**Smith+Nephew**

CATALYSTEM<sup>◇</sup>  
Primary Hip System

Design Rationale

# CATALYSTEM<sup>◇</sup> Primary Hip System

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## Reference

1. Smith+Nephew 2024. CATALYSTEM Comparison to Legacy Stem: POLARSTEM. Internal Report. 10143423 REVA. 2. Smith+Nephew 2024. CATALYSTEM comparison to legacy stem: SYNERGY. Internal Report. 10143364 REV A. 3. Smith+Nephew 2024. Surgeon Feedback on the CATALYSTEM Total Hip System. Internal Report. CSD.REC.24.001. 4. Smith+Nephew 2024. Three-dimensional assessment stem fit and restoration of natural biomechanics in a globally diverse population of femoral CT scans. Internal Report. TM-24-034. 5. Smith+Nephew 2024. CATALYSTEM summary of global evaluations. Internal Report. 10143591. 6. Smith+Nephew 2024. Cadaveric study of the repeatability of CATALYSTEM broach and femoral stem seating level. Internal Report. 10144794. 7. Smith+Nephew 2023. Cadaver validation of the CATALYSTEM Total Hip System. Internal Report. OR-23-106. 8. Smith+Nephew 2024. CATALYSTEM claims tray efficiency memo. Internal Report. 10143458 REV A. 9. 2005 ASM International Engineering Materials Achievement Award. 10. Hunter G, Dickinson J, Herb B, et al. Creation of oxidized zirconium orthopaedic implants. Journal of ASTM International. 2005;2:1-14. 11. Parikh A, Hill P, Pawar V, Sprague J. Long-term simulator wear performance of an advanced bearing technology for THA. Poster presented at: 2013 Annual Meeting of the Orthopaedic Research Society. Poster no. 1028. 12. Papannagari R, Hines G, Sprague J, Morrison M. Long-term wear performance of an advanced bearing technology for TKA. Poster presented at: 2011 Annual Meeting of the Orthopaedic Research Society. Poster no. 1141. 13. Dalal A, Pawar V, McAllister K, Weaver C, Hallab NJ. Orthopedic implant cobalt-alloy particles produce greater toxicity and inflammatory cytokines than titanium alloy and zirconium alloy-based particles in vitro, in human osteoblasts, fibroblasts, and macrophages. J Biomed Mater Res Part A. 2012;100A:2147-2158. 14. Pawar V, Jones B, Sprague J, Salehi A, Hunter G. Acidic fretting tests of oxidized Zr-2.5Nb, CoCr and SS femoral heads. Paper presented at: Materials & Processes for Medical Devices Conference, ASMI, St. Paul, Minnesota, 2004. 15. Davis ET, Pagkalos J, Kopjar B. Bearing surface and survival of cementless and hybrid total hip arthroplasty in the National Joint Registry of England, Wales, Northern Ireland and the Isle of Man. Journal of Bone Joint Surgery. 2020;5(2):pe0075. 16. Innocenti M, Matassi F, Carulli C, Nistri L, Civinini C. Oxidized zirconium femoral component for TKA: A follow-up note of a previous report at a minimum of 10 years. The Knee. 2014;21:858-861. 17. Mayman DJ, Patel AR, Carroll KM. Hospital related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patients. Poster presented at: ISPOR Symposium; May 19-23, 2018; Baltimore, Maryland, USA. 18. Patrick C, Delhougne G, Patel A. Retrospective analysis of oxidized zirconium bearing surface in hip replacement 90-day episode claims. ISPOR EU Poster, 2019. 19. Smith+Nephew 2024. Finite element analysis of the CATALYSTEM hip stem design. Internal Report. OR-24-025. 20. Smith+Nephew 2024. CATALYSTEM distal taper and cross section comparisons of competitive stems. Internal Report. 10142796 Rev B. 21. Smith+Nephew 2024. Morphological and mechanical properties of a dual layered porous plasma spray coating. Internal Report. 10145064. 22. Smith+Nephew 2024. CATALYSTEM ACCUBROACH femur to bone interaction design rationale. Internal Report. 10142827. 23. Smith+Nephew 2024. CATALYSTEM comparison to competitive neck geometries. Internal Report. 10143689. 24. Nerys-Figueroa J, Parsa A, Curley A et al. Slightly reduced early subsidence with similar outcomes and complications rate in collared stems – a systematic review of randomized clinical trials. Journal of Orthopaedics, 2024; 50; 170-176. 25. Smith+Nephew 2024. CATALYSTEM foam femur torsion stability. Internal Report. OR-24-082. 26. Qiu, J. et al. Risk factors for iliopectoral impingement after total hip arthroplasty using a collared femoral prosthesis. Journal of Orthopaedic Surgery and Research (2020) 15:267. <https://doi.org/10.1186/s13018-020-01787-3>. 27. National Joint Registry for England, Wales and Northern Ireland. 20th Annual Report, 2023. Available at: <https://reports.njrcentre.org.uk>. Accessed May 24, 2024. 28. Australian Orthopaedic Association National Joint Replacement Registry Hip, Knee & Shoulder Arthroplasty: 2023 Annual Report. Available at: <https://aoanjrr.sahmri.com/annual-reports-2023>. Accessed May 24, 2024. 29. Sheth NP, Lementowski P, Hunter G, Garino JP. Clinical applications of oxidized zirconium. J Surg Orthop Adv. 2008;17(1):17-26. 30. Davidson JA, et al. Friction and UHMWPE wear of cobalt alloy, zirconia, titanium nitride, and amorphous diamond-like carbon implant bearing surfaces. Poster presented at: 4th World Biomaterials Con1992; Berlin, FRG. 31. Hobbs L, et al. Int J Appl Ceram Technol. 2005. 32. Long M, et al, Hunter G. Nano-Hardness Measurements of Oxidized Zr.2.5Nb and Various Orthopaedic Materials. 1998. 33. Cartner J, et al. Hss j. 2017;13(1):35-41. 34. Cartner J, et al. Journal of Orthopaedic Research. 2016;34. 35. Peters RM, et al. Acta Orthopaedica. 2018;89(2):163-169. 36. Atrey A, et al. Poster presented at: Canadian Orthopaedic Association; June 20-23, 2018; Victoria, British Columbia, Canada. 37. Clavé A, et al. In Orthopaedics & Traumatology, Surgery & Research. 2015;101:647-653. 38. Renkawitz T, et al. J Arthroplasty. 2014;29(5):1021-1025. 39. Ulivi M, et al. J Arthroplasty. 2014;29(5):1026-1029. 40. Naudie D, et al. J Arthroplasty. 2013;28(8 Suppl):48-52. 41. Bourne R, et al. Orthopedics. 2008;31(12 Suppl 2). 42. Boby J, et al. Clin Orthop Relat Res. 1980(150):263-270. 43. Kienapfel H, et al. J Arthroplasty. 1999;14(3):355-368.

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©2025 Smith+Nephew
®Trademark of Smith+Nephew