+ Evidence in Focus

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Referencias

1. Gwam CU, Mistry JB, Mohamed NS, *et al.* Current Epidemiology of Revision Total Hip Arthroplasty in the United States: National Inpatient Sample 2009 to 2013. *J Arthroplasty*. 2017;32:2088–2092.
2. Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR) Hip, Knee & Shoulder Arthroplasty: 2019 Annual Report. Available at: <https://aoanjrr.sahmri.com/annual-reports-2019>. Accessed April 30, 2020.
3. Moore C, Orlandini L. The dual mobility concept in total hip arthroplasty. *Bone & Joint Science*. 2012;4:1–4.
4. Blakeney WG, EpinetteJ-A, Vendittoli P-A. Dual mobility total hip arthroplasty: should everyone get one? *EFORT Open Reviews*. 2019;4:541–547.
5. Jones CW, De Martino I, D'Apolito R, Nocon AA, Sculco PK, Sculco TP. The use of dual-mobility bearings in patients at high risk of dislocation. *Bone Joint J*. 2019;101-B(1 Suppl A):41–5.
6. Abdel MP. Simplifying the hip-spine relationship for total hip arthroplasty: when do I use dual-mobility and why does it work? *J Arthroplasty*. 2019;34:S74–S75.
7. Epinette J-A, Lafuma A, Robert J, Doz M. Cost-effectiveness model comparing dual-mobility to fixed-bearing designs for total hip replacement in France. *Orthopaedics & Traumatology: Surgery & Research*. 2016;102:143–148.
8. De Martino I, Triantafyllopoulos GK, Sculco PK, Sculco TP. Dual mobility cups in total hip arthroplasty. *World J Orthop*. 2014;18:180–187.
9. Neri T, Philippot R, Klasan A, *et al.* Dual mobility acetabular cups for total hip arthroplasty: advantages and drawbacks. *Expert Rev Med Devices*. 2018;15:835–845.
10. Orthopaedic Data Evaluation Panel (ODEP). Available at: <http://www.odep.org.uk>. Accessed March 27, 2020.
11. Cypres A, Fiquet A, Girardin P, *et al.* Long-term outcomes of a dual-mobility cup and cementless triple-taper femoral stem combination in total hip replacement: a multicenter retrospective analysis. *J Orthop Surg Res*. 2019;14(1):376.
12. Bauchu P, Bonnard O, Cypres A, Fiquet A, Girardin P, Noyer D. The dual-mobility POLARCUP: first results from a multicenter study. *Orthopedics*. 2008;31(12 Suppl 2):97–99.
13. Fiquet A, Noyer D. "Polarsystem" dual mobility hip prosthesis and "minimally invasive surgery" (MIS). *Interactive Surgery*. 2006;22:1–5.
14. Cartner J, Aldinger P, Li C, Collins D. Characterization of femoral head taper corrosion features using a 22-year retrieval database. *HSSJ*. 2017;13:35–41.
15. Hunter G, Dickinson J, Herb B, Graham R. Creation of oxidized zirconium orthopaedic implants. *J ATSM Int*. 2005;2:7.
16. ASTM F75-01. Standard specification for cobalt-28 chromium-6 molybdenum alloy castings and casting alloy for surgical implants (UNS R30075), ASTM International, West Conshohocken, PA, 2012.
17. Morrison M, Jani S. Comparison of sequential and single-dose irradiation effects on the mechanical and physical properties of UHMWPE. Poster 1786 presented at Orthopaedic Research Society Congress; February 14–18, 2007; CA, USA.

Referencias

18. Morrison ML, Jani S. Evaluation of sequentially crosslinked ultra-high molecular weight polyethylene. *J Biomed Mater Res B Appl Biomater*. 2009;90:87–100.
19. Davis ET, Pagkalos J, Kopjar B. Polyethylene manufacturing characteristics have a major effect on the risk of revision surgery in cementless and hybrid total hip arthroplasties: an analysis of the National Joint Registry Of England, Wales, Northern Ireland and The Isle of Man. *Bone Joint J*. 2020;102-B(1):90–101.
20. Jassim S, Patel S, Wardle N, et al. Five-year comparison of wear using oxidised zirconium and cobalt–chrome femoral heads in total hip arthroplasty. *Bone Joint J*. 2015;97:883–889.
21. Parikh A, Hill P, Pawar V. The effect of femoral head diameter and material on the wear of highly crosslinked polyethylene liners. Abstract #042 Presented at Society for Biomaterials Annual Meeting; April 13–16, 2011; FL, USA.
22. Pawar V, Weaver C, Parikh A, Jani S. Evaluation of Diffusion Hardened Oxidized Zr-2.5wt%Nb for Hard-on-Hard Articulation in Total Hip Arthroplasty. Medical Device Materials VI; Proceedings from the Materials & Processes for Medical Devices Conference 2011 (ASM International). 2013;149–152.
23. Parikh A, Weaver C, Sprague J, Pawar V. Abrasion resistance, wear performance, and taper corrosion of a new composition of ceramicised metal. Poster 1910 presented at Orthopaedic Research Society Congress; March 5–8, 2016; FL, USA.
24. Parikh A, Weaver C. Wear testing of a new modular dual mobility system under subluxation conditions. Poster 1012 presented at Orthopaedic Research Society Congress; February 2–5, 2019; TX, USA.
25. Parikh A, Weaver C. Wear Testing of a New Modular Dual Mobility System. Poster 1041 presented at ORS Congress; February 2–5, 2019; TX, USA.
26. Lee Y-K, Kim K-C, Jo W-L, Ha Y-C, Parvizi J, Koo K-H. Effect of inner taper angle of acetabular metal shell on the malseating and dissociation force of ceramic liner. *J Arthroplasty*. 2017;32:P1360–1362.
27. Bingenheimer H. Deformation and roundness analysis of diffusion hardened oxidized zirconium dual mobility liners in r3 acetabular shells. Poster 1905 presented at Orthopaedic Research Society Congress; February 2–5, 2019; TX, USA.
28. Bingenheimer H, Li C. Fatigue test of a cross-linked polyethylene dual mobility insert against diffusion hardened oxidized zirconium liner. Poster 1056 presented at Orthopaedic Research Society Congress; February 2–5, 2019; TX, USA.
29. US Department of Health and Human Services. 510(k) Summary Smith & Nephew OR30 Dual Mobility System (K191002). Available at: https://www.accessdata.fda.gov/cdrh_docs/pdf19/K191002.pdf. Accessed April 30, 2020.

Smith & Nephew, Inc.
1450 Brooks Road,
Memphis, TN 38116,
USA

Smith & Nephew S.A.U
Av. Baix Llobregat, 30
08970 Esplugues de Llobregat
Barcelona
Telf: + 34 933 737 301
Fax: + 34 933 737 453
www.smith-nephew.es

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© 2020 Smith & Nephew