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References

1. Latarjet M. A propos du traitement des luxations récidivantes de l'épaule. [Treatment of recurrent dislocations of the shoulder]. *Lyon Chir* 1954;49:994-7. **2.** Walch G. La luxation récidivante antérieure d'épaule. [Recurrent anterior shoulder instability]. *Rev Chir Orthop Reparatrice Appar Mot* 1991;77(Suppl 1):177-91. **3.** Walch G, Boileau P. Latarjet-Bristow procedure for recurrent anterior instability. *Tech Shoulder Elbow Surg* 2000;1:256-61. **4.** Balg F, Boileau P. The instability severity index score. A simple preoperative score to select patients for arthroscopic or open shoulder stabilisation. *J Bone Joint Surg Br* 2007;89:1470-7. **5.** Beran MC, Donaldson CT, Bishop JY. Treatment of chronic glenoid defects in the setting of recurrent anterior shoulder instability: a systematic review. *J Shoulder Elbow Surg* 2010;19:769-80. <http://dx.doi.org/10.1016/j.jse.2010.01.011>. **6.** Bhatia S, Frank RM, Ghodadra NS, Hsu AR, Romeo AA, Bach BRJ, et al. The outcomes and surgical techniques of the Latarjet procedure. *Arthroscopy* 2014;30:227-35. <http://dx.doi.org/10.1016/j.arthro.2013.10.013>. **7.** Burkhart SS, De Beer JF, Barth JRH, Cresswell T, Criswell T, Roberts C, et al. Results of modified Latarjet reconstruction in patients with anteroinferior instability and significant bone loss. *Arthroscopy* 2007;23:1033-41. <http://dx.doi.org/10.1016/j.arthro.2007.08.009>. **8.** Provencher MT, Bhatia S, Ghodadra NS, Grumet RC, Bach BR, Dewing CB, et al. Recurrent shoulder instability: current concepts for evaluation and management of glenoid bone loss. *J Bone Joint Surg Am* 2010;92(Suppl 2):133-51. <http://dx.doi.org/10.2106/JBJS.J.00906>. **9.** Schmid SL, Farshad M, Catanzaro S, Gerber C. The Latarjet procedure for the treatment of recurrence of anterior instability of the shoulder after operative repair: a retrospective case series of forty-nine consecutive patients. *J Bone Joint Surg Am* 2012;94:e75. <http://dx.doi.org/10.2106/JBJS.K.00380>. **10.** Shah AA, Butler RB, Romanowski J, Goel D, Karadagli D, Warner JJP. Short-term complications of the Latarjet procedure. *J Bone Joint Surg Am* 2012;94:495-501. <http://dx.doi.org/10.2106/JBJS.J.01830>. **11.** Giles JW, Boons HW, Elkinson I, Faber KJ, Ferreira LM, Johnson JA, et al. Does the dynamic sling effect of the Latarjet procedure improve shoulder stability? A biomechanical evaluation. *J Shoulder Elbow Surg* 2013;22:821-7. <http://dx.doi.org/10.1016/j.jse.2012.08.002>. **12.** Patte D, Bernageau J, Bancel P. The anteroinferior vulnerable point of the glenoid rim. In: Bateman JE, Welsh, editors. *Surgery of the shoulder*. New York: Marcel Dekker; 1985. p. 94-9. **13.** Boileau P, Thelu CE, Mercier N, Ohl X, Houghton-Clemmey R, Carles M, et al. Arthroscopic Bristow-Latarjet combined with Bankart repair restores shoulder stability in patients with glenoid bone loss. *Clin Orthop Relat Res* 2014;472:2413-24. <http://dx.doi.org/10.1007/s11999-014-3691-x>. **14.** Mizuno N, Denard PJ, Raiss P, Melis B, Walch G. Long-term results of the Latarjet procedure for anterior instability of the shoulder. *J Shoulder Elbow Surg* 2014;23:1691-9. <http://dx.doi.org/10.1016/j.jse.2014.02.015>. **15.** Allain, J, Goutallier, D, Glorion, C. Long-term results of the Latarjet procedure for the treatment of anterior instability of the shoulder. *J Bone Joint Surg Am* 1998;80:841-52. **16.** Bessiere C, Trojani C, Carles M, Mehta SS, Boileau P. The open Latarjet procedure is more reliable in terms of shoulder stability than arthroscopic Bankart repair. *Clin Orthop Relat Res* 2014;472:2345-51. <http://dx.doi.org/10.1007/s11999-014-3550-9>. **17.** Bhatia DN, De Beer JF, du Toit DF. Coracoid process anatomy: implications in radiographic imaging and surgery. *Clin Anat* 2007;20:774-84. <http://dx.doi.org/10.1002/ca.20525>. **18.** Collin P, Rochcongar P, Thomazeau H. Résultat de la butée coracoïdienne type Latarjet pour instabilité antérieure chronique de l'épaule. [Treatment of chronic anterior shoulder instability using a coracoid bone block (Latarjet procedure): 74 cases]. *Rev Chir Orthop Reparatrice Appar Mot* 2007;93:126-32. **19.** Hovelius L, Sandstrom B, Olofsson A, Svensson O, Rahme H. The effect of capsular repair, bone block healing, and position on the results of the Bristow-Latarjet procedure (study III): long-term followup in 319 shoulders. *J Shoulder Elbow Surg* 2012;21:647-60. <http://dx.doi.org/10.1016/j.jse.2011.03.020>. **20.** Boileau, P, et al. A guided surgical approach and novel fixation method for arthroscopic Latarjet. *J Shoulder Elbow Surg* 2015;1-13. <http://dx.doi.org/10.1016/j.jse.2015.02.019>.

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©Trademark of Smith+Nephew