

# + Redefining healing potential

An integral part of  
Advanced Healing Solutions

**Smith+Nephew**



**REGENETEN**   
Bioinductive Implant

# Changing the course of rotator cuff disease

Rotator cuff disease is a significant and costly problem<sup>2-4</sup> that causes ongoing pain and limits patients' mobility.<sup>5</sup> Progressive in nature, small tears tend to grow in size and severity over time, eventually requiring surgery.<sup>1-3</sup>

- Up to 80% of partial-thickness tears increase in size within two years<sup>6</sup>
- Untreated rotator cuff tendinosis can progress to a partial- or full-thickness tear<sup>7</sup>
- Larger tears requiring surgery tend to re-tear over 40% of the time<sup>8-10</sup>

## Now you can disrupt rotator cuff disease progression biologically<sup>1</sup>


The REGENETEN Bioinductive Implant stimulates the body's natural healing response to support new tendon growth and disrupt disease progression.<sup>1,2</sup> Derived from highly purified bovine Achilles tendon, it creates an environment that is conducive to healing.<sup>1,2</sup>

### Biologically improve healing

- Proprietary, highly porous implant design facilitates the formation of new tendon-like tissue<sup>1,2</sup>
- New tissue reduces the peak strain at the site of the tear<sup>11</sup>
- Gradually absorbs within 6 months and leaves a layer of new tendon-like tissue to biologically augment the existing tendon<sup>12</sup>



Arthroscopic view of rotator cuff tear



Implant in situ



12 months post-op



### Demonstrated clinical efficacy

- Induction of new tendon-like tissue in all patients (N=33)
- Mean increase of tendon thickness of 2.2 mm ( $P < 0.0001$ ) at 3 months
- Potentially reduce re-tears<sup>13</sup>



### Excellent safety profile


- No foreign body/ inflammatory reaction
- No implant-related complications



### Impressive patient outcomes

- High patient satisfaction (94%) after 1 year
- Rapid recovery: average 23 days of sling time
- Significantly improved ASES pain score at 1 year ( $P < 0.0001$ )<sup>‡</sup>

## Natural progression of rotator cuff disease



\*Results from a prospective multi-center study of patients with partial-thickness tears. Patients had chronic, degenerative, intermediate grade (n=12) or high grade (n=21) partial-thickness tears of the supraspinatus tendon. The REGENETEN Bioinductive Implant was attached following arthroscopic subacromial decompression without repair. Clinical outcomes were assessed pre-op and at 3 and 12 months post-op using American Shoulder and Elbow Surgeons (ASES) and Constant-Murley (CM) scores. Post-op tendon healing and thickness was assessed with MRI.

‡ASES pain score improved from  $4.2 \pm 0.4$  standard error of mean (SEM) at baseline to  $0.6 \pm 0.2$  (SEM) at 1 year.



# Ordering information

Implants	
Reference #	Description
4565	Medium Bioinductive Implant with Arthroscopic Delivery (1)
4566	Large Bioinductive Implant with Arthroscopic Delivery (1)
Anchors	
Reference #	Description
4403	Bone Anchors (3) with Advanced Delivery System
2504-1	Tendon Anchors (8)
Accessory Devices	
Reference #	Description
4173-1	Tendon Marker (2)
4402	Tendon Stabilizing Guide (1)
2503-S	Bone Anchor (1)

Learn more at [smith-nephew.com](http://smith-nephew.com)

Sports Medicine  
Smith & Nephew, Inc.  
150 Minuteman Road  
Andover, MA 01810

[www.smith-nephew.com](http://www.smith-nephew.com)  
T +978 749 1000  
US Customer Service:  
+1 800 343 5717

\*Trademark of Smith+Nephew.  
©2020 Smith+Nephew. All rights  
reserved. All trademarks acknowledged.  
Printed in USA. 13508 V4 03/20

## References

**1.** Schlegel TF, Abrams JS, Bushnell BD, Brock JL, Ho CP. Radiologic and clinical evaluation of a bioabsorbable collagen implant to treat partial-thickness tears: a prospective multicenter study. J Shoulder Elbow Surg. 2017. doi: <http://dx.doi.org/10.1016/j.jse.2017.08.023>. **2.** Bokor DJ, Sonnabend D, Deady L et al. Evidence of healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a collagen implant: a 2-year MRI follow-up. MLTJ. 2016;6(1):16-25. **3.** Washburn R, Anderson TM, Tokish JM. Arthroscopic rotator cuff augmentation: Surgical technique using bovine collagen bioinductive implant. Arthroscopy Techniques. 2017;6(2):e297-e301. **4.** Mather RC, Koenig L, Acevedo D et al. The societal and economic value of rotator cuff repair. J Bone Joint Surg Am. 2013;95:1993-2000. **5.** Lin JC, Weintraub N, Aragaki DR. Nonsurgical treatment for rotator cuff injury in the elderly. Am Med Dir Assoc. 2008;9(9):626-32. doi: 10.1016/j.jamda.2008.05.003. **6.** Yamanaka K and Matsumoto T. The joint side tear of the rotator cuff: A followup study by arthrography. Clinical Orthopaedics and Related Research. 1994: 304,68-73. **7.** Keener JD, Galatz LM, Teefey SA et al. A prospective evaluation of survivorship of asymptomatic degenerative rotator cuff tears. J Bone Joint Surg Am. 2015;97:89-98. **8.** Bishop J, Klepps S, Lo IK, Bird J, Gladstone JN, Flatow EL. Cuff integrity after arthroscopic versus open rotator cuff repair: A prospective study. J Shoulder Elbow Surg. 2006;15(3):290-299. **9.** Heuberger PR, Smolen D, Pauzenberger L et al. Longitudinal long-term magnetic resonance imaging and clinical follow-up after single-row arthroscopic rotator cuff repair. Am J Sports Med. 2017;45(6):1283-1288. **10.** Henry P, Wasserstein D, Park S, et al. Arthroscopic repair for chronic massive rotator cuff tears: A systematic review. Arthroscopy. 2015;31(12):2472-80. **11.** Chen Q. Proof-of-concept finite element modelling of effect of tissue induction on rotator cuff tears. Material and Structural Testing Core, Mayo Clinic, Rochester, MN, 2011. **12.** Van Kampen C, et al. Tissue-engineered augmentation of a rotator cuff tendon using a reconstituted collagen scaffold: A histological evaluation in sheep. MLTJ. 2013;3:229-235. **13.** Bokor DJ, Sonnabend D, Deady L et al. Preliminary investigation of a biological augmentation of rotator cuff repairs using a collagen implant: a 2-year MRI follow-up. MLTJ. 2015;5(3):144-150.

FOR DISTRIBUTION IN THE US AND ITS TERRITORIES ONLY