

+ Evidence in focus

# Compendium of peer-reviewed clinical evidence

JOURNEY<sup>II</sup>  
Total Knee Arthroplasty

July 2023

Smith+Nephew



# JOURNEY<sup>◇</sup> II TKA Compendium of Evidence

Key evidence

All studies

Key outcomes

Function

Recovery

Patient  
satisfaction

Survivorship

## PDF navigation tips

Interactivity has been included throughout the compendium to aid navigation



will take you to this slide



will take you to a summary of each outcome



will take you to a grid summary of the 51 peer-reviewed studies highlighted in this compendium. Key evidence is represented in **blue** and supporting evidence in **turquoise**.

Full summary

will take you to the Evidence in Focus full summary of the relevant study


Published paper

will take you to the official online journal website on which the original paper is published

Arrows at the bottom of each page will take you to previous or subsequent pages, as indicated

# What are the issues with conventional TKA?

A “Forgotten Joint” is the desired goal for both patients and surgeons following TKA. However, the majority of patients do not obtain a normal feeling knee post-TKA with high numbers of patients experiencing reduced function and dissatisfaction:



**80%** patients feel that their joint is artificial<sup>1</sup>



**>50%** patients report a degree of functional limitation<sup>1,2</sup>



**20%** patients are dissatisfied<sup>3</sup>

Studies show that following conventional TKA, patients demonstrate abnormalities in their gait and functional kinematics of the knee:<sup>4-8</sup>

Velocity



Stride length



Flexion during walking



**Quadricep muscle strength following conventional TKA fails to recover to the same level as those of healthy subjects, resulting in reduced physical functioning of the knee<sup>9,10</sup>**




Conventional TKAs have been successful in their aim to relieve pain and provide long term survivorship.<sup>11,12</sup> However, conventional TKAs do not live up to all patients functional expectations,<sup>2</sup> leaving a high percentage feeling dissatisfied.<sup>3</sup>

# Kinematics and patient satisfaction — how do they relate?

Fluoroscopic analysis of weight-bearing motion has led to a better understanding of TKA kinematics. During flexion, TKA kinematic patterns can vary considerably from the normal knee.<sup>13</sup>

Paradoxical motion (Figure), where there is a sudden anterior translation of the femur relative to the tibia, has been observed in many studies and is related to reduced knee flexion and quadriceps efficiency.<sup>13</sup>



**An implant designed to better replicate the shape and position of the normal knee could improve knee function and subsequently patient satisfaction.**

For the first time, a relationship between kinematic patterns of TKA and patient satisfaction has been demonstrated.<sup>14,15</sup>

Kinematic patterns for one weight-bearing and one non-weight-bearing motion were analysed for JOURNEY II BCS and two posterior-stabilised TKAs. Whilst no differences were observed during non-weight bearing activity, during weight-bearing motion, poor patient-reported outcomes were associated with:<sup>14</sup>

Pronounced paradoxical anterior motion (medial side)

Less stable medial compartment in midflexion

Less posterior translation in deep flexion (lateral side)

**“Reproduction of optimal kinematic patterns during TKA could be instrumental in improving patient satisfaction.”<sup>14</sup>**



# Why choose JOURNEY II TKA for your patients?

## Shape



## Position

## Motion

Unlike conventional TKA design, the shape and position of the JOURNEY II TKA have been designed to replicate the normal knee


### Shape<sup>16-19</sup>

Replication of anatomic asymmetric femoral and tibial profiles



### Position<sup>16,19,20</sup>


Restoration of native anterior/posterior (A/P) starting position and the anatomic 3° varus joint line




### Motion<sup>16-18,20-23</sup>

Replication of native femoral rollback and axial rotation


#### Normal knee



#### JOURNEY II TKA knee



# Why choose JOURNEY II TKA for your patients?



\*Compared to non-JOURNEY II knees. †As observed in JOURNEY II BCS and JOURNEY II CR.



## Key outcome: Improved function

38 studies

Postoperative ROM is one of the most important factors influencing patient satisfaction following TKA, as a full ROM is critical for a patient to perform daily activities.<sup>41</sup>

JOURNEY II TKA has demonstrated better kinematics,<sup>17,24,25</sup> ROM<sup>26-29</sup> and gait<sup>24,30</sup> compared to conventional TKAs.

### Improved knee flexion<sup>24,27-29</sup>

Post-TKA:



3 months

**Significantly improved knee flexion** during walking ( $p < 0.01$ ) with JOURNEY II CR compared to Attune™ CR\* TKA.<sup>24</sup>

**A significantly increased mean ROM** with JOURNEY II BCS compared to conventional PS TKAs:

+6° ( $p = 0.04$ )<sup>27</sup>

+6° ( $p = 0.002$ )<sup>28</sup>

+23° ( $p < 0.001$ )<sup>29</sup>

1 year



### Improved gait<sup>24,30</sup> and stair climbing<sup>42</sup>

Post-TKA:

3 months


**Significantly greater increase in external rotation whilst walking** with JOURNEY II CR compared to Attune™ CR\* TKA ( $p < 0.01$ ).<sup>24</sup>

6 months

**Significant improvements in gait** with JOURNEY II BCS compared to LEGION<sup>®</sup> PS TKA ( $p = 0.03$ ).<sup>30</sup>

1 year

**Significantly more patients able to climb stairs** with JOURNEY II BCS compared to LEGION PS TKA ( $p < 0.05$ ).<sup>42</sup>



### Better kinematic outcomes<sup>16,17</sup>

Compared to the normal knee after 1 year post-TKA:



**A similar axial rotation** in early and late flexion with JOURNEY II BCS<sup>16</sup>







## Key outcome: Smoother recovery

 **13 studies**

Pain relief, function and return to activities are some of the most common patient expectations post-TKA.<sup>2</sup>

JOURNEY II TKA\* has been shown to result in speedier recovery<sup>29,31,32</sup> and improved muscle activation and muscle strength in the early recovery period,<sup>24,30</sup> compared to conventional implants.

### Improved quadricep activation and strength<sup>24,30</sup>

Improved muscle activation and muscle strength in the early recovery period with JOURNEY II BCS and JOURNEY II CR, compared to LEGION<sup>®</sup> PS TKA and Attune<sup>™†</sup> CR, respectively.<sup>24,30</sup>

Improving quadriceps function is important for limiting post-TKA functional deficits<sup>9</sup>



### Speedier recovery<sup>31,32</sup>


Compared with patients receiving conventional TKA, JOURNEY II TKA patients are:



More likely to be discharged home ( $p < 0.001$ )<sup>31</sup>




Discharged from hospital sooner ( $p < 0.0001$ )<sup>32</sup>



### Higher function outcome measures<sup>26,29</sup>

Compared to conventional TKA patients, JOURNEY II TKA patients have demonstrated **higher scores** in functional outcome measures at 6 and 12 weeks post-TKA.<sup>26,29</sup>



Similar levels of patient satisfaction and patient-reported outcomes between JOURNEY II BCS and THA patients at 3 months post-op<sup>33</sup>




\*As observed in JOURNEY II BCS and JOURNEY II CR. †Manufactured by DePuy Synthes Orthopaedics, Warsaw, IN, USA.







## Key outcome: Higher patient satisfaction

 19 studies

The primary determinant of patient satisfaction is the fulfilment of patient expectations, of which pain relief, improved knee function and return to sports are the most common.<sup>2</sup>

JOURNEY II TKA has demonstrated improvements in patient-reported outcome measures post-TKA,<sup>26,29,33-37</sup> compared to pre-TKA.

### Improved WOMAC scores<sup>26</sup>

Significant reductions in WOMAC scores with JOURNEY II CR compared to P.F.C. Sigma<sup>™</sup> CR\* at 6 (p=0.018) and 12 months (p=0.008) post-TKA.<sup>26</sup>



WOMAC scores<sup>43</sup>, KSS<sup>75</sup>, and KOOS<sup>76</sup> directly correlate to patient satisfaction

### Better KSS<sup>25,26,29,36</sup> and KOOS<sup>35</sup>

Significantly better KSS outcomes with:

JOURNEY II CR compared to P.F.C. Sigma<sup>™</sup> CR\* at 3, 6, and 24 months post-TKA (p<0.05) and at 12 months post-TKA (p<0.001)<sup>26</sup>

JOURNEY II BCS compared to Persona<sup>™†</sup> PS TKA (p<0.05)<sup>25</sup> and to LEGION<sup>®</sup> PS TKA (p<0.001)<sup>29</sup> at 12 months post-TKA

JOURNEY II BCS at 24 months post-TKA compared to pre-TKA (p<0.01)<sup>36</sup>



### Quality of life<sup>33</sup>

Patients receiving JOURNEY II BCS have reported **similar levels of satisfaction as those receiving THA** at 3 months and 1 year post-op, with no significant differences in patient quality of life.<sup>33</sup>

### Return to work and sports<sup>33</sup>

One study reported a  
**2 month**

median time to return to work, activities of daily living and sporting activities<sup>33</sup>



Improving patient satisfaction is important for maintaining quality of care and patient loyalty which is linked to reimbursement to payors and healthcare providers<sup>44,45</sup>



\*Manufactured by DePuy Synthes Orthopaedics, Warsaw, IN, USA. †Manufactured by Zimmer Biomet, Warsaw, Indiana, USA.





## Key outcome: High survivorship

5 studies

Revision TKA is a technically challenging procedure and is associated with a high risk of complications.<sup>46</sup>

Studies on survivorship of JOURNEY II BCS and JOURNEY II CR have shown favourable survivorship outcomes, with high early<sup>35,38</sup> and mid-term<sup>39,40</sup> survivorship.

### High early survivorship

Early results from prospective studies on JOURNEY II CR (n=174)<sup>38</sup> have demonstrated a low risk of revision:

**<1.2%**

Risk of revision at 1 year<sup>35,38</sup>



Revision TKA poses significant economic costs to the healthcare industry.<sup>46</sup>


**\$75,000**  
Average estimated US cost of revision TKA<sup>47</sup>




### High mid-term survivorship

JOURNEY II BCS has demonstrated high mid-term survivorship\*, comparable to AOANJRR cemented PS TKA control.<sup>39,40</sup>

Retrospective, international study of **2,059 JOURNEY II BCS TKAs**<sup>39</sup>



Prospective, multicentre, US study of **209 JOURNEY II BCS TKAs**<sup>40</sup>



JOURNEY II BCS resulted in:

**<1% major revisions, lower than registry control<sup>39,40</sup> and a lower revision rate in <55 year olds compared to registry control (3 vs 7%)<sup>39</sup>**

\*Kaplan-Meier (K-M) approach was used to obtain estimates of implant survival.





# Does a uni “feel better” than a total knee? Not necessarily, when using modern implant designs<sup>77</sup>

Ast MP, Kolin DA, Carroll KM, Davis D, Pearle AD, Mayman DJ, Ong AC. *HSS Journal: The Musculoskeletal Journal of Hospital for Special Surgery*. 2022;19(1),32–36

## Retrospective, multi-centre review of:

200 JOURNEY II BCS TKA    200 Restoris™ UKA\*    200 Vanguard™ TKA†

Length of study: **3 years**

## Assessed:

Forgotten Joint Score (FJS)

Mean 1.6 year follow-up



## Results:

JOURNEY II TKA demonstrated:



Significantly higher FJS compared with Vanguard™ posterior stabilised (PS) TKA†  
(p=0.007)



Similar FJS to Restoris™ UKA\*



Higher FJS in patients with increased BMI compared with their UKA and other TKA counterparts

## Conclusion

JOURNEY II BCS as a modern design TKA had no significant difference in FJS when compared to UKAs. However had a significantly higher FJS scores over “traditional” knee designs.

\*Manufactured by Stryker Orthopedics, Mahwah, NJ, USA.

†Manufactured by Zimmer Biomet, Warsaw, IN, USA.



# Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study<sup>24</sup>

Di Benedetto P, Vidi D, Colombo A, Buttironi MM, Cainero V, Causero A. *Acta Biomed.* 2019;90:91–97

## Independent pilot study of:

12 JOURNEY II CR  
Mean age: 70.3 years

14 Attune™ CR\*  
Mean age: 71.8 years

Follow up: **3 months**

## Assessed:

Kinematics were assessed at flexion-extension internal-external rotation and abduction-adduction

Muscle strength and activation using electromyography

KSS

KOOS




## Results

Compared to Attune™ CR\* at 3 months post-TKA, JOURNEY II CR patients demonstrated:



Significantly improved knee flexion during walking compared to pre-TKA ( $p < 0.01$ )



Significantly improved external rotation during walking compared to pre-TKA ( $p < 0.01$ )



Better muscle activation timing, comparable to the physiological knee



**37.5**



**35.4**

A numerically higher KSS satisfaction score (maximum score: 40)



**84.9**



**76.2**

A numerically higher KOOS quality of life score (maximum score: 100)

JOURNEY II CR

Attune CR

## Conclusion

JOURNEY II CR patients experienced significantly better flexion, external rotation and muscle strength during walking at 3 months post-TKA compared to Attune™ CR\* TKA. In addition, JOURNEY II CR patients demonstrated physiological activation timing of most muscles.



# In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy<sup>16</sup>

Grieco TF, Sharma A, Dessinger GM, Cates HE, Komistek RD. *J Arthroplasty*. 2018;33:565–571

## Retrospective, single-surgeon study of:

40 JOURNEY II BCS  
Mean age: 69.8 years

10 normal knees  
Mean age: 57.4 years

Mean follow up: **14.3 months**

## Assessed:

Kinematics assessed during a weight-bearing deep knee bend at 30° increments from full extension to 120° of knee flexion to determine whether the dual cam-post mechanism is able to replicate the cruciate ligament (ACL and PCL) function



## Results



**0–30°**

JOURNEY II BCS subjects exhibited similar patterns of femoral rollback and axial rotation compared with normal knee subjects



**30–60°**

JOURNEY II BCS subjects experienced minimal anterior-posterior motions and axial rotation, whereas normal knees continued to roll back and externally rotate



**60–90°**

JOURNEY II BCS resumed posterior motion with axial rotation increasing in a normal-like fashion after 90°

## Conclusion

JOURNEY II BCS exhibited normal-like kinematic patterns and moved as designed under in vivo observation. Similarities in early and late kinematic patterns between the two groups suggest the dual cam-post design and asymmetric articular geometries of the JOURNEY II BCS adequately replicate ACL and PCL function of the normal knee.



# Midterm clinical results of bi-cruciate stabilized total knee arthroplasty compared with posterior stabilized total knee arthroplasty with single radius design<sup>78</sup>

Inui H, Yamagami R, Kono K, Kawaguchi K, Taketomi S, Saita K, Tanaka S. *Journal of Joint Surgery and Research*. 2023;1(1),117–122.

## Retrospective, single-centre review of:

61 JOURNEY II BCS TKA

81 Scorpio™ NRG & Triathlon™\*

## Assessed at 6 months, 2 years and 5 years post TKA followup:

Maximum Flexion Angle (MFA)



Knee injury and Osteoarthritis Outcome Score (KOOS)



Activities of Daily Living (ADL)



Quality of Life (QOL)



## Results:

JOURNEY II TKA demonstrated:



6 months: MFA of the JOURNEYII cohort was larger than the single radii PS cohort



2 years: MFA, KOOS pain and KOOS ADL scores were significantly better for the JOURNEY II cohort compared to the single radii PS cohort



5 years: MFA, KOOS pain and KOOS QOL scores were significantly better for the JOURNEY II cohort compared to the single radii PS cohort

## Conclusion

JOURNEY II BCS showed significantly better midterm clinical results compared to single radii PS designed knees for MFA and KOOS pain scores.

\*Manufactured by Stryker Orthopedics, Mahwah, NJ, USA.



# The bicruciate substituting knee design and initial experience<sup>29</sup>

Nodzo SR, Carroll KM, Mayman DJ. *Tech Orthop.* 2018;33:37–41

## Retrospective, single-surgeon, study of:

100 JOURNEY II BCS      100 LEGION<sup>®</sup> Total Knee System  
Mean age: 51.0 years

Mean follow up: **1 year**

## Assessed at pre-TKA and 6 weeks and 1 year post-TKA:


ROM

KSS



## Results

JOURNEY II BCS patients, compared to PS TKA patients, demonstrated:



Significantly greater mean ROM at  
1 year post-TKA (119° vs 96°;  $p < 0.0001$ )



Numerically higher mean KSS score  
at 6 weeks post-TKA (88 vs 73)



Significantly higher mean KSS scores  
at 1 year post-TKA (89 vs 81;  $p < 0.001$ )

## Conclusion

JOURNEY II BCS led to significant improvements in ROM and patient-reported outcomes at 1 year post-TKA, compared with standard PS TKA. The authors noted that the results suggest that the more anatomic design of the implant, which is intended to replicate a more normal knee position and kinematic patterns, may be responsible for the improved flexion and patient satisfaction in JOURNEY II BCS patients.





# A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY II bi-cruciate stabilizing knee system and total hip arthroplasty patients<sup>33</sup>

Snyder MA, Sympton A, Gregg J, Levit A. *Orthop Trauma Prosth.* 2018;3:5–10

## Retrospective review of clinically matched data from a regional US total joint registry (Cincinnati, Ohio):

48 JOURNEY II BCS

Mean age: 58.3 years

48 THA

Mean age: 55.9 years

Mean follow up: **1 year**

## Assessed at 3 months and 1 year post-op:

Patient satisfaction


UCLA activity scores

EQ-5D-5L




## Results


JOURNEY II BCS patients, compared to THA patients, reported:



No significant difference in overall satisfaction at 3 months ( $p=0.398$ ) or 1 year ( $p=0.590$ )



Significantly greater median UCLA activity scores at 3 months (8 vs 7;  $p=0.028$ ) and 1 year (8 vs 7;  $p<0.001$ ; maximum score: 10)



Significantly greater median EQ-5D-5L scores at 3 months (90 vs 80;  $p<0.001$ ; maximum score: 100); no significant difference at 1 year

## Conclusion

JOURNEY II BCS TKA patients reported similar levels of satisfaction and patient reported outcomes as THA patients at 3 months and 1 year post-op. The kinematic design advancements of JOURNEY II BCS demonstrate improvements in patient satisfaction following TKA compared to past TKA procedures, via comparison with THAs.



# Midterm performance of a guided-motion bicruciate-stabilized total knee system: results from the international study of over 2000 consecutive primary total knee arthroplasties<sup>39</sup>

Harris AI, Christen B, Malcorps JJ, O'Grady CP, Kopjar B, Sensiba PR, Vandenneucker H, Huang BK, Cates HE, Hur J, Marra DA. *J Arthroplasty*. 2019;34:S201–S208

## Retrospective, multicentre, international case series study of:

2,059 JOURNEY II BCS

Mean age, 64.3 years

Mean follow up: **3.87 years**

## Assessed at last follow-up:

Survivorship compared to all other PS cemented TKA in the AOANJRR




## Results


At 5 years, JOURNEY II BCS, compared to the AOANJRR control, resulted in:



Similar high survivorship rate (96.4 vs 95.9%)



Lower percentage of major revisions, accounting for <1% of JOURNEY II BCS TKAs and defined as all revisions involving tibial and/or femoral component removal (29.9 vs 41.6%)



Lower revision rate in patients <55 year olds (3 vs 7%)

## Conclusion

The JOURNEY II BCS knee system performs favourably when compared to PS TKA control from the AOANJRR, particularly in patients <55 years of age. At up to 6.1 years post TKA, less than 1% of all JOURNEY II BCS TKAs required major revision.



# Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty<sup>30</sup>

Hyodo K, Kanamori A, Kadone H, Takahashi T, Kajiwara M, Yamazaki M. *Arthroplast Today*. 2020;6:338–342

## A retrospective analysis of:

12 JOURNEY II BCS TKA

Mean age: 69.4 years

12 LEGION<sup>®</sup> TKA

Mean age: 70.0 years

Mean length of study: **6 months**

## Assessed at 6 months post-TKA:

Gait and motion capture analysis with a force plate and electromyogram of the lower limb muscles



## Results


JOURNEY II BCS TKA, compared to conventional TKA, demonstrated:



Significantly longer step length ( $p=0.03$ )




Significantly faster gait speed ( $p=0.03$ )



Significantly greater maximum knee extension moment ( $p=0.04$ )


■ JOURNEY II TKA ■ Conventional TKA



Significantly greater maximum knee flexion angle during initial stance phase ( $p=0.04$ )

## Conclusion

JOURNEY II BCS patients demonstrate faster gait speed, longer step length and greater knee extension moment, compared to conventional PS TKA. This suggests that the quadriceps muscles are more effectively activated and that anterior stability function of the anterior cruciate ligament can be achieved with JOURNEY II BCS TKA compared to conventional PS TKA.



# A comparison of rollback ratio between bicruciate substituting total knee arthroplasty and Oxford unicompartmental knee arthroplasty<sup>21</sup>

Iriuchishima T and Ryu K. *J Knee Surg.* 2018;31:568–572

## Retrospective, single-surgeon study of:

64 JOURNEY II BCS Mean age: 71.3 years	50 Oxford™ Partial Knee* Mean age: 73.8 years	50 normal knees N/A
---	--	------------------------



Follow up: **6–9 months**

## Assessed:


When subjects had recovered their range of knee flexion (6–9 months post-op), lateral radiographs in active flexion were taken to measure rollback ratio and flexion angle




## Results



JOURNEY II BCS



Oxford™ UKA\*



Normal knee

No significant difference in rollback ratio and no significant difference in knee flexion angle among the three groups

A significant correlation between rollback ratio and knee flexion angle among the three groups ( $p=0.002$ )

## Conclusion

JOURNEY II BCS showed no significant difference in rollback ratio when compared with UKA or asymptomatic control knees. The results suggest that JOURNEY II BCS reproduces native anterior and posterior cruciate function and native knee rollback.

\*Manufactured by Zimmer Biomet, Warsaw, Indiana, USA.



# Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics<sup>26</sup>

**Lutes W and Fitch D. Presented at: 39<sup>th</sup> SICOT Orthopaedic World Congress; 2018; Montreal, Canada**

Retrospective, single-surgeon study of:		Assessed at 3, 6, 12 and 24 months post-TKA:	
52 JOURNEY II CR	60 P.F.C. Sigma™ CR*	KSS	✓
Mean age: 67.3 years	Mean age: 70.2 years	WOMAC scores	✓
Mean follow up: <b>2 years</b>		ROM (only assessed up to 12 months)	✓

## Results:

JOURNEY II CR patients, compared to P.F.C. Sigma™ CR\* patients, reported:



Significantly greater KSS scores at 3 (69.5 vs 63.0), 6 (84.4 vs 70.1), 12 (93.0 vs 86.1) and 24 (96.4 vs 91.7) months post-TKA ( $p < 0.05$ )



Significantly lower WOMAC scores at 6 (17.8 vs 24.6) and 12 (12.4 vs 18.5) months post-TKA ( $p < 0.05$ )



Significantly greater change in ROM from baseline at 3 (-4.4 vs -10.1), 6 (5.8 vs -1.8) and 12 (11.4 vs 4.0) months post-TKA ( $p < 0.05$ )

## Conclusion

JOURNEY II CR patients reported significant improvements in functional outcomes compared to P.F.C. Sigma™ CR\* TKA.

\*Manufactured by DePuy Synthes Orthopaedics, Warsaw, IN, USA.



# Hospital-related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patient<sup>31</sup>

Mayman DJ, Patel AR, Carroll KM. Poster presented at: ISPOR Symposium; 2018; Baltimore, Maryland, USA

## A retrospective cohort study with real world evidence of:

1,692 JOURNEY II BCS TKA  
Mean age: 64.1 years

1,692 other TKA devices  
Mean age: 63.9 years

## Assessed:


Hospital-related clinical and economic outcomes  
(reported as 2016 US dollars)

1:1 propensity score matching to control for patients and  
provider characteristics






## Results:

JOURNEY II BCS, compared to other TKA knees, was associated with:



Patients receiving JOURNEY II BCS were also:



## Conclusion

JOURNEY II BCS led to a significantly lower total hospital cost and significantly shorter hospital stay, and patients were less likely to be readmitted within 30 days. Patients receiving this device were also more likely to be discharged to home when compared with patients undergoing primary TKA with other posterior-stabilised TKA systems.



# Evaluation of anteroposterior accelerometric change after bi-cruciate stabilized total knee arthroplasty and posterior stabilized total knee arthroplasty<sup>25</sup>

Tomite T, Saito H, Kijima H, Ishikawa N, Hatakeyama Y, Tazawa H, Miyakoshi N, Shimada Y. *Knee*. 2021;32:121–130

## Independent, single-surgeon, prospective study of:

30 JOURNEY II BCS

Mean age: 77.9 years

30 PERSONA<sup>™</sup> PS\*

Mean age: 77.3 years

Mean follow up: **1 year**

## Assessed at pre-TKA and 12 months post-TKA:

ROM

KSS

FTA

Lateral X-rays of the standing extended knee

Accelerometer data



## Results


JOURNEY II BCS patients, compared to PS TKA patients, demonstrated:



Significantly lower posterior offset ratio (2.1 vs 17.9%, respectively), with the AP positioning of the femur and tibia close to that of the normal knee



Lower anteroposterior (AP) acceleration on the femoral side of the knee



## Conclusion






JOURNEY II BCS resulted in better functional kinematics, closer positioning to that of the normal knee on lateral X-ray, and lower anteroposterior acceleration on the femoral side compared to PS TKA at 12 months post-TKA.

\*Manufactured by Zimmer Biomet, Warsaw, IN, USA.








# All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Ast MP, et al. 2022 <sup>77</sup>	Does a uni “feel better” than a total knee? Not necessarily, when using modern implant designs	HSS Journal		BCS			✓		
Biały M, et al. 2021 <sup>51</sup>	A comparison of the JOURNEY II <sup>®</sup> bi-cruciate stabilized total knee system and Genesis II cruciate-retaining implant	Physiotherapy Review		BCS	✓	✓	✓		
Christen B, et al. 2018 <sup>52</sup>	Second-generation bi-cruciate stabilized total knee system has a lower reoperation and revision rate than its predecessor	Arch Orthop Trauma Surg		BCS		✓			
Coomer S, et al. 2021 <sup>53</sup>	Determining patella function in non-implanted knees having functional cruciate ligaments and subjects having a bi-cruciate stabilized total knee arthroplasty	Orthopaedic Proceedings		BCS	✓				
Di Benedetto P, et al. 2019 <sup>24</sup>	Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study	Acta Biomed		CR	✓	✓	✓		






# All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Grieco TF, et al. 2018 <sup>16</sup>	In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy	J Arthroplasty		BCS	✓				
Harris AI, et al. 2018 <sup>37</sup>	Short-term safety and effectiveness of a second-generation motion-guided total knee system	Arthro Today		BCS	✓	✓	✓		
Harris AI, et al. 2019 <sup>39</sup>	Midterm performance of a guided-motion bicruciate-stabilized total knee system: results from the international study of over 2000 consecutive primary total knee arthroplasties	J Arthroplasty		BCS				✓	
Harris AI, et al. 2019 <sup>40</sup>	Guided motion total knee arthroplasty system: five-year outcomes of the prospective multicenter US study	EFFORT Congress		BCS				✓	
Harris AI, et al. 2019 <sup>54</sup>	Guided motion total knee arthroplasty (TKA) system in younger patients has a lower revision rate than registry controls: results from the international multicenter study with up to 6 Years follow-up	EKS Congress		BCS				✓	






# All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Harris AI, et al. 2019 <sup>55</sup>	Guided motion total knee arthroplasty (TKA) in patients with BMI of 40kg/m <sup>2</sup> or more: results from the international multicentre study of 2,059 primary TKAs with up to 6 years follow-up	EKS Congress		BCS				✓	
Heir S, et al. 2019 <sup>36</sup>	Clinical and functional outcomes of a second-generation guided motion total knee arthroplasty system: two-year results of a prospective multicentre study	EKS Congress		BCS	✓	✓	✓		
Hino K, et al. 2018 <sup>56</sup>	Bi-cruciate substituting total knee arthroplasty provides varus–valgus stability throughout the midflexion range	Knee		BCS	✓				
Hommel H, et al. 2017 <sup>57</sup>	Good early results obtained with a guided-motion implant for total knee arthroplasty: A consecutive case series	Open Orthop J		BCS	✓				
Hyodo K, et al. 2020 <sup>30</sup>	Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty	Arthroplasty Today		BCS	✓	✓			



# All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Inui H, et al. 2018 <sup>58</sup>	The relationship between anteroposterior stability and medial-lateral stability of the bi-cruciate stabilized total knee arthroplasty	Knee		BCS	✓				
Inui H, et al. 2019 <sup>59</sup>	The relationship between soft-tissue balance and intraoperative kinematics of guided motion total knee arthroplasty	J Knee Surg		BCS	✓				
Inui H, et al. 2023 <sup>78</sup>	Midterm clinical results of bi-cruciate stabilized total knee arthroplasty compared with posterior stabilized total knee arthroplasty with single radius design	Journal of Joint Surgery and Research		BCS	✓		✓		
Inoue A, et al. 2023 <sup>79</sup>	The clinical results of bi-cruciate vs posterior stabilized total knee arthroplasty for flexion contracture in osteoarthritic knee	Journal of Orthopaedic Surgery		BCS			✓		
Iriuchishima T, et al. 2018 <sup>21</sup>	A comparison of rollback ratio between bicruciate substituting total knee arthroplasty and Oxford unicompartmental knee arthroplasty	J Knee Surg		BCS	✓				






# All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Iriuchishima T, et al. 2019 <sup>60</sup>	Bicruciate substituting total knee arthroplasty improves stair climbing ability when compared with cruciate-retain or posterior stabilizing total knee arthroplasty	Indian J Orthop		BCS	✓				
Ishibashi T, et al. 2020 <sup>42</sup>	Kinematics of bicruciate and posterior stabilized total knee arthroplasty during deep knee flexion and stair climbing	J Orthop Res		BCS	✓				
Ishida K, et al. 2017 <sup>61</sup>	Comparison of intra-operative navigation-based kinematics between bi-cruciate-stabilised total knee arthroplasty (TKA) and conventional posterior-stabilised TKA	Orthop Proceedings		BCS	✓				
Itou J, et al. 2021 <sup>62</sup>	Anterior prominence of the femoral condyle varies among prosthesis designs and surgical techniques in total knee arthroplasty	BMC Musculoskeletal Disorders		BCS	✓				
Kage T, et al. 2021 <sup>15</sup>	The association between in vivo knee kinematics and patient-reported outcomes during squatting in bicruciate-stabilized total knee arthroplasty	J Knee Surg		BCS	✓		✓		






# All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Kaneko T, et al. 2017 <sup>19</sup>	Bi-cruciate substituting total knee arthroplasty improved medio-lateral instability in midflexion range	J Orthop		BCS	✓		✓		
Kaneko T, et al. 2018 <sup>63</sup>	The influence of compressive forces across the patello-femoral joint on patient reported outcome after bi-cruciate stabilized total knee arthroplasty	Bone Joint J		BCS			✓		
Kaneko T, et al. 2020 <sup>64</sup>	The influence of tibiofemoral joint forces on patient-reported outcome measurements after bicruciate stabilized total knee arthroplasty	J Orthop Surg		BCS	✓		✓		
Kiyohara M, et al. 2021 <sup>65</sup>	Comparison of in vivo knee kinematics before and after bicruciate-stabilized total knee arthroplasty during squatting	BMC Musculoskeletal Disorders		BCS	✓	✓			
Kono K, et al. 2019 <sup>66</sup>	Bicruciate-stabilised total knee arthroplasty provides good functional stability during high-flexion weight-bearing activities	Knee Surg Sports Traumatol Arthosc		BCS	✓				

# All evidence






Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Kopjar B, et al. 2019 <sup>38</sup>	Clinical and functional outcomes of JOURNEY <sup>®</sup> II CR total knee system. Interim results of an ongoing, prospective, multicenter study	ISTA Congress		CR	✓		✓	✓	
Kosse NM, et al. 2018 <sup>27</sup>	Minor adaptations in implant design bicruciate-substituted total knee system improve maximal flexion	EFFORT Congress		BCS	✓				
Lutes W, et al. 2018 <sup>26</sup>	Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics	SICOT OW Congress		CR	✓	✓	✓		
Mayman DJ, et al. 2018 <sup>31</sup>	Hospital related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patients	ISPOR Symposium		BCS		✓			✓
Murakami K, et al. 2018 <sup>18</sup>	In vivo kinematics of gait in posterior-stabilized and bicruciate-stabilized total knee arthroplasties using image-matching techniques	Int Orthop		BCS	✓				








# All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Murakami K, et al. 2018 <sup>22</sup>	Knee kinematics in bi-cruciate stabilized total knee arthroplasty during squatting and stair-climbing activities	J Orthop		BCS	✓				
Murakami K, et al. 2018 <sup>68</sup>	Preoperative tibial mechanical axis orientation and articular surface design influence on the coronal joint line orientation relative to the ground during gait after total knee arthroplasties	Knee Surg Sports Traumatol Arthrosc		BCS	✓				
Moewis P, et al. 2020 <sup>69</sup>	Retention of posterior cruciate ligament alone may not achieve physiological knee joint kinematics after total knee arthroplasty: a retrospective study	J Bone Joint Surg Am		BCS, CR	✓		✓		
Nodzo SR, et al. 2018 <sup>29</sup>	The bicruciate substituting knee design and initial experience	Tech Orthop		BCS	✓	✓	✓		
Oikonomidis L, et al. 2020 <sup>70</sup>	The Journey bicruciate knee replacement: design modifications yield better early functional results and reduce complications	J Knee Surg		BCS	✓		✓		






# All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Patel AR, et al. 2019 <sup>32</sup>	Hospital related clinical and economical outcomes of two premium knee systems in total knee arthroplasty (TKA) patients	ISPOR Symposium		BCS		✓			✓
Salzman M, et al. 2017 <sup>71</sup>	Does postoperative mechanical axis alignment have an effect on clinical outcome of primary total knee arthroplasty? A retrospective cohort study	Open Orthop J		BCS	✓		✓		
Snyder MA, et al. 2018 <sup>33</sup>	A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY <sup>®</sup> II bi-cruciate stabilizing knee system and total hip arthroplasty patients	Orthop Trauma Prosth		BCS	✓	✓	✓		
Takubo A, et al. 2017 <sup>28</sup>	Comparison of muscle recovery following bicruciate substituting versus posterior stabilized total knee arthroplasty in an Asian population	J Knee Surg		BCS	✓				
Tomite T, et al. 2016 <sup>72</sup>	Gait analysis of conventional total knee arthroplasty and bicruciate stabilized total knee arthroplasty using a triaxial accelerometer	Case Report Orthop		BCS	✓				

# All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Tomite T, et al. 2021 <sup>25</sup>	Evaluation of anteroposterior accelerometric change after bi-cruciate stabilized total knee arthroplasty and posterior stabilized total knee arthroplasty	Knee		BCS	✓				
West JA, et al. 2019 <sup>34</sup>	Clinical outcomes and patient satisfaction after total knee arthroplasty: a follow-up of the first 50 cases by a single surgeon	J Int Med Res		BCS	✓	✓	✓		
Zambianchi F, et al. 2018 <sup>74</sup>	Changes in total knee arthroplasty design affect in vivo kinematics in a redesigned total knee system: A fluoroscopy study	Clin Biomech		BCS	✓				

Smith & Nephew Pty Ltd  
Australia  
T +61 2 9857 3999  
F +61 2 9857 3900  
[smith-nephew.com/en-au](http://smith-nephew.com/en-au)

Smith & Nephew Ltd  
New Zealand  
T +64 9 820 2840  
F +64 9 820 2841  
[smith-nephew.com/en-nz](http://smith-nephew.com/en-nz)

♦Trademark of Smith+Nephew  
All trademarks acknowledged  
© February 2025 Smith+Nephew  
17058-anz V4 REVA 02/25

This material is intended for healthcare professionals.  
For detailed product information, including indications  
for use, contraindications, precautions and warnings,  
please consult the product's applicable Instructions  
for Use (IFU) prior to use.

Developed by Evidence Communications,  
Global Clinical & Medical Affairs

**Smith+Nephew**

## References

- Collins M, Lavigne M, Girard J, Vendittoli PA. Joint perception after hip or knee replacement surgery. *Orthop Traumatol Surg Res.* 2012;98:275–280.
- Noble PC, Gordon MJ, Weiss JM, Reddix RN, Conditt MA, Mathis KB. Does total knee replacement restore normal knee function? *Clin Orthop Relat Res.* 2005;431:157–165.
- Scott CEH, Howie CR, MacDonald D, Biant LC. Predicting dissatisfaction following total knee replacement. *J Bone Joint Surg Am.* 2010;92-B(9):1253–1258.
- Lee A, Park J, Lee S. Gait analysis of elderly women after total knee arthroplasty. *J Phys Ther Sci.* 2015;27:591–595.
- McClelland JA, Webster KE, Feller JA, Menz HB. Knee kinematics during walking at different speeds in people who have undergone total knee replacement. *The Knee.* 2010;18:151–5.
- Dorr LD, Ochsner JL, Gronley J, Perry J. Functional Comparison of Posterior Cruciate-Retained versus Cruciate-Sacrificed Total Knee Arthroplasty. *Clin Orthop Relat Res.* 1988;236:36–43.
- Kramers-de Quervain IA, Stussi E, Muller R, Drobny T, Munzinger U, Gschwend N. Quantitative Gait Analysis After Bilateral Total Knee Arthroplasty With Two Different Systems Within Each Subject. *J Arthroplasty.* 1997;12(2):168–179.
- Saari T, Tranberg R, Zügner R, Uvehammer J, Karrholm J. Changed gait pattern in patients with total knee arthroplasty but minimal influence of tibial insert design. *Acta Orthop.* 2005;76(2):253–260.
- Thomas AC, Stevens-Lapsley JE. Importance of attenuating quadriceps activation deficits after total knee arthroplasty. *Exerc Sport Sci Rev.* 2012;40:95–101.
- Lester DK, Shantharam R, Zhang K. Dynamic electromyography after cruciate-retaining total knee arthroplasty revealed a threefold quadriceps demand compared with the contralateral normal knee. *J Arthroplasty.* 2013;28(4):557–562.
- Birdsall PD, Hayes JH, Cleary R, Pinder IM, Moran CG, Sher JL. Health outcome after total knee replacement in the very elderly. *J Bone Joint Surg Br.* 1999;81:660–662.
- Evans JT, Evans JP, Walker RB, Blom AW, Sayers A, Whitehouse MR. How long does a hip replacement last? A systematic review and meta-analysis of case series and national registry reports with more than 15 years of follow-up. *Lancet.* 2019;393:647–654.
- Angerame MR, Holst DC, Jennings JM, Komistek RD, Dennis DA. Total knee arthroplasty kinematics. *J Arthroplasty.* 2019;34:2502–2510.
- Van Onsem S, Verstraete M, Van Eenoo W, Van Der Straten C, Victor J. Are TKA kinematics during closed kinetic chain exercises associated with patient-reported outcomes? A preliminary analysis. *Clin Orthop Relat Res.* 2020;478:255–263.
- Kage T, Inui, H, Tomita T, et al. The association between in vivo knee kinematics and patient-reported outcomes during squatting in bicruciate-stabilized total knee arthroplasty. *J Knee Surgery.* 2021. doi:10.1055/s-0041-1723763.
- Grieco TF, Sharma A, Dessinger GM, Cates HE, Komistek RD. In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy. *J Arthroplasty.* 2018;33:565–571.
- Smith LA, Nachtrab J, LaCour M, Cates H, Freeman MG, Komistek RD. In vivo knee kinematics: how important are the roles of femoral geometry and the cruciate ligaments? *J Arthroplasty.* 2020; doi: <https://doi.org/10.1016/j.arth.2020.10.020>.
- Murakami K, Hamai S, Okazaki K, et al. In vivo kinematics of gait in posterior-stabilized and bicruciate-stabilized total knee arthroplasties using image matching techniques. *Int Orthop.* 2018;42:2573–2581.
- Kaneko T, Kono N, Mochizuki Y, Hada M, Toyoda S, Musha Y. Bi-cruciate substituting total knee arthroplasty improved medio-lateral instability in mid-flexion range. *J Orthop.* 2017;14:201–206.
- Catani F, Ensini A, Belvedere C, et al. In vivo kinematics and kinetics of a bi-cruciate substituting total knee arthroplasty: a combined fluoroscopic and gait analysis study. *J Orthop Res.* 2009;27(12):1569–1575.
- Iriuchishima T, Ryu K. A Comparison of rollback ratio between bicruciate substituting total knee arthroplasty and Oxford unicompartmental knee arthroplasty. *J Knee Surg.* 2018;31:568–572.
- Murakami K, Hamai S, Okazaki K, et al. Knee kinematics in bicruciate stabilized total knee arthroplasty during squatting and stair-climbing activities. *J Orthop.* 2018;15:650–654.
- Carpenter RD, Brilhault J, Majumdar S, Ries MD. Magnetic resonance imaging of in vivo patellofemoral kinematics after total knee arthroplasty. *Knee.* 2009;16(5):332–336.
- Di Benedetto P, Vidi D, Colombo, Buttironi MM, Cainero V, Causero A. Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study. *Acta Biomed.* 2019;90:91–97.
- Tomite, T, Saito, H, Kijima H, et al. Evaluation of anteroposterior accelerometric change after bi-cruciate stabilized total knee arthroplasty and posterior stabilized total knee arthroplasty. *Knee.* 2021;32:121–130.
- Lutes W, Fitch D. Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics. Presented at: 39th SICOT Orthopaedic World Congress; October 10–13, 2018; Montreal, Canada.
- Kosse NM, Heesterbeek PJC, Defoort KC, Wymenga AB, van Hellemond G. Minor adaptations in implant design bicruciate-substituted total knee system improve maximal flexion. Poster presented at: 2nd World Arthroplasty Congress; 19–21 April, 2018; Rome, Italy.
- Takubo A, Ryu K, Iriuchishima T, Tokuhashi Y. Comparison of muscle recovery following bicruciate substituting versus posterior stabilized total knee arthroplasty in an Asian population. *J Knee Surg.* 2017;30:725–729.
- Nodzo SR, Carroll KM, Mayman DJ. The bicruciate substituting knee design and initial experience. *Tech Orthop.* 2018;33:37–41.
- Hyodo K, Kanamori A, Kadone H, Takahashi T, Kajiwara M, Yamazaki M. Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty. *Arthroplast Today.* 2020;6:338–342.
- Mayman DJ, Patel AR, Carroll KM. Hospital related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patients. Poster presented at: ISPOR Symposium; May 19–23, 2018; Baltimore, Maryland, USA.
- Patel AR, Delhougne G. Hospital related clinical and economical outcomes of two premium knee systems in total knee arthroplasty (TKA) patients. Poster presented at: ISPOR Symposium; May 18–22, 2019; New Orleans, Louisiana, USA.
- Snyder MA, Simpson A, Gregg J, Levit A. A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY II bi-24. cruciate stabilizing knee system and total hip arthroplasty patients. *Orthop Trauma Prost.* 2018; 2018; 3:5–10.
- West JA, Scudday T, Anderson S, Amin NH. Clinical outcomes and patient satisfaction after total knee arthroplasty; a follow-up of the first 50 cases by a single surgeon. *J Int Med Res.* 2019;47:1667–1676.
- Boese K, MacDonald J, Huang W, et al. Early clinical and patient-reported results of a bi-cruciate retaining total knee implant: six-month results of a prospective multicentre study of 149 primary TKAs. Poster presented at: European Orthopaedic Research Society (EORS); October 2–4, 2019; Maastricht, the Netherlands.
- Heir S, Catani F, van Hellemond G, et al. Clinical and functional outcomes of a second generation guided motion total knee arthroplasty system: two-year results of a prospective multicentre study. Abstract number P37 presented at: European Knee Society; May 2–3, 2019; Valencia, Spain.
- Harris AI, Luo TD, Lang JE, Kopjar B. Short-term safety and effectiveness of a second generation motionguided total knee system. *Arthroplast Today.* 2018;4:240–243.
- Kopjar B, Archan P, MacDonald J, et al. Clinical and functional outcomes of JOURNEY II CR total knee system. Interim results of an ongoing, prospective, multicenter study. Abstract presented at: International Society for Technology in Arthroplasty Congress (ISTA); October 2–5, 2019; Toronto, Canada.

## References (cont.)


39. Harris AI, Christen B, Malcorps JJ, et al. Mid-term performance of a guided motion bicruciate stabilized total knee system; results from the international study of over 2,000 consecutive primary total knee arthroplasties. *J Arthroplasty*. 2019;34:S201–S208.
40. Harris A, Luo TD, Lang JE, et al. Guided motion total knee arthroplasty system: five-year outcomes of the prospective multicentre US study. Presentation number 2761 presented at: EFORT; June 5–7, 2019; Lisbon, Portugal.
41. Matsuda S, Kawahara S, Okazaki K, Tashiro Y, Iwamoto Y. Postoperative alignment and ROM affect patient satisfaction after TKA. *Clin Orthop Relat Res*. 2013;471:127–133.
42. Ishibashi T, Tomita T, Yamazaki T, Tsuji S, Yoshikawa H, Sugamoto K. Kinematics of bicruciate and posterior stabilized total knee arthroplasty during deep knee flexion and stair climbing. *J Orthop Res*. 2020; doi: 10.1002/jor.24773.
43. Walker LC, Clement ND, Bardgett M, et al. The WOMAC score can be reliably used to classify patient satisfaction after total knee arthroplasty. *Knee Surg Sports Traumatol Arthrosc*. 2018;26:3333–3341.
44. Lyu H, Wick E, Housman M, Freischlag JA, Makary MA. Patient satisfaction as a possible indicator of quality surgical care. *JAMA Surg*. 2013;148:362–367.
45. Otani K, Waterman B, Faulkner K, et al. Patient satisfaction: focusing on “excellent”. *J Healthc Manag*. 2009;54:93–102.
46. Hamilton DF, Howie CR, Burnett R, Simpson AHRW, Patton JT. Dealing with the predicted increase in demand for revision total knee arthroplasty; challenges, risks and opportunities. *Bone Joint J*. 2015;97-b:723–728.
47. Delanois RE, Mistry JB, Gwam CU, Mohamed NS, Choksi US, Mont MA. Current epidemiology of revision total knee arthroplasty in the United States. *J Arthroplasty*. 2017;32:2663–2668.
48. Amemiya K, Kaneko T, Omata M, et al. Anatomical bi-cruciate retaining TKA improves gait ability earlier than bi-cruciate stabilized TKA based on triaxial accelerometry data: a prospective cohort study. *AP-SMART*. 2021;25:35–41.
49. Kono K, Inui H, Tomita T, Yamazaki T, Taketomi S, Tanaka S. In vivo kinematics of bicruciate-retaining total knee arthroplasty with anatomical articular surface under high-flexion conditions. *J Knee Surg*. 2021; 34 4:452–459.
50. Husain A, Lee GC. Establishing realistic patient expectations following total knee arthroplasty. *J Am Acad Orthop Surg*. 2015;23:707–713.
51. Biały M, Corte RD, Dec J, Pierzchała A, Gnat R. Evaluation of functional status in patients after total knee arthroplasty: a comparison of the JOURNEY II bi-cruciate stabilized total knee system and Genesis II cruciate-retaining implant. *Physiotherapy review*. 2021;25;3:24–34. doi:10.5114/phr.2021.109028.
52. Christen B, Kopjar B. Second-generation bi-cruciate stabilized total knee system has a lower reoperation and revision rate than its predecessor. *Arch Orthop Trauma Surg*. 2018;138:1591–1599.
53. Coomer S, LaCour M, Khasian M, Cates H, Komistek R. Determining patella function in non-implanted knees having functional cruciate ligaments and subjects having a bicruciate-stabilised total knee arthroplasty. *Orthopaedic Proceedings*. 2021;103:4–4.
54. Harris A, O’Grady C, Sensiba PR, et al. Guided motion total knee arthroplasty (TKA) system in younger patients has a lower revision rate than registry controls: results from the international multicenter study with up to 6 years follow-up. Abstract number P35 presented at: European Knee Society; May 2–3, 2019; Valencia, Spain.
55. Harris A, O’Grady C, Sensiba PR, et al. Guided motion total knee arthroplasty (TKA) in patients with BMI of 40kg/m<sup>2</sup> or more: results from the international multicenter study of 2,059 primary TKAs with up to 6 years follow-up. Abstract number P36 presented at: European Knee Society; May 2–3, 2019; Valencia, Spain.
56. Hino K, Kutsuna T, Watamori K, et al. Bi-cruciate substituting total knee arthroplasty provides varus–valgus stability throughout the midflexion range. *Knee*. 2018;25:897–902.
57. Hommel H, Wilke K. Good early results obtained with a guided-motion implant for total knee arthroplasty: A consecutive case series. *Open Orthop J*. 2017;11:51–56.
58. Inui H, Taketomi S, Yamagami R, Kawaguchi K, Nakazato K, Tanaka S. The relationship between anteroposterior stability and medial-lateral stability of the bi-cruciate stabilized total knee arthroplasty. *Knee*. 2018;25:1247–1253.
59. Inui H, Taketomi S, Yamagami R, Shirakawa N, Kawaguchi K, Tanaka S. The relationship between soft-tissue balance and intraoperative kinematics of guided motion total knee arthroplasty. *J Knee Surg*. 2019;32:91–96.
60. Iriuchishima T, Ryu K. Bicruciate substituting total knee arthroplasty improves stair climbing ability when compared with cruciate-retain or posterior stabilizing total knee arthroplasty. *Indian J Orthop*. 2019. DOI:10.4103/ortho. IJOOrtho\_392\_18.
61. Ishida K, Shibamura N, Toda A, et al. Comparison of intra-operative navigation-based kinematics between bi-cruciate-stabilised total knee arthroplasty (TKA) and conventional posterior stabilised TKA. *Orthopaedic proceedings*. 2017;99B:Suppl 4.
62. Itou J, Kuwashima U, Itoh M, Okazaki K. Anterior prominence of the femoral condyle varies among prosthesis designs and surgical techniques in total knee arthroplasty. *BMC musculoskeletal disorders*. 2021;22:1–8.
63. Kaneko T, Kono N, Mochizuki Y, et al. The influence of compressive forces across the patellofemoral joint on patient reported outcome after bi-cruciate stabilized total knee arthroplasty. *Bone Joint J*. 2018; 100-B:1585–1591.
64. Kaneko T, Kono N, Mochizuki Y, et al. The influence of tibiofemoral joint forces on patient-reported outcome measurements after bicruciate stabilized total knee arthroplasty. *J Orthop Surg*. 2020;28:1–2.
65. Kiyohara M, Hamai S, Gondo H, et al. Comparison of in vivo knee kinematics before and after bicruciate-stabilized total knee arthroplasty during squatting. *BMC Musculoskeletal Disorders*. 2021;22:1–10.
66. Kono K, Inui H, Tomita T, et al. Bicruciate-stabilised total knee arthroplasty provides good functional stability during high-flexion weight-bearing activities. *Knee Surg Sports Traumatol Arthrosc*. 2019;27:2096–2103.
67. Kono K, Inui H, Tomita T, et al. In vivo kinematics and cruciate ligament forces in bicruciate-retaining total knee arthroplasty. *Scientific Reports*. 2021;11:1–10.
68. Murakami K, Hamai S, Okazaki K, et al. Preoperative tibial mechanical axis orientation and articular surface design influence on the coronal joint line orientation relative to the ground during gait after total knee arthroplasties. *Knee Surg Sports Traumatol Arthrosc*. 2018;26:3368–3376.
69. Moewis P, Duda GN, Trepczynski A, et al. Retention of posterior cruciate ligament alone may not achieve physiological knee joint kinematics after total knee arthroplasty: a retrospective study. *J Bone Joint Surg Am*. 2020; doi: 10.2106/ JBJS.20.00024.
70. Oikonomidis L, Santini AJA, Davidson JS, Banks JV, Phillipson A, Pope J. The JOURNEY bicruciate knee replacement: design modifications yield better early functional results and reduce complications. *J Knee Surg*. 2020; doi: 10.1055/s-0040-1718599.
71. Salzmann M, Fennema P, Becker R, Hommel H. Does postoperative mechanical axis alignment have an effect on clinical outcome of primary total knee arthroplasty? A retrospective cohort study. *Open Orthop J*. 2017;11:1330–1336.
72. Tomite T, Saito H, Aizawa T, Kijima H, Miyakoshi N, Shimada Y. Gait analysis of conventional total knee arthroplasty and bicruciate stabilized total knee arthroplasty using a triaxial accelerometer. *Case Report Orthop*. 2016;6875821.
73. Tria Jr AJ, Bateman DK, Preston JS, Diamond KB. Bicruciate retaining total knee arthroplasty: short-term clinical outcomes using a novel prosthesis. *J. long-term eff med implants*. 2021;31:3.
74. Zambianchi F, Fiacchi F, Lombardi V et al. Changes in total knee arthroplasty design affect in vivo kinematics in a redesigned total knee system: A fluoroscopy study. *Clin Biomech*. 2018;54:92–102.
75. Goh GS, Bin Abd Razak HR, Tay DKJ, Lo N, Yeo S. Early post-operative oxford knee score and knee society score predict patient satisfaction 2 years after total knee arthroplasty. *Arch Orthop Trauma Surg*. 2021;141:129–137. doi: 10.1007/s00402-020-03612-2.
76. Wendelspiess S, Kaelin R, Vogel N, Rychen T, Arnold MP. No difference in patient-reported satisfaction after 12 months between customised individually made and off-the-shelf total knee arthroplasty. *Knee Surg Sports Traumatol Arthrosc*. 2022. doi: 10.1007/s00167-022-06900-z.
77. Ast MP, Kolin DA, Carroll KM, Davis D, Pearle AD, Mayman DJ, Ong AC. Does a uni “feel better” than a total knee? Not necessarily, when using modern implant designs. *HSS Journal®: The Musculoskeletal Journal of Hospital for Special Surgery*. 2022;19(1), 32–36. doi:10.1177/15563316221131251.
78. Inui H, Yamagami R, Kono K, Kawaguchi K, Taketomi S, Saita K, Tanaka S. Midterm clinical results of bi-cruciate stabilized total knee arthroplasty compared with posterior stabilized total knee arthroplasty with single radius design. *Journal of Joint Surgery and Research*. 2023;1(1), 117–122. doi:10.1016/j.jjoirs.2023.04.001.
79. Inoue A, Arai Y, Nakagawa S, Yoshihara Y, Kobayashi M, Takahashi K. The clinical results of bi-cruciate vs posterior stabilized total knee arthroplasty for flexion contracture in osteoarthritic knee. *Journal of Orthopaedic Surgery*. 2023;31(3). doi:10.1177/10225536231190524.

# Abbreviations

<b>ACL</b>	Anterior cruciate ligament	<b>PFR</b>	Posterior femoral rollback
<b>AOANJRR</b>	Australian Orthopaedic Association National Joint Replacement Registry	<b>PS</b>	Posterior stabilised
<b>BCS</b>	Bicruciate stabilised	<b>ROM</b>	Range of motion
<b>CR</b>	Cruciate retaining	<b>THA</b>	Total hip arthroplasty
<b>FJS-12</b>	Forgotten joint score-12	<b>TKA</b>	Total knee arthroplasty
<b>ITB</b>	Iliotibial band	<b>TKS</b>	Total knee system
<b>KOOS</b>	Knee injury and osteoarthritis outcome score	<b>UCLA</b>	University of California Los Angeles
<b>KSS</b>	Knee society score	<b>UKA</b>	Unicompartmental knee arthroplasty
<b>LOS</b>	Length of stay	<b>WOMAC</b>	The Western Ontario and McMaster Universities Arthritis Index

# JOURNEY II TKA: further information

## Learn more:



Access additional  
JOURNEY II TKA  
resources and training at  
**+ Education Unlimited**

## Explore similar resources:



Robotics compendium  
of clinical evidence



Contact your sales representative for further information