

# Evidence in focus

## REGENESORB Absorbable Biocomposite Material

A novel formulation of materials with long histories of clinical use




### Summary

REGENESORB is a bioabsorbable biocomposite material designed to remain mechanically stable for a minimum of six months before being absorbed and replaced by bone within 24 months.<sup>1</sup>

### Material overview

REGENESORB is comprised of the copolymer PLGA (poly [L-lactide co-glycolide]; 65%) combined with two fillers, beta-tricalcium phosphate (β-TCP; 15%) and calcium sulfate (20%), which is specific to REGENESORB. Each of these components have long histories of clinical use in implanted surgical devices.<sup>2-4</sup>

65%




**PLGA**

**Purpose**  
Responsible for providing overall structural integrity of the implant

**Action**  
PLGA is composed of two polymers, poly L-lactic acid (PLLA) and polyglycolic acid (PGA). The absorption characteristics of PLGA are dependent on the ratio of PLLA:PGA,<sup>5</sup> for REGENESORB the ratio is 85:15.<sup>1</sup>

**Time to absorbance**  
24 months<sup>1</sup>

20%



**Calcium sulfate**


**Purpose**  
Promotes bone growth in the early stage of healing<sup>6</sup>

**Action**  
Calcium sulfate is osteoconductive.<sup>4,7,8</sup> When calcium sulphate dissolves it chemically activates the cycle of bone growth, reacting with platelets to promote bone growth. The release of calcium ions activates platelets to produce local growth factors which stimulate osteogenic differentiation of mesenchymal stem cells.<sup>9-11</sup>

**Time to absorbance**  
4-8 weeks<sup>6</sup>

**Specific to REGENESORB**

15%



**β-TCP**

**Purpose**  
Acts as a scaffold for new bone ingrowth

**Action**  
β-TCP is a calcium phosphate filler which is osteoconductive.<sup>7,8</sup> Exposed calcium phosphate on the implant supports local osteoblast adhesion and new bone build up. β-TCP has a similar compressive and tensile strength to cancellous bone.<sup>7</sup>

**Time to absorbance**  
6-18 months<sup>6</sup>

# Evidence in focus

## Preclinical testing

The purpose of biocomposite, bioabsorbable surgical implants is to provide fixation during the initial healing period and for some time after until the material is completely absorbed and replaced by new native tissue. The aim is to reduce the time of absorption and replacement, providing fast restoration of new native tissue.

**REGENESORB hypothesis:** The unique composition of REGENESORB will be an improvement on currently available products. REGENESORB will be replaced by bone within 24 months and work faster than implants composed of PLLA-hydroxyapatite (HA).

### In preclinical studies, REGENESORB is replaced by bone within 24 months<sup>1</sup>



#### Study design

9 x 10mm interference screws moulded in REGENESORB were implanted into the distal femur or proximal tibia of sheep

Histological and micro-computer tomography analysis was conducted at 6, 12, 18 and 24 months



#### Results

6 months: REGENESORB screws largely intact

12 months: REGENESORB screws had begun to degrade and were partially fragmented with indications of new bone growth

18 months: REGENESORB screws degrading and evidence of new bone growth

**At 24 months, REGENESORB material was absorbed and replaced by bone**




Figure 1. Absorption profile of the REGENESORB implant

### In preclinical studies, REGENESORB is absorbed and replaced by bone faster than PLLA-HA<sup>12</sup>



#### Study design

REGENESORB or PLLA-HA suture anchors were implanted into the distal femur or proximal tibia of sheep

Histological and micro-CT analysis was conducted at 12 and 18 months



#### Results

12 months: Varying levels of degradation of REGENESORB, no degradation of PLLA-HA implant

**At 18 months, REGENESORB demonstrated significant new bone fill and PLLA-HA showed no visible signs of absorption**




Figure 2. REGENESORB material vs PLLA-HA at 18 months

## References

1. Al-Beik J. Evaluation of new bioresorbable screws in an ovine bone defect model. Smith & Nephew internal report: 15000897. 2. Suchenski M, McCarthy MB, Chowanec D, Hansen D, KcKinnon W, Apostolakis J, Arciero R, Mazzocca AD. Material properties and composition of soft-tissue fixation. *Arthroscopy*. 2010;26:821–831. 3. Albee FH. Studies in bone growth: triple calcium phosphate as a stimulus to osteogenesis. *Ann Surg*. 1920;71:32–36. 4. Peltier LF, Bickel EY, Lillo R, Thein MS. The use of plaster of paris to fill defects in bone. *Ann Surg*. 1957;146:61–69. 5. Chu C. Biodegradable polymeric biomaterials. In: Bronzino JD, ed. *The Biomedical Engineering Handbook*. Hartford: CRC Press; 1995. 6. Costantino PD & Friedman CD. Synthetic bone graft substitutes. *Otolaryngol Clin North Am*. 1994;27:1037–74. 7. Hak DJ. The use of osteoconductive bone graft substitutes in orthopaedic trauma. *J Am Acad Orthop Surg*. 2007;15:525–536. 8. Allison DC, Lindberg AW, Samimi B, Mirzayan R, Menendez LR. A comparison of mineral bone graft substitutes for bone defects. *US Oncology and Hematology*. 2011;7:38–49. 9. Strocchi, R., Orsini, G., Iezzi, G., Scarano, A., Rubini, C., Pecora, G., Piattelli, A., Bone regeneration with calcium sulfate: evidence for increased angiogenesis in rabbits. *J Oral Implant*. 2002;28:273–278. 10. Intini G, Andreana S, Margarone JE, 3rd, Bush PJ, Dziak R. Engineering a bioactive matrix by modifications of calcium sulfate. *Tissue Eng*. 2002;8:997–1008. 11. Groeneveld EH, Burger EH. Bone morphogenetic proteins in human bone regeneration. *Eur J Endocrinol*. 2000;142:9–21. 12. Al-Beik J. Evaluation of bioresorbable suture anchors in an ovine bone defect model. Smith & Nephew internal report: 15001194. 13. Prior A. Histological evaluation of specimens from resorbable screw study (RS-II) 6-month post-implantation. Smith & Nephew internal report: WRP-TE045-700-05, 2009. 14. Butcher K. Micro-CT analysis and histological evaluation of specimens from resorbable screw study (RS-II/OMI-08) 12-month post-implantation. Smith & Nephew internal report: WRP-TE045-700-06, 2009. 15. Prior A. Micro-CT analysis and histological evaluation of specimens from resorbable screw study (RS-II/OMI-08) 18-month post-implantation. Smith & Nephew internal report: WRP-TE045-700-07, 2009. 16. Butcher K. Micro-CT analysis and histological evaluation of specimens from resorbable screw study (RS-II/OMI-08) 24-month post-implantation. Smith & Nephew internal report: WRP-TE045-700-08, 2010. 17. Anderson M & Butcher K. Histological evaluation of specimens from study IVB/TE045/900/01 (18-month time-point). Smith & Nephew internal report: WRP TE045/900/03, 2011.