Publication summary: Schlegel TF, et al. J Shoulder Elbow Surg (2020)\*

# Smith-Nephew

REGENETEN<sup>o</sup> Bioinductive Implant reduced the size of partial-thickness rotator cuff tears by over 50% in all tear locations

## + Plus points



### Overview

- Prospective, multicentre case series of isolated bioinductive repair with the REGENETEN Implant for intermediate- (n=12) or high-grade (n=21) partial-thickness rotator cuff tears
  - 11 articular-sided, 10 bursal-sided, 4 intrasubstance and 8 hybrid tears were enrolled
- 31/33 patients were available at 24-month follow-up
- Outcomes assessed at 3, 12 and 24 months included:
  - MRI assessment of defect fill, tendon thickness and progression to a full-thickness tear
  - Clinical outcomes: American Shoulder and Elbow Surgeons (ASES) and Constant-Murley Shoulder (CMS) scores

### Results

- At 24 months, 87% of tears (26/30) had reduced in size by >50% from preoperative measurements (Figure)
- No compliant patient progressed to a full-thickness tear; one non-compliant patient progressed ~1 month after surgery
- Significant increase in mean tendon thickness from preoperative values to 24 months of 1.2 and 1.8mm in intermediate- and high-grade tears, respectively (p≤0.012)
- Improvements in CMS score met or exceeded minimal clinically important differences (MCIDs) for all patients with intermediate-grade tears and 94% of patients with high-grade tears
- No significant difference in clinical outcomes, tear healing or tendon thickness based on tear location: articular-sided, bursal-sided, intrasubstance or hybrid
- All 28 surveyed patients would recommend the procedure to a friend



Figure: Degree of reduction in tear size from baseline to 24 months after isolated bioinductive repair with the REGENETEN Implant. <sup>†</sup>Tear size could not be determined for one patient with a high-grade tear.

www.smith-nephew.com/education

### Conclusions

In a prospective, multicentre case series of isolated bioinductive repair with the REGENETEN Implant, 87% of partial-thickness tears reduced in size by >50% at 24 months. Improvements in clinical outcomes, tear healing and tendon thickness were consistent in all tear locations including articular-sided and intrasubstance tears.

#### Citation

\*Schlegel TF, Abrams JS, Angelo RL, et al. Isolated bioinductive repair of rotator cuff tears. J Shoulder Elbow Surg. 2020;Nov 18 [ePub ahead of print]

#### Available at: Journal of Shoulder and Elbow Surgery

Sports Medicine, Smith & Nephew, Inc. 150 Minuteman Road, Andover, MA 01810 USA. 27635-en V1 1120. Published November 2020. ©2020 Smith+Nephew. All rights reserved. °Trademark of Smith+Nephew. All Trademarks acknowledged.