

## + Evidence in focus

### Advanced healing solutions

Redefining healing potential for rotator cuff repair

[View the evidence >>](#)




HEALICOIL<sup>®</sup>  
REGENESORB<sup>®</sup>  
Suture Anchor




HEALICOIL<sup>®</sup>  
KNOTLESS<sup>®</sup>  
Suture Anchor



REGENETEN<sup>®</sup>  
Bioinductive Implant



REGENESORB<sup>®</sup>  
Material



# Advanced healing solutions: evidence collection

This Evidence Collection summarises the clinical evidence on the S+N **Advanced Healing Solutions** portfolio:

**Redefining  
healing potential**  
for rotator cuff repair

## Key



Imaging  
outcomes



Patient  
outcomes



Biopsy



Biomechanical

## Functionality



Link to  
full summary



View and share on  
Education + Evidence



Return to menu




Reference list

## Abbreviations


3D-FEM = 3-dimensional finite element method  
ASES = American Shoulder and Elbow Surgeons  
FT = full-thickness  
MCID = minimal clinically important difference  
MRI = magnetic resonance imaging  
PT = partial-thickness

RCR = rotator cuff repair  
SANE = single assessment numeric evaluation  
SS = supraspinatus  
VR-12 PCS = Veterans RAND 12 Item Health Survey (VR-12)  
Physical Component Score  
WORC = Western Ontario Rotator Cuff Index


Click to view the evidence by product:



**REGENETEN<sup>®</sup>**  
Bioinductive Implant




**HEALICOIL<sup>®</sup>**  
Suture Anchor




**REGENESORB<sup>®</sup>**  
Material

# Advanced healing solutions: evidence collection


---



**REGENETEN<sup>®</sup>**  
Bioinductive Implant




**HEALICOIL<sup>®</sup>**  
Suture Anchor



**REGENESORB<sup>®</sup>**  
Material

---


# Advanced healing solutions: evidence collection



**REGENETEN<sup>®</sup>**  
Bioinductive Implant







**HEALICOIL<sup>®</sup>**  
Suture Anchor



**REGENESORB<sup>®</sup>**  
Material

## Changing the course of rotator cuff disease


The REGENETEN Bioinductive Implant stimulates the body's natural healing response to support new tendon growth and disrupt disease progression.<sup>1,2</sup>

Author	Key finding	Tear type				
<b>Key evidence</b>						
Arnoczky SP, et al.	<a href="#">REGENETEN<sup>®</sup> Bioinductive Implant rapidly absorbed and replaced by tendon-like tissue within six months<sup>3</sup></a>	FT			✓	
Bokor DJ, et al.	<a href="#">Improved tendon thickness and integrity was sustained five years after treatment with the REGENETEN<sup>®</sup> Bioinductive Implant<sup>4</sup></a>	PT	✓	✓		
Schlegel TF, et al.	<a href="#">REGENETEN<sup>®</sup> Bioinductive Implant induces new tendon-like tissue in patients with partial-thickness tears of the supraspinatus tendon<sup>2</sup></a>	PT	✓	✓		
Thon SG, et al.	<a href="#">REGENETEN<sup>®</sup> Bioinductive Implant leads to tissue induction and high rates of tendon healing in patients with large and massive rotator cuff tears<sup>5</sup></a>	FT	✓	✓		
<b>Supporting evidence</b>						
Bokor DJ, et al.	<a href="#">REGENETEN<sup>®</sup> Bioinductive Implant maintains repair integrity in full-thickness rotator cuff tears<sup>6</sup></a>	FT	✓	✓		
Mcintyre LF, et al	<a href="#">Clinically meaningful improvements in pain and function in a multicentre registry study of the REGENETEN<sup>®</sup> Bioinductive Implant in rotator cuff repair<sup>7</sup></a>	Both		✓		


# Advanced healing solutions: evidence collection



**REGENETEN<sup>®</sup>**  
Bioinductive Implant







**HEALICOIL<sup>®</sup>**  
Suture Anchor



**REGENESORB<sup>®</sup>**  
Material

## Open-architecture: Now open to even greater possibilities


The unique open-architecture design of HEALICOIL anchors reduces the amount of implanted material in the shoulder compared to traditional, solid-core anchors and may provide a biologic healing advantage.<sup>8</sup>

Author	Key finding	Tear type				
<b>Key evidence</b>						
Chahla J, et al.	<a href="#">Significantly greater bone density surrounding HEALICOIL<sup>®</sup> PK Suture Anchors versus TWINFIX<sup>®</sup> PK Suture Anchors six months after RCR<sup>9</sup></a>	FT	✓	✓		
Clark TR, et al.	<a href="#">Rotator cuff thickness significantly greater in patients who received vented compared with non-vented suture anchors at six weeks<sup>8</sup></a>	FT	✓			
<b>Supporting evidence</b>						
Kim J-H, et al.	<a href="#">Significantly improved bone ingrowth with HEALICOIL<sup>®</sup> PK Suture Anchors versus TWINFIX<sup>®</sup> HA Suture Anchors six months after rotator cuff repair<sup>10</sup></a>	Both	✓	✓		
Sano H, et al.	<a href="#">Fixation properties, stress distribution and failure patterns differ between coil-type and screw-type suture anchors for rotator cuff repair<sup>11</sup></a>	N/A				✓


# Advanced healing solutions: evidence collection



**REGENETEN<sup>◇</sup>**  
Bioinductive Implant







**HEALICOIL<sup>◇</sup>**  
Suture Anchor



**REGENESORB<sup>◇</sup>**  
Material

## A unique formulation of proven materials

REGENESORB Material is designed to provide a jump start in bone healing and formation, and is absorbed and replaced by bone within 24 months.<sup>12</sup>

Author	Key finding	Tear type				
<b>Key evidence</b>						
Vonhoegen J, et al.	<a href="#">HEALICOIL<sup>◇</sup> REGENESORB Suture Anchor mostly resorbed and replaced by new bone material within 21 months of arthroscopic rotator cuff repair<sup>13</sup></a>	FT	✓			
<b>Competitor comparison</b>						
Vonhoegen J, Sgroi M.	<a href="#">Resorption comparison of HEALICOIL REGENESORB anchors and Arthrex Biocomposite anchors</a>	FT	✓			

# REGENETEN<sup>®</sup> Bioinductive Implant rapidly absorbed and replaced by tendon-like tissue within six months

Arnoczky SP, et al. *Arthroscopy* (2017)<sup>3</sup>

## Overview

- Retrospective study of biopsies collected from 7 patients (6 FT tears, 1 PT tear) between 5 weeks and 6 months following RCR augmented with the REGENETEN Bioinductive Implant

## Results

- Rapid host cell ingrowth and early collagen formation seen at 5 weeks
- Increased collagen formation, maturation and organisation on the surface of the implant at 3 months (Figure 1)
- Implant no longer visible at 6 months, with new tendon-like tissue and highly oriented collagen indicative of functional loading (Figure 2)

## Conclusions

In a series of second-look biopsies, the REGENETEN Bioinductive Implant was rapidly resorbed and replaced by new tissue. At six months postoperatively, new tissue was tendon-like with oriented collagen indicative of functional loading.




Figure 1.  
Photomicrograph  
of implant  
surface at  
3 months.  
Remnants of  
implant still  
present (→)  
(haematoxylin  
and eosin x 100)




Figure 2.  
Photomicrograph  
of implant  
surface at  
6 months.  
No evidence  
of implant  
remaining  
(haematoxylin  
and eosin x 100)





# REGENETEN<sup>®</sup> Bioinductive Implant maintains repair integrity in FT rotator cuff tears

Bokor DJ, et al. *Muscles Ligaments Tendons J* (2015)<sup>6</sup>

## Overview

- Prospective study of 9 patients with supraspinatus tendon tears (8 medium-sized FT tears, 1 PT tear converted to a FT tear during surgery) who received the REGENETEN Bioinductive Implant in conjunction with standard repair

## Results

- No MRI evidence of re-tear or gap formation at 24 months
- Significantly greater mean tendon thickness on MRI in study patients versus published values from young, healthy adults, sustained to 24 months ( $p < 0.01$ ; Figure)
- New tissue was indistinguishable from the native tendon on MRI at 12 months
- Significant improvements in ASES score and Constant-Murley score from preoperative values at 24 months (both  $p < 0.001$ )

## Conclusions

REGENETEN Bioinductive Implant facilitates restoration of the normal tendon footprint and ultimately maintains repair integrity of FT tears at 24 months.




Figure. Mean ( $\pm$  standard error) tendon thickness in male ( $n=6$ ) and female ( $n=3$ ) patients compared to published average thicknesses (horizontal lines) for healthy males and females;  $p < 0.01$  at all time points





# REGENETEN<sup>®</sup> Bioinductive Implant induces tendon-like tissue formation in patients with PT tears of the supraspinatus tendon

Schlegel TF, et al. *J Shoulder Elbow Surg* (2017)<sup>2</sup>

## Overview

- Prospective, multicentre study of 33 patients with PT tears of the supraspinatus tendon, who received the REGENETEN Bioinductive Implant following subacromial decompression without repair

## Results

- Reduction in tear size of  $\geq 1$  grade from baseline in 31/33 patients (94%) at 1 year
  - 8/33 patients (24%) had no visible defect
- Significant increase in mean tendon thickness on 1-year MRI versus preoperative values ( $p < 0.0001$ ; Figure)
- Significant improvements in ASES pain, ASES shoulder index and Constant-Murley scores from baseline to 1 year ( $p < 0.0001$ ) were more than twice the respective MCIDs
- 30/33 patients (94%) were satisfied with the results of their procedure at 1 year

## Conclusions

REGENETEN Bioinductive Implant biologically augments healing, increasing tendon thickness and reducing tear size. Patients reported improved function, high satisfaction and rapid recovery.




Figure. Mean ( $\pm$  standard error) tendon thickness on MRI in intermediate- and high-grade tears.

At month 12,  $p = 0.003$  and  $p < 0.0001$  for intermediate- and high-grade tears, respectively, versus preoperative measures





# REGENETEN<sup>®</sup> Bioinductive Implant leads to tissue induction and high rates of tendon healing in large and massive rotator cuff tears

Thon SG, et al. *Am J Sports Med* (2019)<sup>5</sup>

## Overview

- Prospective study of 23 patients (mean age, 57.9 years) receiving a REGENETEN Bioinductive Implant in conjunction with double-row repair of large (n=11) or massive (n=12) rotator cuff tears

## Results

- No implant-related adverse events
- Tendon healing in 22/23 patients (96%) on ultrasound at 24 months
- Treatment success in 21/23 patients (91%) at 24 months; one additional clinical failure due to progression of glenohumeral osteoarthritis (Figure)
- Mean tendon thickness on ultrasound increased from 6.29mm at 3 months to 7.72mm at 12 months, decreasing to 7.28mm at 24 months

## Conclusions

In conjunction with repair of large and massive rotator cuff tears, REGENETEN Bioinductive Implant was safe, induced tissue formation and led to a high tendon healing rate in both primary and revision settings.




Figure. Percentage of patients achieving treatment success in primary and revision repairs





# Improved tendon thickness and integrity was sustained five years after treatment with the REGENETEN<sup>®</sup> Bioinductive Implant

Bokor DJ, et al. *Muscles Ligaments Tendons J* (2019)<sup>4</sup>

## Overview


- Five-year follow-up of a prospective, single-arm study evaluating REGENETEN Bioinductive Implant in lieu of repair in patients with PT rotator cuff tears
- 11/13 enrolled patients were available for assessment at 5 years

## Results

- 8/11 patients (73%) had no decline in tendon integrity between 2 and 5 years (Figure)
- Mean tendon thickness significantly decreased between 2 and 5 years (5.9 vs 5.2mm;  $p=0.0012$ ), but remained significantly greater than preoperative values (4.3mm;  $p<0.0001$ )
- Significant improvements from baseline in pain and function were sustained to 5 years ( $p\leq 0.01$ ) and were not significantly different to 2-year values

## Conclusions

REGENETEN Bioinductive Implant demonstrated sustained effectiveness in improving tendon thickness and integrity, with no decline in pain and function between 2 and 5 years.



2 years



5 years





# Clinically meaningful improvements in pain and function in a multicentre registry study of REGENETEN<sup>®</sup> Bioinductive Implant in RCR

McIntyre L, et al. *Arthroscopy* (2019)<sup>7</sup>

## Overview

- A multicentre registry study of 173 patients who received a REGENETEN Bioinductive Implant in lieu of PT repair (n=90) or in conjunction with standard repair of FT tears (n=83)

## Results

- In PT tears, mean improvements exceeded MCIDs in VAS pain from 2 weeks and ASES score from 6 weeks postoperatively (p<0.001)
- In FT tears, mean improvements exceeded MCIDs in VAS pain from 2 weeks and ASES score from 3 months postoperatively (p<0.001)
- SANE score, WORC and VR-12 PCS also improved significantly from baseline during the study period in both tear types (p<0.001)
- Postoperative recovery was rapid in PT tears (Table)

## Conclusions

REGENETEN Bioinductive Implant led to clinically meaningful improvements in pain and function for patients with PT and FT tears, with rapid postoperative recovery in PT tears.

Table. Duration of postoperative recovery in patients with partial-thickness tears

Measure	REGENETEN Bioinductive Implant	
Time in sling	10.6 days	27.7 days
	No biceps surgery	Concomitant tenodesis
Return to driving	14.6 days	
Return to work	9.4 days	72.9 days
	Sedentary work	Physical work
Return to athletics	65.6 days	117.9 days
	Overall	Overhead athletics
Duration of opioid use	18.3 days	





# HEALICOIL<sup>◇</sup> REGENESORB<sup>◇</sup> Suture Anchor mostly resorbed and replaced by new bone material within 21 months of arthroscopic RCR

Vonhoegen J, et al. *J Orthop Surg Res* (2019)<sup>13</sup>

## Overview

- Retrospective study evaluating resorption of HEALICOIL REGENESORB Suture Anchors in 48 patients (82 anchors) at a mean follow-up of 21 months after RCR

## Results

- At 21 months, 65/82 anchors (79%) could not be distinguished from adjacent bone material on MRI (Figure)
- Osteolysis was detected in only 2/82 anchors (2.4%); no reaction exceeded the diameter of the former suture anchor (5.5mm) and no peri-anchor cysts were formed
- Complete healing was achieved in 46/48 (96%) patients and no anchor pull-out complications were detected

## Conclusions

HEALICOIL REGENESORB anchors provide strong primary stability, reliable degradation and maintains bone quality of the rotator cuff footprint.




Figure. HEALICOIL REGENESORB anchor resorption at 21 months postoperatively

Study summary



Arthrex  
Biocomposite  
comparison



# Significantly greater bone density surrounding HEALICOIL<sup>◇</sup> PK Suture Anchors versus TWINFIX<sup>◇</sup> PK Suture Anchors six months after RCR

Chahla J, et al. *Arthroscopy* (2020)<sup>9</sup>

## Overview

- Single-centre randomised controlled trial comparing outcomes after double-row repair of FT rotator cuff tears, with the medial row fixated with either:
  - Coil-type HEALICOIL PK anchors (n=21)
  - Screw-type TWINFIX PK anchors (n=19)

## Results

- HEALICOIL PK anchors had significantly greater bone density at and up to 1.50mm from the anchor surface compared to TWINFIX PK anchors at 6 months ( $p<0.05$ ; Figure)
- HEALICOIL PK anchors had significantly more total bone mass within the anchor site compared to TWINFIX PK anchors at 6 months ( $p<0.01$ ); there was no significant difference in density
- Pain and shoulder function improved significantly from baseline to 12 months in both groups ( $p<0.05$ )

## Conclusions

Coil-type HEALICOIL PK anchors had significantly greater bone density surrounding the anchor compared to screw-type TWINFIX PK anchors 6 months after rotator cuff repair. The authors suggested that increased bone density around the anchor may contribute to a stronger construct during postoperative rehabilitation.




Figure. Bone density at and up to 2.50mm away from the surface of HEALICOIL PK anchors and TWINFIX PK anchors on 6-month CT scans

\*  $p<0.01$ , †  $p<0.05$





# Rotator cuff thickness significantly greater in patients who received vented compared with non-vented suture anchors at six weeks

Clark TR, et al. *Am Sport Med Res* (2016)<sup>8</sup>

## Overview


- Retrospective study comparing tendon healing following primary RCR in patients who received either vented (HEALICOIL<sup>®</sup> REGENESORB<sup>®</sup> Suture Anchor; n=40) or non-vented suture anchors (Healix Advance<sup>™</sup> anchor, DePuy Synthes, Raynham, MA, USA; n=30)

## Results

- Mean rotator cuff thickness was significantly greater in patients who received HEALICOIL REGENESORB anchors versus Healix Advance anchors, on 6-week ultrasound (Figure)
- In addition to anchor type, mean rotator cuff thickness was also significantly related to gender (p=0.022), age (p<0.001) and days since surgery (p=0.004)

## Conclusions

Patients treated with vented HEALICOIL REGENESORB anchors had a significantly greater rotator cuff thickness at six weeks versus those treated with non-vented anchors.





# Significantly improved bone ingrowth with HEALICOIL<sup>◇</sup> PK Suture Anchors versus TWINFIX<sup>◇</sup> HA Suture Anchors six months after RCR

Kim J-H, et al. *Arthroscopy* (2020)<sup>10</sup>

## Overview

- Randomised controlled trial in two South Korean hospitals comparing bone ingrowth and clinical outcomes six months after RCR with open-architecture HEALICOIL PK anchors (n=36) or non-vented TWINFIX HA anchors (n=33)

## Results

- Significantly more HEALICOIL PK anchors had good or excellent bone ingrowth compared to TWINFIX HA anchors (69.5 vs 36.3%;  $p<0.001$ ; Figure)
- No significant differences between HEALICOIL PK anchor and TWINFIX HA anchor groups in rates of cyst formation (14 vs 12%) and re-tear (5 vs 5%)
- Significant improvements in shoulder function and pain relief from baseline with both anchors ( $p<0.001$ ); no significant differences between groups

## Conclusions

Open-architecture HEALICOIL PK anchors had significantly improved bone ingrowth compared to non-vented TWINFIX HA anchors 6 months after rotator cuff repair.




Figure. Bone ingrowth (Modified Barber's ossification scale) with HEALICOIL PK anchors and TWINFIX HA anchors on 6-month CT scans





# Fixation properties, stress distribution and failure patterns differ between coil-type and screw-type suture anchors for RCR

Sano H, et al. *J Orthop Sci* (2016)<sup>11</sup>

## Overview

- Virtual pull-out testing of screw-type (TWINFIX<sup>®</sup> Ti Suture Anchor) and coil-type (HEALICOIL<sup>®</sup> PK Suture Anchor and HEALICOIL REGENESORB<sup>®</sup> Suture Anchor) suture anchors using 3D-FEM
- Computer models of each anchor were inserted into a model of cancellous bone prior to simulation of a traction force

## Results

- Site of highest stress distribution and element failure differed by anchor type:
  - Around the proximal threads of the screw-type TWINFIX Ti anchor
  - Near the distal tip and site of suture thread attachment in both coil-type HEALICOIL anchors
- HEALICOIL REGENESORB anchors showed the least displacement of the three anchors (Figure)

## Conclusions

Failure of the screw-type TWINFIX Ti anchor occurred closer to the bone surface than with the coil-type HEALICOIL anchors in virtual pull-out testing. As proximal bony tissue is often damaged during repair, this may indicate a greater risk of pull-out with screw-type anchors.




Figure. Mean value of displacement for coil-type and screw-type suture anchors



# References

---

1. Bokor DJ, Sonnabend D, Deady L et al. Evidence of healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a collagen implant: a 2-year MRI follow-up. *Muscles, Ligaments Tendons J.* 2016;6(1):16-25.
2. Schlegel TF, Abrams JS, Bushnell BD, Brock JL, Ho CP. Radiologic and clinical evaluation of a bioabsorbable collagen implant to treat partial-thickness tears: a prospective multicenter study. *J Shoulder Elbow Surg.* 2018;27(2):242-251.
3. Arnoczky SP, Bishai SK, Schofield B, et al. Histologic evaluation of biopsy specimens obtained after rotator cuff repair augmented with a highly porous collagen implant. *Arthroscopy.* 2017;33(2):278-283.
4. Bokor DJ, Sonnabend DH, Deady L, et al. Healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a highly-porous collagen implant: a 5-year clinical and MRI follow-up. *Muscles Ligaments Tendons J.* 2019;9(3):338-347.
5. Thon SG, O'Malley L, O'Brien MJ, Savoie FH. Evaluation of healing rates and safety with a bioinductive collagen patch for large and massive rotator cuff tears. *Am J Sports Med.* 2019;47(8):1901-1908.
6. Bokor DJ, Sonnabend D, Deady L, et al. Preliminary investigation of a biological augmentation of rotator cuff repairs using a collagen implant: a 2-year MRI follow-up. *Muscles Ligaments Tendons J.* 2015;5(3):144-150.
7. McIntyre LF, Bishai SK, Brown PB, Bushnell BD, Trenhaile SW. Patient-reported outcomes after use of a bioabsorbable collagen implant to treat partial and full-thickness rotator cuff tears. *Arthroscopy.* 2019;35(8):2262-2271.
8. Clark TR, Guerrero EM, Song A, O'Brien MJ and Savoie FH. Do vented suture anchors make a difference in rotator cuff healing. *Ann Sport Med Res.* 2016;3(3):1068.
9. Chahla J, Liu JN, Manderle B, et al. Bony ingrowth of coil-type open-architecture anchors compared with screw-type PEEK anchors for the medial row in rotator cuff repair: A randomized controlled trial. *Arthroscopy.* 2020;36(4):952-961.
10. Kim J-H, Kim Y-S, Park I, et al. A comparison of open-construct PEEK suture anchor and non-vented biocomposite suture anchor in arthroscopic rotator cuff repair: A prospective randomized clinical trial. *Arthroscopy.* 2020;36(2):389-396.
11. Sano H, Tokunaga M, Noguchi M, et al. Comparison of fixation properties between coil-type and screw-type anchors for rotator cuff repair: A virtual pullout testing using 3-dimensional finite element method. *J Orthop Sci.* 2016;21(4):452-457.
12. Data on file at Smith & Nephew, report 15000897, 2012.
13. Vonhoegen J, John D, Hagermann C. Osteoconductive resorption characteristics of a novel biocomposite suture anchor material in rotator cuff repair. *J Orthop Surg Res.* 2019;14(1):12.

## REGENETEN® Bioinductive Implant is rapidly absorbed and replaced by tendon-like tissue within 6 months

Cell ingrowth by 5 weeks, with progressive maturation to tendon-like tissue



### Study overview

- A retrospective study of biopsies taken between 5 weeks and 6 months following arthroscopic rotator cuff repair augmented with the REGENETEN Bioinductive Implant
- Biopsies were collected from 7 patients (6 full-thickness tears and 1 partial-thickness) requiring a second procedure
- Specimens were examined for host-tissue ingrowth, host-tissue maturation and host-implant biocompatibility



### Key results

- At the earliest time period (5 weeks), the biopsy showed rapid host cell ingrowth and early collagen formation
- At 3 months, there was increased collagen formation, maturation and organization on the surface of the implant (Figure 1)
- By 6 months, the implant was no longer visible, with new tendon-like tissue and oriented collagen indicative of functional loading (Figure 2)
- No evidence of foreign body or inflammatory reactions at any time point




Figure 1. Photomicrograph of surface of implant at 3 months. Remnants of implant still present (→)




Figure 2. Photomicrograph of surface of implant at 6 months. No evidence of implant remaining

Images included with permission from Dr Craig L. Van Kampen.



### Conclusion

The REGENETEN Bioinductive implant is the first to clinically demonstrate regeneration of tendon tissue. Biopsies demonstrated rapid host cell ingrowth and collagen formation, leading to progressive maturation and functional loading of new tissue.



### Considerations

- Due to the relatively small sample size and the associated inter-patient variability to healing and graft uptake, the progress of tendon-like tissue has not been directly measured



### Study citation

\*Arnoczky SP, Bishai SK, Schofield B, et al. Histologic evaluation of biopsy specimens obtained after rotator cuff repair augmented with a highly porous collagen implant. *Arthroscopy* 2017;33(2):278-283.  
Available at: [Arthroscopy: The Journal of Arthroscopy & Related Surgery](#)

## REGENETEN<sup>®</sup> Bioinductive Implant maintains repair integrity in full-thickness (FT) rotator cuff tears

Repair integrity supported by rapid induction of new tendon-like tissue which matures and integrates with native tendon



### Study overview

- A preliminary analysis of a prospective study of 9 patients (mean age, 56.4 years) with tears of the supraspinatus tendon (8 medium-sized FT tears; 1 high-grade partial-thickness tear converted to a FT tear during surgery)
- All patients received a REGENETEN Bioinductive Implant over the bursal surface of the tendon following standard repair
- MRI and clinical outcome assessments were conducted preoperatively and at 3, 6, 12 and 24 months postoperatively
- Tendon thickness measurements were compared to published values from young healthy adults to determine the relative amount of tissue generation



### Key results

- No MRI evidence of re-tear or gap formation, with the integrity of all repaired tendons intact at 24 months
- Significant increase in mean tendon thickness versus published values at 3, 6, 12 and 24 months ( $p < 0.01$ ), with an average of 2mm new tissue over the bursal surface (Figure)
- New tissue rapidly matured, improved in quality and was indistinguishable from the native tendon by 12 months
- From 12 to 24 months, tendon thickness slightly decreased, likely reflecting continued functional remodeling
- Significant improvement in clinical scores at 24 months versus preoperative measures
  - Constant-Murley score and Constant-Murley pain score (both  $p < 0.001$ )
  - American Shoulder and Elbow Surgeons (ASES) score and ASES pain score (both  $p < 0.001$ )
- Outcomes were satisfactory for 8/9 patients (89%) at 24 months




Figure. Mean ( $\pm$  standard error of mean) tendon thickness in male ( $n=6$ ) and female ( $n=3$ ) patients compared to published normal average thicknesses (horizontal lines) for healthy males and females;  $p < 0.01$  at all time points



### Conclusion

Through the generation of rapidly maturing new tendon-like tissue, REGENETEN Bioinductive Implant facilitates restoration of the normal tendon footprint and ultimately maintains repair integrity of full-thickness tears over 24 months. These findings are consistent with previous pre-clinical research and a finite element analysis.



### Study citation

\*Bokor DJ, Sonnabend D, Deady L, et al. Preliminary investigation of a biological augmentation of rotator cuff repairs using a collagen implant: a 2-year MRI follow-up. *Muscles Ligaments Tendons J*. 2015;5(3):144-150.  
Available at: [Muscle, Ligaments and Tendons Journal](#)

## REGENETEN<sup>®</sup> Bioinductive Implant induces tendon-like tissue formation in patients with partial-thickness (PT) tears of the supraspinatus (SS) tendon

Significant and clinically meaningful benefits were demonstrated in validated assessments of pain and function



### Study overview

- A prospective, multi-center, open-label trial in 33 patients (mean age, 54.6 years) with PT tears of the SS tendon
- All patients received REGENETEN Bioinductive Implant over the bursal surface of the tendon following arthroscopic subacromial decompression without repair
- Tendon thickness and tear size were assessed by MRI preoperatively and at 3 months and 1 year following surgery
- Clinical outcomes were measured using American Shoulder and Elbow Surgeons (ASES) and Constant-Murley assessments, preoperatively and at 3 months and 1 year postoperatively



### Key results

#### MRI outcomes

- At 1 year, 23 patients (70%) had a reduction in tear size of at least one grade from baseline
- Additionally, 8 patients (24%) had no visible defect at 1 year
- Tear progression only occurred in one patient, who did not follow the rehabilitation protocol
- No patients underwent revision surgery
- Significant increase in mean tendon thickness in both intermediate and high grade tears at 1 year ( $p < 0.01$ ; Figure 1)
- No significant differences in tendon thickness between:
  - Intermediate and high-grade tears
  - Articular surface and bursal-sided defects




Figure 1. Mean change in tendon thickness ( $\pm$  standard error) across intermediate and high grade tears. At month 12,  $p = 0.003$  and  $p < 0.0001$  for intermediate and high-grade tears respectively, versus preoperative measures



## Key results (continued)

### Patient outcomes

- Significant improvements in ASES shoulder index, pain and shoulder function scores at 1 year (all  $p < 0.0001$ )
  - Improvements in ASES pain and ASES shoulder index scores were approximately twice the minimal clinically important differences (MCIDs; Figures 2 and 3)
- Significant improvement in Constant-Murley shoulder score from 57.1 at baseline to 81.4 at 1 year ( $p < 0.0001$ ), greater than twice the MCID of 10.4
- At 1 year, 30 patients (94%) agreed or strongly agreed that they were satisfied with the results of their procedure
- Recovery was considered rapid by the investigators when compared with patients undergoing tear conversion and repair:
  - Mean sling time:  $23.3 \pm 2.4$  days
  - Mean return to work:  $30.5 \pm 12.0$  days
  - Mean duration of physical therapy:  $18 \pm 1.6$  visits




Figure 2. Mean ASES pain score before and after treatment with the REGENETEN Bioinductive Implant; MCID: 1.4




Figure 3. Mean ASES shoulder index score before and after treatment with the REGENETEN Bioinductive Implant; MCID: 12.01-16.92



## Conclusion

The REGENETEN Bioinductive Implant biologically augments healing, increasing tendon thickness and creating an environment conducive to healing. Therefore, REGENETEN Bioinductive Implant represents a promising treatment for patients with intermediate and high-grade PT tears of the SS tendon.



## Study citation

\*Schlegel TF, Abrams JS, Bushnell BD, Logan Brock J, Ho CP. Radiologic and clinical evaluation of a bioabsorbable collagen implant to treat partial thickness tears: a prospective multicenter study. *J Shoulder Elbow Surg* 2017;27(2):242-251.  
Available at: [Journal of Shoulder and Elbow Surgery](#)

## REGENETEN<sup>®</sup> Bioinductive Implant leads to tissue induction and high rates of tendon healing in patients with large and massive rotator cuff tears

Consistent treatment success in primary repairs and revisions with no implant-related adverse events



### Study overview

- A prospective study of 23 patients (mean age, 57.9 years) with large (two tendon, n=11) and massive (three tendon, n=12) rotator cuff tears receiving primary (n=7) or revision (n=16) repairs
- Following a double-row repair, a REGENETEN Bioinductive Implant was applied over the repaired supraspinatus and infraspinatus tendons
- Primary outcome was safety. Secondary outcomes included tendon thickness as an assessment of tissue induction on each ultrasound (US) examination (3, 6, 12 and 24 months) and on a single MRI (mean follow-up, 13 months), and American Shoulder and Elbow Surgeons (ASES) score at 24 months
- Standard postoperative rehabilitation protocol for large/massive rotator cuff tears was followed



### Key results

- No implant-related adverse events were reported
- Complete tendon healing in 22/23 patients (96%) on both imaging modalities (US and MRI)
- Treatment success in 21/23 patients (91%) at 24 months; 1 healing failure and 1 clinical failure due to progression of glenohumeral osteoarthritis
- Mean tendon thickness increased from 6.29mm at 3 months to 7.72mm at 12 months, decreasing to 7.28mm at 24 months
- Mean ASES score was 82.87 at 24 months
- No significant difference in treatment success (Figure), tendon thickness or ASES score between primary and revision repair groups or between large and massive tear groups




Figure. Percentage of patients achieving treatment success in primary and revision repairs



### Conclusion

Tendon healing is often unsuccessful following repair of large and massive rotator cuff tears, especially in the revision setting. In conjunction with repair of large and massive tears, REGENETEN Bioinductive Implant was safe, induced tissue formation and resulted in a relatively high tendon healing rate in both primary and revision settings. These findings are consistent with the healing response seen in partial-thickness tears.



### Considerations

- Tendon thickness was measured at the lateral edge of the articular cartilage and slightly posterior to the bicipital groove; the authors suggested a usual thickness in normal rotator cuffs of 8mm with this measurement technique



### Study citation

Thon SG, O'Malley L, O'Brien MJ, Savoie FH. Evaluation of healing rates and safety with a bioinductive collagen patch for large and massive rotator cuff tears: 2-year safety and clinical outcomes. *Am J Sports Med*. 2019;47(8):1901-1908.

## REGENETEN® Bioinductive Implant promotes rapid and sustained healing of partial-thickness (PT) rotator cuff tears

Tear healing is linked to significant improvement in function and reduction in pain compared to preoperative values over 24 months



### Study overview

- A prospective study of 13 patients (mean age, 53.8 years) with various grades and locations of PT tears of the supraspinatus tendon
- All patients received a REGENETEN Bioinductive Implant over the bursal surface of the tendon following arthroscopic subacromial decompression without repair
- MRI and clinical outcome assessments were conducted preoperatively and at 3, 6, 12 and 24 months postoperatively



### Key results

- Significant mean increase in tendon thickness of 2.2mm at 3 months versus preoperative values ( $p < 0.0001$ )
- At 12 months, new tissue was indistinguishable from underlying tissue in 12/13 patients (92%)
- Tendon thickness at 24 months was significantly thicker ( $p < 0.0001$ ) than preoperative values
- At 12 months, all assessable patients had a reduction in tear size of  $\geq 1$  grade, with complete tear disappearance in 7 of 10 patients with measurable tear size (70%) (Figure)
- Significant improvement in clinical scores throughout 24-month follow-up period
  - Constant-Murley score ( $p \leq 0.01$ ) and Constant-Murley pain score ( $p \leq 0.001$ )
  - American Shoulder and Elbow Surgeons (ASES) total score and ASES pain score (both  $p \leq 0.001$ )
- Outcomes were satisfactory for 12/13 patients (92%) at 24 months, suggesting a benefit over acromioplasty alone




Figure. Sequential improvement and healing of cuff defects over 24 months (n=10)



### Conclusion

When treated with REGENETEN Bioinductive Implant, partial-thickness rotator cuff tears can decrease in size and in most cases completely heal. Tear healing is associated with the formation of load-bearing, tendon-like tissue, ultimately leading to improved clinical outcomes.



### Study citation

\*Bokor DJ, Sonnabend D, Deady L, et al. Evidence of healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a collagen implant: a 2-year MRI follow-up. *Muscles Ligaments Tendons J*. 2016;6(1):16-25.

Available at: [Muscle, Ligaments and Tendons Journal](#)

**+ Evidence in focus**Publication summary: Bokor DJ, et al. *MLTJ*. (2019)\***Improved tendon thickness and integrity was sustained five years after treatment with the REGENETEN<sup>®</sup> Bioinductive Implant****+ Plus points**

Comparison of 5-year findings with 2-year follow-up:

**73%**


Most patients had no decline in tendon integrity

**Significant improvements**  
in pain and function  
sustained to 5 years**No additional complications****Overview**

- Five-year follow-up of a prospective, single-arm study evaluating REGENETEN Bioinductive Implant in lieu of repair in patients with partial-thickness rotator cuff tears
- 11/13 enrolled patients (mean age, 54.0 years) were available for assessment
- Outcomes were assessed 5 years postoperatively and were compared to findings at 2 years (previously reported):
  - MRI assessment of tendon integrity and thickness
  - Patient-reported outcomes (American Shoulder and Elbow Surgeons [ASES] shoulder scale and Constant-Murley shoulder score)

**Results**

- 8/11 patients (73%) had no decline in tendon integrity between 2 and 5 years (Figure)
- Tendon thickness significantly decreased between 2 and 5 years (5.9 vs 5.2mm;  $p=0.0012$ ), but remained significantly greater than preoperative values (4.2mm;  $p<0.0001$ ; Figure)
- Significant improvements from baseline in pain and function were sustained to 5 years ( $p\leq 0.01$ ) and were not significantly different to 2-year values
- No additional complications between 2 and 5 years

**11 patients available for 5-year follow-up****n=8**  
No decline in  
tendon integrity**n=2**  
New low grade  
tears at original site**n=1**  
New low grade  
tear distinct from  
original site*Figure. Change in tendon integrity (MRI) between 2- and 5-year follow-up***Conclusions**

REGENETEN Bioinductive Implant demonstrated sustained effectiveness in improving tendon thickness and integrity, with no reduction in patient-reported outcomes between 2 and 5 years.

**Citation**

\*Bokor DJ, Sonnabend DH, Deady L, et al. Healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a highly-porous collagen implant: a 5-year clinical and MRI follow-up. *MLTJ*. 2019;9(3):338-347.

Available at: [Muscles, Ligaments and Tendons Journal](https://www.smith-nephew.com/education)

## Clinically meaningful improvements in pain and function at one year in a multicentre registry study of REGENETEN® Bioinductive Implant for the treatment of partial- and full-thickness rotator cuff tears

Time to return to driving and sport compared favourably to published results from patients treated with standard surgical techniques



### Study overview

- A multicentre registry study of 173 patients who received a REGENETEN Bioinductive Implant in lieu of partial-thickness repair (n=90) or in conjunction with standard repair of full-thickness tears (n=83)
- Patients available at follow-up had a mean age of 54.2 years, were 57% male and included diabetics (10.4%), smokers (14.4%), patients involved in a worker's compensation claim (12.1%) and chronic opioid users (7.5%)
- Outcomes included postoperative recovery and patient-reported outcome measures of pain and function (ASES, SANE, VAS pain, VR-12 and WORC scores) at baseline and at regular intervals to 1 year of study follow-up



### Key results

#### Partial-thickness tears

- MCIDs were achieved in VAS pain from 2 weeks (5.3 at baseline to 3.3 at 2 weeks and 1.1 at 1 year;  $p<0.001$ ) and ASES score from 6 weeks (47.0 at baseline to 60.6 at 6 weeks and 85.6 at 1 year;  $p<0.001$ )
  - Mean VAS pain and ASES score at 3 months were significantly improved compared to published results from patients undergoing transtendon or takedown repair without a REGENETEN Bioinductive Implant ( $p<0.05$ )
- Significant improvements from 6 weeks in SANE score (42.5 at baseline to 59.4 at 6 weeks and 86.0 at 1 year;  $p<0.001$ ), VR-12 PCS (35.8 at baseline to 39.1 at 6 weeks and 49.7 at 1 year;  $p<0.002$ ) and WORC index (38.2 at baseline to 53.5 at 6 weeks and 84.4 at 1 year;  $p<0.001$ )
- Time in sling and time to return to driving and sport was rapid compared to published results for standard surgical techniques (Table 1)






Measure		REGENETEN Bioinductive Implant	
	Time in sling	10.6 days No biceps surgery	27.7 days Concomitant tenodesis
		14.6 days	
	Return to driving	14.6 days	
		14.6 days	
	Return to work	9.4 days Sedentary work	72.9 days Physical work
		18.3 days	
	Return to athletics	65.6 days Overall	117.9 days Overhead athletics
		18.3 days	
	Duration of opioid use	18.3 days	

Table 1. Duration of postoperative recovery in patients with partial-thickness tears



## Key results (continued)

### Full-thickness tears

- MCIDs were achieved in VAS pain from 2 weeks (5.2 at baseline to 3.7 at 2 weeks and 1.2 at 1 year;  $p < 0.001$ ) and ASES score from 3 months (45.5 at baseline to 68.4 at 3 months and 83.8 at 1 year;  $p < 0.001$ )
- Significant improvements from 6 weeks in WORC index (35.0 at baseline to 41.1 at 6 weeks and 80.1 at 1 year;  $p < 0.001$ ), and from 3 months in SANE score (39.2 at baseline to 63.3 at 3 months and 80.7 at 1 year;  $p < 0.001$ ) and VR-12 PCS (34.5 at baseline to 40.8 at 3 months and 45.7 at 1 year;  $p < 0.001$ )
- Time to return to driving and sport was rapid compared to published results for standard surgical techniques (Table 2)






Measure		REGENETEN® Bioinductive Implant	
	Time in sling	34.0 days No biceps surgery	39.4 days Concomitant tenodesis
	Return to driving	24.5 days	
	Return to work	21.8 days Sedentary work	62.5 days Physical work
	Return to athletics	119.2 days Overall	143.7 days Overhead athletics
	Duration of opioid use	26.9 days	

Table 2. Duration of postoperative recovery in patients with full-thickness tears

### Complications

- Complications requiring revision surgery were failure of cuff healing ( $n=4$ ), postoperative infection ( $n=1$ ), deep vein thrombosis and adhesive capsulitis ( $n=1$ ), postoperative stiffness ( $n=1$ ) and recurrent effusions ( $n=1$ )



## Conclusion

**REGENETEN Bioinductive Implant led to clinically meaningful improvements in pain and function for patients with partial- and full-thickness tears. Return to driving and sport compared favourably to standard surgical techniques, suggesting that REGENETEN Bioinductive Implant may also enhance postoperative recovery.**




## Considerations

- Only comparisons to MCIDs for ASES score and VAS pain were made by the authors. However, the improvements from preoperative values in SANE score, VR-12 PCS and WORC index also exceeded the MCIDs reported in existing literature for rotator cuff repair<sup>1-3</sup>



## Study citation

\*McIntyre LF, Bishai SK, Brown PB, Bushnell BD, Trenhaile SW. Patient-reported outcomes following use of a bioabsorbable collagen implant to treat partial and full-thickness rotator cuff tears. *Arthroscopy*. 2019 July 23. [Epub ahead of print]  
Available at: [Arthroscopy](#) 

## References

1. Kirkley A, Griffin S, Dainty K. Scoring systems for the functional assessment of the shoulder. *Arthroscopy*. 2003; 19(10):1109-1120.
2. Cvetanovich GL, Gowd AK, Liu JN, et al. Establishing clinically significant outcome after arthroscopic rotator cuff repair. *J Shoulder Elbow Surg*. 2019; 28:939-948.
3. Zhou L, Natarajan M, Miller BS, Gagnier JJ. Establishing minimal important differences for the VR-12 and SANE scores in patients following treatment of rotator cuff tears. *Orthop J Sports Med*. 2018;6(7):2325967118782159.

## Abbreviations

ASES = American Shoulder and Elbow Surgeons, MCID = minimal clinically important difference, NR = not reported, SANE = Single Assessment Numeric Evaluation, VAS = visual analogue scale, VR-12 = Veterans RAND 12-item Health Survey, VR-12 PCS = Veterans RAND 12-item Health Survey Physical Component Score, WORC = Western Ontario Rotator Cuff

# HEALICOIL<sup>®</sup> REGENESORB biocomposite suture anchor mostly resorbed and replaced by new bone material within 21 months of arthroscopic rotator cuff repair (RCR)

No severe osteolysis or cyst formation observed at any anchor site



## Study overview

- Retrospective, single-centre study assessing the resorption and osteoconductive properties of a novel biocomposite material REGENESORB, comprising 65% polylactic-co-glycolic acid (PLGA), 15% beta-tricalcium phosphate ( $\beta$ -TCP) and 20% calcium sulfate (CS)
  - 48 patients underwent arthroscopic single-row RCR with 5.5mm HEALICOIL REGENESORB Suture Anchor (82 suture anchors, average 1.71 anchors per patient)
- Outcomes included MRI evaluation of implant resorption, osteolysis and re-tear rate at a mean follow-up of 21.2 months



## Key results

- At 21 months, 79% of implants (75% patients) could not be distinguished from adjacent bone material (Figure)
  - No significant correlation between anchor resorption and age, re-tear rate, defect size, gender, number of anchors, and grade of retraction
- Osteolysis was detected in only 2/82 anchors (2.4%), with no reaction exceeding the diameter of the former suture anchor (5.5mm) and no peri-anchor cyst formation
  - No significant correlation between osteolysis and patient age, gender, re-tear rate, or size of the defect
- Complete healing was achieved in 46/48 (96%) patients and no anchor pull-out complications were detected




Figure. Level of resorption of 82 implanted suture anchors



## Conclusion

HEALICOIL REGENESORB biocomposite suture anchor provides strong primary stability, reliable degradation and maintains bone quality of the rotator cuff footprint. Preserving bone quality aids the clinical situation when revision surgery is required. Resorption characteristics and osteolysis occurrence appeared superior compared to existing evidence of commonly used anchor materials containing PLLA (poly-L-lactide) and PDLDA (poly-D-L-lactide).



## Study citation

\*Vonhoegen J, John D, Hägermann C. Osteoconductive resorption characteristics of a novel biocomposite suture anchor material in rotator cuff repair. *J Orthop Surg Res*. 2019;14(1):12.

Available at: [Journal of Orthopaedic Surgery and Research](#)

## Focus on biocomposite suture anchor resorption

## Introduction

Biocomposite suture anchors composed of biodegradable polymers and osteoconductive materials are commonly used in rotator cuff repair.<sup>1</sup> Osteoconductive materials were introduced in response to the poor bone replacement and bone-derived complications associated with polymer-only anchors.<sup>1-5</sup> However, despite their widespread use, the resorption process of biocomposite suture anchors for rotator cuff repair has only recently been studied.<sup>1,2</sup>

HEALICOIL<sup>®</sup> REGENESORB<sup>®</sup> Suture Anchors<sup>1</sup>

## Overview

Retrospective study evaluating resorption of 18.5mm HEALICOIL REGENESORB Suture Anchors in 48 patients (82 anchors).

Anchor resorption was evaluated by MRI 21 months after rotator cuff repair.

## Anchor material composition

65%

PLGA

15%


 $\beta$ -TCP

20%

Calcium sulphate

## Results

- Resorption at 21 months (Figure 1):
- Complete or near complete resorption of 65/82 suture anchors (79%)
- Only 17/82 suture anchors (21%) were visible or clearly visible
- No severe osteolysis

Arthrex BioComposite<sup>™</sup> Suture Anchors<sup>2</sup>

## Overview

Retrospective study evaluating resorption of 14.7mm Biocomposite Corkscrew<sup>®</sup> and 19.1mm Biocomposite SwiveLock<sup>®</sup> anchors (Arthrex, Naples, FL, USA) in 25 patients (84 anchors).

Anchor resorption was evaluated by MRI 28 months after rotator cuff repair.

## Anchor material composition

85%


PLLA

15%

 $\beta$ -TCP

## Results

- Resorption at 28 months (Figure 2):
- Complete or near complete resorption of only 28/84 suture anchors (33%)
- Two-thirds of suture anchors (56/84, 67%) remained visible or clearly visible
- No severe osteolysis



## Commentary

- The material composition and physical structure of biocomposite suture anchors affects their resorption<sup>1,2</sup>
- In separate clinical studies, resorption of more than twice as many HEALICOIL REGENESORB Suture Anchors was complete or near complete at 21 months compared to 85% PLLA and 15%  $\beta$ -TCP suture anchors at 28 months<sup>1,2</sup>
- Successful anchor resorption and bone replacement offers benefits for patients requiring revision surgery<sup>1,2</sup>

## Abbreviations

$\beta$ -TCP =  $\beta$ -tricalcium phosphate, MRI = magnetic resonance imaging, PLGA = poly-L-lactic co-glycolic acid, PLLA = poly L-lactic acid

## References

- Vonhoegen J, John D, Hägermann C. *J Orthop Surg Res.* 2019;14(1):12.
- Sgroi M, Friesz T, Schocke M, Reichel H, Kappe T. *Clin Orthop Relat Res.* 2019;477(6):1469-1478.
- Dhawan A, Ghodadra N, Karas V, et al. *Am J Sports Med.* 2012;40(6):1424-1430.
- Duralde XA. *Clin Orthop Relat Res.* 2019;477:1479-1482.
- Milewski MD, Diduch DR, Hart JM, et al. *Am J Sports Med.* 2012;40(6):1392-1401.

## + Evidence in focus

Publication summary: Chahla J, et al. *Arthroscopy* (2020)\*

## Significantly greater bone density surrounding coil-type HEALICOIL<sup>®</sup> PK Suture Anchors versus screw-type TWINFIX<sup>®</sup> PK Suture Anchors 6 months after rotator cuff repair

### + Plus points



**Significantly greater bone density** at and around the surface of HEALICOIL PK anchors versus TWINFIX PK anchors ( $p < 0.05$ )



**Significantly more total bone mass** within the anchor sites of HEALICOIL PK anchors versus TWINFIX PK anchors ( $p < 0.01$ )

### Overview

- Single-centre randomised controlled trial comparing outcomes after double-row repair of full-thickness rotator cuff tears, with the medial row fixated with either:
  - Coil-type HEALICOIL PK anchors ( $n=21$ )
  - Screw-type TWINFIX PK anchors ( $n=19$ )
- Outcomes included:
  - Bone density at and around the anchor site at 6 months
  - Composition of synovial fluid-bone marrow aspirate collected from the anchor site during surgery
  - Pain and shoulder function at 6 and 12 months

### Results

- HEALICOIL PK anchors had significantly greater bone density at and up to 1.50mm from the anchor surface compared to TWINFIX PK anchors ( $p < 0.05$ ; Figure)
- HEALICOIL PK anchors had significantly more total bone mass within the anchor site compared to TWINFIX PK anchors ( $p < 0.01$ ); there was no significant difference in density
- Pain and shoulder function improved significantly from baseline to 12 months in both groups ( $p < 0.05$ )
- Growth factors and stem cells detected in synovial fluid-bone marrow aspirate during surgery were similar in both groups




Figure. Bone density at and up to 2.50mm away from the surface of HEALICOIL PK anchors and TWINFIX PK anchors  
\* $p < 0.01$ , † $p < 0.05$

### Conclusions

Coil-type HEALICOIL PK anchors had significantly greater bone density surrounding the anchor compared to screw-type TWINFIX PK anchors 6 months after rotator cuff repair. The authors suggested that increased bone density around the anchor may contribute to a stronger construct during postoperative rehabilitation.

### Citation

\*Chahla J, Liu JN, Manderle B, et al. Bony ingrowth of coil-type open-architecture anchors compared with screw-type PEEK anchors for the medial row in rotator cuff repair: a randomized controlled trial. *Arthroscopy*. 2020;36(4):952-961.

Available at: [Arthroscopy](https://www.arthroscopy.com)

## Mean rotator cuff thickness significantly greater in patients who received vented compared with non-vented suture anchors at 6 weeks

Vented suture anchors for rotator cuff repair may provide greater healing potential at 6 weeks postoperatively



### Study overview

- Retrospective study comparing tendon healing following primary rotator cuff repair in patients who received either vented suture anchors (HEALICOIL® REGENESORB suture anchor; n=40) or non-vented suture anchors (Healix Advance™ suture anchor, DePuy Synthes, Raynham, MA, USA; n=30)
- Patients in the vented group were younger than those in the non-vented group (mean age: 55.0±10.1 years vs 62.5±10.7 years, respectively), but there was no difference between groups in terms of days post-op, handedness or gender
- To assess healing, rotator cuff thickness was measured at the medial anchor site by ultrasound at 6 week follow-up
- Both patients and physicians conducting the ultrasound were blinded to the type of anchor used



### Key results

- Mean rotator cuff thickness was significantly greater in patients who received HEALICOIL REGENESORB suture anchor vs Healix Advance suture anchor (0.59cm vs 0.48cm;  $p=0.0074$ ; Figure)
- Mean rotator cuff thickness was significantly greater in:
  - Males vs females, 0.59±0.17cm and 0.49±0.17cm, respectively ( $p=0.022$ )
  - Younger vs older patients ( $p<0.001$ )
  - Patients with more days post-op vs those with fewer ( $p=0.004$ )




Figure. Mean ( $\pm$  standard deviation) rotator cuff thickness at 6 weeks



### Conclusion

This is the first study comparing tendon healing following rotator cuff repair with vented and non-vented suture anchors. Patients treated with vented suture anchors had a significant increase in rotator cuff thickness at six weeks versus those treated with non-vented anchors. These findings support the theory that vented suture anchors provide a biologic healing advantage.



### Considerations

- Two surgeons operated on all patients. Most patients in the non-vented group received surgery from a single surgeon, whereas the vented group had a smaller ratio between operating surgeons
- Results were only collected at 6-week follow-up, other studies have denoted rotator cuff and muscle healing can take up to six months to see full effect



### Study citation

\*Clark TR, Guerrero EM, Song A, O'Brien MJ and Savoie FH. Do vented suture anchors make a difference in rotator cuff healing. *Am Sport Med Res*. 2016;3(3):1068.

Available at: [Annals of Sports Medicine and Research](#) 

## + Evidence in focus

Publication summary: Kim J-H, et al. *Arthroscopy* (2020)\*

## Significantly improved bone ingrowth with HEALICOIL<sup>◇</sup> PK Suture Anchors versus TWINFIX<sup>◇</sup> HA Suture Anchors 6 months after rotator cuff repair

### + Plus points



Significantly more HEALICOIL PK anchors had good or excellent bone ingrowth compared to TWINFIX HA anchors ( $p < 0.001$ )



No significant difference in cyst formation between HEALICOIL PK anchor and TWINFIX HA anchor groups (14 vs 12%)

### Overview

- Randomised controlled trial in two South Korean hospitals comparing bone ingrowth and clinical outcomes after rotator cuff repair with open-architecture HEALICOIL PK anchors ( $n=36$ ) or non-vented TWINFIX HA anchors ( $n=33$ ) for medial row fixation
- Postoperative outcomes included:
  - CT evaluation of bone ingrowth and cyst formation at 6 months
  - MRI or ultrasound evaluation of re-tear at 12 months
  - Shoulder function and pain at 3, 6 and 12 months, and final follow-up (mean, 25.2 months)

### Results

- Significantly more HEALICOIL PK anchors had good or excellent bone ingrowth compared to TWINFIX HA anchors (69.5 vs 36.3%;  $p < 0.001$ ; Figure)
  - Excellent bone ingrowth was observed in 11 patients (30.6%) treated with HEALICOIL PK anchors compared to only 1 patient (3.0%) treated with TWINFIX HA anchors
- No significant difference in cyst formation between HEALICOIL PK anchor and TWINFIX HA anchor groups (14 vs 12%)
- No significant difference in re-tear rate between HEALICOIL PK anchor and TWINFIX HA anchor groups (5 vs 5%)
- Significant improvements in shoulder function and pain relief from baseline with both anchors ( $p < 0.001$ ); no significant differences between groups




Figure. Bone ingrowth (Modified Barber's ossification scale) with HEALICOIL PK anchors and TWINFIX HA anchors at 6 months

### Conclusions

Open-architecture HEALICOIL PK Suture Anchors had significantly improved bone ingrowth compared to non-vented TWINFIX HA Suture Anchors 6 months after rotator cuff repair. Rates of cyst formation were similar and the authors suggested that the open-architecture of HEALICOL PK anchors may allow early and sufficient bone ingrowth before cyst formation occurs.

### Citation

\* Kim J-H, Kim Y-S, Park I, et al. A comparison of open-construct PEEK suture anchor and non-vented biocomposite suture anchor in arthroscopic rotator cuff repair: a prospective randomized clinical trial. *Arthroscopy*. 2020;36(2):389-396.  
Available at: [Arthroscopy](https://www.arthroscopyjournal.org/)

## Fixation properties, stress distribution and failure patterns differ between coil-type and screw-type suture anchors for rotator cuff repair

HEALICOIL<sup>®</sup> REGENESORB suture anchor had the best initial fixation properties of anchors tested



### Study overview

- An independent study conducting virtual pull-out testing using 3-dimensional finite element method (3D-FEM)
- Computer models of three anchors; one screw-type anchor (TWINFIX<sup>®</sup> Ti suture anchor), and two coil-type anchors (HEALICOIL PK suture anchor and HEALICOIL REGENESORB suture anchor) were inserted into the isotropic cube model that simulated cancellous bone
- A tensile load (500 N) along the long axis of the inserted anchor was applied to the site of suture thread attachment to simulate a traction force



### Key results

- With TWINFIX Ti screw-type suture anchor, the highest stress and element failure occurred around the anchor threads, closest to the surface of the cube
- Conversely, the highest stress and element failure with both coil-type anchors occurred deeper, near the anchor tip and site of suture thread attachment
- HEALICOIL REGENESORB suture anchor showed the least displacement of any anchor tested, with less than 0.1mm displacement at a load of 500N, vs 0.1mm displacement at 400N and 370N for TWINFIX Ti suture anchor and HEALICOIL PK suture anchor, respectively (Figure)




Figure. Mean value of displacement



### Conclusion

In virtual pullout testing of the screw-type anchor, stress distribution and element failure occurred around the proximal threads, whereas in coil-type anchors, stress and element failure occurred nearer the distal anchor tip. As proximal bony tissue is often damaged during repair, this may lead to a greater risk of pull-out with screw-type anchors. HEALICOIL REGENESORB suture anchor had the best initial fixation properties of all anchors tested.



### Considerations

- 3D-FEM is a computer aided engineering tool, which has been validated in the prediction of femoral and vertebrae fractures and has been used to predict the failure risk of inserted implants



### Study citation

\*Sano H, Tokunaga M, Noguchi M, et al. Comparison of fixation properties between coil-type and screw-type anchors for rotator cuff repair: A virtual pullout testing using 3-dimensional finite element method. *J Orthop Sci*. 2016;21(4):452-457.