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ADVANCED HEALING SOLUTIONS

Biologisches Portfolio

- REGENETEN[◇] bioinduktives Implantat
- HEALICOIL[◇] REGENESORB[◇] Fadenanker
- HEALICOIL KNOTLESS knotenloser Fadenanker
- REGENESORB Material

Februar 2024

Smith+Nephew

Übersicht

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Histologic evaluation of biopsy specimens obtained after rotator cuff repair augmented with a highly porous collagen implant¹

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Evidence of healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a collagen implant: a 2-year MRI follow-up²

Bokor DJ, Sonnabend D, Deady L, et al. *Muscles Ligaments Tendons J.* 2016;6:16–25.

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Partial-Risse

Healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a highly-porous collagen implant: a 5-year clinical and MRI follow-up³

Bokor DJ, Sonnabend DH, Deady L, et al. *Muscles Ligaments Tendons J.* 2019;9:338–347.

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
Treatment of partial-thickness rotator cuff tears with a resorbable bioinductive bovine collagen implant: 1-year results from a prospective multicenter registry⁴

Bushnell BD, Bishai SK, Krupp RJ, et al. Orthop J Sports Med. 2021;9:23259671211027850.

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Bioinductive collagen implants facilitate tendon regeneration in rotator cuff tears⁵

Camacho-Chacon JA, Cuenca-Espierrez J, Roda-Rojo V, et al. *J Exp Orthop.* 2022;9:53.

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Augmentation of a transosseous equivalent repair in posterosuperior non-acute rotator cuff tears with a bioinductive collagen implant decreases the re-tear rate at one year

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Isolated bioinductive repair of partial-thickness rotator cuff tears using a resorbable bovine collagen implant: two-year radiologic and clinical outcomes from a prospective multicenter study⁷

Schlegel TF, Abrams JS, Angelo RL, Getelman MH, Ho CP, Bushnell BD. *J Shoulder Elbow Surg.* 2021;30:1938–1948.

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Evaluation of healing rates and safety with a bioinductive collagen patch for large and massive rotator cuff tears: 2-year safety and clinical outcomes⁸

Thon SG, O'Malley L 2nd, O'Brien MJ, Savoie FH 3rd. *Am J Sports Med.* 2019;47:1901–1908.

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Preliminary investigation of a biological augmentation of rotator cuff repairs using a collagen implant: a 2-year MRI follow-up⁹

Bokor DJ, Sonnabend D, Deady L, et al. *Muscles Ligaments Tendons J.* 2015;5:144–150.

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Two-year outcomes with a bioinductive collagen implant used in augmentation of arthroscopic repair of full-thickness rotator cuff tears: final results of a prospective multi-center study¹⁰

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Collagen-based bioinductive implant for treatment of partial thickness rotator cuff tears¹¹

Dai A, Campbell A, Bloom D, Baron S, Begly J, Meislin R. *Bull Hosp Jt Dis.* 2020;78:195–201.

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Patient-reported outcomes after use of a bioabsorbable collagen implant to treat partial and full-thickness rotator cuff tears¹²

McIntyre LF, Bishai SK, Brown PB 3rd, Bushnell BD, Trenhaile SW. *Arthroscopy*. 2019;35:2262–2271.

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Full-thickness rotator cuff tears can be safely treated with a resorbable bioinductive bovine collagen implant: one-year results of a prospective, multicenter registry¹³

McIntyre LF, McMillan S, Trenhaile SW, Bishai SK, Bushnell BD. *Arthrosc Sports Med Rehabil.* 2021;3:e1473–e1479.

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Bony ingrowth of coil-type open-architecture anchors compared with screw-type PEEK anchors for the medial row in rotator cuff repair: a randomized controlled trial¹⁴

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Do vented suture anchors make a difference in rotator cuff healing¹⁵

Clark TR, Guerrero EM, Song A, O'Brien MJ, Savoie FH. *Ann Sports Med Res.* 2016;3:1068.

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
Figure. Mean±SD rotator cuff thickness at 6 weeks post-operatively

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Biomechanical analysis of bioabsorbable suture anchors for rotator cuff repair using osteoporotic and normal bone models¹⁶

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A comparison of open-construct PEEK suture anchor and non-vented biocomposite suture anchor in arthroscopic rotator cuff repair: a prospective randomized clinical trial¹⁷

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Comparison of fixation properties between coil-type and screw-type anchors for rotator cuff repair: a virtual pullout testing using 3-dimensional finite element method¹⁸

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Osteointegration of a biocomposite suture anchor after arthroscopic shoulder labral repair¹⁹

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Osteoconductive resorption characteristics of a novel biocomposite suture anchor material in rotator cuff repair²⁰

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Literaturangaben

1. Arnoczky SP, Bishai SK, Schofield B, et al. Histologic evaluation of biopsy specimens obtained after rotator cuff repair augmented with a highly porous collagen implant. *Arthroscopy*. 2017;33:278–283.
2. Bokor DJ, Sonnabend D, Deady L, et al. Evidence of healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a collagen implant: a 2-year MRI follow-up. *Muscles Ligaments Tendons J*. 2016;6:16–25.
3. Bokor DJ, Sonnabend DH, Deady L, et al. Healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a highly-porous collagen implant: a 5-year clinical and MRI follow-up. *Muscles Ligaments Tendons J*. 2019;9:338–347.
4. Bushnell BD, Bishai SK, Krupp RJ, et al. Treatment of partial-thickness rotator cuff tears with a resorbable bioinductive bovine collagen implant: 1-year results from a prospective multicenter registry. *Orthop J Sports Med*. 2021;9:23259671211027850.
5. Camacho-Chacon JA, Cuenca-Espierrez J, Roda-Rojo V, et al. Bioinductive collagen implants facilitate tendon regeneration in rotator cuff tears. *J Exp Orthop*. 2022;9:53.
6. Ruiz Ibán MA, Navlet MG, Marco SM, et al. Augmentation of a transosseous equivalent repair in posterosuperior non-acute rotator cuff tears with a bioinductive collagen implant decreases the re-tear rate at one year. A randomised controlled trial. *Arthroscopy*. 2023 Dec 28. [Epub ahead of print]
7. Schlegel TF, Abrams JS, Angelo RL, Getelman MH, Ho CP, Bushnell BD. Isolated bioinductive repair of partial-thickness rotator cuff tears using a resorbable bovine collagen implant: two-year radiologic and clinical outcomes from a prospective multicenter study. *J Shoulder Elbow Surg*. 2021;30:1938–1948.
8. Thon SG, O'Malley L 2nd, O'Brien MJ, Savoie FH 3rd. Evaluation of healing rates and safety with a bioinductive collagen patch for large and massive rotator cuff tears: 2-year safety and clinical outcomes. *Am J Sports Med*. 2019;47:1901–1908.
9. Bokor DJ, Sonnabend D, Deady L, et al. Preliminary investigation of a biological augmentation of rotator cuff repairs using a collagen implant: a 2-year MRI follow-up. *Muscles Ligaments Tendons J*. 2015;5:144–150.
10. Bushnell BD, Connor P, Harris HW, Ho CP, Trenhaile SW, Abrams JS. Two-year outcomes with a bioinductive collagen implant used in augmentation of arthroscopic repair of full-thickness rotator cuff tears: final results of a prospective multi-center study. *J Shoulder Elbow Surg*. 2022 Jul 1. [Epub ahead of print]
11. Dai A, Campbell A, Bloom D, Baron S, Begly J, Meislin R. Collagen-based bioinductive implant for treatment of partial thickness rotator cuff tears. *Bull Hosp Jt Dis*. 2020;78:195–201.
12. McIntyre LF, Bishai SK, Brown PB 3rd, Bushnell BD, Trenhaile SW. Patient-reported outcomes after use of a bioabsorbable collagen implant to treat partial and full-thickness rotator cuff tears. *Arthroscopy*. 2019;35:2262–2271.
13. McIntyre LF, McMillan S, Trenhaile SW, Bishai SK, Bushnell BD. Full-thickness rotator cuff tears can be safely treated with a resorbable bioinductive bovine collagen implant: one-year results of a prospective, multicenter registry. *Arthrosc Sports Med Rehabil*. 2021;3:e1473–e1479.
14. Chahla J, Liu JN, Manderle B, et al. Bony ingrowth of coil-type open-architecture anchors compared with screw-type PEEK anchors for the medial row in rotator cuff repair: a randomized controlled trial. *Arthroscopy*. 2020;36:952–961.
15. Clark TR, Guerrero EM, Song A, O'Brien MJ, Savoie FH. Do vented suture anchors make a difference in rotator cuff healing. *Ann Sports Med Res*. 2016;3:1068.
16. Yamauchi S, Tsukada H, Sasaki E, et al. Biomechanical analysis of bioabsorbable suture anchors for rotator cuff repair using osteoporotic and normal bone models. *J Orthop Sci*. 2022;27:115–121.
17. Kim JH, Kim YS, Park I, et al. A comparison of open-construct PEEK suture anchor and non-vented biocomposite suture anchor in arthroscopic rotator cuff repair: a prospective randomized clinical trial. *Arthroscopy*. 2020;36:389–396.
18. Sano H, Tokunaga M, Noguchi M, et al. Comparison of fixation properties between coil-type and screw-type anchors for rotator cuff repair: a virtual pullout testing using 3-dimensional finite element method. *J Orthop Sci*. 2016;21:452–457.
19. Sugaya H, Suzuki K, Yoshimura H, et al. Osteointegration of a biocomposite suture anchor after arthroscopic shoulder labral repair. *Arthroscopy*. 2019;35:3173–3178.
20. Vonhoegen J, John D, Hägermann C. Osteoconductive resorption characteristics of a novel biocomposite suture anchor material in rotator cuff repair. *J Orthop Surg Res*. 2019;14:12.
21. Bushnell BD, Connor PM, Harris HW, Ho CP, Trenhaile SW, Abrams JS. Retear rates and clinical outcomes at 1 year after repair of full-thickness rotator cuff tears augmented with a bioinductive collagen implant: a prospective multicenter study. *JSES Int*. 2021;5:228–237.

Literaturangaben

22. Micheloni GM, Salmaso G, Zecchinato G, Giaretta S, Barison E, Momoli A. Bio-inductive implant for rotator cuff repair: our experience and technical notes. *Acta Biomed.* 2020;91:e2020004.
23. Schlegel TF, Abrams JS, Bushnell BD, Brock JL, Ho CP. Radiologic and clinical evaluation of a bioabsorbable collagen implant to treat partial-thickness tears: a prospective multicenter study. *J Shoulder Elbow Surg.* 2018;27:242–251.
24. Van Kampen C, Arnoczky S, Parks P, et al. Tissue-engineered augmentation of a rotator cuff tendon using a reconstituted collagen scaffold: a histological evaluation in sheep. *Muscles Ligaments Tendons J.* 2013;3:229–235.
25. Rashid MS, Cooper C, Cook J, et al. Increasing age and tear size reduce rotator cuff repair healing rate at 1 year. *Acta Orthop.* 2017;88:606–611.
26. Lo IK, Burkhart SS. Transtendon arthroscopic repair of partial-thickness, articular surface tears of the rotator cuff. *Arthroscopy.* 2004;20:214–220.
27. Millett PJ, Wilcox RB, O'Holleran JD, Warner JJP. Rehabilitation of the rotator cuff: an evaluation-based approach. *J Am Acad Orthop Surg.* 2006;14:599–609.
28. Arai E, Nakashima H, Tsukushi S, et al. Regenerating the fibula with beta-tricalcium phosphate minimizes morbidity after fibula resection. *Clin Orthop Relat Res.* 2005:233–237.
29. Chu C-C. Biodegradable polymeric biomaterials: an updated overview. In: Bronzino JD, ed. *The biomedical engineering handbook.* CRC Press; 1995:Section IV:44.
30. Calori GM, Mazza E, Colombo M, Ripamonti C. The use of bone-graft substitutes in large bone defects: any specific needs? *Injury.* 2011;42(SUPP 2):S56–S63.
31. Costantino PD, Friedman CD. Synthetic bone graft substitutes. *Otolaryngol Clin North Am.* 1994;27:1037–1074.
32. Gaasbeek RD, Toonen HG, van Heerwaarden RJ, Buma P. Mechanism of bone incorporation of beta-TCP bone substitute in open wedge tibial osteotomy in patients. *Biomaterials.* 2005;26:6713–6719.
33. Park K, Skidmore S, Hadar J, et al. Injectable, long-acting PLGA formulations: Analyzing PLGA and understanding microparticle formation. *J Control Release.* 2019;304:125–134.
34. Walsh WR, Morberg P, Yu Y, et al. Response of a calcium sulfate bone graft substitute in a confined cancellous defect. *Clin Orthop Relat Res.* 2003:228–236.
35. Smith+Nephew 2010. Micro-CT and histological evaluation of specimens from resorbable screw study (RS-II / OM1-08) 24-month post-implantation. Internal Report WRP-TE045-700-08.
36. Smith+Nephew 2016. Healicoil Regenesorb Suture Anchor – a study to assess implant replacement by bone over a 2 year period. NCS248.
37. Milewski MD, Diduch DR, Hart JM, Tompkins M, Ma S-Y, Gaskin CM. Bone replacement of fast-absorbing biocomposite anchors in arthroscopic shoulder labral repairs. *Am J Sports Med.* 2012;40:1392–1401.
38. Sgroi M, Friesz T, Schocke M, Reichel H, Kappe T. Biocomposite suture anchors remain visible two years after rotator cuff repair. *Clin Orthop Relat Res.* 2019;477:1469–1478.
39. Duralde XA. CORR Insights®: Biocomposite suture anchors remain visible two years after rotator cuff repair. *Clin Orthop Relat Res.* 2019;477:1479–1481.

Abkürzungen

| | | | |
|----------------|--|------------------|---|
| 3D-FEM | dreidimensionale Finite-Elemente-Methode | PT | Partial-Thickness (Partial-Riss) |
| ASES | American Shoulder and Elbow Surgeons | SANE | Single-Assessment Numeric Evaluation (Numerische Bewertung einer einzelnen Beurteilung) |
| KI | Konfidenzintervall | SD | Standard Deviation (Standardabweichung) |
| CMS | Constant-Murley-Punktzahl | SECEC | Société Européenne pour la Chirurgie de l'Epaule et du Coude (Europäische Gesellschaft für Schulter- und Ellenbogenchirurgie) |
| CT | Computertomographie | SEM | Standard error of the mean (Standardfehler des Mittelwerts) |
| FT | Full-Thickness (kompletter Riss) | Ti | Titan |
| HRQoL | Health-Related Quality of Life (gesundheitsbezogene Lebensqualität) | VAS | Visuelle Analogskala |
| JSS-SIS | Shoulder Instability Score der Japanese Shoulder Society | VR-12 MCS | Veterans RAND 12-Item Health Survey Mental Component Score |
| MCID | Minimal Clinically Important Difference (minimale klinisch bedeutsame Differenz) | VR-12 PCS | Veterans RAND 12-Item Health Survey Physical Component Score |
| MRT | Magnetresonanztomografie | WORC | Western Ontario Rotator Cuff Index |
| N | Newton | β-TCP | β-Trikalziumphosphat |
| OQS | Ossifikations-Qualitäts-Punktzahl | | |
| PLGA | Poly-L-Lactid-co-Glycolsäure | | |
| PROM | Patient-Reported Outcome Measure (von Patienten angegebene Ergebnisse) | | |

Sports Medicine,
Smith & Nephew, Inc.
150 Minuteman
Road, Andover, MA
01810, USA.

www.smith-nephew.com

◊Marke von Smith+Nephew

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© 2024 Smith+Nephew

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Smith & Nephew GmbH
Friesenweg 30
22763 Hamburg
T +49 (0)40 87 97 44-0
F +49 (0)40 87 97 44-375
info@smith-nephew.com
www.smith-nephew.de

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Smith & Nephew GmbH
Concorde Business Park 1/C/3
2320 Schwechat
Österreich
T +43 1 70 79102
F +43 1 70 79101
Info.austria@smith-nephew.com
www.smith-nephew.com

Kontakt Schweiz

Smith & Nephew Schweiz AG
Theilerstrasse 1A
CH-6300 Zug
Schweiz
T +41 41 766 22 66
F +41 41 766 39 93
CustomerService.CH@smith-nephew.com
www.smith-nephew.com

Erstellt von Evidence Communications,
Global Clinical & Medical Affairs

Smith+Nephew