

A route to more effective infection management

Improve patient outcomes¹ with accurate decision making, a fast response and effective treatment choices

T.I.M.E.
Clinical decision support tool


A
B
C
D
E

Start with following steps to undertake a comprehensive assessment²

- A Assess** patient, wellbeing and wound
- B Bring in** a multi-disciplinary team and informal carers to promote holistic patient assessment
- C Control** and treat the underlying causes and barriers to wound healing
- D Decide** appropriate treatment
- E Evaluate** and reassess the treatment and wound management outcomes

What clinical signs and symptoms of infection are present?

<div style="display: flex; justify-content: space-around;"> </div> <div style="display: flex; justify-content: space-around; margin-top: 5px;"> <div style="border: 1px solid black; padding: 5px; width: 45%;"> <p style="text-align: center;">Biofilm^{1,3-5}</p> <ul style="list-style-type: none"> Antibiotic/antimicrobial treatment failure Recurrence of delayed healing on cessation of antibiotic treatment Delayed healing despite optimal wound/patient management Low level chronic inflammation Low level erythema Friable granulation Covert (subtle) signs of infection </div> <div style="border: 1px solid black; padding: 5px; width: 45%;"> <p style="text-align: center;">Covert (subtle)^{1,3}</p> <ul style="list-style-type: none"> Delayed wound healing Serous drainage with concurrent inflammation Hypergranulation Bleeding, friable granulation Epithelial bridging and pocketing in granulation tissue Wound breakdown & enlargement New or increasing pain Increasing malodour </div> </div>	<div style="display: flex; justify-content: space-around;"> </div> <div style="display: flex; justify-content: space-around; margin-top: 5px;"> <div style="border: 1px solid black; padding: 5px; width: 45%;"> <p style="text-align: center;">Overt (classic)^{1,3}</p> <ul style="list-style-type: none"> Erythema Warmth Oedema/swelling Purulent discharge Pain Increasing malodour Delayed wound healing </div> <div style="border: 1px solid black; padding: 5px; width: 45%;"> <p style="text-align: center;">Spreading or systemic infection^{1,3}</p> <ul style="list-style-type: none"> Spreading erythema, warmth May include cellulitis, crepitus Wound breakdown/dehiscence with or without satellite lesions Malaise/lethargy Loss of appetite Systemic inflammatory response Sepsis Organ dysfunction </div> </div>
--	--



TWO-WEEK CHALLENGE^{1,6,90}

Antimicrobial dressings are recommended to be used for a minimum of two weeks' duration. After two weeks, re-evaluate and either:

- discontinue if signs and symptoms of infection have resolved,
- continue with antimicrobial if wound is progressing but there are still signs and symptoms, or
- consider an alternative antimicrobial and refer to an appropriate specialist if no improvement.

* No one sign or symptom can reliably confirm the presence of infection, and those with immunosuppression may not exhibit signs and symptoms of clinical infection.
 † Cleanse wound and periwound skin thoroughly. Should an antiseptic cleanser be selected, the product's Instructions for Use (IFU) and soak time should be followed.
 ‡ Consider the use of **DURAFIBER[®] Ag Silver Gelling Fibre Dressing** for deep infected wounds.
 § Unless iodine contraindicated.
 ¶ For very-high risk patients and wounds (e.g. osteomyelitis), it may be appropriate to use antimicrobial treatment for longer than the two-week challenge.
 For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's Instructions for Use (IFU).

References 1. International Wound Infection Institute (IWII) Wound infection in clinical practice. Wounds International (2016). 2. Moore Z, et al. J Wound Care 28(3):154-161 (2019). 3. Weir D, Schultz G. Assessment and Management of Wound-Related Infections. In Doughty D & McNichol L (Eds.). Wound, Ostomy and Continence Nurses Society Core Curriculum: Wound Management (p. 156-180). 2016. Philadelphia: Wolters-Kluwer. 4. Wolcott RD, et al. J Wound Care 19(2): 45-53 (2010). 5. Schultz G, et al. Wound Repair Regen 25(5): 744-757 (2017). 6. Ayello EA, et al. Wounds Int 1-24 (2012). 7. Roche ED, et al. Int Wound J 1-10 (2019). 8. Malone M, et al. J Antimicrob Chemother 72, 2093-2101 (2017). 9. Schwarzer S, et al. J Infect 80(3):261-270 (2020). 10. Gago M, Garcia F, Gaztelu V, Verdu J, Lopez P, Nolasco A. A comparison of three silver-containing dressings in the treatment of infected, chronic wounds. Wounds. 2008;20(10):273-278.
 Photographs (from left to right) courtesy of Kerlyn Carville, Kevin Woo, and Henri Post.
 ©Trademark of Smith+Nephew. All Trademarks acknowledged. ©June 2020 Smith+Nephew. 25993 | CA55374 10/20