

Quality Improvement Audit for hard-to-heal wounds using sNPWT

Aim

To evaluate the clinical effectiveness, patient-related benefits and economic impact of a -80mmHg single-use negative pressure wound therapy device (sNPWT)* in routine clinical practice for hard-to-heal wounds.

Background

Hard-to-heal wounds¹ are associated with delayed healing, increased healthcare resource utilisation and reduced patient quality of life². sNPWT has been introduced as a practical solution to improve healing outcomes and reduce treatment burden in routine care³.

Methods

- Quality Improvement Audit conducted in a specialised wound care centre in Germany
- Data collected between May and November 2025
- n = 29 patients with hard-to-heal wounds
- Data collection tool: SnapSurveys (Bristol, UK)⁴
- Wounds were assessed before initiation of sNPWT and weekly during therapy
- Data collection was limited to the treatment period
- sNPWT was discontinued when sufficient wound improvement or closure was achieved

Outcomes assessed:

- Wound healing progression
- Wound area reduction
- Dressing change frequency
- Patient comfort and mobility
- Staff satisfaction



*PICO^o sNPWT. ^oTrademark of Smith+Nephew, Hull, UK. All Trademarks are acknowledged.

BREMEN, GERMANY 6 - 8 MAY 2026

The 36th Conference of the European
Wound Management Association
EWMA 2026

DE WU
DEUTSCHER
WUNDKONGRESS


 **WundD-A-CH**

Quality Improvement Audit for hard-to-heal wounds using sNPWT

A total of 29 hard-to-heal wounds were included. Before sNPWT, wounds were predominantly stagnant or delayed in healing, with some showing clinical deterioration. Most wounds had been present for several weeks to months, indicating a predominantly chronic wound population [Fig.1].

Wound types: surgical wound dehiscence, diabetic foot ulcers, leg ulcers, haematomas.

Fig.1 Wound duration before sNPWT




Results

Clinical outcomes

- 48% of wounds achieved complete healing
 - 48% demonstrated wound area reduction
 - 3% showed no response [Fig.2]
- ✓ 4 wounds healed within 4 weeks
 - ✓ 6 additional wounds healed by 8 weeks
 - ✓ A further 4 wounds healed by 12 weeks

Fig 2. Healing outcomes




Resource utilisation

Dressing change frequency decreased from 4.9/week before sNPWT to 1.7/week during therapy with sNPWT (65% reduction); it then increased to 2.4/week after return to standard care [Fig.3]

Patient and staff experience
 96.4% rated the therapy as good or very good.

Fig. 3 Dressing change frequency



Quality Improvement Audit for hard-to-heal wounds using sNPWT

Conclusion

- In this real-world quality improvement audit, the use of sNPWT was associated with a clear shift from previously stagnant or delayed wounds towards healing.
- Clinical improvement was observed across a range of hard-to-heal wound types, supporting its use in routine practice. In addition, a reduction in dressing change frequency indicates a meaningful decrease in treatment burden in daily care.
- Overall, these findings highlight sNPWT as a practical and effective option for the management of hard-to-heal wounds in routine clinical settings.
- This initiative demonstrated the practical value of sNPWT in supporting healing and reducing treatment burden in everyday clinical practice.

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use prior to use. The information presented is not, intended to serve as, medical advice. It is the responsibility of the treating physician to determine and utilize the appropriate products and techniques according to their own clinical judgment for each of their patients.

References

1. Dowsett C, Hampton J, Myers D, Styche T. Use of PICO™ to improve clinical and economic outcomes in hard-to-heal wounds. *Wounds International*. 2017; 8(2):52-58
2. Probst S, Atkin L, Kapp S, et al. Clinical inertia in the care of patients with chronic wounds. *Wounds International*. 2025;16(1):48-50.
3. Kirsner RS, Dove C, Reyzelman A, Vayser D, Jaimes H. A prospective, randomized, controlled clinical trial on the efficacy of a single-use negative pressure wound therapy system, compared to traditional negative pressure wound therapy in the treatment of chronic ulcers of the lower extremities. *Wound Repair Regen*. 2019;27(5):519-529. <https://doi.org/10.1111/wrr.12727>
4. <https://www.snapsurveys.com/> (last accessed 13th April 2026)