



# Improved patient recovery and fewer 90-day readmissions after RI.KNEE total knee arthroplasty (TKA) on CORI Surgical System versus conventional methods

Cochrane NH, Kim BI, Leal J, Hallows RK, Seyler TM. Comparing a robotic imageless second-generation system to traditional instrumentation in total knee arthroplasty: a matched cohort analysis. *J Orthop.* 2024;57:1–7.

Available at: [Journal of Orthopaedics](#) 

## Key points

Compared with conventional TKA, RI.KNEE TKA on CORI Surgical System resulted in:




## Overview

- Multi-surgeon (two), multi-centre retrospective review of 430 TKAs conducted with RI.KNEE TKA on CORI Surgical System ( $n=215$ ) or conventional methods ( $n=215$ )
- Patients were propensity score matched between groups (by age, sex, BMI, American Society of Anesthesiologists [ASA] score) on a surgeon-specific basis
- Intra-, peri- and post-operative outcomes were assessed including:
  - Procedure time, blood loss, wound drain use
  - Length of stay, discharge location, readmissions (30- and 90-day)
  - Pain scores (pre-operatively, in the care unit, and at 3, 6 and 12 months post-TKA) and PROMs (KOOS JR and PROMIS Global 10; pre-operatively, at 3 months and final follow-up)
  - Survivorship to all-cause revision

## Results

- Compared with conventional TKA, use of RI.KNEE TKA on CORI Surgical System resulted in:
  - Significantly lower pain scores in the post-operative care unit (3 vs 2;  $p=0.016$ ) and at 1-year post-TKA (2 vs 1;  $p=0.022$ )
  - Significantly fewer patients required a post-operative drain (34.4 vs 12.6%;  $p<0.001$ )
  - Significantly earlier discharge by 12 hours on average (mean length of stay: 54 vs 42 hours;  $p=0.022$ )
  - Significantly more patients were discharged to home or self care (60.0 vs 86.5%;  $p<0.001$ )
  - Significantly fewer 90-day readmissions (Figure;  $p=0.001$ ) and numerically fewer 30-day readmissions (6.5 vs 2.3%)
  - No difference in procedure time, blood loss, 3-year survivorship, improvements in PROMs and pain scores at 3 and 6 months



## Conclusions

Use of RI.KNEE TKA on CORI Surgical System resulted in significantly quicker recovery (reduced pain in the care unit and earlier discharge) with significantly fewer 90-day readmissions, compared with conventional TKA. These benefits were seen without an increase in procedure time.

**Abbreviations:** KOOS JR = Knee Injury and Osteoarthritis Outcome Score for Joint Replacement; PROMIS = Patient-Reported Outcomes Measure Information System; PROMs = patient-reported outcome measures.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area. For detailed product information, including indications for use, contraindications, warnings and precautions, please consult the product's Instructions for Use (IFU) prior to use.