



Smith+Nephew

TSF[◇] ALLY

Deformity Correction

Powered by Real Intelligence

Examples of TSF[◇] Ally Reports

Through the collaboration of surgeon and expert engineer, TSF[◊] Ally serves to optimize TSF treatment for your patient and practice, by using dedicated software to conduct the x-ray analysis and accurately plan the hardware assembly

Pre-operative Report — Tibial Malunion with Knee Arthroplasty	4
Post-operative Report — Equinus treated with TSF	14
Pre-operative Report — Multi-apical Tibial Deformity	20
Post-operative Report — Multi-apical Tibial Deformity treated with TSF	34
References and Glossary of Terms	40
Terms of Use	41

TSF Ally is hosted at www.TSFALLY.com

Example of TSF[◊] Ally Report

Pre-operative Planning Report

Malunion with Knee Arthroplasty

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Clinical Details

Surgery Date	XX/XX/XXXX
Length Discrepancy	—
Axial Rotation Deformity	None
Preferred Proximal Ring Size and Type	155 mm
Preferred Distal Ring Size and Type	155 mm
Preferred Reference Ring	Proximal
Preferred Strut Type	Fast FX
Additional Notes	Proximal ring block, with middle ring being the reference ring

Initial Frontal Plane Radiograph







Initial Sagittal Plane Radiograph





Images courtesy of Jeffrey Holmes MD

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia


Deformity Analysis: Frontal Plane														
Pre-Correction Joint Angle Measurements		Joint Line Placement												
Angle/Length	Measured	Normal Range¹												
mLPFA	—	85°-95°												
mLDFA	—	85°-90°												
mMPTA	85°	85°-90°												
mLDTA	99°	86°-92°												
JLCA	0°	0° +/- 2°												
MAD	—	+/- 3 mm												
<p>Pre-Correction Length Measurements</p> <table border="1"> <thead> <tr> <th>Length</th> <th>Measured</th> <th>Discrepancy†</th> </tr> </thead> <tbody> <tr> <td>Femur</td> <td>—</td> <td>—</td> </tr> <tr> <td>Tibia</td> <td>360 mm</td> <td>—</td> </tr> <tr> <td>Total Length</td> <td>—</td> <td>—</td> </tr> </tbody> </table> <p>†Discrepancies measured from bilateral radiograph.</p> <p>Scaling method: 25.4 mm Calibration Sphere</p> 			Length	Measured	Discrepancy†	Femur	—	—	Tibia	360 mm	—	Total Length	—	—
Length	Measured	Discrepancy†												
Femur	—	—												
Tibia	360 mm	—												
Total Length	—	—												
														
														
														

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia


Deformity Analysis: Sagittal Plane		
Pre-Correction		Joint Line Placement
Joint Angle Measurements		
Angle	Measured	Normal Range ¹
PDFA	—	79°-87°
PPTA	89°*	77°-84°
ADTA	78°	78°-82°
ACL	—	Within 5° of total bony deformity
<p>*Please note: Knee implant and rotation of lateral image may lead to inaccurate joint measurements.</p>		
Soft-tissue Contribution †		
STC= ACL—PDFA—PPTA		
N/A		
<p>†Please note that only bony deformity is corrected within this report. Soft-tissue contribution is provided for reference only. Scaling method: 25.4 mm Calibration Sphere</p>		
		



PPTA = 89°



ADTA = 78°




PPTA = 89°

ADTA = 78°

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Deformity Analysis: Apex			
Frontal Plane		Sagittal Plane	
Proximal Axis	Normal mMPTA of 90°	Proximal Axis	Mid-diaphyseal line
Distal Axis	Normal mLDTA of 89°	Distal Axis	Mid-diaphyseal line
Apex Angulation	15° Varus	Apex Angulation	4.7° Recurvatum
Apex Location (to respective joint line)	125 mm proximal to tibiotalar joint line	Apex Location (to respective joint line)	137.9 mm proximal to tibiotalar joint line




Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Correction Analysis

Planned Osteotomy Level	125 mm proximal to tibiotalar joint line
Frontal Plane Angulation Correction	15° Varus
Sagittal Plane Angulation Correction	4.7° Recurvatum
Length Correction (along the axis)	4.9 mm
Anticipated Post-Correction Bony Translation	The distal segment will translate 1 mm posterior to the proximal segment.

Frontal Plane Correction	Sagittal Plane Correction
---------------------------------	----------------------------------



- | | |
|---|---|
| — Proximal Segment | — Ruler |
| — Distal Segment | — Transverse Bisector |

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia



Correction Analysis			Frontal Plane	Sagittal Plane		
Post-Correction Frontal Joint Angle Measurements						
Angle/Length	Measured	Normal				
mLPFA	—	85°-95°				
mLDFA	—	85°-90°				
mMPTA	90°	85°-90°				
mLDTA	89°	86°-92°				
JLCA	0°	0°+/-2°				
MAD	—	+/- 3mm				
Post-Correction Sagittal Joint Angle Measurements						
Angle	Measured	Normal Range¹				
PDFA	—	79°-87°				
PPTA	87°*	77°-84°				
ADTA	81°	78°-82°				
*Please note: Knee implant and rotation of lateral image may lead to inaccurate joint measurements.						
Post-Correction Length Measurements						
Length	Measured	Residual Discrepancy†				
Femur	—	—				
Tibia	367 mm	—				
Total Length	—	—				
†Discrepancies measured from bilateral radiograph.						

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

TSF Planning



Length Correction (along the axis)	4.9 mm
Osteotomy Level (to respective joint line)	125 mm proximal to tibiotalar joint line
Reference Ring Placement	130 mm distal to knee joint line
Neutral Frame Height	190 mm
Neutral Strut Length	175 mm
Ring Sizes	155 mm Full (both proximal), 155 mm Full

Frontal Plane Correction	Sagittal Plane Correction
--------------------------	---------------------------



- | | |
|---|---|
| <ul style="list-style-type: none"> — Reference Ring(s) — Moving Ring(s) — Joint Lines | <ul style="list-style-type: none"> — Ruler — Osteotomy ● Center of Ring |
|---|---|

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

TSF Planning					
Deformity Parameters					
(Method of Planning: Apex = Corresponding Point)					
AP View Angulation	AP View Translation	Lateral View Angulation	Lateral View Translation	Axial View Angulation	Axial Translation
15° Varus	0 mm	4.7° Apex Posterior	1.3 mm Anterior	0°	4.9 mm Short
Mounting Parameters					
Referencing			Proximal		
AP View Frame Offset			0 mm		
Lateral View Frame Offset			35 mm posterior to origin		
Rotary Frame Angle			0°		
Axial Frame Offset			119.7 mm proximal to (extrinsic) origin		
The distance from the reference ring to the origin is increased when extrinsic length is added.					
Initial Strut Settings					
Strut 1 ■	Strut 2 ■	Strut 3 ■	Strut 4 ■	Strut 5 ■	Strut 6 ■
178	154	151	170	186	193
Medium Fast Fx [®] Strut	Medium Fast Fx [®] Strut	Medium Fast Fx [®] Strut	Medium Fast Fx [®] Strut	Medium Fast Fx [®] Strut	Medium Fast Fx [®] Strut
Initial Frame in Frontal Plane			Initial Frame in Sagittal Plane		
					

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Case Summary					
Deformity Analysis					
Frontal Plane Angulation			15° Varus		
Sagittal Plane Angulation			4.7° Recurvatum		
Leg Length Discrepancy Pre-Correction			—		
Correction Analysis					
Planned Osteotomy Level			125 mm proximal to tibiotalar joint line		
Total Anticipated Length Correction			4.9 mm along the axis		
Leg Length Discrepancy Post-Correction			—		
TSF Planning					
Ring Sizes/Types			Proximal: 155 mm Full (both proximal rings)		
Reference Ring (Must be placed orthogonal to bone segment.)			2 nd Ring of proximal ring block (Proximal referencing)		
Reference Ring Placement			130 mm distal to knee joint line		
Number of Strut Change-outs			0		
Initial Strut Settings					
Strut 1 ■	Strut 2 ■	Strut 3 ■	Strut 4 ■	Strut 5 ■	Strut 6 ■
178	154	151	170	186	193
Medium Fast Fx ® Strut	Medium Fast Fx ® Strut	Medium Fast Fx ® Strut	Medium Fast Fx ® Strut	Medium Fast Fx ® Strut	Medium Fast Fx ® Strut
Parts List					
Part (Catalog Number)			Quantity		
155mm Ring (7107-0114)			2		
Fast Fx ® Identification Band Kit (7107-0340)			1		
Medium Fast Fx ® Strut (7107-0720)			6		
Notes					
<p>Knee implant and rotation of lateral image may lead to inaccurate joint measurements.</p> <p>We selected a mid-diaphyseal line for the proximal axis in the sagittal view due to the implantation of the tibia baseplate and the rotation of the lateral x-ray. The remaining abnormal PPTA may be inaccurate based on these factors.</p> <p>Planned TSF using a 60 mm proximal ring block, with 2nd ring of block as the reference ring (proximal referencing). 60 mm Threaded Sockets may be used to assemble the proximal ring block as shown.</p> <p>Planned TSF using 155 mm distal ring based on preferences; however, a 130 mm ring would also fit based on soft-tissue measurements.</p>					

Example of TSF[◊] Ally Report

Post-operative Planning Report

Equinus treated with TSF


Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Left Foot

Post-Op Clinical Details

Surgery Date	XX/XX/XXXX
Prescription Start Date	XX/XX/XXXX
Length Discrepancy	Distract 5 mm
Axial Rotation Deformity	None
Proximal Ring Size and Type	155 mm Full
Distal Ring Size and Type	155 Foot Plate
Reference Ring	Distal (Foot Plate)
Structure at Risk (SAR)	Posterior Ankle Capsule
Maximum Safe Distraction Rate	1 mm/day
Additional Notes	Soft tissue correction





Initial Strut Settings

Strut 1 ■	Strut 2 ■	Strut 3 ■	Strut 4 ■	Strut 5 ■	Strut 6 ■
248	230	200	98	101	190
Long Fast FX	Long Fast FX	Medium Fast FX	XS Fast FX	XS Fast FX	Medium Fast FX



Images courtesy of Eugene Stautberg MD

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Left Foot


Deformity Analysis			
Method: TraumaCad			
Proximal Axis	Normal mLDTA of 89°	Proximal Axis	Mid-diaphyseal line of tibia
Distal Axis	Center of talar dome	Distal Axis	Sagittal axis of foot
AP View Angulation	3° Varus	Lateral View Angulation	48° Equinus
AP View Translation	0 mm	Lateral View Translation	0 mm
Frontal Plane Deformity Analysis		Sagittal Plane Deformity Analysis	
			
 Reference Ring/Segment		 Moving Segment	

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Left Foot


TSF Mounting Parameters

AP View Frame Offset	7.5 mm Medial to Origin
Lateral View Frame Offset	9.4 mm Anterior to Origin
Axial Frame Offset	58.1 mm Distal to Origin (Extrinsic)
Rotary Frame Angle	8° Externally Rotated

Frontal Plane Reference Ring



Sagittal Plane Reference Ring



— Ruler


— Reference Axis (Center of Talus)
 — Center of Reference Ring

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Left Foot


Structure at Risk Parameters





Structure at Risk	Posterior Ankle Capsule
AP View SAR Offset	0 mm
Lateral View SAR Offset	50 mm Posterior to Origin
Axial SAR Offset	9.2 mm Distal to Origin (Extrinsic)

**Frontal Plane
Reference Ring and SAR Placement**



**Sagittal Plane
Reference Ring and SAR Placement**



	Structure at Risk		Reference Axis (Center of Talus)
	Ruler		Center of Reference Ring

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Left Foot

Case Summary					
Deformity Parameters					
AP View Angulation	3° Varus		AP View Translation	0 mm	
LAT View Angulation	58° Plantar Flexion		LAT View Translation	0 mm	
Axial View Angulation	0°		Axial View Translation	5 mm Short	
Mounting Parameters					
AP View Frame Offset			7.5 mm Medial to Origin		
LAT View Frame Offset			9.4 mm Anterior to Origin		
Axial View Frame Offset			58.1 mm Distal to Origin (Extrinsic)		
Rotary Frame Angle			8° Externally Rotated		
Structure at Risk Parameters					
Structure at Risk			Posterior Ankle Capsule		
AP View SAR Offset			0 mm		
Lateral View SAR Offset			50 mm Posterior to Origin		
Axial SAR Offset			9.2 mm Distal to Origin (Extrinsic)		
Correction Analysis					
Length Correction			5 mm along the axis		
Anticipated Post-Correction Bony Translation			0 mm		
Maximum Safe Distraction Rate			1 mm/day		
Correction Time			54 days		
Number of Strut Change-Outs			6		
Strut Change-Outs					
Change-Out	Strut	Overlap Start Date	Overlap End Date	Strut Change From	Strut Change To
a	2 (Orange)	38 (03/05/20)	50 (03/17/20)	7107-0730 Long FAST FX [®]	7107-0720 Medium FAST FX [®]
b	3 (Yellow)	0 (01/27/20)	25 (02/21/20)	7107-0720 Medium FAST FX [®]	7107-0730 Long FAST FX [®]
c	4 (Green)	14 (02/10/20)	17 (02/13/20)	7107-0705 Extra Short FAST FX [®]	7107-0710 Short FAST FX [®]
d	4 (Green)	29 (02/25/20)	34 (03/01/20)	7107-0710 Short FAST FX [®]	7107-0720 Medium FAST FX [®]
e	5 (Blue)	13 (02/09/20)	16 (02/12/20)	7107-0705 Extra Short FAST FX [®]	7107-0710 Short FAST FX [®]
f	5 (Blue)	30 (02/26/20)	36 (03/03/20)	7107-0710 Short FAST FX [®]	7107-0720 Medium FAST FX [®]

Example of TSF[◊] Ally Report

Pre-operative Planning Report

Multi-apical Tibial Deformity

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Clinical Details

Surgery Date	XX/XX/XXXX
Length Discrepancy	Right Leg 5 cm shorter than left
Axial Rotation Deformity	0° Rotational Deformity
Preferred Proximal Ring Size and Type	180 mm Full Ring
Preferred Distal Ring Size and Type	180 mm Full Ring
Preferred Reference Ring	Middle Ring
Preferred Strut Type	Fast FX
Additional Notes	Patient has a tibia malunion and was treated with a previous frame.

Initial Frontal Plane Radiograph	Initial Sagittal Plane Radiograph
---	--



Images courtesy of Evgeny Dyskin MD





Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Deformity Analysis: Frontal Plane

Pre-Correction Joint Angle Measurements			Joint Line Placement	
Angle/Length	Measured	Normal Range ¹		
mLPFA	85°	85°-95°		
mLDFA	85°	85°-90°		
mMPTA	64°	85°-90°		
mLDTA	84°	86°-92°		
JLCA	1°	0° +/- 2°		
MAD	69 mm	+/- 3 mm		
Pre-Correction Length Measurements				
Length	Measured	Discrepancy [†]		
Femur	484 mm	2 mm		
Tibia	344 mm	32 mm		
Total Length	819 mm	48 mm short		
[†] Discrepancies measured from bilateral radiograph. Scaling method: 25.4 mm Calibration Sphere				

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Deformity Analysis: Sagittal Plane


Pre-Correction			Joint Line Placement	
Joint Angle Measurements				
Angle	Measured	Normal Range¹		
PDFA	—	79°-87°		
PPTA	82°	77°-84°		
ADTA	80°	78°-82°		
ACL	—	Within 5° of total bony		
Soft-tissue Contribution †				
STC= ACL—PDFA—PPTA				
N/A				
<p>†Please note that only bony deformity is corrected within this report. Soft-tissue contribution is provided for reference only.</p> <p>Scaling method: 25.4 mm Calibration Sphere</p> 				

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Deformity Analysis: Frontal Plane Apices

Proximal Axis	Match normal contralateral mMPTA of 85°
Middle Axes	Mid-diaphyseal lines
Distal Axis	Match normal contralateral mLDTA of 86°
Apex Angulations	Proximal: 37° Varus Middle: 14° Valgus Distal: 3° Valgus
Apex Locations (to respective joint line)	Proximal: 77.8 mm distal to proximal joint line Middle: 167 mm proximal to distal joint line

Proximal Apex	Middle Apex	Distal Apex
----------------------	--------------------	--------------------



— Proximal Segment
— Proximal Middle Segment


— Distal Middle Segment
— Distal Segment

— Ruler
⋯ Transverse Bisector

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Deformity Analysis: Sagittal Plane Apices

Proximal Axis	Normal PPTA of 84°
Middle Axes	Mid-diaphyseal lines
Distal Axis	Normal ADTA of 82°
Apex Angulations	Proximal: 5° Procurvatum Middle: 3° Recurvatum Distal: 3° Recurvatum
Apex Locations (to respective joint line)	Proximal: 89.8 mm distal to proximal joint line Middle: 233.9 mm proximal to distal joint line




Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Correction Analysis

	Proximal Frame	Middle Frame	Distal Frame
Planned Osteotomy Level	65 mm distal to	215 mm proximal to tibiotalar joint line	81 mm proximal to tibiotalar joint line
Frontal Plane Angulation Correction	36.9° Varus	14.5° Valgus	3° Valgus
Sagittal Plane Angulation Correction	5.3° Procurvatum	3° Recurvatum	3.2° Recurvatum
Length Correction (along the axis)	12.6 mm	2.1 mm	0.7 mm
Anticipated Post-Correction Bony Translation	Proximal Frame: The proximal middle segment is 7.8 mm medial and 2.2 mm posterior to the proximal segment.		
	Middle Frame: The proximal middle segment is 12 mm medial and 1 mm anterior to the distal segment.		
	Distal Frame: The distal segment is 1.9 mm medial and 1.9 mm anterior to the distal middle segment.		

Frontal Plane Correction	Sagittal Plane Correction
---------------------------------	----------------------------------




- Proximal Segment
- Distal Middle Segment
- Proximal Middle Segment
- Distal Segment
- Ruler
- Osteotomy

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Correction Analysis			Frontal Plane	Sagittal Plane
Post-Correction Frontal Joint Angle Measurements				
Angle/Length	Measured	Normal Range¹		
mLPFA	85°	85°-95°		
mLDFA	85°	85°-90°		
mMPTA	85°	85°-90°		
mLDTA	86°	86°-92°		
JLCA	1°	0°+/-2°		
MAD	3 mm	+/- 3mm		
Post-Correction Sagittal Joint Angle Measurements				
Angle	Measured	Normal Range¹		
PDFA	—	79°-87°		
PPTA	84°	77°-84°		
ADTA	82°	78°-82°		
Post-Correction Length Measurements				
Length	Measured	Residual Discrepancy†		
Femur	484 mm	2 mm		
Tibia	369 mm	7 mm		
Total Length	856 mm	11 mm short		
†Discrepancies measured from bilateral radiograph.				

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

TSF Planning			
	Proximal Frame	Middle Frame	Distal Frame
Length Correction (along the axis)	15.6 mm	4.1 mm	6.7 mm
Osteotomy Level (to respective joint line)	65 mm distal	215 mm proximal	81 mm proximal
Reference Ring Placement	135.5 mm distal		31.1 mm proximal
Neutral Frame Height	118 mm	135 mm	77.7 mm
Neutral Strut Length	122 mm	135 mm	98 mm
Ring Sizes	180 mm, 180 mm	180 mm, 180 mm	180 mm, 155 mm



Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

TSF Planning: Proximal Frame

Proximal Frame Deformity Parameters
 (Method of Planning: Apex = Corresponding Point)

AP View Angulation	AP View Translation	Lateral View Angulation	Lateral View Translation	Axial View Angulation	Axial Translation
36.9° Varus	0 mm	5.3° Apex Anterior	1.1 mm Posterior	0°	15.6 mm short

Proximal Frame Mounting Parameters

Referencing	Distal
AP View Frame Offset	0 mm
Lateral View Frame Offset	40 mm posterior to origin
Rotary Frame Angle	0°
Axial Frame Offset	70.6 mm proximal to (extrinsic) origin
The distance from the reference ring to the origin is increased when extrinsic length is added.	

Proximal Frame Initial Strut Settings

Strut 1 ■	Strut 2 ■	Strut 3 ■	Strut 4 ■	Strut 5 ■	Strut 6 ■
130	95	84	89	119	155
Medium Strut	Short Strut	Extra Short Strut	Extra Short Strut	Medium Strut	Medium Strut

Initial Frame in Frontal Plane **Initial Frame in Sagittal Plane**



Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

TSF Planning: Middle Frame

Middle Frame Deformity Parameters
 (Method of Planning: Apex = Corresponding Point)

AP View Angulation	AP View Translation	Lateral View Angulation	Lateral View Translation	Axial View Angulation	Axial Translation
14.5° Valgus	0 mm	3° Apex Posterior	3.6 mm Anterior	0°	4.1 mm Short


Middle Frame Mounting Parameters

Referencing	Proximal
AP View Frame Offset	0 mm
Lateral View Frame Offset	40 mm posterior to origin
Rotary Frame Angle	0°
Axial Frame Offset	21 mm proximal to (extrinsic) origin
The distance from the reference ring to the origin is increased when extrinsic length is added.	



Middle Frame Initial Strut Settings

Strut 1 ■	Strut 2 ■	Strut 3 ■	Strut 4 ■	Strut 5 ■	Strut 6 ■
126	131	157	151	118	120
Short Fast Fx	Short Fast Fx	Medium Fast Fx	Medium Fast Fx	Short Fast Fx	Short Fast Fx

Initial Frame in Frontal Plane **Initial Frame in Sagittal Plane**



Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

TSF Planning: Distal Frame					
Distal Frame Deformity Parameters (Method of Planning: Apex = Corresponding Point)					
AP View Angulation	AP View Translation	Lateral View Angulation	Lateral View Translation	Axial View Angulation	Axial Translation
3° Valgus	0 mm	3.2° Apex Posterior	3.9 mm Anterior	0°	6.7 mm Short
Distal Frame Mounting Parameters					
Referencing			Distal		
AP View Frame Offset			0 mm		
Lateral View Frame Offset			20 mm posterior to origin		
Rotary Frame Angle			0°		
Axial Frame Offset			56.7 mm distal to (extrinsic) origin		
The distance from the reference ring to the origin is increased when extrinsic length is added.					
Distal Frame Initial Strut Settings					
Strut 1 ■	Strut 2 ■	Strut 3 ■	Strut 4 ■	Strut 5 ■	Strut 6 ■
94	95	95	101	97	91
Extra Short Fast Fx	Extra Short Fast Fx	Extra Short Fast Fx	Extra Short Fast Fx	Extra Short Fast Fx	Extra Short Fast Fx
Initial Frame in Frontal Plane			Initial Frame in Sagittal Plane		
					

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Parts List

<u>Part (Catalog Number)</u>	<u>Quantity</u>
180 mm Ring (7107-0115)	3
155 mm Ring (7107-0114)	1
Standard Identification Band Kit (7107-0320)	1
Fast Fx Identification Band Kit (7107-0340)	2
Medium Strut (7107-0220)	3
Short Strut (7107-0210)	3
Extra Short Strut (7107-0205)	2
Medium Fast Fx Strut (7107-0720)	2
Short Fast Fx Strut (7107-0710)	6
Extra Short Fast Fx Strut (7107-0705)	6

Notes

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia

Case Summary						
Deformity Analysis						
Proximal Deformity Angulations			36.9° Varus, 5.3° Procurvatum			
Middle Deformity Angulations			14.5° Valgus, 3° Recurvatum			
Distal Deformity Angulations			3° Valgus, 3.2° Recurvatum			
Leg Length Discrepancy Pre-Correction			Right Lower Limb is 48 mm shorter than Left			
Correction Analysis						
Planned Osteotomy Levels			Proximal: 65 mm distal to knee joint line Middle: 215 mm proximal to tibiotalar joint line Distal: 81 mm proximal to tibiotalar joint line			
Total Anticipated Length Correction			26.4 mm along the axis			
Leg Length Discrepancy Post-Correction			0 mm			
TSF Planning						
Ring Sizes/Types			Proximal: 180 mm Full Middle: 180 mm Full and 180 mm Full Distal: 155 mm Full			
Reference Ring			Proximal Middle Ring and Distal			
Reference Ring Placement			Proximal Middle: 135.5 mm distal to knee joint line			
Proximal Frame Initial Strut Settings						
Strut 1 ■	Strut 2 ■	Strut 3 ■	Strut 4 ■	Strut 5 ■	Strut 6 ■	
130	95	84	89	119	155	
Medium Strut	Short Strut	Extra Short Strut	Extra Short Strut	Medium Strut	Medium Strut	
Middle Frame Initial Strut Settings						
Strut 1 ■	Strut 2 ■	Strut 3 ■	Strut 4 ■	Strut 5 ■	Strut 6 ■	
126	131	157	151	118	120	
Short Fast Fx	Short Fast Fx	Medium Fast Fx	Medium Fast Fx	Short Fast Fx	Short Fast Fx	
Distal Frame Initial Strut Settings						
Strut 1 ■	Strut 2 ■	Strut 3 ■	Strut 4 ■	Strut 5 ■	Strut 6 ■	
94	95	95	101	97	91	
Extra Short Fast Fx	Extra Short Fast Fx	Extra Short Fast Fx	Extra Short Fast Fx	Extra Short Fast Fx	Extra Short Fast Fx	

Example of TSF[◊] Ally Report

Post-operative Planning Report

Multi-apical Tibial Deformity treated with
TSF

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia (Proximal Frame)


Post-Op Clinical Details

Surgery Date	XX/XX/XXXX
Prescription Start Date	XX/XX/XXXX
Length Discrepancy	Per pre-op planning
Axial Rotation Deformity	0°
Proximal Ring Size and Type	180 mm, 2/3rd (Open between Struts 4 & 5)
Distal Ring Size and Type	180mm Full
Reference Ring	Distal
Structure at Risk (SAR)	Concavity of the deformity
Maximum Safe Distraction Rate	1 mm/day
Additional Notes	

Initial Strut Settings

Strut 1 ■	Strut 2 ■	Strut 3 ■	Strut 4 ■	Strut 5 ■	Strut 6 ■
143	134	91	127	157	188
Short Fast FX (7107-0710)	Short Fast FX (7107-0710)	X-Short Fast FX (7107-0705)	Short Fast FX (7107-0710)	Medium Fast FX (7107-0720)	Medium Fast FX (7107-0720)

Initial Frontal Plane Radiograph



Initial Sagittal Plane Radiograph



Images courtesy of Evgeny Dyskin MD

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia (Proximal Frame)

TSF Mounting Parameters



AP View Frame Offset	4.8 mm Medial to Origin
Lateral View Frame Offset	32.9 mm Posterior to Origin
Axial Frame Offset	127.6 mm Distal to Origin (Extrinsic)
Rotary Frame Angle	0° - based on TraumaCad Measurements


Structure at Risk Parameters


Structure at Risk	Concavity of the deformity
AP View SAR Offset	20.8 mm Medial to Origin
Lateral View SAR Offset	19.2 mm Posterior to Origin
Axial SAR Offset	0 mm Proximal to Origin

**Frontal Plane
Reference Ring and SAR Placement**


**Sagittal Plane
Reference Ring and SAR Placement**



 Structure at Risk

 Master Tab

 Reference Ring

 Osteotomy


Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia (Proximal Frame)

Deformity Analysis


Method: TraumaCad

Proximal Axis	Normal mMPTA of 87°	Proximal Axis	Normal PPTA of 84°
Distal Axis	Mid-diaphyseal line	Distal Axis	Mid-diaphyseal line
AP View Angulation	26.3° Varus	Lateral View Angulation	3.7° Apex Anterior
AP View Translation	0 mm	Lateral View Translation	0.6 mm Posterior

Frontal Plane Deformity Analysis



Sagittal Plane Deformity Analysis




Structure at Risk	Master Tab	Reference Ring/Segment	Osteotomy
		Moving Segment	Transverse Bisector

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia (Proximal Frame)


TSF Correction Analysis







Frontal Plane Angulation Correction	26.3° Varus
Sagittal Plane Angulation Correction	3.7° Procurvatum
Length Correction	15.6 mm about the axis
Anticipated Post-Correction Bony Translation	The proximal segment will translate 9.5 mm lateral to the distal segment.

Frontal Plane Correction Analysis



Sagittal Plane Correction Analysis



 Structure at Risk	 Master Tab	 Reference Ring/Segment	 Osteotomy
		 Moving Segment	 Transverse Bisector

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia (Proximal Frame)

Case Summary					
Deformity Parameters					
AP View Angulation	26.3° Varus		AP View Translation	0 mm	
LAT View Angulation	3.7° Apex Anterior		LAT View Translation	0.6 mm Posterior	
Axial View Angulation	0°		Axial View Translation	15.6 mm Short	
Mounting Parameters					
AP View Frame Offset			4.8 mm Medial to Origin		
LAT View Frame Offset			32.9 mm Posterior to Origin		
Axial View Frame Offset			127.6 mm Distal to Origin		
Rotary Frame Angle			0		
Structure at Risk Parameters					
Structure at Risk			Concavity of Deformity		
AP View SAR Offset			20.8 mm Medial to Origin		
Lateral View SAR Offset			19.2 mm Posterior to Origin		
Axial SAR Offset (to Extrinsic Origin)			0		
Correction Analysis					
Length Correction			15.6 mm about the axis		
Anticipated Post-Correction Bony Translation			The proximal segment will translate 9.5 mm medial to the distal segment.		
Maximum Safe Distraction Rate			1 mm/day		
Correction Time			27 Days		
Number of Strut Change-Outs			2		
Strut Change-Outs					
Change-Out	Strut	Overlap Start Date	Overlap End Date	Strut Change From	Strut Change To
a	3 (Yellow)	12 (01/10/20)	13 (01/11/20)	7107-0705 Extra Short FAST FX ®	7107-0710 Short FAST FX ®
b	4 (Green)	9 (01/07/20)	14 (01/12/20)	7107-0710 Short FAST FX ®	7107-0720 Medium FAST FX ®
Notes					
Case will be sent via spatialframe.com once plan is approved.					

Patient ID	XXXXXXXX
Patient	Patient Name
Surgeon	Dr. Smith
Planning Date	XX/XX/XXXX
Anatomy	Right Tibia (Proximal Frame)

References

1. Standard SC, Herzenberg JE, Conway JD, Siddiqui NA, McClure PK. The Art of Limb Alignment. 8th ed. Baltimore: Rubin Institute for Advanced Orthopedics, Sinai Hospital of Baltimore, 2019.

Glossary

Joint Angle Names

Angle (°)	Complete Name of Joint Angle
mLPFA	mechanical Lateral Proximal Femoral Angle
mLDFA	mechanical Lateral Distal Femoral Angle
mMPTA	mechanical Medial Proximal Tibial Angle
mLDTA	mechanical Lateral Distal Tibial Angle
JLCA	Joint Line Convergence Angle
MAD	Mechanical Axis Deviation
PDFA	Posterior Distal Femoral Angle
PPTA	Posterior Proximal Tibial Angle
ADTA	Anterior Distal Tibial Angle
ACL	Anterior Cortical Lines

Terms of Use

PERMITTED USE

These Terms of Use permit you to use TAYLOR SPATIAL FRAME[◇] Ally (TSF[◇] Ally). TSF Ally enables you to provide clinical assessment data, imaging, and other context in order to provide support to assist in surgical planning of the patients you treat (e.g., limb alignment analysis and deformity correction reports, case examples, case discussions and conference calls, measurements and correction of the deformity being treated, TSF configuration being utilized, the initial/final settings of the TSF fixation device, inputs/outputs from spatialframe.com (the Website) and any other necessary information), referred to in this Terms of Use as "Support". You undertake full responsibility for non-compliance with these Terms of Use.

INDEPENDENT MEDICAL JUDGMENT

You are solely responsible for the Support generated through the use of TSF Ally. TSF Ally, or the Support generated, shall not be construed as a substitute for the independent medical judgment of a physician relative to diagnostic and treatment options of a specific patient's medical condition. You will use your independent medical judgment in treating patients, and you will not use or rely upon any information provided through TSF Ally Support as a diagnosis or treatment of your patients. You understand and agree that information provided through TSF Ally is for planning and scheduling purposes only, and you will use it only in that capacity. Any decisions made following use of TSF Ally Support are the sole responsibility of the user, who should determine whether or not to use any information provided through TSF Ally Support exercising his/her own independent medical judgment, and how to use such information to make his or her own clinical decisions.

NO RETENTION OF RECORDS

TSF Ally is not a part of your medical practice files or patient treatment records, including your electronic health records. You should download or print a copy of the Support materials or other information for each case. We do not represent or warrant that we will maintain a copy of your cases and disclaim all responsibility and liability for the loss of data you have input into TSF Ally. You are solely responsible for complying with any state, federal or international medical records retention laws applicable to you or your medical practice. You further agree that the Company (Smith+Nephew) has no obligation to store any content you may enter, upload or share with TSF Ally after termination of your account, and we shall have the right to delete your content permanently.