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Advanced Wound Management  
Smith & Nephew Medical Ltd  
101 Hesse Road  
Hull HU3 2BN, UK  
T +44 (0) 1482 225181  
F +44 (0) 1482 328326

[www.smith-nephew.com](http://www.smith-nephew.com)

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**REFERENCIAS:** 1. Karlakki SL, Hamad AK, Whittall C, Graham NM, Banerjee RD, Kuiper JH. Incisional negative pressure wound therapy dressings (NPWTd) in routine primary hip and knee arthroplasties: a randomised controlled trial. *Bone & joint research*. 2016 Aug;5(8):328-37. 2. Mhaske AS, Budithi S, Karlakki S. Incisional Negative Pressure Wound Therapy Can Reduce Wound Complications And Hospital Stay In Revision Hip And Knee Arthroplasty. Paper presented at: EFORT; 2019; Lisbon. 3. Nherera LM, Trueman P, Karlakki SL. Cost-effectiveness analysis of single-use negative pressure wound therapy dressings (sNPWT) to reduce surgical site complications (SSC) in routine primary hip and knee replacements. *Wound Repair and Regeneration*. 2017 May;25(3):474-82. 4. Patel VP, Walsh M, Sehgal B, Preston C, DeWal H, Di Cesare PE. Factors associated with prolonged wound drainage after primary total hip and knee arthroplasty. *JBJS*. 2007 Jan 1;89(1):33-8. 5. Holm B, Kristensen MT, Bencke J, Husted H, Kehlet H, Bandholm T. Loss of knee-extension strength is related to knee swelling after total knee arthroplasty. *Archives of physical medicine and rehabilitation*. 2010 Nov 1;91(11):1770-6. 6. Tucker A, Walls A, Leckey B, Hill JC, Phair G, Bennett DB, O'Brien S, Beverland DE. Postdischarge unscheduled care burden after lower limb arthroplasty. *The Journal of arthroplasty*. 2018 Sep 1;33(9):2745-51. 7. Avram V, Petruccielli D, Winemaker M, de Beer J. Total joint arthroplasty readmission rates and reasons for 30-day hospital readmission. *The Journal of arthroplasty*. 2014 Mar 1;29(3):465-8. 8. Poultsides LA, Ma Y, Della Valle AG, Chiu YL, Sculco TP, Memtsoudis SG. In-hospital surgical site infections after primary hip and knee arthroplasty—incidence and risk factors. *The Journal of arthroplasty*. 2013 Mar 1;28(3):385-9. 9. Anis HK, Sodhi N, Klika AK, Mont MA, Barsoum WK, Higuera CA, Molloy RM. Is operative time a predictor for post-operative infection in primary total knee arthroplasty? *The Journal of arthroplasty*. 2019 Jul 1;34(7):S331-6. 10. Leekha S, Sampathkumar P, Berry DJ, Thompson RL. Should national standards for reporting surgical site infections distinguish between primary and revision orthopedic surgeries? *Infection Control & Hospital Epidemiology*. 2010 May;31(5):503-8. eek 11. Drew JM, Griffin WL, Odum SM, Van Doren B, Weston BT, Stryker LS. Survivorship after periprosthetic femur fracture: factors affecting outcome. *The Journal of arthroplasty*. 2016 Jun 1;31(6):1283-8. 12. Matharu GS, Pynsent PB, Dunlop DJ, Revell MP. Clinical outcome following surgical intervention for periprosthetic hip fractures at a tertiary referral centre. *Hip International*. 2012 Sep;22(5):494-9. 13. Everhart JS, Andridge RR, Scharschmidt TJ, Mayerson JL, Glassman AH, Lemeshow S. Development and validation of a preoperative surgical site infection risk score for primary or revision knee and hip arthroplasty. *JBJS*. 2016 Sep 21;98(18):1522-32. 14. Sugrue M, Ciprandi G, Djohan R, et al. World Union of Wound Healing Societies (WUWHS) Consensus Document. Closed surgical incision management: Understanding the role of NPWT. *Wounds Int* [Internet]. 2016. [www.woundsinternational.com/wuwhs/view/consensus-document-closed-surgical-incision-management-understanding-the-roleof-npwt](http://www.woundsinternational.com/wuwhs/view/consensus-document-closed-surgical-incision-management-understanding-the-roleof-npwt). (Last accessed July 29, 2017). 15. The World Health Organisation (2016) Global Guidelines for the Prevention of Surgical Site Infection. [online] Available from: <https://apps.who.int/iris/bitstream/handle/10665/250680/9789241549882-eng.pdf?sequence=8>. Last Accessed June 2020 16. Saunders C, Buzza K, Nherera L. 2019. A single use negative pressure system reduces surgical site complications Compared with conventional dressings in closed surgical incisions: a systematic literature review with meta-analysis. Poster presented at the European Wound Management Association annual meeting, June 5-7, 2019, Gothenburg, Sweden. 17. NICE (2019) PICO negative pressure wound dressings for closed surgical incisions [online] accessible from: <https://www.nice.org.uk/guidance/mtg43>. Last accessed June 2020 18. Canonico S, Campitiello F, Della Corte A. Therapeutic possibilities of portable NPWT. *Acta Vulnologica* 10 (2012): 57-64. 19. Hyldig N, Birke-Sorensen H, Kruse M, Vinter C, Joergensen JS, Sorensen JA, Mogensen O, Lamont RF, Bille C. Meta-analysis of negative-pressure wound therapy for closed surgical incisions. *British Journal of Surgery*. 2016 Apr;103(5):477-86. 20. Loveluck J, Copeland T, Hill J, Hunt A, Martin R. Biomechanical modeling of the forces applied to closed incisions during single-use negative pressure wound therapy. *ePlasty*. 2016;16. 21. Malmjö M, Huddleston E, Martin R. Biological effects of a disposable, canisterless negative pressure wound therapy system. *ePlasty*. 2014;14. 22. Pellino G, Sciaudone G, Candilio G, Campitiello F, Selvaggi F, Canonico S. Effects of a new pocket device for negative pressure wound therapy on surgical wounds of patients affected with Crohn's disease: a pilot trial. *Surgical Innovation*. 2014 Apr;21(2):204-12. 23. Data on file reference 1102010 – Bacterial Barrier Testing (wet-wet) of PICO Dressing with a 7 day test Duration against *S. marcescens*. 24. Smith & Nephew October 2017. Project Opal PICO 7 System Stability Testing, Initial Time Point. Internal report. DS/17/253/R. 25. Data on file, report DS.17/666/R2. Comparison of PICO 1.6 and 2.1 Device Air Leak Tolerance. January 2018. 26. Smith & Nephew January 2019. Acoustic Testing Report: Comparison of PICO v2 (PICO 7 and PICO 14) Devices to PICO v1.6 (PICO) devices. Internal report. RD/19/005. 27. Smith & Nephew December 2018. Summary of PICO v2 (PICO 7 and PICO 14) Human Factors Testing. Internal report. RD/18/136. 28. Birke-Sorensen H, Malmjö M, Rome P, et al. Evidence-based recommendations for negative pressure wound therapy: treatment variables (pressure levels, wound filler and contact layer)—steps towards an international consensus. *J Plast Reconstr Aesthet Surg*. 2011;64 Suppl:S1-16. 29. Scalise A, Calamita R, Tartaglione C, et al. Improving wound healing and preventing surgical site complications of closed surgical incisions: a possible role of Incisional Negative Pressure Wound Therapy. A systematic review of the literature. *Int Wound J*. 2016;13(6):1260-1281. 30. Shim HS, Choi JS, Kim SW. A Role for Postoperative Negative Pressure Wound Therapy in Multitissue Hand Injuries. *Biomed Res Int*. 2018;2018. 31. Kilpaal DV, Cunningham MR. Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system. *Wound Repair Regen*. 2011;19(5):588-596. 32. Ma Z, Shou K, Li Z, et al. Negative pressure wound therapy promotes vessel destabilization and maturation at various stages of wound healing and thus influences wound prognosis. *Exp Ther Med*. 2016;11(4):1307-1317. 33. Xia CY, Yu AX, Qi B, et al. Analysis of blood flow and local expression of angiogenesis associated growth factors in infected wounds treated with negative pressure wound therapy. *Mol Med Rep*. 2014;9(5):1749-1754. 34. Smith & Nephew 2018. PICO 14 Service Life Testing: 14 Day Device Lifespan. Internal report. RD/18/132. 35. Smith & Nephew 2018. The review of evidence supporting the use of PICO in wounds ≥2cm in depth. Internal report EO AWM. PCS230.001 v2 36. Smith & Nephew 2019. Research & Development Report. PICO 14 and PICO 7 Initial Pump Down and Maintenance Pump Down Time Outs RD/19/084. 37. Hudson DA, Adams KG, Van Huyssteen A, Martin R, Huddleston EM. Simplified negative pressure wound therapy: clinical evaluation of an ultraportable, no-canister system. *Int Wound J*. 2015;12(2):195-201. 38. Smith & Nephew 2018. Summary of routine QA testing on MVP of PICO dressings. 2018. Internal report. DS/18/153/R. 39. Payne C, Edwards D. 2014. Application of the single use negative pressure wound therapy device (PICO) on a heterogeneous group of surgical and traumatic wounds. *ePlasty*. 14. 40. Rossington A. 2015. A prospective, open, non-comparative, multicentre study to evaluate the functionality and dressing performance of a new negative pressure enhanced dressing (NPEd) in acute wounds. CT09/02, May. 41. Sharp E. Single use NPWT for the treatment of complex orthopaedic surgical and trauma wounds. *Journal of Wound Care*. 2013;22(10):55-59. 42. Stryja J, Staffa R, Řiřa D, Stryjová K, Nicelínková K. Cost-effectiveness of negative pressure wound therapy in outpatient setting. *Prolekare*. 2015. 43. Data on file 0810010 – The antimicrobial activity of ACTICOAT<sup>®</sup> and ACTICOAT FLEX 3 while under negative pressure; Lumb, H; 2011. 44. Smith & Nephew 2019. Summary report of *in vitro* Wound Model and Negative Pressure Delivery (Nominal -80mmHg) testing for PICO v2 (PICO 7 and PICO 14) System. Internal report. RD/18/134 V2 45. Helito CP, Sobrado MF, Giglio PN, Bonadio MB, Pecora JR, Demange MK, Gobbi RG. The use of negative-pressure wound therapy after total knee arthroplasty is effective for reducing complications and the need for reintervention. 46. Villeco JP. Edema: a silent but important factor. *Journal of Hand Therapy*. 2012 Apr 1;25(2):153-62. 47. Smith & Nephew 2018. Wound Model Testing: Project Opal PICO 15x20cm dressings with ACTICOAT FLEX 7, Moderate Flow – 4 days. Internal report. DS/17/448/R V2. 48. Smith & Nephew January 2018. Wound model testing of PICO 7 Project Opal, low flow with ACTICOAT FLEX 7 wound contact layer – 7 days. Internal report. DS/17/627/R V2. 49. Sharpe A, Myers D, Searle R. 2018. Using single use negative pressure wound therapy for patients with complicated diabetic foot ulcers: an economic perspective. *Wounds UK*, 14(4). 50. Dowsett C, Hampton K, Myers D, Styché T. Use of PICO to improve clinical and economic outcomes in hard-to-heal wounds. *Wounds International*. 2017;8(2):52-58. 51. Bozic KJ, Katz P, Cisternas M, Ono L, Ries MD, Showstack J. Hospital resource utilization for primary and revision total hip arthroplasty. *JBJS*. 2005 Mar 1;87(3):570-6.

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©Octubre de 2020 Smith+Nephew