

+ ALLEVYN[®] THIN Foam Dressing is a self-adhesive,^{1,2} flexible and conformable^{1,3,4} dressing designed for the management of low to moderate exuding wounds.

Smith+Nephew

ALLEVYN[®] THIN
Foam Dressing

Application



1

- Remove the release paper



2

- Apply the self-adhesive side of the dressing to the wound



3

- Ensure good contact
- Smooth the dressing into place, ensuring that the edges of the dressing are not wrinkled

Removal

To remove the ALLEVYN THIN Dressing, lift one corner and pull the dressing back on itself.

Frequency of change

The frequency of dressing change depends on the level of wound exudate. Change dressing when the secondary dressing rises or strikethrough occurs. The interval between dressing changes is one to several days.

A dressing change is required at least once a week.

The diagram below is for guidance only.

ALLEVYN THIN Dressing should be inspected frequently.



ALLEVYN THIN Dressing ordering codes

Size	Code	Carton
5cm x 6cm / 2in x 2 ³ / ₈ in	66047576	10
10cm x 10cm / 4in x 4in	66047578	5
15cm x 20cm / 6in x 8in	66047579	3

Advanced Wound Management
Smith & Nephew Medical Ltd
Croxley Park, Building 5,
Hatters Lane, Watford,
Hertfordshire WD18 8YE, UK

www.smith-nephew.com/allevyn
T +44 (0) 1923 477100
F +44 (0) 1923 477101

◇Trademark of Smith+Nephew
All Trademarks acknowledged
© January 2021 Smith+Nephew
AWM-AWC-26659 | GMC1249

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

References

1. Beiersdorf AG. Clinical evaluation of Cutinova™ Thin. Internal report. 2. Schmidt A, Gassmueller J, Hughes-Formella B, Bielfeldt S. Treating hypertrophic scars for 12 or 24 hours with a self-adhesive hydroactive polyurethane dressing. *J Wound Care*. 2001;10(5):149–153. 3. Smith + Nephew 2009. ALLEVYN THIN Dressing physical properties. Internal report. DS/09/013/R6. 4. Smith+Nephew 2007. A prospective, multi-centre, randomised and comparative evaluation to assess the relative cost effectiveness of treatment of Stage II pressure ulcers with ALLEVYN THIN and Saline soaked gauze. Internal report. SR/SR/CE/011/ALT.