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## **RI.HIP NAVIGATION**

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## **+ Evidence in focus**

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## Studio fondamentale

### Intraoperative validation of navigated limb measurements in THA using a pinless femoral array<sup>4</sup>

Ulivi M, Orlandini L, Pascale W, Consonni O, Sansone V. *J Arthroplasty*. 2014;29:1026–1029

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## Studio fondamentale

A new method of registration in navigated hip arthroplasty without the need to register the anterior pelvic plane<sup>18</sup>

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## Studio fondamentale

Impingement-free range of movement, acetabular component cover and early clinical results comparing 'femur-first' navigation and 'conventional' minimally invasive total hip arthroplasty: a randomised controlled trial<sup>23</sup>

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## Studio fondamentale

Visual intraoperative estimation of range of motion is misleading in minimally invasive total hip arthroplasty<sup>31</sup>

Woerner M, Weber M, Sendtner E, et al. *Arch Orthop Trauma Surg.* 2016;136:1015–1020

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## Studio fondamentale

A new system of computer-assisted navigation leading to reduction in operating time in uncemented total hip replacement in a matched population<sup>24</sup>

Chaudhry FA, Ismail SZ, Davis ET. *Eur J Orthop Surg Traumatol.* 2018;28:645–648

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## Studio fondamentale

Computer guided total hip arthroplasty is associated with a reduced risk of revision and increased patient satisfaction. An analysis of a single manufacturer acetabular components from the National Joint Registry of England, Wales, Northern Ireland and the Isle of Man<sup>26</sup>

**Davis ET, McKinney KD, Kamali A, Kuljaca S, Pagkalos J. Poster presented at: World Arthroplasty Congress (WAC) Virtual Meeting; April 22–24, 2021**

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
Sendtner  
et al.  
2011<sup>20</sup>

Accuracy of acetabular cup placement in computer-assisted, minimally-invasive THR in a lateral decubitus position<sup>20</sup>  
Sendtner E, Schuster T, Wörner M, Kalteis T, Grifka J, Renkawitz T. *Int Orthop.* 2011;**35**:809–815

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Rizzi et al.  
2012<sup>17</sup>

Acetabular cup positioning using computer navigation through direct anterior approach<sup>16</sup>  
Rizzi L, Gotti V, Castelli CC. *Hip Int.* 2012;22:431–432

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Renkawitz  
et al.  
2014<sup>20</sup>

Femoral pinless length and offset measurements during computer-assisted, minimally invasive total hip arthroplasty<sup>20</sup>  
**Renkawitz T, Sendtner E, Schuster T, Weber M, Grifka J, Woerner M. *J Arthroplasty*. 2014;29:1021–1025**

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Clavé et al.  
2015<sup>19</sup>

Comparison of the reliability of leg length and offset data generated by three hip replacement CAOS systems using EOS™ imaging<sup>19</sup>

Clavé A, Fazilleau F, Cheval D, Williams T, Lefèvre C, Stindel E. *Orthop Traumatol Surg Res.* 2015;101:647–653

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Femur first navigation can reduce impingement severity compared to traditional free hand total hip arthroplasty<sup>21</sup>

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## Bibliografia

1. National Joint Registry for England, Wales and Northern Ireland: 17th Annual Report. Available at: <http://www.njrcentre.org.uk>. Accessed 15 February 2021.
2. Beswick AD, Wylde V, Goberman-Hill R, Blom A, Dieppe P. What proportion of patients report long-term pain after total hip or knee replacement for osteoarthritis? A systematic review of prospective studies in unselected patients. *BMJ Open*. 2012;2(1):e000435.
3. NHS Digital. Finalised patient reported outcome measures (PROMs) in England for hip and knee replacement procedures (April 2019 to March 2020). Available at: <http://digital.nhs.uk/pubs/promsapr19mar20fin>. Accessed 18 May 2021.
4. Ulivi M, Orlandini L, Pascale W, Consonni O, Sansone V. Intraoperative validation of navigated limb measurements in THA using a pinless femoral array. *J Arthroplasty*. 2014;29:1026–1029.
5. Callanan MC, Jarrett B, Bragdon CR, et al. The John Charnley Award: risk factors for cup malpositioning: quality improvement through a joint registry at a tertiary hospital. *Clin Orthop Relat Res*. 2011;469:319–329.
6. Bosker BH, Verheyen CC, Horstmann WG, Tulp NJ. Poor accuracy of freehand cup positioning during total hip arthroplasty. *Arch Orthop Trauma Surg*. 2007;127:375–379.
7. Dorr LD, Malik A, Wan Z, Long WT, Harris M. Precision and bias of imageless computer navigation and surgeon estimates for acetabular component position. *Clin Orthop Relat Res*. 2007;465:92–99.
8. Feng JE, Anoushiravani AA, Eftekhary N, Wiznia D, Schwarzkopf R, Vigdorichik JM. Techniques for optimizing acetabular component positioning in total hip arthroplasty: defining a patient-specific functional safe zone. *JBJS Rev*. 2019;7:e5.
9. Maloney WJ, Keeney JA. Leg length discrepancy after total hip arthroplasty. *J Arthroplasty*. 2004;19:108–110.
10. Parvizi J, Sharkey PF, Bissett GA, Rothman RH, Hozack WJ. Surgical treatment of limb-length discrepancy following total hip arthroplasty. *J Bone Joint Surg Am*. 2003;85:2310–2317.
11. Renkawitz T, Weber T, Dullien S, et al. Leg length and offset differences above 5mm after total hip arthroplasty are associated with altered gait kinematics. *Gait Posture*. 2016;49:196–201.
12. Lewinnek GE, Lewis JL, Tarr R, Compere CL, Zimmerman JR. Dislocations after total hip-replacement arthroplasties. *J Bone Joint Surg Am*. 1978;60:217–220.
13. Abdel MP, von Roth P, Jennings MT, Hanssen AD, Pagnano MW. What safe zone? The vast majority of dislocated THAs are within the lewinnek safe zone for acetabular component position. *Clin Orthop Relat Res*. 2016;474:386–391.
14. Schwarz T, Benditz A, Springorum HR, et al. Assessment of pelvic tilt in anteroposterior radiographs by means of tilt ratios. *Arch Orthop Trauma Surg*. 2018;138:1045–1052.
15. Zhu J, Wan Z, Dorr LD. Quantification of pelvic tilt in total hip arthroplasty. *Clin Orthop Relat Res*. 2010;468:571–575.
16. Sendtner E, Schuster T, Worner M, Kalteis T, Grifka J, Renkawitz T. Accuracy of acetabular cup placement in computer-assisted, minimally-invasive THR in a lateral decubitus position. *Int Orthop*. 2011;35:809–815.
17. Rizzi L, Gotti V, Castelli CC. Acetabular cup positioning using computer navigation through direct anterior approach. *Hip Int*. 2012;22:431–432.
18. Davis ET, Schubert M, Wegner M, Haimerl M. A new method of registration in navigated hip arthroplasty without the need to register the anterior pelvic plane. *J Arthroplasty*. 2015;30:55–60.
19. Clave A, Fazilleau F, Cheval D, Williams T, Lefevre C, Stindel E. Comparison of the reliability of leg length and offset data generated by three hip replacement CAOS systems using EOS™ imaging. *Orthop Traumatol Surg Res*. 2015;101:647–653.
20. Renkawitz T, Sendtner E, Schuster T, Weber M, Grifka J, Woerner M. Femoral pinless length and offset measurements during computer-assisted, minimally invasive total hip arthroplasty. *J Arthroplasty*. 2014;29:1021–1025.
21. Palit A, Williams MA, Turley GA, Renkawitz T, Weber M. Femur first navigation can reduce impingement severity compared to traditional free hand total hip arthroplasty. *Sci Rep*. 2017;7:7238.
22. Weber M, Woerner ML, Sendtner E, Vollner F, Grifka J, Renkawitz TF. Even the intraoperative knowledge of femoral stem anteversion cannot prevent impingement in total hip arthroplasty. *J Arthroplasty*. 2016;31(11):2514–2519.
23. Renkawitz T, Weber M, Springorum H, et al. Impingement-free range of movement, acetabular component cover and early clinical results comparing 'femur-first' navigation and 'conventional' minimally invasive total hip arthroplasty: a randomised controlled trial. *Bone Joint J*. 2015;97:890–898.
24. Chaudhry FA, Ismail SZ, Davis ET. A new system of computer-assisted navigation leading to reduction in operating time in uncemented total hip replacement in a matched population. *Eur J Orthop Surg Traumatol*. 2018;28:645–648.
25. Sugano N, Takao M, Sakai T, Nishii T, Miki H. Does CT-Based navigation improve the long-term survival in ceramic-on-ceramic THA? *Clinical Orthopaedics and Related Research*. 2012;470:3054–3059.
26. Davis ET, McKinney KD, Kamali A, Kuljaca S, Pagkalos J. Computer guided total hip arthroplasty is associated with a reduced risk of revision and increased patient satisfaction. An analysis of a single manufacturer acetabular components from the National Joint Registry of England, Wales, Northern Ireland and the Isle of Man. Poster presented at: World Arthroplasty Congress (WAC) Virtual Meeting; April 22–24, 2021.
27. National Joint Registry for England, Wales and Northern Ireland: POLARSTEM Cementless (Oxinium/XLPE/R3 cup) implant summary report. 14 August 2019. Data on file. [Link: [http://bit.ly/POLAR3\\_Aug2019](http://bit.ly/POLAR3_Aug2019)]
28. Orthopaedic Data Evaluation Panel (ODEP). Available at <http://www.odep.org.uk>. Accessed 17 May 2021.
29. Smith+Nephew 2020. Navigation assisted THA vs conventional THA: frequently reported study outcomes and key differences in outcomes between the techniques – EA\_RECON\_NAVIO\_004\_v2.
30. Keshmiri A, Schroter C, Weber M, Craiovan B, Grifka J, Renkawitz T. No difference in clinical outcome, bone density and polyethylene wear 5-7 years after standard navigated vs conventional cement free total hip arthroplasty. *Arch Orthop Trauma Surg*. 2015;135(5):723–730.
31. Woerner M, Weber M, Sendtner E, et al. Visual intraoperative estimation of range of motion is misleading in minimally invasive total hip arthroplasty. *Arch Orthop Trauma Surg*. 2016;136:1015–1020.
32. Parratte S, Ollivier M, Lunebourg A, Flecher X, Argenson JN. No Benefit after THA performed with computer-assisted cup placement: 10-year results of a randomized controlled study. *Clin Orthop Relat Res*. 2016;474:2085–2093.

**Fabbricante**  
**Smith & Nephew Inc.**  
1450 Brooks Road  
Memphis, TN 38116  
U.S.A.

[www.smith-nephew.com](http://www.smith-nephew.com)

**Contatto**  
**Smith & Nephew S.r.l.**  
Via De Capitani 2A  
20864 Agrate Brianza MB  
Italia

[www.smith-nephew.it](http://www.smith-nephew.it)  
T +39 039 60941  
F +39 039 651535

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©Luglio 2021 Smith+Nephew