


JOURNEY[®] II Bi-Cruciate Stabilised (BCS) OXINIUM[®] Total Knee System (TKS) demonstrates high long-term survivorship with lower revision risk versus other TKA implants in the AOANJRR¹

- The AOANJRR annual report includes the revision rate for JOURNEY II BCS OXINIUM TKS both with and without patella resurfacing (6.3% at 5 years)²
- Data from the AOANJRR (up to January 2025) were obtained via the Automated Industry Report System (AIRS) to compare JOURNEY II BCS OXINIUM TKS versus all other TKAs when the patella was resurfaced¹
- This independent analysis includes 10 years of data from 3,430 JOURNEY II BCS OXINIUM TKS procedures¹
- Revision risk, expressed as a hazard ratio (HR), has been adjusted for variables such as age and gender¹

10 years Risk of revision with JOURNEY II BCS OXINIUM TKS is **50% lower than other TKAs in the registry** when the patella is resurfaced (Adjusted HR: 0.50 [95% CI: 0.33–0.76]; p=0.001)¹



Number at Risk	0 Yr	1 Yr	2 Yrs	3 Yrs	4 Yrs	5 Yrs	6 Yrs	7 Yrs
Journey II Oxinium/Journey	3430	2207	956	174	110	107	106	103
Other Total Knee	627948	575477	514006	463179	411597	364915	318165	274587
Number at Risk	8 Yrs	9 Yrs	10 Yrs	11 Yrs	12 Yrs	13 Yrs	14 Yrs	15 Yrs
Journey II Oxinium/Journey	98	98	46	7	0	0	0	0
Other Total Knee	234516	197973	165085	136280	111237	89540	70690	54955
Number at Risk	16 Yrs	17 Yrs	18 Yrs	19 Yrs	20 Yrs	21 Yrs	22 Yrs	23 Yrs
Journey II Oxinium/Journey	0	0	0	0	0	0	0	0
Other Total Knee	42346	32025	23486	16841	11557	7197	3965	1568

Figure 1. Cumulative percentage revision of primary TKA by model (all diagnoses). (Figure 1, AOANJRR AIRS report 08/01/2025)

10 years **98% survivorship** for JOURNEY II BCS OXINIUM TKS when the patella is resurfaced¹

Table 1. Yearly cumulative percent revision of primary TKA by model (all diagnoses). (Table 12, AOANJRR AIRS report 08/01/2025)

CPR	1 Yr	2 Yrs	3 Yrs	4 Yrs	5 Yrs	
Journey II Oxinium/Journey	0.4 (0.2, 0.7)	0.8 (0.5, 1.3)	1.1 (0.7, 1.9)	2.0 (0.8, 4.9)	2.0 (0.8, 4.9)	
Other Total Knee	0.9 (0.9, 0.9)	1.6 (1.5, 1.6)	2.1 (2.0, 2.1)	2.4 (2.4, 2.5)	2.7 (2.7, 2.8)	
CPR	6 Yrs	7 Yrs	8 Yrs	9 Yrs	10 Yrs	11 Yrs
Journey II Oxinium/Journey	2.0 (0.8, 4.9)	2.0 (0.8, 4.9)	2.0 (0.8, 4.9)	2.0 (0.8, 4.9)	2.0 (0.8, 4.9)	
Other Total Knee	3.0 (3.0, 3.0)	3.3 (3.2, 3.3)	3.5 (3.5, 3.6)	3.8 (3.7, 3.9)	4.1 (4.0, 4.1)	4.4 (4.3, 4.4)
CPR	12 Yrs	13 Yrs	14 Yrs	15 Yrs	16 Yrs	17 Yrs
Journey II Oxinium/Journey						
Other Total Knee	4.7 (4.6, 4.7)	5.0 (4.9, 5.1)	5.3 (5.2, 5.4)	5.6 (5.5, 5.7)	6.0 (5.9, 6.1)	6.3 (6.2, 6.4)
CPR	18 Yrs	19 Yrs	20 Yrs	21 Yrs	22 Yrs	23 Yrs
Journey II Oxinium/Journey						
Other Total Knee	6.6 (6.5, 6.8)	6.9 (6.8, 7.1)	7.2 (7.1, 7.4)	7.5 (7.3, 7.7)	7.7 (7.5, 8.0)	7.8 (7.6, 8.1)

Conclusions

At 10 years follow-up, JOURNEY II BCS OXINIUM TKS results in high survivorship and lower risk of revision than the average for other TKAs in the AOANJRR (both groups had patella resurfaced).¹

This material is intended for healthcare professionals. For detailed product information, including indications for use, contraindications, warnings and precautions, please consult the product's Instructions for Use (IFU) prior to use.

Abbreviations: AIRS = automated industry report system; AOANJRR = Australian Orthopaedic Association National Joint Replacement Registry; BCS = bi-cruciate stabilised; HR = hazard ratio.

References: 1. Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR), Automated Industry Report System (AIRS), ID No.16955* for Smith & Nephew Australia, JOURNEY II OXINIUM Total Knee (Procedures from 1 September 1999 – 7 January 2025), Accessed 31 January 2025, report generated January 8, 2025, AOA, Adelaide: 1–16. 2. Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR) Hip, Knee and Shoulder Arthroplasty: 2024 Annual Report, AOA, Adelaide; 2024:1–629. Available at: <https://aoanjrr.sahmri.com/annual-reports-2024>. Accessed 03 February 2025.

*AOANJRR is confident in the accuracy of the data included in this report, at the time it was provided. However, it was generated using an automated reporting system and has not been reviewed by the AOANJRR personnel.