
Smith+Nephew

OR30[◇]

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

Designed with safety in mind – Corrosion avoidance

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OXINIUM femoral heads are associated with decreased corrosion damage compared to CoCrMo femoral heads.

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OR3O[◇] Dual Mobility uses clinically Proven design elements

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Wettability and friction – Reducing wear potential in the joint

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Wettability test results Oxidized Zirconium (OxZr) vs CoCrMo

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OXINIUM[◇] DH – A next-generation Advanced bearing material for hip arthroplasty

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<https://aoanjrr.sahmri.com/>

** 2021 Australian Orthopaedic Association National Joint Replacement Registry

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
Notes

References

1. Parikh, A; Hill, P; Sprague, J. "Long-term Simulator Wear Performance of an Advanced Bearing Technology for THA." Poster Presented at: ORS 2013 Annual Meeting; Poster No: 1028 2. Parikh, A; Weaver, C. "Wear Testing of a New Modular Dual Mobility System." Poster Presented at: ORS 2019 Annual Meeting; Poster No: 1041 3. Boyer, B; Philippot, R; Geringer J; Farizon F. "Primary total hip arthroplasty with dual mobility socket to prevent dislocation: a 22-year follow-up of 240 hips." International Orthopaedics (SICOT) , 36, 2012, pp. 511-518 4. Bauchu, P; Bonnard, O; Cypres, A; Fiquet A; Girardin, P; Noyer, D. "The Dual -Mobility POLARCUP: First Results From a Multicenter Study." ORTHOPEDICS, 31, 2008, pp. 97-99 5. Reina, N; Pareek, A; Krych, AJ; Pagnano, MW; Berry, DJ; Abdel MP. "Dual-Mobility Constructs in Primary and Revision Total Hip Arthroplasty: A Systematic Review of Comparative Studies." The Journal of Arthroplasty, 34, 2019, pp. 594-603 6. Heffernan, C; Bhimji, S; Macintyre, J; Coustance, A; Ron, M; Markel, D; Cabanela, M; Nevelos, J. "Development and Validation of a Novel Modular Dual Mobility Hip Bearing." Poster Presented at: ORS 2011; Poster No: 1165 7. Fabry, C; Langlois, J; Hamadouche, M; Bader, R. "Intra-prosthetic dislocation of dual-mobility cups after total hip arthroplasty: potential causes from a clinical and biomechanical perspective." International Orthopaedics, 40, 2016, pp. 901-906 8. Addona, JL; Gu, A; De Martino, I; Malahias, MA; Sculco, TP; Sculco PK. "High Rate of Early Intraprosthetic Dislocations of Dual Mobility Implants: A Single Surgeon Series of Primary and Revision Total Hip Replacements." The Journal of Arthroplasty, 34, 2019, pp. 2793-2798 9. Nam, D; Salih, R; Brown, KM; Nunley, RM; Barrack RL. "Metal Ion Levels in Young, Active Patients Receiving a Modular, Dual Mobility Total Hip Arthroplasty." The Journal of Arthroplasty, 32, 2017, pp. 1581-1585 10. Matsen Ko, LJ; Pollag, KE; Yoo, JY; Sharkey, PF. "Serum Metal Ion Levels Following Total Hip Arthroplasty With Modular Dual Mobility Components." The Journal of Arthroplasty, 31, 2016, pp. 186-189 11. Zachwieja, E; Sharkey, PF. "Severe Corrosion of a Modular Dual Mobility Acetabular Component." International Congress For Joint Reconstruction, Abstract Presented at: Rothman Institute for Orthopedics Grand Rounds, January 29, 2020 12. Padgett, DE; Romero, J; Wach, A; Wright, TM. "Incidence of Liner Malseating in Dual Mobility Implants." Abstract Presented at, The Hip Society Summer Meeting, 2019 13. Di Laura, A; Hothi, H; Battisti, C; Cerquiglini, A; Henckel, J; Skinner, J; Hart, A. "Wear of dual-mobility cups: a review article." International Orthopaedics, 41, 2017, pp. 625-633 14. Fabry, C; Kaehler, M; Herrmann, S; Woernle, C; Bader, R. "Dynamic behavior of tripolar hip endoprostheses under physiological conditions and their effect on stability." 15. Philippot, R; Boyer, B; Farizon, F. "Intraprosthetic Dislocation: A Specific Complication of the Dual-mobility System." Clinical Orthopaedics and Related Research, 471, 2013, pp. 965-970 16. Salehi, A; Hunter, G. "Laboratory and clinical performance of oxidized zirconium alloy." Clinical Cases in Mineral and Bone Metabolism, 7(3), 2010, pp. 163-175 17. Brown, TD; Elkins, JM; Pedersen, DR; Callaghan, JJ. "Impingement and dislocation in total hip arthroplasty: mechanisms and consequences." The Iowa Orthopaedic Journal, 34, 2012, 1-15 18. Range of motion and jump distance memos 19. Parikh, A; Weaver, C; Pawar, V. "Electrochemical Evaluation of Zr2.5Nb, Diffusion Hardened Oxidized Zr2.5Nb, and CoCrMo." Poster Presented at: ORS 2020, Poster No.: 0528 20. Cartner, J; Aldinger, P; Li, C; Collins, D. "Characterization of Femoral Head Taper Corrosion Features Using a 22-Year Retrieval Database." The Musculoskeletal Journal of Hospital for Special Surgery, 13, 2017, pp. 35-41 21. Hampton, C; Weitzler, L; Baral, E; Wright TM; Bostrom MPG. "Do oxidized zirconium heads decrease tribocorrosion in total hip arthroplasty." Bone Joint J 101-B, 2019, 22. Parikh, A; Pawar, V; Sprague, J. "Long-term Simulator Wear Performance of an Advanced Bearing Technology for THA." Poster Presented at: ORS 2013, Poster No.:1028 23. Parikh, A; Weaver, C. "Wear testing of a new modular dual mobility system." Poster Presented at: ORS 2019, Poster No.: 24. Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR). Hip, Knee & Shoulder Arthroplasty: 2019 Annual Report. Adelaide: AOA, 2019. 25. Peters RM, Van Steenberg LN, Stevens M, Rijk PC, Bulstra SK, Zijlstra WP. The effect of bearing type on the outcome of total hip arthroplasty. Acta Orthop. 2018;89:163-169. 26. Atrey A, Ancarani C, Fitch D, Bordini B. Impact of bearing couple on long-term component survivorship for primary cementless total hip replacement in a large arthroplasty registry. Poster presented at: Canadian Orthopaedic Association; June 20-23, 2018; Victoria, British Columbia, Canada. 27. Parikh, A; Weaver, C; Pawar, V. "Evaluation of Dislocation Damage Resistance using a Novel Test Method." Poster Presented at: ORS 2017, Poster No.: 1070 28. Pawar, V; Weaver, C; Jani, S. "Physical characterization of a new composition of oxidized zirconium - 2.5wt% niobium produced using a two step process for biomedical applications." Applied Surface Science, 257, 2011, pp.6118-6124 29. Hunter, G; Marc, L. "Abrasive Wear of Oxidized Zr-2.5Nb, CoCrMo, and Ti-6AL-4V Against Bone Cement." Transactions of the World Biomaterials Conference. Kamuela, HI : Society for Biomaterials, 2000. 835. 30. Bingenheimer, H; Li, C. "Fatigue Test of a Cross-linked Polyethylene Dual Mobility Insert against Diffusion Hardened Oxidized Zirconium Liner." Poster Presented at: ORS 2019, Poster No.: 1036 31. Rister, D; Parikh, A. "Neck Impingement Testing of a Novel Hard on Hard Hip Bearing." Poster Presented at: 2015 Society for Biomaterials, Poster No.: 527 32. Bingenheimer, H. "Deformation and Roundness Analysis of Diffusion Hardened Oxidized Zirconium Dual Mobility Liners in R3 Acetabular Shells." Poster Presented at: ORS 2019, Poster No.: 1905 33. Davis et al. SAFETY AND EFFICACY OF CERAMIC ON CERAMIC TOTAL HIP REPLACEMENT, World Arthroplasty Congress 2018. 34. Data on file. 35. Lee, YK; Kim, KC; Jo, WL; Ha, YC; Parvizi, J; Koo, KH. "Effect of Inner Taper Angle of Acetabular Metal Shell on the Malseating and Dissociation Force of Ceramic Liner." The Journal of Arthroplasty, 32, 2017, pp. 1360-1362 36. Liner removal memo 37. Latest ratings can be found at www.odep.org.uk (Accessed on April 1 2020). 38. Naudie, DDR; Somerville, L; Korczak, A; Yuan, X; McCalden, RW; Holdsworth, D; Bourne, RB. "A Randomized Trial Comparing Acetabular Component Fixation of Two Porous Ingrowth Surfaces Using RSA." The Journal of Arthroplasty, 28, 2013, pp. 48-52 39. Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR). Hip, Knee & Shoulder Arthroplasty: 2019 Annual Report. Adelaide: AOA, 2020. 40. Serum Metal Ions in Contemporary Monoblock and Modular Dual Mobility Articulations Alexander Greenberg, MD , *, Allina Nocon, PhD, MPH a, Ivan De Martino, MD David J. Mayman, MD , Thomas P. Sculco, MD , Peter K. Sculco, MD 41. Fretting and Corrosion at the Backside of Modular Cobalt Chromium Acetabular Inserts: A Retrieval Analysis T. David Tarity, MD , Chelsea N. Koch, BS , *, Jayme C. Burket, PhD, Timothy M. Wright, PhD , Geoffrey H. Westrich, MD 42. Adverse reaction to metal debris due to fretting corrosion between the acetabular components of modular dual-mobility constructs in total hip replacement: a systematic review and meta-analysis Jonathan M. R. French Paul Bramley Sean Scattergood Nemandra A. Sandiford

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© 2022 Smith+Nephew
®Trademark of Smith+Nephew