

Evidence in focus

A systematic literature review and meta-analysis

Adults

High treatment success rate
with the TAYLOR SPATIAL FRAME
in adults

 **smith&nephew**
**TAYLOR SPATIAL
FRAME[®]**

External Fixator

Supporting healthcare professionals




Purpose

To systematically evaluate the available evidence to determine the overall treatment success rate of TAYLOR SPATIAL FRAME® (TSF) in adults with acute trauma, non-unions/mal-unions and deformities.

Background

TSF is an external device for limb correction, lengthening and/or straightening, with a long history of clinical use:



20

More than 20 years
of clinical use



200


More than 200 publications
detailing the use of TSF in
adults and children




Methods

Literature search


A search for clinically relevant results was conducted using Embase and PubMed across three indications (September 6, 2018):



Acute trauma



Non-unions/mal-unions



Deformity correction

Study suitability

Abstracts were analysed to determine study relevance. Additional studies were identified from other sources, such as by reviewing reference lists. To be considered eligible, a study had to fulfil the following criteria:

Inclusion criteria:

- Published from 2008 onwards in a peer-reviewed journal
- English language publication
- Adult population
- Proportion of successful cases identifiable in study

Exclusion criteria:

- Single case report
- Off-label product use

Only studies with >10 patients in the TAYLOR SPATIAL FRAME® treatment group were included in the meta-analyses (Figure 1). Studies with 2-10 patients are reported in the Appendices.




Figure 1. Search strategy

Meta-analyses

The results of each suitable study were analysed to determine the proportion of patients treated with TSF who successfully achieved the treatment goal.

The goals used to indicate treatment success were:

- **Consolidation** in patients with acute trauma
- **Bony union** in patients with mal-unions/non-unions
- **Deformity correction**

Meta-analyses were then conducted to determine the overall success rate of TSF per indication.



Results

Total number of studies meeting the inclusion criteria with ≥ 10 patients



15 studies

Combined treatment success

Across all three indications in adult patients, the meta-analyses demonstrated consistently high treatment success rates (Figure 2).





Figure 2. Combined treatment success in adult patients treated with TSF.

Full details of studies included in the meta-analysis are included in the Appendices.



Conclusion

The TAYLOR SPATIAL FRAME[®] has a long history of clinical use and has been reported in more than 200 peer-reviewed publications. These meta-analyses demonstrate consistently high success rates in adult patients for the treatment of acute trauma, mal-unions/non-unions and deformities.



Appendix 1. Literature review and meta-analysis in acute trauma

Table 1. Characteristics of relevant studies.

Study, year	Evidence Level				n	Age (years)	Indication
	Level I: Randomised controlled trials	Level II: Prospective, comparative	Level III: Retrospective, comparative	Level IV: Case series			
n ≥ 10; included in meta-analysis							
Ahearn et al, 2014 ¹					21	Mean: 44* Range: 17-78*	Unstable bicondylar tibial plateau fractures
Menakaya et al, 2014 ²					37	Mean: 45* Range: NR*	Various
O'Neill et al, 2016 ³					15 fractures	Mean: 39* Range: 16-79*	Various
Rampurada et al, 2008 ⁴					26	Mean: 40 Range: 22-59	Tibial plateau and pilon fractures
Sala et al, 2017 ⁵					20	Mean: 37 Range: 11-72	Open supracondylar-intracondylar femoral fractures
n < 10; not included in meta-analysis							
Lahoti et al, 2013 ⁶					7	Mean: 38 Range: 15-70 [†]	Various
Sharma and Nunn, 2013 ⁷					2	Mean: 54 Range: 48-60	Open tibial fractures

Forest plot for consolidation in acute trauma




Figure 3. Proportional meta-analysis of studies (with ≥10 patients) assessing the use of TSF for acute trauma


* Mean age and range of overall patient population.

[†] Data for adults and children not separated but mean age >18 years.

[‡] Data shown is for number of fractures successfully treated.

Abbreviations

CI = confidence interval; NR = not reported



Appendix 2. Literature review and meta-analysis in non-unions/mal-unions

Table 2. Characteristics of relevant studies.

Study, year	Level I: Randomised controlled trials	Level II: Prospective, comparative	Level III: Retrospective, comparative	Level IV: Case series	n	Age (years)	Indication
n≥10; included in meta-analysis							
Arvensen et al, 2017 ⁸					37	Mean: 50 Range: NR	Deformities associated with distal tibial non-unions
Khunda et al, 2016 ⁹					40	Mean: 40 Range: 9-69*	Non-unions
Napora et al, 2017 ¹⁰					75	Mean: 46 Range: NR†	Infected posttraumatic non-unions of the tibia
Rozbruch et al, 2008 ¹¹					38	Mean: 43 Range: 8-72*	Tibial non-unions
n<10; not included in meta-analysis							
Schoenleber & Hutson, 2015 ¹²					5	Mean: 38† Range 22-51†	Non-unions/mal-unions

Forest plot for bony union in mal-unions/non-unions





Figure 4. Proportional meta-analysis of studies (with ≥10 patients) assessing the use of TSF for non-unions or mal-unions

* Data for adults and children not separated but mean age >18 years.

† Mean age of overall patient population.

Abbreviations

CI = confidence interval; NR = not reported



Appendix 3. Literature review and meta-analysis in deformity correction

Table 3. Characteristics of relevant studies.

Study, year	Level I: Randomised controlled trials	Level II: Prospective, comparative	Level III: Retrospective, comparative	Level IV: Case series	n	Age (years)	Indication
n≥10; included in meta-analysis							
Ashfaq et al, 2012 ¹³					57*	Mean: 39 Range: 21-72	Proximal tibia varus
Horn et al, 2011 ¹⁴					52	Mean: 44 Range 18-79	Various
Nakase et al, 2009 ¹⁵					10	Mean: 29 Range: 10-71 [†]	Various
Rozbruch et al, 2010 ¹⁶					102 (122 tibiae)	Mean: 39 Range: 5-72 [†]	Various
Sokucu et al, 2013 ¹⁷					37 (50 limbs)	Mean: 23 Range: 10-58 [†]	Various deformities around the knee
Thiryayi et al, 2010 ¹⁸					10	Mean: 59 Range: 48-71	Various
n<10; not included in meta-analysis							
Manggala et al, 2017 ¹⁹					7	Mean: 43 Range: 18-63	Foot and ankle deformities
Baumgartner and Weber, 2017 ²⁰					6	Mean: NR Range: NR	Post-infectious bony malalignments
Docquier et al, 2008 ²¹					2	Mean: 27 Range 24-31	Various
Robinson et al, 2011 ²²					9	Mean: 49 Range: 37-59	OA deformities
Tawari et al, 2018 ²³					2	Mean: 64 Range: 63-64	Tibial deformities

Forest plot for deformity correction




Figure 5. Proportional meta-analysis of studies (with ≥10 patients) assessing the use of TSF for deformity correction in adults

* Two n values are given in this study. n=55 is described in the Materials & Methods. n=57 is calculated from treatment success in the results section and has therefore been used in the meta-analysis.

[†] Data for adults and children not separated but mean age >18 years.

[‡] Data reported as number of limbs rather than patients.

Abbreviations

CI = confidence interval; NR = not reported; OA = osteoarthritis



References

1. Ahearn N, Oppy A, Halliday R, et al. The outcome following fixation of bicondylar tibial plateau fractures. *Bone Joint J.* 2014;96-B(7):956-962.
2. Menakaya C, Rigby A, Hadland Y, Barron E, Sharma H. Fracture healing following high energy tibial trauma: Ilizarov versus Taylor Spatial Frame. *Ann R Coll Surg Engl.* 2014;96(2):106-110.
3. O'Neill B, Fox C, Molloy A, O'hEireamhoin S, Moore D. The use of circular external fixators in the management of lower limb trauma in Dublin: a single surgeon's 20-year experience. *Ir J Med Sci.* 2016;185(1):133-138.
4. Rampurada A, Madan S, Tadross T. Treatment of complex tibial plateau and distal tibial fractures with Taylor spatial frame: experience in a district general hospital. *Eur J Orthop Surg Traumatol.* 2008;18:521-524.
5. Sala F, Thabet A, Capitani P, Bove F, Abdelgawad A, Lovisetti G. Open supracondylar-intercondylar fractures of the femur treatment with Taylor spatial frame. *J Orthop Trauma.* 2017;31(10):546-553.
6. Lahoti O, Findlay I, Shetty S, Abhishetty N. Intentional deformation and closure of soft tissue defect in open tibial fractures with Taylor spatial frame: 'sat nav' technique. *J Orthop Trauma.* 2013;27(8):451-456.
7. Sharma H, Nunn T. Conversion of open tibial IIIb to IIIa fractures using intentional temporary deformation and the Taylor Spatial Frame. *Strategies Trauma Limb Reconstr.* 2013;8(2):133-140.
8. Arvesen J, Watson J, Israel H. Effectiveness of treatment for distal tibial nonunions with associated complex deformities using a hexapod external fixator. *J Orthop Trauma.* 2017;31(2):e43-e48.
9. Khunda A, Al-Maiyah M, Eardley W, Montgomery R. The management of tibial fracture non-union using the Taylor Spatial Frame. *J Orthop.* 2016;13(4):360-363.
10. Napora J, Weinberg D, Eagle B, Kaufman B, Sontich J. Hexapod frame stacked transport for tibial infected nonunions with bone loss: analysis of use of adjunctive stability. *J Orthop Trauma.* 2017;31(7):393-399.
11. Rozbruch S, Pugsley J, Fragomen A, Ilizarov S. Repair of tibial nonunions and bone defects with the Taylor Spatial Frame. *J Orthop Trauma.* 2008;22(2):88-95.
12. Schoenleber S, Hutson J. Treatment of hypertrophic distal tibia nonunion and early malunion with callus distraction. *Foot Ankle Int.* 2015;36(4):400-407.
13. Ashfaq K, Fragomen A, Nguyen J, Rozbruch S. Correction of proximal tibia varus with external fixation. *J Knee Surg.* 2012;25(5):375-384.
14. Horn D, Fragomen A, Rozbruch S. Supramalleolar osteotomy using circular external fixation with six-axis deformity correction of the distal tibia. *Foot Ankle Int.* 2011;32(10):986-993.
15. Nakase T, Kitano M, Kawai H, et al. Distraction osteogenesis for correction of three-dimensional deformities with shortening of lower limbs by Taylor Spatial Frame. *Arch Orthop Trauma Surg.* 2009;129(9):1197-1201.
16. Rozbruch S, Segal K, Ilizarov S, Fragomen A, Ilizarov G. Does the Taylor Spatial Frame accurately correct tibial deformities? *Clin Orthop Relat Res.* 2010;468(5):1352-1361.
17. Sokucu S, Karakoyun O, Arikan Y, Kucukkaya M, Kabukcuoglu Y. Efficacy of the Taylor spatial frame in the treatment of deformities around the knee. *Acta Orthop Traumatol Turc.* 2013;47(2):86-90.
18. Thiryayi W, Naqui Z, Khan S. Use of the Taylor spatial frame in compression arthrodesis of the ankle: a study of 10 cases. *J Foot Ankle Surg.* 2010;49(2):182-187.
19. Manggala Y, Angthong C, Primadhi A, Kungwan S. The deformity correction and fixator-assisted treatment using Ilizarov versus Taylor spatial frame in the foot and ankle. *Orthop Rev (Pavia).* 2018;9(4):114-117.
20. Baumgartner H, Weber A. Segmental bone transport and deformity correction in patients with osteomyelitis using the Taylor spatial frame (TSF). *Archives of Orthopaedic and Trauma Surgery.* 2017;137(11):1602.
21. Docquier P, Rodriguez D, Mousny M. Three-dimensional correction of complex leg deformities using a software assisted external fixator. *Acta Orthop Belg.* 2008;74(6):816-822.
22. Robinson P, Papanna M, Somanchi B, Khan S. High tibial osteotomy in medial compartment osteoarthritis and varus deformity using the Taylor spatial frame: early results. *Strat Traum Limb Recon.* 2011;6(3):137-145.
23. Tawari G, Maheshwari R, Madan S. Extra-articular deformity correction using Taylor spatial frame prior to total knee arthroplasty. *Strategies Trauma Limb Reconstr.* 2018;13(3):179-184.