

✚ Evidence in focus

Interactive collection of evidence

OXINIUM[®] Oxidized Zirconium -
Materials science

June 2022

Smith+Nephew



OXINIUM Oxidized Zirconium - An Overview


What is OXINIUM Technology?

Oxidized zirconium (**OXI**dized zirco**NIUM**) is a patented, advanced bearing material for hip and knee arthroplasty

- Alloy of 97.5% zirconium and 2.5% niobium
- Virtually no nickel, cobalt or chromium^{1,2}

How is OXINIUM Technology manufactured?

The zirconium alloy core undergoes a transformation at its surface, upon heating, to form a ~5µm thick ceramicized surface³



Key properties of OXINIUM Technology

Wear resistance*



OXINIUM Technology demonstrates reduced polyethylene wear compared with cobalt chrome (CoCr),⁴ and similar wear to ceramic at 45 million cycles⁵

Compared with CoCr, OXINIUM Technology is also:

- harder^{6†} and more abrasion resistant^{7†}
- more wettable^{8†} and has lower friction^{9†}

Strength and stability




- OXINIUM Technology behaves with the strength of a metal¹⁰ and has been shown to withstand fracture under application of strong force¹¹
- Zero recorded instances of brittle fracture in total joint arthroplasty in the published literature¹²

Biocompatibility



OXINIUM Technology (zirconium oxide and zirconium alloy) demonstrates reduced cytotoxicity and cellular inflammatory response[‡] in vitro compared with CoCr and other materials¹³

Corrosion and fretting resistance



OXINIUM Technology shows lower corrosion damage compared with CoCr^{14,15} and ceramic^{16,17}

OXINIUM Technology combines the strength of metal with the wear resistance and biocompatibility of ceramic and corrosion resistance better than both⁴⁻¹⁷

*The results of in vitro wear simulation testing have not been proven to quantitatively predict clinical wear performance. †Material property refers specifically to the oxide surface of OXINIUM.

‡The results of in vitro cytokine expression analyses have not been proven to quantitatively predict clinical cytokine expression.

Key evidence

● Knee
 ● Hip
 ● Material

Click on the **arrow** by the study name to view the study summary, and **click on the property icons** to the right to view the corresponding property summary

Ezzet K, et al. 2012 ¹⁸	✓			
Good V, et al. 2003 ¹⁹	✓			
Morrison M, et al. 2015 ²⁰	✓			
Papannagari R, et al. 2011 ²¹	✓			
Parikh A, et al. 2009 ²²	✓			
Parikh A, et al. 2013 ⁴	✓			
Ries MD, et al. 2002 ²³	✓			
Hobbs LW, et al. 2005 ²⁴		✓		
Sprague J, et al. 2004 ¹¹		✓		
Tsai S, et al. 2001 ¹⁰		✓		
Caicedo M, et al. 2014 ²⁵			✓	
Rose SF, et al. 2012 ²⁶			✓	
Cartner J, et al. 2017 ¹⁴				✓
Hampton C, et al. 2019 ¹⁵				✓
Pawar V, et al. 2005 ²⁷				✓

*Particles tested in these studies were derived from diffusion-hardened OXINIUM Oxidized Zirconium material.

Further studies

Interactivity has been included throughout the collection of evidence to aid navigation

Menu

will bring you back to this menu page

Banners at the bottom of each page will take you to previous or subsequent pages, as indicated

Key evidence is represented in **turquoise** and further evidence in **brown**



Property summary:

Wear resistance

Wear resistance refers to the susceptibility of implant components to wear processes such as abrasion, scratching and fracture. Component and liner wear-related complications have a significant impact on implant longevity and are shown to be a major cause of revision and reoperation^{4,19,28}



OXINIUM Oxidized Zirconium shows reduced in vitro wear compared with CoCr, as evidenced by:^{4,7,18,19,21-23,29-31*}

Lower wear rates following articulation against polyethylene liners^{4,18,19,22,23,29,31}




Figure 1. Wear rate (mm³/million cycles) of OXINIUM and CoCr femoral heads articulated against highly cross-linked polyethylene (XLPE) liners

OXINIUM femoral heads exhibit a lower wear rate than ceramic heads in vivo:³²




Figure 2. Linear wear rate of OXINIUM and ceramic femoral heads articulated against XLPE liners

OXINIUM Technology demonstrates superior wear resistance and reduced surface damage compared with CoCr in laboratory and retrieval studies^{4,7,18,19,21-23,29-32*}

*The results of in vitro wear simulation testing have not been proven to quantitatively predict clinical wear performance.



Property summary:

Strength and stability

Strength refers to a material's toughness and, as a result, its ability to avoid chipping or fracture under application of force. The chemical stability of materials contributes to their susceptibility to mechanical and chemical wear during implantation



OXINIUM Oxidized Zirconium exhibits high fracture resistance in laboratory strength testing^{10,11}

OXINIUM	Zirconia
<p>20,000lbf</p>	<p>Burst after unloading</p>
<p>OXINIUM femoral heads withstood strong force without fracture, whilst zirconia heads withstood the maximum load but burst after unloading¹¹</p>	

OXINIUM Technology behaves like metal under application of strong force¹⁰

OXINIUM knee components withstood cyclic loads estimated to be 13 times body weight without fracture¹⁰

Under application of substantial loads when fractures did occur, OXINIUM fractures demonstrated:¹⁰

- A clean break with a metal-like fracture mechanism
- No evidence of oxide delamination

The surface stability of OXINIUM Technology reduces the impact of ageing and mechanical stress compared with ceramic, including:¹¹

Less surface roughening¹¹

Reduced in vitro and in vivo phase transformation¹¹

OXINIUM Technology behaves with the strength of a metal, resisting fracture under high mechanical stress in laboratory testing,^{10,11} and its surface stability reduces the impact of ageing and mechanical stress versus ceramic¹¹



Property summary: Biocompatibility

Nickel (Ni), cobalt (Co) and chromium (Cr) are sensitizers in 10–15% of the general population.³³ Inflammatory responses to metal debris released by orthopaedic implants may contribute to aseptic loosening and periprosthetic osteolysis, potentially accelerating implant failure^{13,25}



OXINIUM Oxidized Zirconium contains almost no Ni, Co and Cr compared with CoCr¹

	OXINIUM	CoCr
Ni	<0.0035%	<0.5%
Co	<0.002%	59–67%
Cr	<0.02%	27–30%

OXINIUM Technology (zirconium oxide and zirconium alloy particles*) leads to improved in vitro cellular viability (ie. reduced cytotoxicity) vs CoCr in:^{13,34}



Osteoblasts



Macrophages



Fibroblasts

OXINIUM Technology (zirconium oxide and zirconium alloy particles*) induces lower inflammatory cytokine expression vs CoCr in vitro^{13,34†}




Figure. Fold-change in inflammatory cytokine expression seen with Co-alloy particles versus zirconium oxide and alloy particles* in vitro

OXINIUM Technology contains very low levels of Ni, Co and Cr,¹ and zirconium oxide and underlying zirconium particles demonstrate improved cellular viability and lower inflammatory response versus CoCr in vitro^{13,25,26,34†}

*Derived from DH-OXINIUM material. †The results of in vitro cytokine expression analyses have not been proven to quantitatively predict clinical cytokine expression.



Property summary:

Corrosion and fretting resistance

Taper corrosion is the degradative wear of the taper junctions of modular implants and may result in the release of metal ions and construct instability. Taper corrosion is a known cause of problems with implant longevity^{35,36}



OXINIUM Oxidized Zirconium shows fewer signs of corrosion damage compared with CoCr, including:^{14-16,27,37*}

Lower mean Goldberg corrosion scores in retrieved hip and knee components^{14,15,27,37}




Figure 1. Mean Goldberg corrosion scores for retrieved OXINIUM and CoCr components across three studies

Reduced incidence of material loss due to MACC in worst-case retrieved femoral heads¹⁴




Figure 2. Proportion of retrieved OXINIUM and CoCr femoral heads with Goldberg corrosion scores >3 that exhibited material loss due to MACC

No corrosion-induced discolouration^{27,37}

Less extensive fretting scars^{27,37}

Less metal transfer²⁷

Significantly lower metal ion concentrations on femoral heads (p<0.05)¹⁶

OXINIUM Technology demonstrates superior corrosion and fretting resistance compared with CoCr and other materials in laboratory and retrieval studies^{14-16,27,37*}

*The results of in vitro wear simulation testing have not been proven to quantitatively predict clinical wear performance.



Wear of polyethylene against oxidized zirconium femoral components: effect of aggressive kinematic conditions and malalignment in total knee arthroplasty¹⁸

Ezzet K, Hermida J, Steklov N, D'Lima D. *J Arthroplasty*. 2012;27:116–121.

Objective

Determine whether OXINIUM Oxidized Zirconium femoral components can reduce polyethylene wear under conditions of varus malalignment (75:25 mediolateral distribution of vertical tibial load) and excessive tibiofemoral rotations (20° tibial axial rotation) in vitro

Results

The mean±SD gravimetric conventional polyethylene (CPE) wear rate was 55% lower with OXINIUM versus CoCr components (17.2±1.2 vs 38.6±1.3mg/million cycles, respectively; p<0.001; Figure)

Volumetric CPE wear rate was also significantly lower with OXINIUM versus CoCr components (p<0.001)

Values reported in the original material have been rounded to 1 d.p. in the present summary.




Figure. Mean gravimetric CPE wear rate with CoCr and OXINIUM femoral components

Conclusion

OXINIUM femoral components substantially reduced CPE wear compared with CoCr under simulated conditions of athletically active patients with modestly malaligned TKA prostheses



Reduced wear with oxidized zirconium femoral heads¹⁹

Good V, Ries M, Barrack RL, Widding K, Hunter G, Heuer D. *J Bone Joint Surg Am.* 2003;85-A Suppl 4:105–110.

Objective

Evaluate the in vitro wear characteristics of CPE and 10-XLPE liners articulating against OXINIUM Oxidized Zirconium femoral heads under both smooth and clinically relevant, roughened conditions using a hip simulator

Results

Wear rate for CPE liners was significantly lower when articulated against OXINIUM heads than CoCr heads under all test conditions ($p < 0.01$; Figure)

Wear rate for 10-XLPE liners was undetectable when articulated against smooth or rough OXINIUM heads, but wear was detected with rough CoCr heads (Figure)

Compared with CoCr heads, smooth and rough OXINIUM heads generated fewer polyethylene particles per cycle against CPE (-30%; $p = 0.01$ and -45%; $p < 0.01$, respectively), and fewer polyethylene particles per cycle against 10-XLPE (-27%; $p = 0.03$ and 63%; $p = 0.03$, respectively)




Figure. Wear rate of CPE and 10-XLPE liners against CoCr and OXINIUM femoral heads under smooth and roughened conditions

*Wear undetectable.

Conclusion

OXINIUM heads produced significantly less wear and fewer polyethylene particles than CoCr heads when articulated against CPE and XLPE liners, under both smooth and clinically relevant, roughened conditions



Design of an advanced bearing system for total knee arthroplasty²⁰

Morrison M, Jani S, Parikh A. *Lubricants*. 2015;3:475–492.

Objective

Evaluate the in vitro wear performance of 5-XLPE and 7.5-XLPE tibial inserts coupled with OXINIUM Oxidized Zirconium and CoCr femoral components under pristine and rough conditions using a knee simulator

Results

Mean±SD wear rates of 5-XLPE and 7.5-XLPE tibial inserts articulating against pristine OXINIUM components were lower (72 and 79%, respectively; $p \leq 0.001$) than with pristine CoCr components (Figure)

Mean±SD wear rates of 5-XLPE and 7.5-XLPE tibial inserts articulating against rough OXINIUM components were also lower (80 and 83%, respectively; $p \leq 0.003$) than with pristine CoCr components (Figure)




Figure. Wear rate of 5-XLPE and 7.5-XLPE tibial inserts against OXINIUM and CoCr femoral components

*†Denotes statistical difference of $p \leq 0.001$ and $p \leq 0.003$, respectively.

Conclusion

OXINIUM femoral components coupled with 5-XLPE and 7.5-XLPE tibial inserts demonstrated superior wear resistance versus CoCr components under both pristine and rough conditions



Long-term wear performance of an advanced bearing technology for TKA²¹

Papannagari R, Hines G, Sprague J, Morrison M. Poster presented at: Orthopaedic Research Society Annual Meeting; 2011; Long Beach, CA.

Objective

Evaluate the long-term in vitro polyethylene wear performance of OXINIUM Oxidized Zirconium femoral components (LEGION[◇] Cruciate Retaining [CR] TKA) coupled with 7.5-XLPE tibial inserts and CoCr femoral components (LEGION[◇] CR TKA) coupled with CPE tibial inserts using a knee simulator

Results

Compared with CoCr/CPE couples after 5 million cycles ($120.4 \pm 12.0 \text{ mm}^3$), the mean \pm SD volumetric wear of OXINIUM/7.5-XLPE couples was 98% lower after 5 million cycles ($2.7 \pm 1.5 \text{ mm}^3$; $p < 0.01$) and 81% lower after 45 million cycles ($22.8 \pm 7.2 \text{ mm}^3$; $p < 0.01$; Figure)

OXINIUM femoral component roughness did not change during wear testing ($p > 0.05$) and there was no measurable loss of oxide surface



Values reported in the original material have been rounded to 1 d.p. in the present summary.




Figure. Mean volumetric wear of CoCr/CPE and OXINIUM/7.5-XLPE couples

*Denotes statistical difference ($p < 0.01$) versus CoCr/CPE.

Conclusion

OXINIUM femoral components coupled with 7.5-XLPE tibial inserts demonstrated superior long-term wear resistance compared with CoCr heads and CPE inserts. OXINIUM components were in excellent condition following 45 million cycles of lab testing



Wear of conventional and highly crosslinked polyethylene liners during simulated fast walking/jogging²²

Parikh A, Hill P, Hines G, Pawar V. Poster presented at: Orthopaedic Research Society Annual Meeting; 2009; Las Vegas; NV.

Objective

Determine if OXINIUM Oxidized Zirconium femoral heads provide an advantage over CoCr femoral heads during simulated fast walking/jogging in vitro

Results

The mean \pm SD wear rate of 10-XLPE liners was significantly lower when articulated against 36mm OXINIUM versus CoCr heads (1.6 \pm 0.2 and 3.5 \pm 0.2mm³/million cycles, respectively; p<0.01)

The cumulative wear of 10-XLPE liners was lower with 36mm OXINIUM versus CoCr heads when tested under fast walking/jogging conditions (Figure)




Figure. Cumulative wear of 10-XLPE liners articulated against 36mm OXINIUM and CoCr heads

Conclusion

OXINIUM femoral heads demonstrated superior wear resistance to CoCr femoral heads during simulated moderate physical activity



Long-term simulator wear performance of an advanced bearing technology for THA⁴

Parikh A, Hill P, Pawar V, Sprague J. Poster presented at: Orthopaedic Research Society Annual Meeting; 2013; San Antonio, TX.

Objective

Assess the long-term in vitro polyethylene wear performance of OXINIUM Oxidized Zirconium femoral heads coupled with 10-XLPE liners using a hip simulator


Results

Mean volumetric liner wear after 45 million cycles with OXINIUM femoral heads and 10-XLPE was 67% lower than with CoCr and 10-XLPE ($p < 0.05$; Figure), and 80% lower than with CoCr and CPE after 7.8 million cycles ($p < 0.05$)*

No significant loss in OXINIUM femoral head surface thickness was observed over 45 million cycles



*Approximately 7.8 and 45 million cycles of testing were completed for CPE and 10-XLPE liners, respectively.



Conclusion

OXINIUM femoral heads combined with 10-XLPE liners demonstrated superior long-term wear resistance versus CoCr femoral heads coupled with either 10-XLPE or CPE liners in hip simulator testing



Polyethylene wear performance of oxidized zirconium and cobalt-chromium knee components under abrasive conditions²³

Ries MD, Salehi A, Widding K, Hunter G. *J Bone Joint Surg.* 2002;84-A (S2):129–135.


Objective

Evaluate the in vitro polyethylene wear performance of OXINIUM Oxidized Zirconium and CoCr femoral components in vitro under abrasive conditions (produced through tumbling)

Results



The aggregate \pm SD wear rate of UHMWPE inserts was 89% lower with OXINIUM versus CoCr components following abrasion (5.3 ± 0.6 vs 46.1 ± 4.7 mm³/million cycles, respectively; $p<0.01$)



OXINIUM components produced 44% fewer submicrometre UHMWPE particles than CoCr ($p=0.02$)



Condyle roughness scores were lower for OXINIUM than CoCr components, suggesting that the former was roughened less by tumbling

Values reported in the original material have been rounded to 1 d.p. in the present summary.

Conclusion

OXINIUM femoral components demonstrate superior resistance to roughening versus CoCr components, which may have contributed to a reduction in wear under abrasive conditions



Oxidation microstructures and interfaces in the oxidized zirconium knee²⁴

Hobbs LW, Rosen VB, Mangin SP, et al. *Int J Appl Ceram Technol.* 2005;2:221–246.

Objective

Characterise the microstructural basis for the surface durability and strong adherence of the oxide layer of OXINIUM Oxidized Zirconium Technology

Results

'Grains' of zirconium oxide arranged in columns were observed at the surface of the oxide layer and may contribute to its surface durability (Figure)

A lack of pores at the oxide-alloy interface and regions of alloy material anchoring the two layers were observed, potentially contributing to the strong oxide-alloy interface adhesion of the OXINIUM oxide layer




Figure. Columnar arrangement of OXINIUM oxide layer and features that may confer its surface durability

✓ Resistance to lateral crack propagation

✓ Maintenance of surface smoothness

✓ Resists extensive spallation

Conclusion

Microstructural features observed in the oxide and alloy layers, as well as at the oxide-alloy interface, may contribute positively to the wear resistance, low friction and strong adhesion of the OXINIUM oxide layer when used in orthopaedic applications



Mechanical behaviour of zirconia, alumina, and oxidized zirconium modular heads¹¹

Sprague J, Salehi A, Tsai S, Pawar V, Thomas R, Hunter G. In: Brown S, Clarke IC, Gustafson A, eds. *International Society for Technology in Arthroplasty*. Birmingham, AL; 2004.

Objective


Compare the in vitro strength and stability of OXINIUM Oxidized Zirconium and ceramic (alumina and zirconia) femoral heads before and after autoclaving that simulated 80 years of in vivo ageing

Results


Pre- (n=3) and post-autoclaved (n=3) OXINIUM femoral heads did not fracture and displayed virtually no signs of damage at 20,000lbs of force during crush testing. By contrast, pre- (n=3) and post-autoclaved (n=2/3) zirconia heads withstood the maximum load but all burst after unloading (Figure), and alumina heads (n=6) all fractured during testing

OXINIUM heads were largely unaffected by autoclaving, whereas zirconia heads displayed phase transformation, surface roughening and microcracks

OXINIUM femoral heads



Zirconia femoral heads



Conclusion

OXINIUM femoral heads did not display brittle fracture during crush testing and exhibited greater stability after autoclaving compared with ceramic heads



Mechanical testing and finite element analysis of oxidized zirconium femoral components¹⁰


Tsai S, Sprague J, Hunter G, Thomas R, Salehi A. *Trans Soc Biomaterials*. 2001;24:163.

Objective

Evaluate the in vitro strength of OXINIUM Oxidized Zirconium femoral components (GENESIS^o II TKA) during cyclic fatigue testing

Results

The medial condyle of OXINIUM knee femoral components (n=24) withstood a maximum fatigue load of 1,000lbf (equivalent to 6.67 times body weight) for 10 million cycles without fracture, and both medial condyles were estimated to be able to withstand the equivalent of 13 times body weight



Conclusion

Medial condyles of OXINIUM knee femoral components behaved like metal under application of strong force, and at higher loads fractured with a metal-like mechanism



Oxidized Zr-alloy particles induce a lower incidence of in vitro lymphocyte metal-sensitivity responses compared to cobalt and titanium implant alloys²⁵

Caicedo M, Pawar V, Hallab NJ. Poster presented at: Orthopaedic Research Society Annual Meeting; 2014; New Orleans, LA.

Objective

Compare the in vitro immune hypersensitivity response to OXINIUM Technology (zirconium oxide, zirconium alloy particles* and zirconium ions) with that of Co-alloy, Ti-alloy and bone cement particles using peripheral blood mononuclear cells from healthy volunteers (n=50)

Results

Zirconium ions did not induce a significant increase in the proliferation of cultured cells (quantified using the stimulation index[†]) in any of the subjects tested, with effectively 0% incidence of zirconium reactivity

The incidence of metal hypersensitivity (at a stimulation index of >4, indicating moderate reactivity) was substantially lower for zirconium oxide and zirconium alloy particles (0%) versus Ti-alloy (5%) and bone cement (5%; Figure)

[†]The stimulation index measures reactivity to different particles, with results that range from 0–2 (non-reactive), 2–4 (mildly reactive), 4–8 (reactive) and >8 (highly reactive).




Figure. Incidence of hypersensitivity at a stimulation index >4 for zirconium oxide and alloy particles, Ti-alloy and bone cement particles

Conclusion

Zirconium oxide and zirconium alloy particles resulted in lower rates of in vitro lymphocyte hypersensitivity compared with Ti-alloy and bone cement



The effect of diffusion hardened oxidized zirconium wear debris on cell viability and inflammation – an in vitro study²⁶

Rose SF, Weaver CL, Fenwick SA, Horner A, Pawar VD. *J Biomed Mater Res Part B*. 2012;100B:1359–1368.

Objective

Evaluate the in vitro biocompatibility of wear debris generated by OXINIUM Technology (zirconium oxide and zirconium alloy particles*) compared with Ti, Co and CoCr

Results

Zirconium oxide and zirconium alloy particles did not induce significant increases in TNF release above the media control at most doses at 4 or 8 hours (except 1mg at 4 hours), in contrast to Ti, Co and CoCr, which dose-dependently stimulated TNF release at all timepoints (Table)

Wear debris generated by zirconium oxide and zirconium alloy particles demonstrated lower cytotoxicity than Co and CoCr in both osteoblast- and fibroblast-like cells

	4hr		8hr		24hr	
	0.1mg	1mg	0.1mg	1mg	0.1mg	1mg
Co	✓	✓	✓	✓	✓	✓
CoCr	✓	✓	✓	✓	✓	✓
Ti		✓	✓	✓	✓	✓
Zirconium oxide and alloy particles		✓			✓	✓

Table. Doses (mg) at which a significant increase in TNF release above the media control was observed following exposure to wear debris particles (tested doses were 0.001mg, 0.01mg, 0.1mg and 1mg)

Conclusion

Wear debris generated by zirconium oxide and zirconium alloy particles was less cytotoxic and elicited a reduced inflammatory response versus Co and CoCr[†]

*Derived from DH-OXINIUM material. †The results of in vitro cytokine expression analyses have not been proven to quantitatively predict clinical cytokine expression.



Characterization of femoral head taper corrosion features using a 22-year retrieval database¹⁴

Cartner J, Aldinger P, Li C, Collins D. *HSS Journal*. 2017;13:35–42.


Objective

Compare taper corrosion severity and extent of material loss in retrieved CoCr and OXINIUM Oxidized Zirconium femoral heads

Results

CoCr femoral heads (n=165) exhibited higher average corrosion scores* compared with OXINIUM heads (n=45; 2.5 ± 1.0 vs 1.9 ± 0.7 ; $p < 0.05$; Figure)

Of femoral heads with a corrosion score ≥ 3 , measurable material loss due to MACC was found on 17/30 CoCr heads, and was not observed on any of the OXINIUM heads (0/8)



*Corrosion was assessed using the semi-quantitative Goldberg scoring system. Scores of 1 (none), 2 (mild), 3 (moderate) or 4 (severe) were assigned to taper surfaces depending on the level of discolouration, fretting and corrosion evidence according to Goldberg et al.³⁸





Figure. Mean corrosion score of retrieved CoCr and OXINIUM femoral heads

Conclusion

OXINIUM femoral heads were associated with less corrosion damage than CoCr heads, with lower corrosion scores and no evidence of material loss due to MACC



Do oxidized zirconium heads decrease tribocorrosion in total hip arthroplasty?¹⁵

Hampton C, Weitzler L, Baral E, Wright TM, Bostrom MPG. *Bone Joint J.* 2019;101-B:386–389.

Objective

Evaluate female taper fretting and corrosion (tribocorrosion) in retrieved OXINIUM Oxidized Zirconium femoral heads and matched* CoCr heads (n=28 for both)

Results

OXINIUM femoral heads showed significantly lower mean [range] corrosion (1.3 [1–2.8] vs 2.2 [1–4]; $p < 0.01$) and fretting scores[†] (1.3 [1–2] vs 1.5 [1–2.3]; $p = 0.02$) compared with CoCr heads (Figure)

Fretting or corrosion was observed in 29% (8/28) of OXINIUM femoral heads, compared with 96% (27/28) of CoCr heads

Values reported in the original material have been rounded to 1 d.p. in the present summary. *Retrieved OXINIUM and CoCr femoral heads were matched according to length of implantation within ± 5 months and patient BMI within $\pm 3 \text{ kg/m}^2$. [†]Corrosion and fretting was assessed using a modified Goldberg scoring system,³⁸ with scores of 1–4 assigned to each femoral head according to the criteria presented in the Figure.




Figure. Mean fretting and corrosion scores of retrieved CoCr and OXINIUM femoral heads

Conclusion

Retrieved OXINIUM femoral heads displayed reduced taper fretting and corrosion compared with matched CoCr heads



Acidic fretting tests of oxidized Zr-2.5Nb, CoCr and SS femoral heads²⁷

Pawar V, Jones B, Sprague J, Salehi A, Hunter G. In: Medical Device Materials II: Proceedings of the Materials and Processes for Medical Devices Conference. Helmus M, Medlin D, eds. ASMI; 2005;403–408.

Objective

Evaluate the relative resistance of OXINIUM Oxidized Zirconium, CoCr and stainless steel (SS) femoral heads combined with Ti or SS stem trunnions to fretting corrosion damage following in vitro cyclic fatigue testing under extreme, supra-clinical conditions of elevated temperature and acidity

Results

Femoral head tapers and trunnions of OXINIUM-Ti/-SS combinations showed the lowest level of chemical attack with no discolouration and fewer signs of damage when compared with other head-trunnion combinations

Couple	Femoral head	Trunnion
CoCr-Ti	Cr-rich deposits and Ti transfer	Co-, Cr- and Mo-rich deposits
OXINIUM-Ti	Few areas of damage to the OXINIUM material and minimal Ti transfer	Fretting scars and weak Zr signal, no chemical etching or pitting
SS-SS	Extensive chemical activity, iron oxide and Cr-rich deposits	Iron oxide deposits
OXINIUM-SS	Few areas of damage to the OXINIUM material and small iron rich areas	Fretting scars, no chemical etching or pitting

Table. Observations on the taper of the femoral head and trunnion from several bearing combinations based on SEM/spectrometry analysis

Conclusion

OXINIUM-Ti/SS femoral-head trunnion combinations demonstrated the lowest levels of corrosion compared with the CoCr-Ti and SS-SS combinations



Further studies

Wear resistance

Study	Design / n=	Key findings
<p>Bourne RB, Barrack R, Rorabeck CH, Salehi A, Good V. Arthroplasty options for the young patient. <i>Clin Orthop Relat Res.</i> 2005;441:159–167.²⁹</p> <p>Available at: Clinical Orthopaedics and Related Research</p>	<p>In vitro hip simulator testing of femoral heads damaged by simulated dislocation and undamaged heads; analysis of retrieved implants</p> <p>Implant tested: hip</p> <p>Six damaged OXINIUM Oxidized Zirconium heads (articulated against CPE liners, n=3; articulated against XLPE liners, n=3); four retrieved femoral heads (CoCr, n=2; OXINIUM, n=2)</p>	<ul style="list-style-type: none"> Hip simulator testing showed equivalent mean±SD wear rates for damaged OXINIUM heads and undamaged CoCr heads that were articulated against conventional UHMWPE liners (37.4±0.7 vs 38.5±0.7mm³/million cycles, respectively; p=0.46)* <ul style="list-style-type: none"> When articulated against XLPE liners, damaged OXINIUM heads demonstrated a low mean±SD wear rate of 2.4±0.3mm³/million cycles* Retrieval analysis revealed scratching of both types of metal femoral heads with in vivo use, though CoCr femoral heads showed more severe scratching compared with OXINIUM heads
<p>Carli AV, Patel AR, Cross MB et al. Long-term performance of oxidized zirconium on conventional and highly cross-linked polyethylene in total hip arthroplasty. <i>J Int Soc Orthop Surg Traum.</i> 2020;6.³⁹</p> <p>Available at: SICOT-J</p>	<p>Single-centre, longitudinal, retrospective study; in vivo radiographic assessment of polyethylene wear rates and osteolysis</p> <p>Implant tested: hip</p> <p>168 patients who underwent THA with OXINIUM Oxidized Zirconium femoral heads (UHMWPE liners, n=80; XLPE liners, n=88)</p>	<ul style="list-style-type: none"> When compared with OXINIUM/UHMWPE couples (average follow-up: 10.5±1.2 years), OXINIUM/XLPE couples (average follow-up: 10.3±1.0 years) were associated with:* <ul style="list-style-type: none"> A significantly lower mean±SD annual wear rate (0.05±0.03 vs 0.21±0.12mm/year, respectively; p<0.001) A significantly lower incidence of acetabular osteolysis (0 vs 21.2%, respectively; p<0.001) and femoral osteolysis (0 vs 46.3%, respectively; p<0.001) No OXINIUM/XLPE patients underwent revision THA, compared with six patients (7.5%) who received OXINIUM/UHMWPE

*Values reported in the original material have been rounded to 1 d.p.



Further studies

Wear resistance

Study	Design / n=	Key findings
<p>Karidakis GK, Karachalios T. Oxidized zirconium head on crosslinked polyethylene liner in total hip arthroplasty: A 7- to 12-year in vivo comparative wear study. <i>Clin Orthop Relat Res.</i> 2015;473:3836–3845.³²</p> <p>Available at: Clinical Orthopaedics and Related Research</p>	<p>Single-centre, prospective, randomised trial; in vivo assessment of polyethylene wear</p> <p>Implant tested: hip</p> <p>85 patients who underwent THA with different femoral head/liner combinations and completed 10 years of follow up (28mm ceramic/CPE, n=20; 28mm ceramic/XLPE, n=21; 28mm OXINIUM Oxidized Zirconium, n=23; 32mm OXINIUM/XLPE, n=22)</p>	<ul style="list-style-type: none"> At 10 years follow-up, compared with both 28mm ceramic groups, the 28mm OXINIUM/XLPE group showed: <ul style="list-style-type: none"> Lower mean±SD volumetric wear (28mm OXINIUM/XLPE: 35.6±6.4; ceramic/CPE: 89.5±13.1; ceramic/XLPE: 162.5±25.3mm³) Lower mean±SD volumetric wear rates (OXINIUM/XLPE: 21.5±5.4; ceramic/CPE: 38.0±8.0; ceramic/XLPE: 48.0±10.1mm³/year) Compared with both 28mm ceramic groups, the 32mm OXINIUM/XLPE group showed: <ul style="list-style-type: none"> Lower mean±SD volumetric wear (OXINIUM/XLPE: 35.4±5.0; ceramic/CPE: 89.5±13.1; ceramic/XLPE: 162.5±25.3mm³) Lower mean±SD volumetric wear rates (OXINIUM/XLPE: 21.0±6.4; ceramic/CPE: 38.0±8.0; ceramic/XLPE: 48.0±10.1mm³/year)
<p>Heyse TJ, Elpers ME, Nawabi DH, Wright TM, Haas SB. Oxidized zirconium versus cobalt-chromium in TKA: profilometry of retrieved femoral components. <i>Clin Orthop Relat Res.</i> 2014;472:277–283.⁷</p> <p>Available at: Clinical Orthopaedics and Related Research</p>	<p>Analysis of retrieved implants</p> <p>Implant tested: knee</p> <p>OXINIUM Oxidized Zirconium femoral components, n=10; matched CoCr femoral components, n=10</p>	<ul style="list-style-type: none"> CoCr femoral components showed more scratching than matched OXINIUM components on visual assessment The average surface roughness for retrieved CoCr components was 83% greater than that of the retrieved OXINIUM components, with CoCr showing significantly higher surface roughness than OXINIUM across all four surface roughness parameters (all p≤0.036) CoCr bearing surfaces showed 39% greater maximum peak-to-peak height roughness and 33% greater 10-point height roughness, when compared with OXINIUM components The increase in surface roughness between retrieved and pristine OXINIUM implants was less than the increase between retrieved and pristine CoCr implants (20 vs 267%, respectively)



Further studies

Wear resistance

Study	Design / n=	Key findings
<p>Morrison ML, Jani S, Parikh A. Development of an advanced bearing couple for total knee arthroplasty. Poster presented at: Orthopaedic Research Society Annual Meeting; 2011; Long Beach, CA.³¹</p> <p>Available at:  Orthopaedic Research Society</p>	<p>In vitro knee simulation of in vivo scratching and wear against different UHMWPE materials</p> <p>Implant tested: knee</p> <p>OXINIUM Oxidized Zirconium femoral components, n=3; CoCr femoral components, n=3</p>	<ul style="list-style-type: none"> OXINIUM femoral components showed lower surface roughness compared with CoCr components following tumbling ($p < 0.001$) When articulated against UHMWPE liners, pristine OXINIUM components showed a significantly lower wear rate than pristine CoCr components ($p \leq 0.001$) OXINIUM femoral components did not show any significant changes in wear rates ($p \geq 0.181$) following tumbling against any UHMWPE liner, in contrast to CoCr components, which showed significantly increased wear rates following tumbling ($p \leq 0.003$)
<p>Anderson FL, Koch CN, Elpers ME, Wright TM, Haas SB, Heyse TJ. Oxidized zirconium versus cobalt alloy bearing surfaces in total knee arthroplasty. <i>Bone Joint J.</i> 2017;99-B:793–798.²⁸</p> <p>Available at:  The Bone & Joint Journal</p>	<p>Analysis of retrieved implants</p> <p>Implant tested: knee</p> <p>Tibial inserts that had been articulated against OXINIUM Oxidized Zirconium femoral components, n=20; matched inserts that had been articulated against CoCr, n=20</p>	<ul style="list-style-type: none"> Compared with retrieved inserts that had been articulated against CoCr components, retrieved inserts articulated against OXINIUM components demonstrated: <ul style="list-style-type: none"> A significantly lower overall mean\pmSD volume loss from the UHMWPE bearing surface (122 ± 87 vs $170 \pm 96 \text{mm}^3$; $p = 0.03$) A smaller mean\pmSD loss of volume in the medial (72 ± 67 vs $92 \pm 60 \text{mm}^3$) and lateral (49 ± 36 vs $79 \pm 61 \text{mm}^3$) compartments ($p = \text{n.s.}$)



Further studies

Wear resistance

Study	Design / n=	Key findings
<p>Brandt JM, Guenther L, O'Brien S, Vecherya A, Turgeon TR, Bohm ER. Performance assessment of femoral knee components made from cobalt–chromium alloy and oxidized zirconium. <i>Knee.</i> 2013;20:388–396.³⁰</p> <p>Available at: Knee</p>	<p>Analysis of retrieved implants</p> <p>Implant tested: knee</p> <p>OXINIUM Oxidized Zirconium femoral components, n=26; matched CoCr knee femoral components, n=26</p>	<ul style="list-style-type: none"> Retrieved CoCr components had a femoral damage score (FDS) 1.8 times greater than retrieved OXINIUM components No significant difference was observed between any surface parameters measured on the condyles of the new, never implanted OXINIUM components and the retrieved OXINIUM components ($p \geq 0.201$) The surface roughness of retrieved CoCr femoral components was found to be significantly higher than that of retrieved OXINIUM femoral components across three parameters ($p \leq 0.031$)



Further studies

Biocompatibility

Study	Design / n=	Key findings
<p>Dalal A, Pawar V, McAllister K, Weaver C, Hallab NJ. Orthopedic implant cobalt-alloy particles produce greater toxicity and inflammatory cytokines than titanium alloy and zirconium alloy-based particles in vitro, in human osteoblasts, fibroblasts, and macrophages. <i>J Biomed Mater Res Part A.</i> 2012;100A:2147–2158.¹³</p> <p>Available at:  Journal of Biomedical Materials Research</p>	<p>In vitro cytotoxicity and inflammatory reactivity testing in human peri-implant cells</p> <p>Implant tested: material only</p> <p>OXINIUM Technology (zirconium oxide and zirconium alloy particles*); particles derived from CoCr-alloy and Ti-alloy</p>	<ul style="list-style-type: none"> Zirconium oxide and zirconium alloy particles led to improved viability of osteoblasts, fibroblasts and macrophages at 48hrs compared with CoCr-alloy particles Osteoblast and fibroblast cell proliferation was inhibited significantly less with zirconium oxide and zirconium alloy particles versus CoCr-alloy particles (p<0.05) Inflammatory cytokine expression was lower following exposure to zirconium oxide and zirconium alloy particles, compared with CoCr- and Ti-alloy particles:[†] <ul style="list-style-type: none"> Significantly lower levels of IL-6 and TNF-α in osteoblasts (p<0.05) Significantly lower levels of IL-8 in macrophages (p<0.05)
<p>Hallab NJ, McAllister H, Jacobs JJ, Pawar V. Zirconium-alloy and zirconium-oxide particles produce less toxicity and inflammatory cytokines than cobalt-alloy and titanium-alloy particles in vitro, in human osteoblasts, fibroblasts and macrophages. Poster presented at: Orthopaedic Research Society Annual Meeting; 2012; San Francisco, CA.³⁴</p> <p>Available at:  Orthopaedic Research Society</p>	<p>In vitro cytotoxicity and inflammatory reactivity testing in human peri-implant cells</p> <p>Implant tested: material only</p> <p>OXINIUM Technology (zirconium oxide and zirconium alloy particles*); Co-alloy and Ti-alloy particles</p>	<ul style="list-style-type: none"> Zirconium oxide, and zirconium alloy and zirconium oxide mix, particles led to improved viability of osteoblasts, fibroblasts and macrophages compared with Co-alloy particles Inflammatory cytokine expression (IL-6 and TNF-α in all cells and IL-8 in macrophages) was significantly lower following exposure to zirconium oxide and zirconium alloy particles, compared with Co- and Ti-alloy particles[†]

*Derived from DH-OXINIUM material. [†]The results of in vitro cytokine expression analyses have not been proven to quantitatively predict clinical cytokine expression.



Further studies

Corrosion and fretting resistance

Study	Design / n=	Key findings
<p>Aldinger P, Williams T, Woodard E. Accelerated fretting corrosion testing of zirconia toughened alumina composite ceramic and a new composition of ceramicised metal femoral heads. Poster presented at: Orthopaedic Research Society Annual Meeting; 2017; San Diego, CA.¹⁶</p> <p>Available at:  Orthopaedic Research Society</p>	<p>In vitro accelerated fretting corrosion testing</p> <p>Implant tested: hip</p> <p>Biolog^o Delta* (ZTA) femoral heads, n=3; CoCr femoral heads, n=6; diffusion-hardened OXINIUM Oxidized Zirconium (DH-OXINIUM) femoral heads, n=6</p>	<ul style="list-style-type: none"> Signs of metal transfer on the femoral heads and fretting on the stem tapers of ZTA couples were observed upon visual examination DH-OXINIUM total metal ion concentration was significantly lower compared with ZTA and CoCr heads (p<0.05) The average wear depth for DH-OXINIUM heads was significantly lower than that of CoCr heads, over 10 million cycles (p<0.05)
<p>Cartner J, Aldinger P, Newman M. Characterization of tapers in TKA revisions from a 16-year retrieval database. Poster presented at: Orthopaedic Research Society Annual Meeting; 2016; Orlando, FL.³⁷</p> <p>Available at:  Orthopaedic Research Society</p>	<p>Analysis of retrieved implants</p> <p>Implant tested: knee</p> <p>70 tapers from 27 retrieved TKAs (Ti6Al4V stems/couplers, n=35; TiAl4V tibial baseplates, n=21; Ti6Al4V femoral stems, n=7; CoCr femoral components, n=3; OXINIUM Oxidized Zirconium femoral components, n=3; CoCr tibial baseplate, n=1)</p>	<ul style="list-style-type: none"> OXINIUM tapers did not show any significant discoloration, wear scars, fretting or corrosion precipitates on their surface, while CoCr and titanium tapers exhibited all of these features, upon SEM analysis Compared with retrieved tapers comprised of CoCr/Ti6Al4V material combinations, OXINIUM/Ti6Al4V couples showed significantly fewer corrosion artifacts upon microscopic analysis and scoring (p<0.05)

References

1. ASTM. Standard specification for wrought zirconium-2.5niobium alloy for surgical implant applications (UNS R60901). Available at: <https://www.astm.org/Standards/F2384.htm>. Accessed June 23, 2022.
2. ASTM. International standard specification for cobalt-28 chromium-6 molybdenum alloy castings and casting alloy for surgical implants (UNS R30075). Available at: <https://www.astm.org/f0075-07.html>. Accessed June 23, 2022.
3. Hunter G, Dickinson J, Herb B, Graham R. Creation of oxidized zirconium orthopaedic implants. *J ASTM Int*. 2005;2:1–14.
4. Parikh A, Hill P, Pawar V, Sprague J. Long-term simulator wear performance of an advanced bearing technology for THA. Poster presented at: Orthopaedic Research Society Annual Meeting; 2013; San Antonio, TX.
5. Smith+Nephew 2010. OR-10-155.
6. Long M, Riester L, Hunter G. Nano-hardness measurements of Oxidized Zr-2.5 Nb and various orthopaedic materials. Poster presented at: Society for Biomaterials Annual Meeting; 1998; San Diego, CA.
7. Heyse TJ, Elpers ME, Nawabi DH, Wright TM, Haas SB. Oxidized zirconium versus cobalt-chromium in TKA: profilometry of retrieved femoral components. *Clin Orthop Relat Res*. 2014;472:277–283.
8. Salehi A, Tsai S, Pawar V, et al. Wettability analysis of orthopaedic materials using optical contact angle methods. *Key Eng Mater*. 2006;309:1199–1202.
9. Poggie RA, Wert JJ, Mishra AK, Davidson JA. Friction and wear characterization of UHMWPE in reciprocating sliding contact with Co-Cr, Ti-6Al-4V, and zirconia implant bearing surfaces. In: *Wear and Friction of Elastomers*. ASTM International; 1992.
10. Tsai S, Sprague J, Hunter G, Thomas R, Salehi A. Mechanical testing and finite element analysis of oxidized zirconium femoral components. *Trans Soc Biomaterials*. 2004;24:163.
11. Sprague J, Salehi A, Tsai S, Pawar V, Thomas R, Hunter G. Mechanical behavior of zirconia, alumina, and oxidized zirconium modular heads. In: Brown S, Clarke IC, Gustafson A, eds. *International Society for Technology in Arthroplasty*. Birmingham, AL; 2004.
12. Smith+Nephew 2020. EA/RECON/POLAR3/007/v1.
13. Dalal A, Pawar V, McAllister K, Weaver C, Hallab NJ. Orthopedic implant cobalt-alloy particles produce greater toxicity and inflammatory cytokines than titanium alloy and zirconium alloy-based particles in vitro, in human osteoblasts, fibroblasts, and macrophages. *J Biomed Mater Res Part A*. 2012;100A:2147–2158.
14. Cartner J, Aldinger P, Li C, Collins D. Characterization of femoral head taper corrosion features using a 22-year retrieval database. *HSS Journal*. 2017;13:35–42.
15. Hampton C, Weitzler L, Baral E, Wright TM, Bostrom MPG. Do Oxidized zirconium heads decrease tribocorrosion in total hip arthroplasty? A study of retrieved components. *Bone Joint J*. 2019;101-B:386–389.
16. Aldinger P, Williams T, Woodard E. Accelerated fretting corrosion testing of zirconia toughened alumina composite ceramic and a new composition of ceramicised metal femoral heads. Poster presented at: Orthopaedic Research Society Annual Meeting; 2017; San Diego, CA.
17. Smith+Nephew 2016. OR-16-127.
18. Ezzet K, Hermida JC, Steklov N, D'Lima D. Wear of polyethylene against oxidized zirconium femoral components: effect of aggressive kinematic conditions and malalignment in total knee arthroplasty. *J Arthroplasty*. 2012;27:116–121.
19. Good V, Ries M, Barrack RL, Widding K, Hunter G, Heuer D. Reduced wear with oxidized zirconium femoral heads. *J Bone Joint Surg Am*. 2003;85-A Suppl 4:105–110.
20. Morrison M, Jani S, Parikh A. Design of an advanced bearing system for total knee arthroplasty. *Lubricants*. 2015;3:475–492.
21. Papannagari R, Hines G, Sprague J, Morrison M. Long-term wear performance of an advanced bearing technology for TKA. Poster presented at: Orthopaedic Research Society Annual Meeting; 2011; Long Beach, CA.
22. Parikh A, Hill P, Hines G, Pawar V. Wear of conventional and highly crosslinked polyethylene liners during simulated fast walking/jogging. Poster presented at: Orthopaedic Research Society Annual Meeting; 2009; Las Vegas; NV.
23. Ries M, Salehi A, Widding K, Hunter G. Polyethylene wear performance of oxidized zirconium and cobalt-chromium knee components under abrasive conditions. *J Bone Joint Surg*. 2002;84-A (S2):129–135.

References

24. Hobbs LW, Rosen VB, Mangin SP, et al. Oxidation microstructures and interfaces in the oxidized zirconium knee. *Int J Appl Ceram Technol*. 2005;2:221–246.
25. Caicedo M, Pawar V, Hallab NJ. Oxidized Zr-alloy particles induce a lower incidence of in vitro lymphocyte metal-sensitivity responses compared to cobalt and titanium implant alloys. Poster presented at: Orthopaedic Research Society Annual Meeting; 2014; New Orleans, LA.
26. Rose SF, Weaver CL, Fenwick SA, Horner A, Pawar VD. The effect of diffusion hardened oxidized zirconium wear debris on cell viability and inflammation – an in vitro study. *J Biomed Mater Res Part B*. 2012;100B:1359–1368.
27. Pawar V, Jones B, Sprague J, Salehi A, Hunter G. Acidic fretting tests of oxidized Zr-2.5Nb, CoCr and SS femoral heads. In: Helmus M, Medlin D, eds. Medical Device Materials II: Proceedings of the Materials and Processes for Medical Devices Conference. ASMI; 2005;403–40.
28. Anderson FL, Koch CN, Elpers ME, Wright TM, Haas SB, Heyse TJ. Oxidized zirconium versus cobalt alloy bearing surfaces in total knee arthroplasty. *Bone Joint J*. 2017;99-B:793–798.
29. Bourne RB, Barrack R, Cecil RH, Salehi A, Good V. Arthroplasty options for the young patient. *Clin Orthop Relat Res*. 2005;441:159–167.
30. Brandt JM, Guenther L, O'Brien S, Vecherya A, Turgeon TR, Bohm ER. Performance assessment of femoral knee components made from cobalt–chromium alloy and oxidized zirconium. *Knee*. 2013;20:388–396.
31. Morrison ML, Jani S, Parikh, A. Development of an advanced bearing couple for total knee arthroplasty. Poster presented at: Orthopaedic Research Society Annual Meeting; 2011; Long Beach, CA.
32. Karidakis GK, Karachalios T. Oxidized zirconium head on crosslinked polyethylene liner in total hip arthroplasty: A 7- to 12-year in vivo comparative wear study. *Clin Orthop Relat Res*. 2015;473:3836–3845.
33. Hallab N, Merritt K, Jacobs JJ. Metal sensitivity in patients with orthopaedic implants. *J Bone Joint Surg*. 2001;83:428–436.
34. Hallab NJ, McAllister K, Jacobs JJ, Pawar V. Zirconium-alloy and zirconium-oxide particles produce less toxicity and inflammatory cytokines than cobalt-alloy and titanium-alloy particles in vitro, in human osteoblasts, fibroblasts and macrophages. Poster presented at: Orthopaedic Research Society Annual Meeting; 2012; San Francisco, CA.
35. Mistry JB, Chughtai M, Elmallah RK, et al. Trunnionosis in total hip arthroplasty: a review. *J Orthopaed Traumatol*. 2016;17:1–6.
36. Pulido R, Kester B, Schwarzkopf R. Trunnionosis in total knee arthroplasty: Is it a clinical problem? *Int J Orth*. 2017;4:837–840.
37. Cartner J, Aldinger P, Newman M. Characterization of tapers in TKA revisions from a 16-year retrieval database. Poster presented at: Orthopaedic Research Society Annual Meeting; 2016; Orlando, FL.
38. Goldberg JR, Gilbert JL, Jacobs JJ, Bauer TW, Paprosky W, Leurgans S. A multicenter retrieval study of the taper interfaces of modular hip prostheses. *Clin Orthop Relat Res*. 2002;401:149–161.
39. Carli AV, Patel AR, Cross MB et al. Long-term performance of oxidized zirconium on conventional and highly cross-linked polyethylene in total hip arthroplasty. *SICOT J*. 2020;6.

Abbreviations

AL	Alabama	IL-6	Interleukin-6
ASMI	American Society for Metals International	IL-8	Interleukin-8
BMI	Body mass index	kg	Kilogram
CA	California	kGy	kilo-Gray (unit of radiation)
Co	Cobalt	LA	Louisiana
CoCr	Cobalt chrome	lbf	Pound-force
CPE	Conventional polyethylene	MACC	Mechanically assisted crevice corrosion
CPM	Counts per minute	mg	Milligram
Cr	Chromium	mm	Millimetre
d.p.	Decimal place	MN	Minnesota
DH-OXINIUM	Diffusion-hardened OXINIUM	Mo	Molybdenum
EtOH	Ethanol	Mrad	Megarad (unit of radiation)
FDS	Femoral damage score	Nb	Niobium
FL	Florida	Ni	Nickel

Abbreviations

NR	Not reported	TNF-α	Tumour necrosis factor alpha
n.s.	Not significant	TX	Texas
NV	Nevada	UHMWPE	Ultra-high molecular weight polyethylene
OxZr	Oxidized Zirconium	5-XLPE/ 7.5-XLPE/ 10-XLPE*	Highly crosslinked polyethylene
PE	Polyethylene	Zr	Zirconium
SD	Standard deviation	Zr-2.5Nb	Zirconium-niobium alloy
SEM	Scanning electron microscopy	ZrO₂	Zirconium oxide
SS	Stainless steel	ZTA	Zirconia toughened alumina
THA	Total hip arthroplasty	μm	Micrometre
Ti	Titanium		
TiAlV/ TiAl4V/ Ti6Al4V	Titanium-aluminium-vanadium alloy		
TKA	Total knee arthroplasty		
TNF	Tumour necrosis factor		

The results of in vitro wear simulation testing have not been proven to quantitatively predict clinical wear performance. The results of in vitro cytokine expression analyses have not been proven to quantitatively predict clinical cytokine expression.

*XLPE is classified as ultra-high molecular weight polyethylene that has been irradiated by high dose (>50kGy) gamma or electron beam radiation. The 5-XLPE, 7.5-XLPE and 10-XLPE acetabular liners described within this collection are manufactured using 5 (50kGy), 7.5 (75kGy) and 10 (100kGy) Mrad doses of gamma radiation, respectively.

Smith & Nephew, Inc
1450 Brooks Road
Memphis,
TN 38116,
USA

www.smith-nephew.com

[◇]Trademark of Smith+Nephew
All Trademarks acknowledged
© 2022 Smith+Nephew
33572 V1 0622. Published June 2022.

Developed by Evidence Communications,
Global Clinical & Medical Affairs

Smith+Nephew