




Smith+Nephew

POLARCUP 
Dual Mobility System

VERILAST 
Oxidized Zirconium with XLPE

Product Information

POLARCUP[®] Dual Mobility System



Insert / cup articulation

Ball head / insert articulation

Third articulation


Neck / insert articulation

POLARCUP combined with a polished neck stem is highly recommended


Disclaimer:

The following product information is for informational and educational purposes only. It is not intended to serve as medical advice. It is the responsibility of treating physicians to determine and utilize the appropriate products and techniques according to their own clinical judgment for each of their patients. For more information on POLARCUP and VERILAST, including indications for use, contraindications, and product safety information, please refer to the label and the Instructions for Use packaged with the product.

Dual Mobility Articulation




Neutral position¹



Low level activity¹


Primary movement occurs in the ball head / insert articulation, allowing the insert to sit in its natural position




High level activity¹

Secondary movement occurs in the insert / cup articulation

Dual Mobility prostheses increase jump distance^{1,2}



Standard hip ø 50/28mm
28 mm head (28 mm diameter head)



POLARCUP® ø 53/28mm
28 mm head (53 mm insert = 47 mm diameter head)

- The Dual Mobility insert acts as a large head intended to increase the jump distance, designed to reduce the risk of dislocation and restoring hip stability.²⁻⁶
- The superior-posterior 6° skirt under the equatorial rim is designed to increase jump distance by ensuring additional stability.

POLARCUP[®] performance

POLARCUP performance demonstrates excellent clinical and functional outcomes^{1,7,8}

- Cementless POLARCUP 7A* ODEP Rating⁷
- Cemented POLARCUP 3A ODEP Rating⁷


POLARCUP is designed to reduce the risk of dislocation and increase joint stability¹⁻⁶

- Dislocation rate of 0.0 and 0.7% reported^{8,9}
- Increased range of motion (ROM) compared to conventional THA and constrained liners²

POLARCUP demonstrates low wear in vitro¹⁰

- Standard Polyethylene (UHMWPE) 8.87 mg/million cycles (StdDev. 0.59) after 5 million cycles
- Cross-linked Polyethylene (XLPE) 0.39 mg/million cycles (StdDev 0.41) between 2 and 5 million cycles

Comparison UHMWPE and XLPE



VERILAST[®] technology for hips

VERILAST Technology for Hips from Smith+Nephew uses the exclusive bearing combination of proprietary OXINIUM[™] and highly cross-linked polyethylene, which provides superior clinical survivorship¹¹ and biocompatibility¹²⁻¹⁴ without sacrificing versatility or introducing the risk of ceramic-like fracture.¹⁵

Wear performance



VERILAST Technology for total hip arthroplasty has been laboratory tested and shown to provide superior wear performance compared to CoCr on highly-crosslinked polyethylene, for up to 45 million cycles. With advanced materials designed to last, VERILAST Technology helps restore patients to their active lifestyles.

Corrosion avoidance

With its biocompatible properties, due to its use of oxidized zirconium, VERILAST Technology has been shown to reduce taper corrosion in total hip arthroplasty, minimizing the concern of trunnionosis.¹²⁻¹³

An article published in the HSS Journal showed that in a 22 year retrieval database, OXINIUM femoral heads are associated with decreased corrosion damage compared to CoCr femoral heads. Furthermore, chromium rich deposits which were present on some CoCr femoral heads were completely absent from OXINIUM.¹³

Cumulative volumetric wear comparison^{16*}



SEM image of worst case CoCrMo taper (at left) and worst case OxZr (at right).¹³

*Cumulative volumetric wear comparison components; XLPE groups consist of R3[®] shells and liners and the CPE group consists of REFLECTION[®] shells and liners.

Product features

Low wear

All versions of the POLARCUP® have a highly polished internal surface designed to minimize wear¹⁰

.....

Self-aligning implant

The self-aligning insert is designed to allow loading following the path of least resistance

.....

Surface

The Ti-Plasma and Ti-Plasma / hydroxyapatite surfaces have a high degree of surface roughness¹⁷

.....

Stability

The equatorial teeth and ribs on the cup are designed to prevent rotation and provide primary stability

.....

Cement

The cemented POLARCUP can be used with previously implanted reinforcement cages and rings¹⁸

.....



Additional fixation with Flanged POLARCUP

One cup, four combinations

The screws and pegs can be used in combination with the Ti Plasma / HA flanged POLARCUP to provide additional fixation



Flanges broken off using the flange cutter



Flanges bent over the acetabular rim without cortical screws




Flanges bent over the acetabular rim fixed with 1 or 2 cortical screws



Flanges bent over the acetabulum fixed with 1 or 2 cortical screws and 2 impacted anchoring pegs

Product overview



Implant Overview

Sizes mm	Ti-plasma/HA cementless cup with flanges	Ti-Plasma cementless cup	Stainless steel cemented cup	PE Insert Ø 22mm head	PE Insert Ø 28mm head	XLPE Insert Ø 22mm head	XLPE Insert Ø 28mm head
43	75100436	75100406	75100451	●		●	
45	75100437	75100407	75100452	●		●	
47	75100438	75100408	75100453		●	●	●
49	75100439	75100409	75100454		●	●	●
51	75100440	75100410	75100455		●	●	●
53	75100441	75100411	75100456		●	●	●
55	75100442	75100412	75100457		●	●	●
57	75100443	75100413	75100458		●	●	●
59	75100444	75100414	75100459		●	●	●
61	75100445	75100415	75100460		●	●	●
63	75100446	75100416	75100461		●	●	●
65	75100447	75100417			●	●	●
67	75100448	75100418			●	●	●

Smith & Nephew Orthopaedics AG
Oberneuhofstrasse 10d
6340 Baar
Switzerland

www.smith-nephew.com

◊Trademark of Smith+Nephew
All Trademarks acknowledged
©2019 Smith+Nephew
01098-en V5 10/19



References

1. Moore, C. and Orlandini, L. (2012). The Dual Mobility Concept in THA. *Bone and Joint Science*, 2004;4(3):1-5. **2.** Stroh, D.A., Naziri, Q., Johnson, A. J. and Mont, M. A. Dual-mobility bearings: a review of the literature. *Expert Rev Med Devices*, 2012; 9(1): 23-31. **3.** Hamadouche, M., Biau, D. J., Hutten, D., Musset, T. and Gaucher, F. The Use of Cemented Dual Mobility Socket to Treat Recurrent Dislocation. *Clin Orthop Relat Res*, 2010; 468:3248-3254. **4.** Guyen, O., Pibarot, V., Wegrzyn, J., Chevillotte, C., Carret, J. P. and Bejui-Hugues, J. The Use of Dual Mobility implants for Primary Total Hip Replacements in Patients at Risk of Dislocation. *Journées Lyonnaises de Chirurgie de la Hanche*, 2008 ; 163-170. **5.** Civinini, R., Carulli, C., Matassi, F., Nistri, L. and Innocenti, M. A Dual-Mobility Cup Reduces the Risk of Dislocation in Isolated Acetabular Revisions. *Clin Orthop Relat Res*, 2012; 470:3542-3548. **6.** Vasukutty, N. L., Middleton, R. G., Matthews, E. C., Young, P. S., Uzoigwe, and Minhas, T. H. A. The Double-Mobility Acetabular Component in Revision Total Hip Replacement. The United Kingdom Experience. *J bone Joint Surg Br*, 2012; 94-B(5): 603-8. **7.** assessed 10/19: <http://www.odep.org.uk/products.aspx> **8.** Bauchu, P., Bonnard, O., Cypres, A., Fiquet, A., Girardin, P. and Noyer, D. The Dual-Mobility POLARCUP: First Results from a Multicenter study. *Orthopaedics*, 2008; 12 (31): 97-99. **9.** Fiquet, A. and Noyer, D. POLARSYSTEM dual mobility hip prosthesis and minimally invasive surgery MIS. *Interactive surgery*, 2006; 22: 1-5. **10.** ENDOLAB Test Report – POLARCUP Dual Mobility System Std UHMWPE v's Crosslinked UHMWPE 2008. Smith & Nephew Internal Data **11.** Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR). Hip, Knee & Shoulder Arthroplasty: 2017 Annual Report. Adelaide: AOA, 2017. **12.** Pawar et al. Acidic Fretting Tests of Oxidized Zr-2.5Nb, CoCr, and SS Femoral Heads, ASMI, 2004. **13.** Cartner, J et al. Characterization of Femoral Head Taper Corrosion Features Using a 22-year Retrieval Database. *HSS Jnl*, 2017; 13: 35. <https://doi.org/10.1007/s11420-016-9517-5>. **14.** Hallab et al. Zirconium and Niobium Affect Human Osteoblasts, Fibroblasts, and Lymphocytes in a Similar Manner to More Traditional Implant Alloy Metals. *Journal of ASTM International*, 2006; 1 (3) **15.** Hunter, G et al. Creation of Oxidized Zirconium Orthopaedic Implants. *Journal of ASTM International*, 2005; 7 (2). **16.** Parikh et al. "Long-term simulator wear performance of an advanced bearing technology for THA," *Orthop Res Soc*, San Antonio, TX, J2013; 1028: 26-29. **17.** Massin, P., Orain, V., Philippot, R., Farizon, F. and Fessy, M. H. Fixation Failures of Dual Mobility Cups. A Mid-term Study of 2601 Hip Replacements. *Clin Orthop Relat Res*, 2012; 470:1932 -1940. **18.** POLARCUP Surgical Technique, 01620-en (1582).