+ Evidence in focus

Publication summary

Smith-Nephew

REGENETEN^o Bioinductive Implant in partial-thickness rotator cuff tears: 12-month results from a prospective multi-centre registry

Bushnell BD, Bishai SK, Krupp RJ, et al. Treatment of partial-thickness rotator cuff tear repairs with a resorbable bioinductive bovine collagen implant: 1-year results from a prospective multi-center registry. Orthop J Sports Med. 2021;9(8):23259671211027850.

Available at: Orthopaedic Journal of Sports Medicine

Key points



I

Early clinical outcomes were significantly improved with isolated bioinductive repair

repair in higher grade tears (p<0.05)

Overview

- An analysis of outcome data from patients with partial-thickness rotator cuff tears enrolled into a prospective registry study (the REBUILD Registry), conducted at 19 sites in the USA
- A total of 272 patients (mean age, 52.1 years) with partialthickness tears received a REGENETEN Implant:
 - 241 in an isolated bioinductive repair (without surgical repair)
 - 31 to augment a takedown and repair, all in grades 2 or 3 tears

Results

- All PROMs were significantly improved at 3, 6 and 12 months from preoperative values, except VR-12 MCS, which was only significant at 12 months (p<0.05; Table)
 - MCIDs for ASES, SANE and WORC scores were met or exceeded by >90% of patients at 12 months
- In higher grade tears (≥grade 2) and compared with the augmented takedown and repair group, isolated bioinductive repair resulted in:
 - Significantly better ASES, SANE and WORC scores at 2 and 6 weeks (p<0.05)
 - No significant differences in 12-month outcomes, except VR-12 PCS, which was significantly improved with isolated bioinductive repair (p=0.0213)
 - Significantly less sling time (19.1 vs 34.3 days; p<0.0001) and faster return to non-overhead sports (72.2 vs 128.9 days; p=0.0192); no other significant differences in recovery

 Preoperative tear size (Ellman classification): 49 grade 1 (<3mm), 101 grade 2 (3-6mm) and 122 grade 3 (>6mm tears)

- Outcomes included postoperative recovery and PROMs, which were assessed preoperatively and postoperatively at 2 and 6 weeks, and 3, 6 and 12 months
- Twelve-month follow-up data were available for 227 patients

Table. Mean PROM scores (± standard deviation) before and after treatment of partial-thickness rotator cuff tears with the REGENETEN Implant.*

	•••••••••••••••••••••••••••••••••••••••		
	Preoperative	3 months	12 months
ASES pain	5.5 ± 2.4	2.1 ± 2.4 [†]	1.1 ± 2.0 [†]
ASES function	14.1 ± 6.2	18.9 ± 6.9 [†]	26.1 ± 5.9 [†]
ASES overall	46.8 ± 18.2	71.9 ± 20.5 [†]	88.1 ± 17.9 [†]
SANE	41.7 ± 19.9	69.9 ± 19.4 [†]	86.2 ± 18.2 ⁺
VR-12 MCS	51.9 ± 13.2	54.6 ± 12.0	55.4 ± 8.9 [†]
VR-12 PCS	35.3 ± 8.4	43.1 ± 8.3 [†]	49.2 ± 9.3 [†]
WORC	36.4 ± 16.6	64.0 ± 22.6 [†]	83.7 ± 21.7 [†]

 * Overall population (n=227): includes both isolated bioinductive repair and augmented takedown and repair groups

⁺ p<0.02 vs preoperative values

Conclusions

In 227 patients with partial-thickness tears receiving a REGENETEN Implant either alone or with takedown and repair, significant improvements in pain, shoulder function and health-related quality of life were reported at 3 and 12 months. Isolated bioinductive repair offered improved early clinical outcomes and comparable 12-month outcomes to augmented takedown and repair.

Abbreviations

ASES = American Shoulder and Elbow Surgeons, MCID = minimal clinically important difference, PROMs = patient-reported outcome measures, SANE = Single Assessment Numeric Evaluation, VR-12 MCS = Veterans RAND 12-Item Health Survey Mental Component Score VR-12 PCS = Veterans RAND 12-Item Health Survey Physical Component Score, WORC = Western Ontario Rotator Cuff Index

Sports Medicine, Smith & Nephew, Inc. 150 Minuteman Road, Andover, MA 01810 USA. 28865 V2 0921 Published September 2021. ©2021 Smith+Nephew. All rights reserved. ®Trademark of Smith+Nephew. All Trademarks acknowledged.