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

The screenshot shows the SPATIALFRAME.COM website. At the top, there are navigation tabs for Home, Product Info, Seminars, and Contact Us. The main content area features a banner for the TAYLOR SPATIAL FRAME with a description and a 'Learn More' link. Below the banner, there is a 'SPATIALFRAME.COM VERSION 5' section with a 'Welcome!' message and several paragraphs of text. On the right side, there are two login options: 'Existing Users - Sign In' with fields for 'User Name' and 'Password', and a 'Sign In' button; and 'New Users' with a 'Request an Account' button. A red arrow points to the 'Request an Account' button with the text 'アカウントを申請するのにクリックしてください。' (Click here to request an account.)

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
The screenshot shows a pop-up form titled 'Request a SpatialFrame.com Account'. The form contains several input fields: 'First Name', 'Last Name', 'E-mail Address', 'Phone Number', 'Smith & Nephew Sales Rep', 'Hospital Affiliation', 'Country', 'Specialty', and 'Classification'. The 'Requested User Name' field is also present. There is a 'Notes' field with a scrollable area. At the bottom right, there are 'Cancel' and 'OK' buttons. A red arrow points to the 'OK' button with the text 'OKをクリックします。48時間以内に確認メールが送信されます。' (Click OK. A confirmation email will be sent within 48 hours.)

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
SPATIALFRAME.com smith&nephew

Home Product Info Seminars Contact Us

Precision alignment. Anatomic perfection. Total control.

TAYLOR SPATIAL FRAME[®]

The Taylor Spatial Frame is a unique external fixation system, providing treatment for a variety of skeletal fractures, malunions and nonunions.

Our easy-to-use software guides you every step of the way.

[Learn More >>](#)

SPATIALFRAME.COM VERSION 5

Welcome!

As the world's most advanced, versatile and clinically proven circular fixator, the TAYLOR SPATIAL FRAME™; system enables uncompromising stability, flexibility and precision. Or, in a word, control. With the TAYLOR SPATIAL FRAME external fixator and SpatialFrame.com, a one-of-a-kind web application that allows you to preoperatively plan or make fine, ongoing corrections postoperatively, you can achieve greater soft tissue control and quality control.

To begin using the TAYLOR SPATIAL FRAME, contact your local Smith & Nephew sales representative or select the "Request an Account" link on the right side of this screen.

For more information related to the product, contact [Drew Barry](#) at +1 (901) 399-5306.

For more information related to the product in Europe or Canada, contact [Marco Motta](#) at +39 039 6094236.

[Click here](#) for surgeon information about the TAYLOR SPATIAL FRAME.

For information on upcoming product training courses, please visit www.orthomeetings.com.

To be included in the TSF Surgeon Locator database, please visit www.tsflocator.com.

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The screenshot shows the SPATIALFRAME.com website. At the top, there are navigation tabs for Home, Product Info, Seminars, and Contact Us. The main banner features the Taylor Spatial Frame product with the text: "Precision alignment. Anatomic perfection. Total control. TAYLOR SPATIAL FRAME. The Taylor Spatial Frame is a unique external fixation system, providing treatment for a variety of skeletal fractures, malunions and nonunions. Our easy-to-use software guides you every step of the way. Learn More >>". Below the banner, there is a section titled "SPATIALFRAME.COM VERSION 5" with a "Welcome!" message and several paragraphs of text providing information about the product and contact details for Drew Barry and Marco Motta. On the right side of the page, there is a "Sign In" section for existing users with fields for "User Name" and "Password", a "Sign In" button, and a link for "Forgot your Username or Password?". Below that is a "New Users" section with a "Request an Account" button.

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The screenshot shows a "Password Reset" pop-up box. It contains the following text: "Please enter you username or email, and a temporary password will be sent to your email account registered with this site." Below this text are two input fields: "User Name:" and "E-mail Address:". Between these fields is the text "- or -". At the bottom of the box, there is a link for "helpdesk". At the very bottom, there are two buttons: "Cancel" with a red 'X' icon and "OK" with a green checkmark icon.

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Case Name:
 Date:
 Case ID:
 Case Notes:

► **Then, Choose Operating Mode:**
 Total Residual Chronic

[Continue to Step 2](#)

Help

Please Note: The images used in this software are graphical representations and are not actual illustrations of your case.

Total Residual
 Frame can be applied without determining a neutral frame position, referred to as "Crooked Frame on a Crooked Bone". Most commonly used. [More info...](#)

Chronic
 Frame matches the patient's deformity and ends in a neutral frame position. [More info...](#)

Notice: Per the Health Insurance Portability and Accountability Act of 1996, **no input fields should include the patient's full name or other patient identifiers.** The user takes full responsibility for non-compliance.

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Ankle

Select 6 x 6 Miter Region:

For 6x6 Miter Hindfoot, select a Tibia Mid Diaphysis and change to Dowel view on the Deformity page

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Preview:

Select 6 + 6 Region:

For 6+6 **Hindfoot**, select a Tibia Mid Diaphysis and change to Dowel view on the Deformity page

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- 6x6 Miter
- 6x6 Butt
- 6+6
- Ankle

Preview:


There is no region selection for an **Ankle Contracture**.

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First, Select Reference Fragment: Case: Test **Help**

Proximal Distal

Next, Enter Deformity Measurements:

AP View
Angulation: (degrees) Translation: (mm)
 Valgus Medial
 Varus Lateral
Med Lat

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Axial View
Angulation: (degrees) Translation: (mm)
 Internal Short
 External Long
Lat Med Lat

Apply Axial Translation First

Help
Reference Fragment - Non-moving fragment from which all deformity measurements are made relative to. In Dowel View, the Reference Fragment is always colored blue. [More...](#)
Origin - Frame can be applied without determining a neutral frame position, referred to as 'Crooked Frame on a Crooked Bone'. Most commonly used. [More...](#)
Corresponding Point - The point on the moving fragment that is coincident with the origin when the bones are reduced. [More...](#)
Views:
AP View Lateral View Axial View
AP View Angulation - Indicated in degrees of varus or valgus, this is the angle between the fragments on an AP radiograph. [More...](#)
AP View Translation - Distance in millimeters that the corresponding point is medial or lateral to the centerline passing through the origin. [More...](#)
[Iconography Key](#)

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[Go Back to Step 2](#) [Continue to Step 4](#)

Help

Selecting Frame Components: In this step, select the rings and struts to construct the patient's frame. Using the drop down boxes, select the type and size of ring. Using the table, select strut family, type, and size.

Ring Types:

Full Ring 2/3 Ring U Ring Foot Ring

Full Ring - Standard spatial ring which can be used in almost any application.

Open Ring Orientation: When using a 2/3 Ring or U-Plate, you can describe the location of the ring opening. This provides a more accurate graphic, and also defines the position of the 7 hole tab, which is used in calculations.

Inner Mount / Outer Mount: This strut mounting option allows for a more customized fit for your corrections. [View a diagram](#) of inner and outer mounts.


Strut Families: The system is comprised of two types of struts, Standard Struts and FAST FX Struts. [More info](#).

Custom Mount: The option to mount your struts in holes other than the default positions is available for some rings. The Custom Mount link will only display for configurations that support this technique.

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
SpatialFrame.com User Account Information

This is your user account information for using the SpatialFrame.com website.

First Name	Test
Last Name	1
E-mail Address	Test_1@smith-nephew.com
Office Phone	###-###-####
Smith & Nephew Sales Representative:	--
Hospital Affiliation:	--
Country:	United States
Specialty:	Limb Length/Deformity ▾
Classification:	S&N Employee ▾
Date Format Preference:	mm/dd/yyyy ▾
Enable alternate strut mounts:	Yes ▾
Change Password	
New Password	
Confirm New Password	

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

▶ First, Select Rings:


Case: TEST

Proximal Ring:	<input checked="" type="radio"/> Full Ring	155 mm, 7107-0114	
	<input type="radio"/> Custom Mount		
Distal Ring:	<input type="radio"/> U Ring	130 mm, 7107-1303 Inner Mount	<input type="radio"/> Between Struts 5 & 6

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► **First, Select Rings:** Case: TEST

Proximal Ring: Full Ring 155 mm, 7107-0114
Custom Mount

Distal Ring: Full Ring 155 mm, 7107-0114
Custom Mount

► **Next, Select Struts:**

	Standard Struts					Fast FX Struts			
	XXS 59-76 mm 7107-0200 Select All	XS 75-96 mm 7107-0205 Select All	Short 90-125 mm 7107-0210 Select All	Medium 116-178 mm 7107-0220 Select All	Long 169-283 mm 7107-0230 Select All	XS 91-121 mm 7107-0705 Select All	Short 116-152 mm 7107-0710 Select All	Medium 143-205 mm 7107-0720 Select All	Long 195-311 mm 7107-0730 Select All
Strut 1	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Strut 2	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Strut 3	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Strut 4	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Strut 5	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Strut 6	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

[Go Back to Step 2](#) [Continue to Step 4](#)

Help

Selecting Frame Components: In this step, select the rings and struts to construct the patient's frame. Using the drop down boxes, select the type and size of ring. Using the table, select strut family, type, and size.

Ring Types:

Full Ring 2/3 Ring U Ring Foot Ring

Full Ring - Standard spatial ring which can be used in almost any application.

Open Ring Orientation: When using a 2/3 Ring or U-Plate, you can describe the location of the ring opening. This provides a more accurate graphic, and also defines the position of the 7 hole tab, which is used in calculations.

Inner Mount / Outer Mount: This strut mounting option allows for a more customized fit for your corrections. [View a diagram](#) of inner and outer mounts.

Strut Families: The system is comprised of two types of struts, Standard Struts and FAST FX Struts. [More info...](#)

Custom Mount: The option to mount your struts in holes other than the default positions is available for some rings. The Custom Mount link will only display for configurations that support this technique.

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
► **First, Specify Neutral Frame Settings:** Case: Test

Neutral Frame Height: (mm) or Neutral Strut Length: (mm)

► **Next, View Final Settings & Frame Previews:**

	Strut 1	Strut 2	Strut 3	Strut 4	Strut 5	Strut 6
Initial Strut Settings:	111	134	160	155	134	106
Final Strut Settings:	136	136	136	136	136	136

Frame Previews View: Bone Dowel

Initial Frame

Final Frame

Help

Neutral Frame Height is the distance from the center of the proximal ring to the center of the distal ring with all struts at their neutral length.

Neutral Strut Length is obtained by measuring the distance from a ring to the interior end of it's fragment, for both proximal and distal fragments

The **Final Strut Settings** will indicate the setting of each strut on the last day of the correction program.

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Change-Out	Strut	Overlap Interval		Strut Change	
		First Day	Last Day	From	To
a	4 (Green)	0 (08/11/10)	2 (08/13/10)	7107-0720 Medium Fast Fx	7107-0710 Short Fast Fx
b	5 (Blue)	0 (08/11/10)	4 (08/15/10)	7107-0720 Medium Fast Fx	7107-0710 Short Fast Fx
c	3 (Yellow)	5 (08/16/10)	11 (08/22/10)	7107-0720 Medium Fast Fx	7107-0710 Short Fast Fx
d	4 (Green)	20 (08/31/10)	24 (09/04/10)	7107-0710 Short Fast Fx	7107-0705 Extra Short Fast Fx

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SPATIALFRAME.com login | smith&nephew

Cases **Test** Resources Account Contact Us

Save Save As Close Create New Case

1. Case Info 2. Deformity 3. Frame 4. Mount 5. Strut Settings 6. Duration / SAR 7. Prescription

Prescription: Case / Patient # Test **Help**

Your Contact Info for this Case: [edit](#) Your Notes for this Prescription: [edit](#)
Click to edit

The Report is intended for the doctor to keep in their records. It is a comprehensive list of all aspects of the correction.

Prescription **Report** Printable PDF Version opens in new window

Spatial Frame Report
The Report is intended to be a part of the Doctor's Records.

General

Anatomy	Bone	Region	Operative Mode
Left	Tibia / Fibula	Mid Diaphysis	Total Residual

Deformity Parameters

AP View Angulation	AP View Translation	Lateral View Angulation	Lateral View Translation	Axial View Angulation	Axial Translation
20.0° Varus	20.0 mm Lateral	10.0° Apex Posterior	10.0 mm Posterior	0.0°	15.0 mm Long

Frame Parameters

Proximal Ring	Distal Ring	Reference
155mm Ring (7107-0114)	155mm Ring (7107-0114)	Proximal

Struts

Strut 1	Strut 2	Strut 3	Strut 4	Strut 5	Strut 6
Medium Fast Fx @ Strut (7107-0720)	Medium Fast Fx @ Strut (7107-0720)	Medium Fast Fx @ Strut (7107-0720)	Medium Fast Fx @ Strut (7107-0720)	Medium Fast Fx @ Strut (7107-0720)	Medium Fast Fx @ Strut (7107-0720)

Mounting Parameters

AP View Frame Offset	Lateral View Frame Offset	Rotary Frame Angle	Axial Frame Offset
10.0 mm Medial to Origin	35.0 mm Posterior to Origin	0.0°	20.0 mm Proximal to Origin

Initial Struts

Strut 1	Strut 2	Strut 3	Strut 4	Strut 5	Strut 6
145.0	160.0	160.0	145.0	145.0	145.0

Final Struts

Strut 1	Strut 2	Strut 3	Strut 4	Strut 5	Strut 6
156.0	145.0	119.0	112.0	129.0	160.0

SAR Data

AP View SAR Offset	Axial SAR Offset	Lateral View SAR Offset
0.0 mm	0.0 mm	0.0 mm

Max Safe Distraction Rate (mm/day)	Correction Time
1.0	27 days

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Strut Change-Outs

Change-Out	Strut	Overlap Interval		Strut Change	
		First Day	Last Day	From	To
a	4 (Green)	0 (08/11/10)	2 (08/13/10)	7107-0720 Medium Fast Fx 	7107-0710 Short Fast Fx 
b	5 (Blue)	0 (08/11/10)	4 (08/15/10)	7107-0720 Medium Fast Fx 	7107-0710 Short Fast Fx 
c	3 (Yellow)	5 (08/16/10)	11 (08/22/10)	7107-0720 Medium Fast Fx 	7107-0710 Short Fast Fx 
d	4 (Green)	20 (08/31/10)	24 (09/04/10)	7107-0710 Short Fast Fx 	7107-0705 Extra Short Fast Fx 

Parts List

Part	Quantity
155mm Ring (7107-0114)	2
Fast Fx  Identification Band Kit (7107-0340)	1
Medium Fast Fx  Strut (7107-0720)	6
Short Fast Fx  Strut (7107-0710)	3
Extra Short Fast Fx  Strut (7107-0705)	1

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
The screenshot shows the SPATIALFRAME.com web application interface. At the top, there is a navigation bar with 'Cases' and 'Test' tabs. Below the tabs, there is a menu bar with 'Save', 'Save As', and 'Close' buttons highlighted in orange. The main content area shows a 'Prescription' form with fields for 'Your Contact Info for this Case' and 'Your Notes for this Prescription'. Below the form, there is a 'Spatial Frame Report' section with a table of 'General' information and a table of 'Deformity Parameters'.

Anatomy	Bone	Region	Operative Mode
Left	Tibia / Fibula	Mid Diaphysis	Total Residual

AP View Angulation	AP View Translation	Lateral View Angulation	Lateral View Translation	Axial View Angulation	Axial Translation
20.0° Varus	20.0 mm Lateral	10.0° Apex Posterior	10.0 mm Posterior	0.0°	15.0 mm Long

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
The screenshot shows a 'Save Case' dialog box. The text inside the dialog box reads: 'Do you want to save the case: Test.' Below the text are two buttons: 'Cancel' with a red 'X' icon and 'OK' with a green checkmark icon. A dotted line points from the 'OK' button to the text on the right.

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The screenshot shows a 'WARNING!' message box. The text inside the box reads: 'ELECTRONICALLY STORED SPATIAL FRAME FILES ARE WORKSPACE ONLY AND NOT PART OF THE PATIENT TREATMENT RECORD. PRIOR TO SAVING THIS DATA ELECTRONICALLY, YOU ARE ADVISED TO PRINT A HARD COPY FOR YOUR PERMANENT RECORDS.' Below the text is an 'OK' button with a green checkmark icon.

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The screenshot shows the SPATIALFRAME.com web application interface. At the top, there is a navigation bar with 'Cases' and 'Test' tabs. The 'Test' tab is active, and within it, the 'Close' button is highlighted with an orange circle. Other buttons in the navigation bar include 'Resources', 'Account', 'Contact Us', 'Save', 'Save as', and 'Create New Case'. Below the navigation bar, there is a breadcrumb trail: '1. Case Info > 2. Deformity > 3. Frame > 4. Mount > 5. Strut Settings > 6. Duration / SAR > 7. Prescription'. The main content area is titled 'Prescription:' and includes fields for 'Your Contact Info for this Case: edit' and 'Your Notes for this Prescription: edit'. A 'Help' section explains that the report is for the doctor's records. Below this, there are tabs for 'Prescription' and 'Report'. The 'Report' tab is active, showing a 'Spatial Frame Report' section with a 'Printable PDF Version' link. The report includes a 'General' section with a table of patient information and a 'Deformity Parameters' section with a table of surgical parameters.

Anatomy	Bone	Region	Operative Mode
Left	Tibia / Fibula	Mid Diaphysis	Total Residual

AP View Angulation	AP View Translation	Lateral View Angulation	Lateral View Translation	Axial View Angulation	Axial Translation
20.0° Varus	20.0 mm Lateral	10.0° Apex Posterior	10.0 mm Posterior	0.0°	15.0 mm Long

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The screenshot shows the 'Prescription' tab in the SPATIALFRAME.com interface. It includes a navigation bar with steps 1 through 7, and a table of strut measurements. The table has columns for Date, Day, and Struts 1 through 6, along with an Actions column. A red circle highlights a 'New Total Residual' icon in the Actions column for the date 08/18/10.

Date	Day	Strut 1	Strut 2	Strut 3	Strut 4	Strut 5	Strut 6	Actions
08/11/10	0	145	160	160	145 ^a	145 ^b	145	[Icons]
08/12/10	1	145	159	158	144 ^a	144 ^b	146	[Icons]
08/13/10	2	146	159	157	143 ^a	144 ^b	146	[Icons]
08/14/10	3	146	158	155	141	143 ^b	147	[Icons]
08/15/10	4	147	158	154	140	143 ^b	147	[Icons]
08/16/10	5	147	157	152 ^c	139	142	148	[Icons]
08/17/10	6	147	157	151 ^c	138	141	148	[Icons]
08/18/10	7	148	156	149 ^c	136	141	149	[Icons]
08/19/10	8	148	156	148 ^c	135	140	149	[Icons]
08/20/10	9	149	155	146 ^c	134	140	150	[Icons]
08/21/10	10	149	154	145 ^c	133	139	151	[Icons]
08/22/10	11	149	154	143 ^c	132	138	151	[Icons]
08/23/10	12	150	153	142	130	136	152	[Icons]

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

The dialog box is titled 'Create New Total Residual' and includes a 'Required fields in bold' note. It shows the 'Existing Case Name' as 'Test' and the 'New Case Name' as 'Test - Second Correction'. There are 'Cancel' and 'Create New Total Residual' buttons at the bottom.

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SPATIALFRAME.com logout

Cases **Test - Second Correction** Resources Account Contact Us

Save Save As Close Create New Case +

1. Case Info 2. Deformity 3. Frame 4. Mount 5. **Strut Settings** 6. Duration / SAR 7. Prescription

► **First, Enter Initial Strut Settings:** Case: Test - Second Correction **Help**

	Strut 1	Strut 2	Strut 3	Strut 4	Strut 5	Strut 6
Initial Strut Lengths:	147 (mm)	158 (mm)	154 (mm)	140 (mm)	143 (mm)	147 (mm)

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SPATIALFRAME.com

logout smith&nephew

Cases Test Resources **Account** Contact Us

Account Info

First Name

Last Name

E-mail Address

Office Phone

Smith & Nephew Sales Representative:

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SPATIALFRAME.com

logout smith&nephew

Cases Test Resources **Account** Contact Us

Resources

Tools

[Rings First Method/Direct Scheduler](#)
Advanced technique that is dependent on accuracy of interoperative ring placement

Links

[OrthoMeetings.com](#)
Smith & Nephew's comprehensive medical education meetings website.

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© 2022 Smith+Nephew
◊Trademark of Smith+Nephew