

Integra®

Integra® XT Revision
Ankle Replacement System

SURGICAL TECHNIQUE



INTEGRA®
LIMIT UNCERTAINTY

Table of Contents

Introduction	3
Indications/Contraindications	3
Instrumentation Concept	4
Preoperative Planning	5
Additional Information	6
Surgical Technique	7
Instruments	20
Implants	24
Disposables	26

Introduction

The Integra® XT Revision Ankle Replacement System is one of the first FDA cleared products indicated for Revision Ankle Arthroplasty only. This represents a strong commitment by Integra LifeSciences to further advance ankle arthroplasty and continue to offer our surgeon customers treatment options to further the continuum of care for ankle arthritis.

Thank you for putting your faith in Integra to provide you and your patients with an implant that replicates the natural ankle anatomy and a surgical procedure that you can trust.

Indications for Use

The Integra XT Revision Ankle Replacement System prosthesis is indicated as a total ankle replacement in revision surgeries only for patients with ankle joints damaged by severe rheumatoid, post-traumatic, or degenerative arthritis.

Components are intended for cemented use only.

Contraindications

Known contraindications to date:

- Sepsis.
- Infection sequelae.
- Systemic infection, fever and/or local inflammation.
- Complete talar necrosis.
- Insufficient quantity of bone stock or poor skin coverage around the ankle joint that would make the procedure unjustifiable.
- Persisting skin lesion.
- Important ligament laxity.
- Severe osteoporosis.
- Ankle arthrodesis with malleolar exeresis.
- Neuromuscular or mental disorders which might jeopardize fixation and postoperative care.
- Neurobiologic diseases.
- Nonfunctional lower limb muscles.
- Complete loss of ankle collateral ligament.
- Charcot's arthropathy.
- Elevation of WBC count.
- Distant foci of infection from genitourinary, pulmonary, skin and other sites, dental focus infection which may cause hematogenous spread to the implant site.
- Bone immaturity.
- Known allergy to one of the materials.
- Patient is pregnant.

Instrumentation Concept

The instrumentation is designed to achieve accurate and reproducible tibiotalar alignment while adapting to various anatomical conditions, depending on the lesions encountered in the ankle or a particular morphotype.

1) Patient Positioning

The patient is placed in a supine position with a bump under ipsilateral hip to reduce external rotation of the extremity. The heel is placed near the end of the table. A bump under the calf should be used throughout the surgery to keep the heel off the table.

2) Removal of Current Implant

The failed implant must be removed with minimal bone loss and with care to avoid fracture of the medial or lateral malleolus. Prophylactic pinning may be needed in the medial malleolus if excessive bone loss is expected. Note: Salto Talaris components are compatible with Integra XT components.

3) Tibial and Talar Preparation

Using a free hand technique or supplied tibial and talar cut guides, minimal bone should be removed to allow for appropriate resections for the selected implants. The tibial and flat cut talar cuts should be parallel to the weight bearing surface and remove just enough bone to maintain good bone stock. The sloped talar dome provides extra posterior augmentation, and this cut should be 14 degrees off the weight bearing surface from anterior to posterior.

4) Final Trial Positioning

Both tibial and talar trials can be inserted together. Lateral fluoroscopy should be utilized to verify correct position in the Sagittal plane. Once correct poly size is chosen the tibial or talar component can be secured first then the ankle can be ranged to locate the corresponding implant based on the patient's own biomechanics.

5) Final Preparation

With the final positioning completed, finish prepping the tibial fin resections and then remove the trials. The XT tibial tray is inserted first followed by the talar component, and lastly the poly insert is installed.

The Integra XT Revision Total Ankle Replacement System is designed to achieve proper positioning of the tibial implant in relation to the talar implant for a successful arthroplasty.

Preoperative Planning:

The preoperative planning for the Integra XT Revision Total Ankle Replacement System prosthesis is carried out using three standard weight-bearing radiological views:

- Anterior view;
- Anterior view with 30° internal rotation to expose the tibial-fibular joint space;
- Straight lateral.

Examination of the healthy side should be used for comparison.

Complementary imaging may be requested to:

- Confirm or reject the indication (CT scan examination for talar necrosis, a relative contraindication for prosthetic replacement);
- Discuss the need for an associated procedure (CT scan of the subtalar joint to determine if fusion is needed for this joint);
- Modify the technical details (assess the anterior deviation of the knee with panoramic x-ray).

Special consideration should be given to two types of pre-existing conditions.

- Malunions responsible for malalignment of the tibia or imbalance of the malleoli, which may require an initial correction.
- Major ligamentous instabilities demonstrated by an examination under stress will require specific intervention (release of the contracted side or possible need for an associated ligamentoplasty on the lengthened side).
- Previous ankle implant removal and subsequent bone loss.

1) Key planning elements determined from the anterior view:

- Choice of an implant size that does not impinge with the lateral malleolus;
- Determination of the ideal joint line level accommodating for articular wear.

Comparative images are often necessary to assess the prosthetic joint line, which should be located at the theoretical anatomic joint line. The thickness of the tibial resection depends on this determination.

2) Key planning elements determined from the lateral view:

- Confirmation of the implant size selected from the anterior view;
- Evaluation of the anterior osteophytic margin and assessment of the proposed bone resection required to expose the roof of the pilon;
- Evaluation of talar dome morphology, particularly its degree of convexity;
- Evaluation of talar positioning, which can be centered or retroplaced beneath the pilon. The relative positioning of the tibial and talar components should take into account a possible off-centered location with the understanding that the prosthesis adapts to this position and does not correct it.

In extreme cases, a pronounced anterior or posterior talar subluxation may preclude implantation of a prosthesis.

General Rules


- The Tibial component size is the same or one size bigger than the Talar component size. Only in the case of the Size O1 Insert, the Size O Tibial Component is smaller to its mating Size 1 Talar Component.
- The polyethylene (PE) insert matches the Talar component. For the size O talar dome use either the size O or size OO PE inserts. For the size 1 talar dome use either the size O1 or size 1 PE inserts.

Additional Information

- The tibial implant comes in 4 symmetrical sizes and 2 thicknesses, that can all be implanted on either the right or the left ankle.
- The Salto Talaris® flat cut talar component and standard keeled tibial tray are compatible with the Integra XT Revision System and surgical techniques for the flat cut talar dome and standard tibial tray are captured in the guide. Separate surgical techniques are available and should be referenced when implanting a Salto Talaris standard keeled Tibial Tray with a Salto Talaris flat or chamfer cut talar dome.

Recommendations about how to associate components with each other are mentioned in the following drawing:

Figure 1. Compatibility of components for Integra XT Revision Ankle Replacement System



Surgical Technique



As the manufacturer of this device, Integra does not practice medicine and does not recommend this or any other surgical technique for use on a specific patient. The surgeon who performs any implant procedure is responsible for determining and using the appropriate techniques for implanting the device in each patient.

Step 1 • Surgical Approach

- 1-1** The ankle is opened with an extensive longitudinal anterior incision lateral to the anterior tibialis tendon. This allows for an anterior release and broad arthrolisis with resection of all the osteophytes. Care should be taken to ensure that all scar tissue surrounding the failed implant has been removed.



Step 2 • Removal of Failed Prosthesis

General instruments may be used to remove failed prosthesis, i.e. an osteotome. If removing Salto components, follow below instructions:

Revision should begin by removing PE Insert. This is disassembled from the tibial base by inserting the Insert extractor blade between the base and the PE insert. A towel clamp holds the PE component for its extraction, after a lever maneuver using the extractor has separated the two components.

If necessary, the tibial base can then be removed as follows:

- To precut the bone around the tibial plug, use the osteotomes provided for this purpose in the instrumentation.
- Hook the posterior aspect of the tibial implant with the tibial component extractor.
- Insert the extractor plug by screwing the slap hammer on the tibial extractor.
- Push and pull vigorously with the slap hammer until the implant is fully removed.

The talar implant is separated from the talus with the thin osteotome.

Care should be taken to avoid any unnecessary bone loss while removing the prosthesis. Attention should be paid to the approach to avoid fracture of the malleoli. Prophylactic pinning may be needed in the medial malleolus if excessive bone loss is expected.

Note:

Must remove all implant components for competitor replacement ankle systems, as the Integra XT System is not compatible with competitor components.

Step 3 • Positioning of Tibial Alignment Guide

3-1a As an option, the 2mm cut guide may be used if the process of removing the prior implant has left the tibial resection relatively flat and at the desired angle.

If 2mm of bone is needed to be resected use Step 3-1b, otherwise follow this step.

The guide should be aligned parallel to the tibia's mechanical axis; this is a determining factor in all the resections performed during the procedure.

If the 7 degree block is desired to maintain the prior Salto Talaris tibial resection use the 7 degree block. It is recommended to use the 3 degree block and must be used for the talar resections later.

First check that all of the guide's set screws are unscrewed.

Drill a 110 mm drill tip pin at the anterior tibial tubercle in the neutral hole of the guide parallel to the tibial crest.

Align the guide to position 0 mm and place at the location of the tibia of the desired resection. In order to minimize bone loss a minimal resection is desired to maintain good quality tibial bone stock. The distal flange of the tibial alignment guide should rest on the anterior cortex of the tibia at the desired level of bone resection. The new resections should not be more than 2mm above the previous resection. Lateral fluoroscopy can assist in determining the desired level of resection.

Once the level of tibial resection is determined, tighten the extension screw located in the middle, then insert a second 110 mm drill tip pin through the guide's medial hole, positioning the guides axis in the center of the inferior metaphysis.

3-1b Place the 2mm Cut Guide (Fig. 3-1b) in the joint space with the platform placed firmly against the resected tibial resection. The joint distractor (MJU928T or MJU929T) should be used to keep the platform flush to the tibial resection. Pin the 2mm Cut Guide to the distal tibia going through the two proximal holes using the 3mm self drilling pins (LJV527T). Insert two 75mm pins (LJV526T) into the distal medial and lateral holes of the guide to protect against saw excursion. Take a lateral fluoroscopic image to ensure the platform is flush to the tibial resection. Perform the additional tibial resection extended carefully to the posterior cortex. A 1.27mm thick saw blade is recommended. Repeat this process as needed until the minimum amount of bone has been resected. Remove the 2mm Cut Guide and pins after the desired bone is resected. Place the tibial alignment guide appropriately and continue to Step 4.

Note: Check final position in ML and AP fluoroscopic views.



Items Used:

MJV553T - Tibial Alignment Guide, No Block
MJU668T - 3 Degree Block, Tibial Alignment
LJV527T - 110mm Self-Drilling Pins
MJV552T - 7 Degree Block, Tibial Alignment
LJV526T - 75mm Pin
MJU928T/MJU929T - Joint Distractors

Step 4 • Adjusting the Alignment Guide

4-1

Frontal plane:

The axis of the tibial resection guide should be made parallel to the tibia's mechanical axis by choosing the proper hole of the proximal pin guide.

Sagittal plane:

With both flanges in contact with the tibia, the resection guide is adjusted parallel to the anterior tibial crest. At this stage, a genu varum or a genu valgum deformation can be corrected by moving the proximal guide medially or laterally over the pin, making it possible to implant the prosthesis strictly perpendicular to the tibial axis to compensate for any axial defect. The timing and degree of this compensation should be discussed for each case as there may be a need for further compensation with additional intervention (i.e. secondary knee surgery or subtalar joint stiffening) to achieve horizontal axis. Once the guide is positioned in the frontal and sagittal planes, the set-up is finalized by tightening the set screw of the superior guide and then the set screw on the distal guide.

Rotational and mediolateral positioning:

The tibial alignment jig, used for mediolateral and rotational adjustment of the implant, is attached to the tibial alignment guide.

Adjustments are made as follows.

- Hold the ankle in a neutral position.
- Rotational adjustment: insert a 110 mm pin into the guide's adjustable arm; the pin in the adjustable arm should also be in line between the 2nd and 3rd metatarsal.

Once the rotational position has been adjusted, the guide's rotational adjustment set screw is tightened.

- Mediolateral adjustment: the tibial implant size planned preoperatively is confirmed through a series of lateral and medial holes on the guide.

The different implant sizes available (0, 1, 2, and 3) are on the guide; hence the size is confirmed by inserting two 75 mm pins in the medial and lateral holes and by positioning them at the M-L joint spaces formed by the malleoli. Once the mediolateral position has been adjusted, the guide's mediolateral adjustment set screw is firmly tightened with the 3.5mm Hex Driver. Check final position in AP fluoroscopic view. Remove the Tibial Alignment Jig from the Tibial Alignment Guide.

Note: This does not commit you to the selected size.

Caution:

Since this guide is not a cutting guide, do not drill through the holes. The pins inserted in the holes are used only to verify that the tibial plate is properly positioned. They are inserted in the holes but not drilled.

4-1



Items Used:

MJU334T - Tibial Alignment Jig

LJV526T - 75mm Pins

MLN113T - Hex Driver, 3.5mm

LJV527T - 110mm Self-Drilling Pin


Step 5 • Placing the Cutting Guide

Note: Skip this step if the 2mm Cut Guide is used.

- 5-1** Depending on the size chosen at the preoperative planning stage and in accordance with the size determined from the tibial alignment jig, tibial cut guide number 0, 1, 2 or 3 is chosen. This unit is attached to the alignment guide by tightening the set screw.

Caution:


Once all the adjustments have been made and before using the oscillating saw, make sure that the guide is sitting on the anterior tibia and all the set screws have been firmly tightened with the screwdriver provided in the instrumentation. An AP fluoroscopy image may be utilized to ensure proper alignment of the cut guide. A lateral fluoroscopy image can also confirm the level of resection of the tibia.



Step 6 • Preparation for Tibial Bone Cut

Note: Skip this step if the 2mm Cut Guide is used.

- 6-1** Drill the pin holes on the guide with 75mm pins. These converging pins protect the sweep of the oscillating saw blade during the horizontal cut as well as the proximal migration of the reciprocating saw blade during vertical cuts.



Step 7 • Tibial Cut

Note: Skip this step if the 2mm Cut Guide is used.

Note: If necessary, a pin pusher tool is available and can be used to further advance the pins through the guide after drilling.

- 7-1** The horizontal tibial resection is performed with an oscillating saw blade extending carefully to the back, just through the posterior cortex. If greater depth is needed, use the wide saw blade.

Care should be taken to avoid over penetration of either oscillating or reciprocating saw blade beyond the posterior cortical bone. Be sure to use only the oscillating saw for the horizontal cut and only the reciprocating saw for the vertical cuts.

If needed, the vertical resection is performed through the slots on the medial and lateral aspects of the guide with the provided end cutting reciprocating saw. Care needs to be taken to avoid unwanted contact to the talus when using the one-sided or two-sided reciprocating saw. If using the one-sided reciprocating saw, point the teeth distally. The cut should extend up to the pins that are seated in the guide.

Once the cuts have been made, remove the guide and remove the anterior half of the distal bone. Care should be taken not to leverage the malleoli while completing the cuts and removing bone. At this stage, the goal is to be able to straighten the foot at a right angle below the tibia.



Caution:

Point the saw blade's teeth distally to prevent wear of the cut guide.

Caution:

Pins may dislodge while cutting and care must be taken to retain pin location.

Items Used:

MJU645T-MJU648T - Tibial Cut Guide, Size 0-3
Saw Blades - Listed on Pg. 24
MLN113T - 3.5mm Hex Driver
LJV526T - 75mm Pins

Step 8 • Preparing the Talar Cut

- 8-1** Select the sloped talar cutting guide or flat talar cutting guide.

First, make sure the cut block is labeled as 3 degrees.

The talar cutting guide is attached to the tibial alignment guide and positioned as close to the talus as possible. Dorsal excision of the talar head with general instrumentation may be necessary prior to attaching the talar cutting guide to achieve proper resection. The set screw is tightened. With the foot in neutral the level of resection can be adjusted based on the desired amount of talus to be removed utilizing the height adjustment of the tibial alignment guide. The resection should be at the level of the anterior border of the talar articular surface, or at the level of the previous chamfer. Place a free saw blade through the guide while taking a lateral fluoroscopic view to properly check the resection level. Make sure the foot maintains a neutral position. Joint distractors may be used to properly position the foot.

8-1



Step 9 • Talar Cut

- 9-1** With the desired resection level determined and the foot in a neutral position, pin the foot using two 75mm or 110mm pins.

Make sure foot is not in varus or valgus before setting pins.

Pin positioning will be variable depending on patient size and resection level. Check final position in ML fluoroscopic view, before cutting.

The resection is completed through the talar cutting guide utilizing the narrow oscillating saw blade.

To protect the malleoli from the sweep of the saw blade, a set of ribbon retractors are provided in the instrumentation.

After cutting the talus remove the tibial alignment guide and the sloped talar cutting guide or flat talar cutting guide and any optional joint distractors. Remove the resected bone.

Ensure the talus has smooth posterior bony surface. General instrumentation may be used as needed.

9-1



Items Used:


MJU985T - Sloped Talar Cutting Guide
LJV527T - 110mm Self-Drilling Pins
MJUo86T - Ribbon Retractors
MJU928T - Joint Distractor, Right
MJU929T - Joint Distractor, Left
MJU910T - Talar Cut Guide, Flat Cut
Saw Blades - Listed on Pg. 24
LJV526T - Pin, 75mm

Step 10 • Trialing and Tibial / Talar Prep

At this step the tibial alignment guide and cutting guides should be removed.

Size Selection

To prevent any risk of medial malleolus fracture, please adhere to the following combinations: do not use thicker tibial baseplate (10 mm) with a poly insert thickness higher than 15 mm (see table below). While trialing select the appropriate thickness poly for an adequate fit. A 10mm Tibial Baseplate is recommended to increase medial and lateral area of titanium plasma coating.


		Polyethelene Insert Thickness (mm)								
		8	9	10	11	13	15	17	19	21
	Tibial baseplate - 4 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Tibial baseplate - 10 mm	✓	✓	✓	✓	✓	✓	NOT ALLOWED		

The plastic trial insert is selected depending on:

- The size and side of the plastic trial insert should match the talar component. For a size O talar dome, use either the size O or size OO insert based on tibial component size. For a size 1 talar dome, use either the size O1 or size 1 insert based on tibial component size. Please reference the chart on pg. 6 for complete component sizing compatibility.
- Thickness: they vary from 8 to 21 mm and correspond to the accumulated thickness of the metallic base and the PE

10-1 When using the sloped or flat cut talar dome:

Insert the selected size of trials for the tibial trial, the PE insert trial, and the sloped or flat cut talar trial. The PE insert trial will snap into place on to the tibial trial. Articulate the foot to find the proper positioning of the trials. Use the handle attached to the talar trial to control rotation and medial/ lateral positioning. For sagittal plane reference use a lateral fluoroscopic view to ensure that the components are centered under the long axis of the tibia. Also check that the handle is pointing between the 2nd and 3rd ray. Once satisfied with the positioning use step **10-1a. for a sloped talar dome** and **10-1b. for flat cut talar dome**.



10-1a When using the sloped talar trial, drill one talar peg hole with a 4.5mm drill, drilling to hard stop, followed immediately by filling the hole with a stabilization peg by hand or with stabilization peg removal tool, then repeat for the second talar peg hole. Do not tamp the stabilization pegs too far (Fig. 10-2). Leave the nail head above the surface of the dome to allow removal. Remove the talar trial handle. **For prepping the standard tibial trial**, drill a 75mm pin into the distal hole to hold the tibial trial base. Next, drill the proximal 3.0mm hole and leave empty. For the larger proximal hole, use the 7.9mm drill bit and drill to its stop.

For prepping the finned tibial trial use the four 3.0 mm holes. Drilling unicortically, drill the two most proximal holes with 45 mm pins, and drill the two distal holes with 75 mm pins. The 45 mm pins should be flush to the tibial trials anterior face to make room for the tapered drill. With the tibial trial secured, drill the large proximal hole with 5 mm tapered drill.

10-1b Remove the stabilization pegs with the stabilization peg remover, and remove all the trials and then skip to Step 11.

When using the Flat talar trial, drill the converging pins with a 75mm pin in the medial hole followed by another 75 mm pin in the lateral hole. Remove the talar trial handle. **For prepping the standard tibial trial**, drill a 75mm pin into the distal hole to hold the tibial trial base. Next, drill the proximal 3.0mm hole and leave empty. For the larger proximal hole, use the 7.9mm drill bit and drill to its stop. **For prepping the finned tibial trial**, use the four 3.0 mm holes. Drilling unicortically, drill the two most proximal holes with 45 mm pins, and the two distal holes with 75 mm pins. The 45 mm pins should be flush to the tibial trials anterior face to make room for the tapered drill. With the tibial trial secured, drill the large proximal hole with 5 mm tapered drill.


Remove only the tibial trial and insert trial. Finish prepping for the talar stem by plantar flexing the foot and inserting the appropriate drill guide for the stem. Next, drill for the stem using the appropriate stem drill to its hard stop. Next, drill for the talar flange using the Talar Flange Drill. Remove the flat cut trial and finish prepping the talar flange with the tibial keel osteotome to the approximate depth of the drill holes. Remove the pins and the talar trial, and then skip to Step 11.

Consideration must be given to possible adjunct soft tissue balancing procedures at this stage (i.e. achilles tendon lengthening, ligament release and repair).

Items Used:

MJU983T - Drill, 4.5mm
 MJV210T-MJV213T - Tibial Trial Jig, Size 0-3, Thick
 MJU981T - Stabilization Peg
 MJU960T-MJU963T - Sloped Talar Trial, Size 0-3, Right
 MJU970T-MJU973T - Sloped Talar Trial, Size 0-3, Left
 MJV200T-MJV203T - Tibial Trial Jig, Size 0-3
 MJU342T - Handle, Lateral Chamfer Guide
 Trial Inserts - Listed on Pg. 19
 MJU939T - Drill, Tapered
 DWD060T / LJV528T - 3mm Drill Bit

MJU380T, MJU384T-MJU386T - Standard Tibial Trials, Size 0-3
 MJV500T-MJV503T - Extended Tibial Trials, Size 0-3
 LJV525T - 45mm Pins
 MWE044T - Remover, Stabilization Peg
 MJU387T - Osteotome, Tibial Keel
 MJV940T, MJU941T-MJU943T - Talar Trial, Flat Cut, Size 0-3, Right
 MJV950T, MJU951T-MJU953T - Talar Trial, Flat Cut, Size 0-3, Left
 LJV526T - 75mm Pins
 MJU362T - Drill Bit, Size 0
 MJU922T - Bell Saw Reamer, Flat Cut
 MJU907T - Tibial Osteotome
 MJU925T - Talar Flange Drill



Note: For size 0 flat cut talar preparation, the size 0 drill bushing (MJU662T) needs to be used along with the size 0 talar drill (MJU362T)

Note: If utilizing a standard tibial tray, skip to Step 15 for instructions on final preparation and implantation.

Step 11 • Final Tibial Preparation - Finned Tibial Tray

11-1 To complete the tibial preparation the tibial window is completed utilizing the thin osteotome to connect the drill holes through cortical bone and finalize distal shape for implant insertion. The tibial trial rasp is then utilized to prep the distal tibia for the final implant. It is inserted into the resection and impacted in and out utilizing the tibial holder and mallet. Care should be made to ensure that the tibial rasp is flush with the inferior and anterior aspect of the tibial resection.


It is recommended to use either the tibial osteotome or the optional Antrum Rasp to remove bone from the intramedullary canal.

Make sure the alignment of the tibial osteotome and the pointed end of the Antrum Rasp are in line with the tibial axis. If necessary, the angle of the Antrum Rasp allows for access to the anterior portion of the intramedullary canal.

To ensure bone removal up to the height of the finned tibial implant, the Antrum Rasp should be inserted until the end of the rasps teeth are at the level of the distal tibial resection.

Note:

The Integra XT is indicated for cemented use only. Cement needs to be applied to the bone contact surfaces of the tibial and talar component: the cement is applied to the top of the tibial tray component and to the bottom surface of the talar dome.



Items Used:

MJU357T - Thin Osteotome
MJV100T-MJV103T - Tibial Trial Rasp, Size 0-3
400236 - Wiener Antrum Rasp
MJU907T - Tibial Osteotome

Step 12 • Tibial Implant Insertion

12-1 The tibial implant must be inserted first.

To implant the tibial component, attach the tibial holder onto the tibial implant.

Position the fin of the tray in the intramedullary canal and rotate the tibial tray into the joint space. At the same time, insert proximally in the intramedullary canal as far as possible.


If a gap still remains between the tibial resection and tibial component use either the Tibial Impactor and Punching Tip Assembly (Fig 12-1) or the Optional Tibial Impactor (Fig 12-2) to impact the tray to its final position.

The tibial holder may be removed at this time or remain in place during impaction.

A hole is positioned in the middle of the tibial keel: it is useful for the surgeon in order to check a proper lateral view with X-rays.

Items Used:

MJU936T - Tibial Holder
MJU937T - Tibial Impactor
MJU938T - Punching Tip
MJU085T - Curette
1129 - Tibial Impactor



Step 13 • Talar Implant Insertion

13-1 When inserting the talar implant, care should be taken to avoid scratching or marking of the talar implant articulating surface. A dedicated instrument, the talar impactor, helps to protect the talar surfaces.

The talar component is then manually inserted into place. For the sloped talar dome, align the 4.5mm pegs into the predrilled holes and impact with the talar impactor.

Lateral fluoroscopy can be utilized to check for proper placement.

Caution:

If contact is made with the surface of the talar dome, apart from the use of the talar impactor, replacement of the talar dome is recommended.

Items Used:

MJU351T - Talar Component Impactor



Step 14 • Poly Insert Insertion

14-1 The PE inserts may be retried after the tibial tray and talar dome have been placed. Refer to the compatibility chart on **Pg. 6** and the sizing guide in **Step 10** to guide selection. The trial will not snap into place on to the tibial tray, as it does with the tibial trial. Once trialing is complete, the holding clamp may be used to remove the trial. If retrying is performed, utilize lateral fluoroscopy to ensure implants have retained proper placement.

Adjust the thickness dial on the in situ poly inserter to the selected PE insert's thickness. Use the attached insert spacer for the 4mm thick tibial tray, but if using the 10mm thick tibial tray, switch this part to the XL insert spacer (to reference spacers, see Fig. 14-1a). (Note only change the spacer when the thickness dial is fully located to one position.) Manually insert the PE insert one third of the way into the tibial tray, making sure the dove tail of the poly is engaged. Attach the in situ poly inserter to the tibial component by threading the attachment screw to the center of the keel. Ensure this screw is tight.

The hex driver may be used to ensure that the attachment screw is tight.

Note:

If the inserter does not attach to the tibial implant, check if additional bone removal is needed in the anterior tibial window.

Caution:

If pressure is not maintained to keep the in situ poly inserter perpendicular while inserting the PE insert, it may bend the inserter's attachment screw or prevent proper poly insertion.

Place the foot against the belly to ensure a neutral position and place a hand on the top of the handle keeping pressure to maintain perpendicular attachment to the tray. Using the T-handle, advance the plastic spacer against the PE insert. Slowly advance the PE insert into the tibial tray while maintaining the pressure on the in situ poly inserter handle towards the foot to prevent the tool from bending on the tibial tray. Stop advancing when the face of the PE insert is flush to the anterior face of the tibial tray. At this point a small gap will exist between the tray and PE insert. Release pressure on the T-handle and remove the attachment screw, and then remove the in situ poly inserter. To close the gap, use thumb pressure or a bone tamp. Be careful to only impact on the anterior chamfer of the PE insert without touching the talar dome. The gap must be fully closed.

Fill in the tibial window with bone graft.

Caution:

If the PE insert does not properly engage the dove tail, remove the PE insert and inspect the dove tail area for damage, and if acceptable, retry the insertion technique. Do not reinsert the PE insert more than twice.

14-1a



Standard Spacer



XL Spacer

14-1b



Items Used:

MJU351T - Talar Impactor
 MJU940T - In Situ Poly Inserter
 MJU984T - XL Insert Spacer
 MJU048T - Holding Clamp
 MLN113T - 3.5mm Hex Driver

Note: Follow Steps 15 and 16 when utilizing a standard tibial tray with the sloped talar dome.

Step 15 • Final Tibial Preparation - Standard Tibial Implant

Note:

When using the tibial keel osteotome or tibial flange rasp, the depth markings are accurate for either standard or XL tibial trays.

15-1 Remove all trials and pins.

The tibial holes are rejoined using a small osteotome; then the thickness and depth of the resection is checked with the engraved lines of the tibial keel osteotome.

The distal part of the anterior groove of the tibia is beveled using the tibial flange rasp, so that the tibial implant lies flush on the resection.


With the different tibial implant sizes (0, 1, 2 and 3) marked on the upper surface of the rasp, the trimming done in this manner perfectly matches the length of the implant selected.

Items Used:

MJU357T - Thin Osteotome

MJU387T - Tibial Keel Osteotome

MJU350T - Tibial Flange Rasp



Step 16 • Placing Final Implants - Standard Tibial Implant

Note:

The Salto Talaris is indicated for cemented use only. Cement must be applied to the bone contact surfaces of the tibial and talar component.

16-1 The size and side selected during the implant trials must be retained.

The talar implant is placed first. When inserting the talar implant, care should be taken to avoid scratching or marking of the talar implant articulating surface. A dedicated instrument, the talar impactor, helps to protect the talar surfaces.

The talar component is then manually inserted into place. For the sloped talar dome, align the 4.5mm pegs into the predrilled holes and impact with the talar impactor. Lateral fluoroscopy can be utilized to check for proper placement.

Caution:

If contact is made with the surface of the talar dome, apart from the use of the talar impactor, replacement of the talar dome is recommended.

16-2 Assembly of the insert on the tibial implant:

The polyethylene insert is assembled to the tibial implant with the insert assembly press.

The tibial implant is mounted onto the press. Adjust the metal bar on the base of the press to the appropriate poly thickness. The PE insert is positioned on the implant and manually pushed to engage approximately 1/3 of the depth. Close the press arm over the assembly and continue to press down until the PE is fully engaged. The tibial assembly must be visually checked by the operator. Generally, an audible “click” is heard.

16-2



16-3 Use of the tibial impactor and method of tibial impact:

Using the tibial impactor, grasp the implant assembly between the metallic plug and the central anterior zone of the tibial tray.

The tibial component is impacted until the position of the tibial trial is reproduced.

During tibial implant impaction, maintain good contact between the superior side of the implant and the tibial resection to prevent any risk of a posterior gap between the tibial cut and the implant.

16-4 Flexion/Extension movements are applied to check the ankle kinematics. Apply bone graft to the tibial window.

Items Used:

- MJU351T - Talar Component Impactor
- MJV549T - Insert Assembly Bench Press
- MJU937T - Tibial Impactor

16-3



Step 17 • Revising or Removing Implants

17-1 If the implant must be revised, revision should begin by removing the PE insert. This is disassembled from the tibial base by inserting the blade of the insert extractor between the base and the PE. The Holding Clamp holds the PE component for its extraction, after a lever maneuver using the extractor has separated the two components.

The talar implant is separated from the talus with the thin osteotome.

For finned tibia:

If necessary, the tibial base can then be removed first by loosening the tibial implant with general instruments. Then attach the Tibial Holder to the tibial implant and extract.

For standard tibia:

To precut the bone around the tibial plug, use the osteotomes provided for this purpose in the instrumentation. Hook the posterior aspect of the tibial implant with the tibial component extractor. Tighten the screw on the tibial extractor to engage the tibial keel. Insert the extractor plug by screwing the slap hammer on the tibial extractor. Push and pull vigorously with the slap hammer until the implant is fully removed.

Items Used:

MJV545T - Insert Extractor
MJU048T - Holding Clamp
MJU357T - Osteotome, Thin
MJU936T - Tibial Holder
MJV556T - Tibial Implant Extractor
MJU356T - Osteotome, tibial barrel, revision
MJU358T - Slap hammer


Instrumentation Case Reference: INTEGRAXTREVBOX1



Caddy



Top Contents



Bottom Contents

INTEGRAXTREVBBOX1 Contents


Top Tray Contents

Item	Reference	Instrument
1	MJU291T	Screw, M4
2	MJU073T	Thumb Screw, M5
3	MJU086T	Ribbon Retractors
4	MJU365T	Pin Pusher
5	MJU981T	Stabilization Peg
6	MJU944T	Handle Screw, In Situ Poly Inserter
7	MJU357T	Osteotome, Thin
8	MLN113T	Hex Driver, 3.5 mm
9	MJU668T	3 Degree Block, Tibial Alignment
10	MJV552T	7 Degree Block, Tibial Alignment
11	MJV553T	Tibial Alignment Guide, No block
12	MJU334T	Tibial Alignment Jig
13	MJU645T	Tibial cut guide, Size 0
13	MJU646T	Tibial cut guide, Size 1
13	MJU647T	Tibial cut guide, Size 2
13	MJU648T	Tibial cut guide, Size 3
14	MJU011T	2 mm Cut Guide
15	MJU085T	Curette
16	MJU928T	Joint Distractor, Right
16	MJU929T	Joint Distractor, Left
17	MJU358T	Slap Hammer
18	MJU359T	Pin Puller
-	CSA-1030-0011	Caddy, Revision
-	CSA-1030-0006	Caddy Lid
-	CSA-1030-0009	Revision Instrument Tray 1


Bottom Tray Contents

Item	Reference	Instrument
19	MJU342T	Handle, Lateral Chamfer Guide
20	MJU362T	Dill Bit, Size 0, Talar Stem
21	MJV555T	Fixation Plug, Size 0
22	MJU662T	Drill Guide, Bell Saw, Size 0
23	MJU351T	Talar component impactor
24	MJU910T	Talar Cut Guide, Flat Cut
25	MJU922T	Bell Saw Reamer, Flat Cut
26	MJU924T	Talar Plug, Flat Cut
27	MJU930T	Extraction adaptor, Flat Cut
28	MJU925T	Drill, Talar Flange, Flat Cut
29	MJV940T	Talar Trial, Flat Cut, Size 0, Rt
29	MJU941T	Talar Trial, Flat Cut, Size 1, Rt
29	MJU942T	Talar Trial, Flat Cut, Size 2, Rt
29	MJU943T	Talar Trial, Flat Cut, Size 3, Rt
29	MJV950T	Talar Trial, Flat Cut, Size 0, Left
29	MJU951T	Talar Trial, Flat Cut, Size 1, Left
29	MJU952T	Talar Trial, Flat Cut, Size 2, Left
29	MJU953T	Talar Trial, Flat Cut, Size 3, Left
30	MJU960T	Sloped talar trial Size 0, Right
30	MJU961T	Sloped talar trial Size 1, Right
30	MJU962T	Sloped talar trial Size 2, Right
30	MJU963T	Sloped talar trial Size 3, Right
30	MJU970T	Sloped talar trial Size 0, Left
30	MJU971T	Sloped talar trial Size 1, Left
30	MJU972T	Sloped talar trial Size 2, Left
30	MJU973T	Sloped talar trial Size 3, Left
31	MJU940T	In situ Poly inserter
32	MJU984T	XL Insert Spacer
33	MJU663T / MJV549T	Insert Assembly Bench Press
34	MJU937T	Tibial impactor
-	CSA-1030-0007	Revision Instrument Case 1

Instrumentation Case Reference: INTEGRAXTREVBBOX2



Top
Contents



Bottom
Contents

INTEGRAXTREVBBOX2 Contents

Top Tray Contents

Item	Reference	Instrument
35	MJU907T	Tibial osteotome
36	MJU983T	Drill, 4.5mm
37	MJU939T	Drill, Tapered
38	MWE044T	Remover, Stabilaztion Peg
39	MJU938T	Punching tip
40	MJU048T	Holding Clamp
41	MJU985T	Sloped talar cutting guide
42	MJU936T	Tibial holder
43	MJV200T	Tibial trial jig Size 0
43	MJV201T	Tibial trial jig Size 1

Item	Reference	Instrument
43	MJV202T	Tibial trial jig Size 2
43	MJV203T	Tibial trial jig Size 3
44	MJV210T	Tibial trial jig Size 0, Thick
44	MJV211T	Tibial trial jig Size 1, Thick
44	MJV212T	Tibial trial jig Size 2, Thick
44	MJV213T	Tibial trial jig Size 3, Thick
45	MJV100T	Tibial trial rasp Size 0
45	MJV101T	Tibial trial rasp Size 1
45	MJV102T	Tibial trial rasp Size 2
45	MJV103T	Tibial trial rasp Size 3
-	CSA-1030-0010	Revision Instrument Tray 2

INTEGRAXTREVBOX2 Contents, continued

Bottom Tray Contents

Item	Reference	Instrument
46	MJU545T	Trial Insert, Size oo, Right, TH8
46	MJU546T	Trial Insert, Size oo, Right, TH9
46	MJU547T	Trial Insert, Size oo, Right, TH10
46	MJU548T	Trial Insert, Size oo, Right, TH11
46	MJV412T	Trial Insert, Size oo, Right, TH13
46	MJV413T	Trial Insert, Size oo, Right, TH15
46	MJV414T	Trial Insert, Size oo, Right, TH17
46	MJV415T	Trial Insert, Size oo, Right, TH19
46	MJV416T	Trial Insert, Size oo, Right, TH21
47	MJU555T	Trial Insert, Size oo, Left, TH8
47	MJU556T	Trial Insert, Size oo, Left, TH9
47	MJU557T	Trial Insert, Size oo, Left, TH10
47	MJU558T	Trial Insert, Size oo, Left, TH11
47	MJV422T	Trial Insert, Size oo, Left, TH13
47	MJV423T	Trial Insert, Size oo, Left, TH15
47	MJV424T	Trial Insert, Size oo, Left, TH17
47	MJV425T	Trial Insert, Size oo, Left, TH19
47	MJV426T	Trial Insert, Size oo, Left, TH21
46	MJV300T	Trial Insert, Size o1, Right, TH8
46	MJV301T	Trial Insert, Size o1, Right, TH9
46	MJV302T	Trial Insert, Size o1, Right, TH10
46	MJV303T	Trial Insert, Size o1, Right, TH11
46	MJV304T	Trial Insert, Size o1, Right, TH13
46	MJV305T	Trial Insert, Size o1, Right, TH15
46	MJV306T	Trial Insert, Size o1, Right, TH17
46	MJV307T	Trial Insert, Size o1, Right, TH19
46	MJV308T	Trial Insert, Size o1, Right, TH21
47	MJV310T	Trial Insert, Size o1, Left, TH8
46	MJV311T	Trial Insert, Size o1, Left, TH9
47	MJV312T	Trial Insert, Size o1, Left, TH10
47	MJV313T	Trial Insert, Size o1, Left, TH11
47	MJV314T	Trial Insert, Size o1, Left, TH13
47	MJV315T	Trial Insert, Size o1, Left, TH15
47	MJV316T	Trial Insert, Size o1, Left, TH17
47	MJV317T	Trial Insert, Size o1, Left, TH19
47	MJV318T	Trial Insert, Size o1, Left, TH21
46	MJU565T	Trial Insert, Size o, Right, TH8
46	MJU566T	Trial Insert, Size o, Right, TH9
46	MJU567T	Trial Insert, Size o, Right, TH10
46	MJU568T	Trial Insert, Size o, Right, TH11
46	MJV219T	Trial Insert, Size o, Right, TH13
46	MJV220T	Trial Insert, Size o, Right, TH15
46	MJV221T	Trial Insert, Size o, Right, TH17
46	MJV222T	Trial Insert, Size o, Right, TH19
46	MJV223T	Trial Insert, Size o, Right, TH21
47	MJU575T	Trial Insert, Size o, Left, TH8
47	MJU576T	Trial Insert, Size o, Left, TH9
47	MJU577T	Trial Insert, Size o, Left, TH10
47	MJU578T	Trial Insert, Size o, Left, TH11
47	MJV229T	Trial Insert, Size o, Left, TH13
47	MJV230T	Trial Insert, Size o, Left, TH15
47	MJV231T	Trial Insert, Size o, Left, TH17
47	MJV232T	Trial Insert, Size o, Left, TH19
47	MJV233T	Trial Insert, Size o, Left, TH21
46	MJU585T	Trial Insert, Size 1, Right, TH8
46	MJU586T	Trial Insert, Size 1, Right, TH9
46	MJU587T	Trial Insert, Size 1, Right, TH10
46	MJU588T	Trial Insert, Size 1, Right, TH11
46	MJU589T	Trial Insert, Size 1, Right, TH13
46	MJU590T	Trial Insert, Size 1, Right, TH15
46	MJU591T	Trial Insert, Size 1, Right, TH17
46	MJV242T	Trial Insert, Size 1, Right, TH19

Item	Reference	Instrument
46	MJV243T	Trial Insert, Size 1, Right, TH21
47	MJU595T	Trial Insert, Size 1, Left, TH8
47	MJU596T	Trial Insert, Size 1, Left, TH9
47	MJU597T	Trial Insert, Size 1, Left, TH10
47	MJU598T	Trial Insert, Size 1, Left, TH11
47	MJU599T	Trial Insert, Size 1, Left, TH13
47	MJU600T	Trial Insert, Size 1, Left, TH15
47	MJU601T	Trial Insert, Size 1, Left, TH17
47	MJV252T	Trial Insert, Size 1, Left, TH19
47	MJV253T	Trial Insert, Size 1, Left, TH21
46	MJU605T	Trial Insert, Size 2, Right, TH8
46	MJU606T	Trial Insert, Size 2, Right, TH9
46	MJU607T	Trial Insert, Size 2, Right, TH10
46	MJU608T	Trial Insert, Size 2, Right, TH11
46	MJU609T	Trial Insert, Size 2, Right, TH13
46	MJU610T	Trial Insert, Size 2, Right, TH15
46	MJU611T	Trial Insert, Size 2, Right, TH17
46	MJV262T	Trial Insert, Size 2, Right, TH19
46	MJV263T	Trial Insert, Size 2, Right, TH21
47	MJU615T	Trial Insert, Size 2, Left, TH8
47	MJU616T	Trial Insert, Size 2, Left, TH9
47	MJU617T	Trial Insert, Size 2, Left, TH10
47	MJU618T	Trial Insert, Size 2, Left, TH11
47	MJU619T	Trial Insert, Size 2, Left, TH13
47	MJU620T	Trial Insert, Size 2, Left, TH15
47	MJU621T	Trial Insert, Size 2, Left, TH17
47	MJV272T	Trial Insert, Size 2, Left, TH19
47	MJV273T	Trial Insert, Size 2, Left, TH21
46	MJU625T	Trial Insert, Size 3, Right, TH8
46	MJU626T	Trial Insert, Size 3, Right, TH9
46	MJU627T	Trial Insert, Size 3, Right, TH10
46	MJU628T	Trial Insert, Size 3, Right, TH11
46	MJU629T	Trial Insert, Size 3, Right, TH13
46	MJU630T	Trial Insert, Size 3, Right, TH15
46	MJU631T	Trial Insert, Size 3, Right, TH17
46	MJV282T	Trial Insert, Size 3, Right, TH19
46	MJV283T	Trial Insert, Size 3, Right, TH21
47	MJU635T	Trial Insert, Size 3, Left, TH8
47	MJU636T	Trial Insert, Size 3, Left, TH9
47	MJU637T	Trial Insert, Size 3, Left, TH10
47	MJU638T	Trial Insert, Size 3, Left, TH11
47	MJU639T	Trial Insert, Size 3, Left, TH13
47	MJU740T	Trial Insert, Size 3, Left, TH15
47	MJU741T	Trial Insert, Size 3, Left, TH17
47	MJV292T	Trial Insert, Size 3, Left, TH19
47	MJV293T	Trial Insert, Size 3, Left, TH21
48	MJU380T	Tibial Trial, Size o
48	MJU384T	Tibial Trial, Size 1
48	MJU385T	Tibial Trial, Size 2
48	MJU386T	Tibial Trial, Size 3
48	MJV500T	Tibial Trial, XL, Size o
48	MJV501T	Tibial Trial, XL, Size 1
48	MJV502T	Tibial Trial, XL, Size 2
48	MJV503T	Tibial Trial, XL, Size 3
49	MJU387T	Osteotome, Tibial Keel
50	MJU350T	Rasp, Tibial Flange
51	MJU353T	Drill Bit, 7.9 mm
52	MJU356T	Osteotome, Tibial Barrel, Revision
53	MJV545T	Insert Extractor
54	MJU361T	Tibial Impactor
55	MJV556T	Tibial Implant Extractor
-	CSA-1030-0008	Revision Instrument Case 2

Implants



TH4

TH10

Tibial Components

Tibial Baseplate	Size	Reference
TH4	Size 0	LJV160T
	Size 1	LJV161T
	Size 2	LJV162T
	Size 3	LJV163T
TH10	Size 0	LJV170T
	Size 1	LJV171T
	Size 2	LJV172T
	Size 3	LJV173T



Sloped

Talar Implants

Side	Type	Size	Reference
Right	Sloped	Size 0	LJV180T
		Size 1	LJV181T
		Size 2	LJV182T
		Size 3	LJV183T
Left	Sloped	Size 0	LJV190T
		Size 1	LJV191T
		Size 2	LJV192T
		Size 3	LJV193T



Flat Cut

Flat Cut Talar Components

Reference - Left	Reference - Right	Size
LJU820T	LJU810T	Size 0
LJU821T	LJU811T	Size 1
LJU822T	LJU812T	Size 2
LJU823T	LJU813T	Size 3

Note:

The Salto Talaris standard implants are ordered separately and implant descriptions are found in accompanying surgical techniques.



Size 00 Fixed Inserts

Insert	Reference
Insert, Size 00, Right, TH8	LJU408T
Insert, Size 00, Right, TH9	LJU409T
Insert, Size 00, Right, TH10	LJU410T
Insert, Size 00, Right, TH11	LJU411T
Insert, Size 00, Right, TH13	LJU412T
Insert, Size 00, Right, TH15	LJU413T
Insert, Size 00, Right, TH17	LJU414T
Insert, Size 00, Right, TH19	LJU415T
Insert, Size 00, Right, TH21	LJU416T
Insert, Size 00, Left, TH8	LJU418T
Insert, Size 00, Left, TH9	LJU419T
Insert, Size 00, Left, TH10	LJU420T
Insert, Size 00, Left, TH11	LJU421T
Insert, Size 00, Left, TH13	LJU422T
Insert, Size 00, Left, TH15	LJU423T
Insert, Size 00, Left, TH17	LJU424T
Insert, Size 00, Left, TH19	LJU425T
Insert, Size 00, Left, TH21	LJU426T

Size 0 Fixed Inserts

Insert	Reference
Insert, Size 0, Left, TH8	LJU225T
Insert, Size 0, Left, TH9	LJU226T
Insert, Size 0, Left, TH10	LJU227T
Insert, Size 0, Left, TH11	LJU228T
Insert, Size 0, Left, TH13	LJU229T
Insert, Size 0, Left, TH15	LJU230T
Insert, Size 0, Left, TH17	LJU231T
Insert, Size 0, Left, TH19	LJU232T
Insert, Size 0, Left, TH21	LJU233T
Insert, Size 0, Right, TH8	LJU215T
Insert, Size 0, Right, TH9	LJU216T
Insert, Size 0, Right, TH10	LJU217T
Insert, Size 0, Right, TH11	LJU218T
Insert, Size 0, Right, TH13	LJV219T
Insert, Size 0, Right, TH15	LJV220T
Insert, Size 0, Right, TH17	LJV221T
Insert, Size 0, Right, TH19	LJV222T
Insert, Size 0, Right, TH21	LJV223T

Size 01 Fixed Inserts

Insert	Reference
Insert, Size 01, Right, TH8	LJV300T
Insert, Size 01, Right, TH9	LJV301T
Insert, Size 01, Right, TH10	LJV302T
Insert, Size 01, Right, TH11	LJV303T
Insert, Size 01, Right, TH13	LJV304T
Insert, Size 01, Right, TH15	LJV305T
Insert, Size 01, Right, TH17	LJV306T
Insert, Size 01, Right, TH19	LJV307T
Insert, Size 01, Right, TH21	LJV308T
Insert, Size 01, Left, TH8	LJV310T
Insert, Size 01, Left, TH9	LJV311T
Insert, Size 01, Left, TH10	LJV312T
Insert, Size 01, Left, TH11	LJV313T
Insert, Size 01, Left, TH13	LJV314T
Insert, Size 01, Left, TH15	LJV315T
Insert, Size 01, Left, TH17	LJV316T
Insert, Size 01, Left, TH19	LJV317T
Insert, Size 01, Left, TH21	LJV318T

Size 1 Fixed Inserts

Insert	Reference
Insert, Size 1, Right, TH8	LJU235T
Insert, Size 1, Right, TH9	LJU236T
Insert, Size 1, Right, TH10	LJU237T
Insert, Size 1, Right, TH11	LJU238T
Insert, Size 1, Right, TH13	LJU239T
Insert, Size 1, Right, TH15	LJU240T
Insert, Size 1, Right, TH17	LJU241T
Insert, Size 1, Right, TH19	LJU242T
Insert, Size 1, Right, TH21	LJU243T
Insert, Size 1, Left, TH8	LJU245T
Insert, Size 1, Left, TH9	LJU246T
Insert, Size 1, Left, TH10	LJU247T
Insert, Size 1, Left, TH11	LJU248T
Insert, Size 1, Left, TH13	LJU249T
Insert, Size 1, Left, TH15	LJU250T
Insert, Size 1, Left, TH17	LJU251T
Insert, Size 1, Left, TH19	LJU252T
Insert, Size 1, Left, TH21	LJU253T









Size 2 Fixed Inserts



Insert	Reference
Insert, Size 2, Right, TH8	LJU255T
Insert, Size 2, Right, TH9	LJU256T
Insert, Size 2, Right, TH10	LJU257T
Insert, Size 2, Right, TH11	LJU258T
Insert, Size 2, Right, TH13	LJU259T
Insert, Size 2, Right, TH15	LJU260T
Insert, Size 2, Right, TH17	LJU261T
Insert, Size 2, Right, TH19	LJU262T
Insert, Size 2, Right, TH21	LJU263T
Insert, Size 2, Left, TH8	LJU265T
Insert, Size 2, Left, TH9	LJU266T
Insert, Size 2, Left, TH10	LJU267T
Insert, Size 2, Left, TH11	LJU268T
Insert, Size 2, Left, TH13	LJU269T
Insert, Size 2, Left, TH15	LJU270T
Insert, Size 2, Left, TH17	LJU271T
Insert, Size 2, Left, TH19	LJU272T
Insert, Size 2, Left, TH21	LJU273T

Size 3 Fixed Inserts

Insert	Reference
Insert, Size 3, Right, TH8	LJU275T
Insert, Size 3, Right, TH9	LJU276T
Insert, Size 3, Right, TH10	LJU277T
Insert, Size 3, Right, TH11	LJU278T
Insert, Size 3, Right, TH13	LJU279T
Insert, Size 3, Right, TH15	LJU280T
Insert, Size 3, Right, TH17	LJU281T
Insert, Size 3, Right, TH19	LJU282T
Insert, Size 3, Right, TH21	LJU283T
Insert, Size 3, Left, TH8	LJU285T
Insert, Size 3, Left, TH9	LJU286T
Insert, Size 3, Left, TH10	LJU287T
Insert, Size 3, Left, TH11	LJU288T
Insert, Size 3, Left, TH13	LJU289T
Insert, Size 3, Left, TH15	LJU290T
Insert, Size 3, Left, TH17	LJU291T
Insert, Size 3, Left, TH19	LJU292T
Insert, Size 3, Left, TH21	LJU293T

Disposables - Sterile Single Use Only

Saw Blades								
Integra No.	SAW5944T/ SAW6944T	SAW5945T/ SAW6945T	SAW5946T/ SAW6946T	SAW5947T/ SAW6947T	SAW5948T/ SAW6948T	SAW5949T/ SAW6949T	SAW5950T/ SAW6950T	SAW6951T
Depth	70.0 mm	85.0 mm	80.0 mm	90.0 mm	75.0 mm	90.0 mm	70.0 mm	70.0 mm
Width	13.0 mm	21.0 mm	13.0 mm	21.0 mm	13.0 mm	21.0 mm	12.5 mm	12.5 mm
Mtl Thk	1.24 mm	1.24 mm	1.24 mm	1.24 mm	1.24 mm	1.24 mm		
Cut Thk	1.27 mm	1.27 mm	1.27 mm	1.27 mm	1.27 mm	1.27 mm	0.94 mm	0.94 mm
Hub	Stryker System 7	Stryker System 7	Hall Versipower	Hall Versipower	Hall Power Pro	Hall Power Pro	Brasseler	Brasseler
								

Pin Pack: LJV529T Includes:		
Description	Reference	
3 x 110 mm Self-drilling Pins	LJV527T	
5 x 75 mm Self-drilling Pins	LJV526T	
3 x 45 mm Self-drilling Pins	LJV525T	
Reamer	LJU097T	

Drill Bit	
Description	Reference
3mm x 220mm	DWD060T
3mm x 135mm	LJV528T

Optional Instruments

Description	Reference
Wiener Antrum Rasp	400236
Tibial Impactor	1129

*Optional instruments are not included in the standard instrument set - special order only.

Page Intentionally Left Blank

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com

Manufacturer:



Ascension Orthopedics Inc
11101 Metric Blvd
Austin, TX 78758 ▪ USA