

Lower risk of prosthetic joint infection (PJI) and revision for infection with OXINIUM[®]/XLPE primary THA versus other bearing combinations^{1,2}

Why is PJI a growing concern?

No. 1
indication
for revision
arthroplasty³



400% projected
increase in demand for
arthroplasty from the
early 2000s to 2030⁴

3.6X higher risk of death
for THA patients after PJI
versus average mortality
risk at 65 years old⁵

\$1.85 billion
projected annual US hospital
expense for PJI by 2030
(THA: \$753.4 million, TKA: \$1.1 billion)⁶

Evidence from two real-world analyses of primary THA patients has shown that:^{1,2}

	Odds of PJI-related infection/ inflammation were: ^{1*}	Risk of revision for infection was: ²
Versus metal/ XLPE	27% lower with OXINIUM/XLPE Medical record-diagnosed PJI-related infection or inflammation validated with procedure code hip imaging: OXINIUM/XLPE 1.29% vs metal/XLPE 1.77% (odds ratio [OR] 0.73; p=0.025)	24% lower with OXINIUM/XLPE (adjusted hazard ratio [HR]: 0.76; p=0.008)
Versus ceramic/ XLPE	34% lower with OXINIUM/XLPE Medical record-diagnosed PJI-related infection or inflammation validated with procedure code hip imaging: OXINIUM/XLPE 1.33% vs ceramic/XLPE 1.99% (OR 0.66; p=0.0027)	20% lower with OXINIUM/XLPE (adjusted HR: 0.80; p=0.016)
Analysis methodology	Included >123,000 patients Propensity score matched primary THA patients with PJI diagnosis from the PREMIER PINC AI Healthcare Database (2017–2022)	Included >784,000 patients Bespoke implant report produced by the NJR (data up to Jan 2026) [†] for primary THA patients with up to 20 years follow-up

Proprietary OXINIUM (OXidized zircoNIUM) Technology:

Contains virtually no
cobalt, chromium
or nickel^{7,8}

Has been shown in laboratory testing to have:
Improved wear resistance^{9–14}
and biocompatibility¹⁵ versus CoCr



LEARN MORE

Has shown the **lowest revision risk of**
all modern bearing combinations
across registries in primary THA (with XLPE)^{16–19§}



LEARN MORE

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area. For detailed product information, including indications for use, contraindications, warnings and precautions, please consult the product's Instructions for Use (IFU) prior to use.

Abbreviations: CoCr = cobalt chrome; HR = hazard ratio; OR = odds ratio; PJI = prosthetic joint infection; THA = total hip arthroplasty; TKA = total knee arthroplasty; XLPE = cross-linked polyethylene.

*As the AJRR reported 96% of THAs used XLPE during the analysis period (2017–2022), it was assumed that all procedures captured used XLPE.¹ †Before validation via cross-referencing with imaging procedure codes (restricting to medical record-diagnosed infection/inflammation), PJI-related infection/inflammation rates for OXINIUM/XLPE vs ceramic/XLPE were significantly lower at 1, 3, and 6 months, but not at 12 or 24 months; versus metal/XLPE, differences were significantly lower at all time points (at 2 years: 7.37% vs 8.37%; OR 0.87 [95% CI: 0.77–0.99], p=0.032). ‡The data used for this analysis was obtained from the NJR Supplier Feedback System. The Healthcare Quality Improvement Partnership (HQIP) and/or the NJR take no responsibility for the accuracy, currency, reliability and correctness of any data used or referred to in this report, nor for the accuracy, currency, reliability and correctness of links or references to other information sources and disclaims all warranties in relation to such data, links and references to the maximum extent permitted by legislation. §See linked registry analysis (QR code) for other modern bearing combinations compared to metal/XLPE. Ceramicised metal/XLPE is equivalent to 'OXINIUM/XLPE'. As in previous years, the Registry urges caution in the interpretation of this result. This bearing is a single company product, used with a small number of femoral stem and acetabular component combinations. See linked registry analysis (QR code) for further details.

References: 1. Schaffler BC, et al. *Hip Int*. Published online February 4, 2026. doi: 10.1177/11207000251410114. 2. National Joint Registry. Bespoke Implant Report for OXINIUM femoral head with XLPE liners. 10 March 2026. Copy available on request. 3. Ayoade F, et al. *Periprosthetic Joint Infection*. Treasure Island (FL): StatPearls Publishing; 2024. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK448131/>. Accessed April 1, 2026. 4. Ahmed SS, Haddad FS, et al. *Bone Joint Res*. 2019;8(11):570–572. 5. Natsuhara KM, et al. *J Arthroplasty*. 2019;34(7S):S337–S342. 6. Premkumar A, et al. *J Arthroplasty*. 2021;36(5):1484–1489e3. 7. ASTM International. ASTM F2384–24 (May 2024). Standard specification for wrought zirconium-2.5niobium alloy for surgical implant applications (UNS R60901). Available at: <https://store.astm.org/f2384-24.html>. Accessed April 2, 2026. 8. Hunter G, et al. *J ASTM Int*. 2005;2(7):1–14. 9. Parikh A, et al. Poster 1028 presented at: Orthopaedic Research Society Annual Meeting; January 26–29, 2013; San Antonio, TX, USA. 10. Long M, et al. In: Transactions: twenty-fourth annual meeting of the Society for Biomaterials. 21:528. 11. Sheth NP, et al. *J Surg Orthop Adv*. 2008;17(1):17–26. 12. Heyse TJ, et al. *Clin Orthop Relat Res*. 2014;472(1):277–283. 13. Salehi A, et al. *Key Eng Mater*. 2006;309–311:1199–1202. 14. Poggie RA, et al. In: Denton R, Keshavan, eds. *Wear and Friction of Elastomers*, ASTM STP1145. ASTM International; 1992:65–81. 15. Dalal A, et al. *J Biomed Mater Res Part A*. 2012;100A:2147–2158. 16. Lewis PL, et al. Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR) Hip, Knee & Shoulder Arthroplasty: 2025 Annual Report Adelaide; AOA, 2025. Available at: <https://aoanjrr.sahmri.com/annual-reports-2024>. Accessed April 2, 2026. 17. Whitehouse MR, et al. *PLoS Med*. 2024;21(11):e1004478. 18. Atrey A, et al. Poster HS2 presented at: Canadian Orthopaedic Association; June 20–23, 2018; Victoria, BC, Canada. 19. Peters RM, et al. *Acta*. 2018;89(2):163–169.