

Smith+Nephew

Back to Basics


Mechanisms of actions



PICO 


Single Use Negative Pressure
Wound Therapy System





Mechanisms of action

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.




*As demonstrated *in vivo*.

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

Summary of the MoAs




Increases the efficiency of functional lymph vessels which helps to **reduce oedema**⁵⁻⁷




Helps to promote **lymphatic drainage**^{4*}




Maintains an **efficient blood supply** to the wound (perfusion), which helps to support the immune response⁸⁻¹⁰



Holds closed incision together, **reducing lateral tension** forces across the incision^{2†}



Provides a **moist wound environment** which helps with wound healing²⁹⁻³²



The combination of macro and micro-deformations (wound contraction and filling of tissue defects with new granulation tissue) leads to **reduction in wound area and wound depth**.¹¹⁻¹³

*As demonstrated *in vivo*. †As demonstrated in biomechanical modelling.

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

Mechanism of action of NPWT



Smith+Nephew Clinical Scientific Medical Affair (CSMA) experts highlight the following key wound effects between **surgically closed incisions and open wounds**.

Closed surgical incisions Bolstering incisions and reduce infection risks	GO TO LINK
Summary	
Reduce lateral tension	
Fluid removal	
Blood flow	
Maintenance of wound homeostasis	


Open wounds Improve granulation tissue and exudate management	GO TO LINK
Summary	
Fluid removal	
Blood flow	
Micro-deformation	
Macro-deformation	
Maintenance of wound homeostasis	

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.



MoA for closed incision wounds

Lateral tension and wound strength



By compressing the wound, PICO[◇] sNPWT helps to reduce lateral tension placed on the closed incision,^{2*} pulling the borders together, facilitating wound healing and increasing breaking strength^{2,8,14}

- Sutures hold wound borders together at the skin surface and in the subcutis; they help reduce the lateral tension that can cause wound borders to move apart¹⁴
- Where wound borders are not closed granulation tissue can form, which increases the risk of undesirable scarring¹⁴
- Application of negative pressure across the incision, to a wider zone of injury, reduces lateral tension on closed incision, pulling the wound borders together^{2,3,8}
- This may result in improvements of quality and appearance of scars.^{15†}

*As demonstrated in biomechanical modelling. †When applied post-operatively compared to care with standard dressings; $p < 0.001$.

Note: PICO is not a scar treatment but rather helps prevent scar formation/ aid in scar appearance, by aiding with the prevention of surgical site complication/ wound healing.



Video

[Back to MoA Summary page](#)


For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

MoA for closed incision wounds



Reduction in lateral tension helps improve the quality and appearance of scars.¹⁵

Visual Analog Scale (VAS) score



Patient Observer Scar Assessment Scale (POSAS)




Fig.1: Within patient difference in VAS score (PICO – SC)¹⁵

Fig.2: Within patient difference in POSAS score (PICO – SC)¹⁵



Video


[Back to MoA Summary page](#)

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

Note: PICO is not a scar treatment but rather helps prevent scar formation/ aid in scar appearance, by aiding with the prevention of surgical site complication/ wound healing.

- **200 patients** undergoing bilateral reduction mammoplasty treated with both **PICO[◇] sNPWT** and standard dressing¹⁵
- **Each patient was their own control;** one breast was treated with **PICO sNPWT** and one breast treated with standard dressing
- Patients treated for up to **14 days** to enable a within patient comparison
- Differences in scar quality and aesthetic appearance assessed using the patient and observer scar assessment scale (POSAS) and The visual analogue scale at **42 and 90 days**¹⁵
- Scar quality as measured by both scoring systems was shown to be **significantly better** on **PICO[◇] sNPWT** than standard dressing, at both **42 and 90 day** assessment ($p < 0.001$)¹⁵

Lateral tension and wound strength




MoA for closed incision wounds

Fluid management | Improve fluid management via lymphatic system

Standard dressings

- Inefficient functioning of the lymphatic system can cause build up of interstitial fluid and potential toxins in the tissues (oedema)⁴⁻⁵
- Consequently, accumulation of serious fluid (seroma) and blood (hematoma) can occur within tissue cavities^{4,6}
- Standard dressings have no known effects on the lymphatic system



PICO[◇] sNPWT

PICO Dressings with AIRLOCK[◇] Technology compress the peri-wound area and remove fluid from wound bed, this simultaneous mechanism may enable fewer dressing changes during the healing process.^{16-18*}

NPWT helps to increase the activity of the lymphatic system,^{4†} in part by reducing oedema and compression of the vessels.⁴⁻⁷

This helps to reduce the incidence of seroma in closed surgical incisions^{19-24‡}



Video

[Back to MoA Summary page](#)

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.


*Compared to care with standard dressings. †As demonstrated *in vivo*.
‡Compared to care with standard dressings.

MoA for closed incision wounds

Fluid management | Effects on oedema and swelling

Standard dressings

- If the lymphatic system is compromised and/or overloaded, fluid can accumulate between the cells (interstitial fluid) causing oedema and swelling²⁵⁻²⁶
- Uncontrolled oedema may delay or compromise wound healing and increase the risk of infection²⁷
- Standard dressings have no known beneficial effects of oedema and swelling



PICO sNPWT

NPWT may help to reduce post-operative oedema⁵⁻⁷

Proposed mechanisms for this possible effect of NPWT may include improvements in the functionality of damaged lymph vessels and increased efficiency of functional lymph vessels^{4,6,14}



Video


[Back to MoA Summary page](#)

MoA for closed incision wounds

Fluid Management | Effects on blood vessels, perfusion and angiogenesis

Standard dressings

- Uncontrolled oedema can cause compression of blood vessels²⁷
- This may lead to decreased blood circulation and impair wound healing²⁷
- Standard dressings have no known beneficial effects on oedema and blood circulation



PICO[◇] sNPWT

May help to reduce build up of interstitial fluid (oedema) and minimise compression of blood vessels²⁵

In vivo, promotes changes in blood flow (perfusion), which may be beneficial to the wound^{8*}


May help to form new blood vessels (angiogenesis) around the incisional margin²⁸

*As demonstrated *in vivo*.




Video

[Back to MoA Summary page](#)



MoA for closed incision wounds

Maintenance of wound homeostasis




Provides a closed, moist wound environment²⁹⁻³²



- Controlled moisture evaporation⁸ (*in-vitro*) and helps to minimise heat loss¹
- Helps to minimise contamination⁵



[Back to MoA Summary page](#)



Summary: NPWT has multiple mechanisms of action which may help promote incisional wound healing and reduce the odds of SSCs.^{1,7,30,31,33}



This pathway is adapted from the WUWHS guidelines document and it shows how NPWT can help reduce SSCs and lateral tension while increasing lymphatic drainage. This effect is likely to contribute to faster and stronger healing, and a reduced risk of infection and dehiscence³⁵

*As demonstrated in biomechanical modelling †As demonstrated *in vivo* [Back to MoA Summary page](#)


For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

MoA for open wounds

Fluid Management | Improve fluid management via lymphatic system

Standard dressings

- Inefficient functioning of the lymphatic system can cause build up of interstitial fluid and potential toxins in the tissues (oedema)^{4,6}
- Uncontrolled oedema may delay or compromise wound healing²⁷
- Standard dressings have no known effects on the lymphatic system



PICO[◇] sNPWT

PICO Dressings with AIRLOCK[◇] Technology compress the peri-wound area and remove fluid from wound bed, this simultaneous mechanism may enable fewer dressing changes during the healing process.^{16-18*}

NPWT helps to increase the activity of the lymphatic system,^{4†} in part by reducing oedema and compression of the vessels.⁴⁻⁷



Video

[Back to MoA Summary page](#)


*Compared to care with standard dressings. †As demonstrated *in vivo*.

MoA for open wounds

Fluid Management | Effects on oedema and swelling

Standard dressings

- If the lymphatic system is compromised and/or overloaded, fluid can accumulate between the cells (interstitial fluid) causing oedema and swelling²⁵⁻²⁶
- Uncontrolled oedema may delay or compromise wound healing and increase the risk of infection²⁷
- Standard dressings have no known beneficial effects of oedema and swelling



PICO sNPWT

NPWT helps been shown to improve the functionality of damaged lymph vessels and increase the efficiency of functional lymph vessels, helping to reduce oedema^{4,6,14}



Video

[Back to MoA Summary page](#)


MoA for open wounds

Blood flow | Increased blood flow helps in the process of wound healing and granulation tissue formation

Standard dressings

- Uncontrolled oedema can cause compression of blood vessels²⁷
- This may lead to decreased blood circulation and impair wound healing²⁷
- Standard dressing have no known beneficial effects on oedema and blood circulation

PICO[◇] sNPWT



Helps to reduce build up of interstitial fluid (oedema) and minimise compression of blood vessels²⁵


Maintains an efficient blood supply to the wound (perfusion), which helps to support the immune response and remove debris^{6,8}

Can promote the formation of new blood vessels (angiogenesis)⁹⁻¹⁰ which facilitates tissue growth and wound healing



Video

[Back to MoA Summary page](#)




Increased blood flow helps in the process of wound healing and granulation tissue formation.

NPWT directly influences blood flow (*in-vivo*)

- Stimulation of blood flow experimental wounds³⁶
- Tissue stretching shown to increase blood flow¹³

NPWT reduces oedema³⁷

- Reduces gap between capillaries and cells
- Improves blood flow to cells




MoA for open wounds

Macro-deformation

The stretch or contraction of tissue when NPWT is applied.

- When NPWT is applied to a wound an immediate reduction in wound area or 'wound contraction' occurs¹² (*in-vivo*)
- This reduction in wound area appears to be permanent as the wounds remained contracted after discontinuation^{11,12}



Video

[Back to MoA Summary page](#)


MoA for open wounds

Micro-deformation

The transfer of an imprint of the surface topography of the compressed wound filler to the tissue surface


- When negative pressure is applied, tissue bulges into the open pores of the wound filler¹¹ (*in-vitro*)
- Cells are stimulated to proliferate leading to new granulation tissue formation^{11,29,38,39} (*in-vitro*)

Note: There is no micro-deformation in closed incision wounds. Micro-deformation happens in open wounds only when filler is being used.




Video

[Back to MoA Summary page](#)



MoA for open wounds

Maintenance of wound homeostasis



NPWT with a wound filler

NPWT without a wound filler

Provides a closed, moist wound environment²⁹⁻³²


- Controlled moisture evaporation⁸ (*in-vitro*) and helps to minimise heat loss¹
- Helps to minimise contamination⁵





Video

[Back to MoA Summary page](#)

Maintenance of wound homeostasis



Summary: NPWT has multiple mechanisms of action that help improve the speed, strength and quality of wound healing,* lead to fewer dressing changes^{17,40,41} and faster healing of chronic wounds^{41†,43}



Schematic representation of the multiple proposed mechanisms of action of sNPWT on acute and chronic wounds

[Back to MoA Summary page](#)

*Compared to care with standard dressing.
†n=9; 44 weeks mean wound duration prior to study
#in partial thickness burns

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.



**Not all NPWT
systems are
the same**

PICO[◇] sNPWT
Product features,
benefits and applications

The PICO[◇] Single-use Negative Pressure Wound Therapy



Features:

Improved device performance

- The PICO 7 pump has a significantly* higher maximum leak rate tolerance than the original PICO pump⁴⁶

Cost savings

- The dressing full indicator on the PICO System led to a low frequency of outpatient clinic visits and home visits, creating a potential saving in treatment costs^{47,48†}

Designed to improve patient quality of life

- The PICO pump is quiet (<30db at 1m)⁴⁹
- The belt clip is to enable portability worn by the patient during use

Increased flexibility

- Multipacks of five dressings available, allowing therapy to be tailored to patients' clinical needs



The PICO System is suitable for use in both hospital and homecare settings.⁵⁰

*p<0.001. †As demonstrated in a 4-patient case study on complex DFUs. Reduction per patient of 1-2 outpatient visits and 1-3 home visits per week; compared to pre-NPWT.

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

Only PICO[◇] sNPWT dressings have AIRLOCK[◇] Technology




Wide delivery. Constant pressure. **Optimal outcomes.**

Top film layer has a high moisture vapour transmission rate and protects the wounds from external contamination^{30,45}

Super-absorbent core locks exudate away from wound^{52*}

Silicone adhesive layer protects the wound environment and helps to minimise pain on removal^{31,33,51}



The PICO System with **AIRLOCK** Technology allows the delivery of negative pressure across the entire dressing to ensure that treatment is delivered to a wider zone beyond the wound itself.^{53,54}


**In vitro testing*

PICO[◇] sNPWT with the proprietary AIRLOCK[◇] Technology




Deliver negative pressure across the entire dressing area,⁵⁴⁻⁵⁶ providing benefits to the wound and surrounding area^{37,57}

PICO[◇] sNPWT



tNPWT




CT imaging of tissue compression


- Delivery of negative pressure across the dressing to surrounding tissue,⁵⁴⁻⁵⁶ which may benefit perfusion⁵⁷ and management of oedema³⁷
- Consistent delivery of negative pressure for the duration of therapy^{31,53}
- **PICO** sNPWT showed less wound edge inflammation, yielding a pro-healing environment^{38,58}
- Fewer changes to skin health and barrier function than with tNPWT³⁸

PICO[◇] Dressing: A pioneering dressing

Closed surgical incisions



Open wounds



Unique design to facilitate delivery of negative pressure at the wound bed.^{8,59}

- Exudate and bacteria are retained in the dressing and away from the wound bed^{52,60}
- Dressing can be disconnected from the pump to enable showering
- Dressing is waterproof⁶¹

Using PICO[◇] sNPWT with wound fillers



Choosing a wound filler

Distribution of negative pressure requires intimate contact with the entire wound surface.


If the wound depth is under 0.5cm the **PICO** Dressing can be applied directly to the wound without the need for a filler

OUS	US
Wound >0.5cm in depth are likely to require a foam or gauze	
Wound >2.0cm in depth must be treated with the use of filler	
	Recommend to treat wound up-to 4.5cm in depth

You may still consider the use of a filler for wounds especially if:

- 1. The wound bed is uneven
- 2. There is undermining
- 3. There is a sinus

Wound types



Did you know?





In a variety of chronic wounds deeper than 2cm, use of PICO sNPWT as part of the integrated care bundle:

- 1. Shown to significantly reduce mean time to wound healing by 11.75 weeks*⁶²
- 2. Shown to have 3.23 days longer mean dressing wear time*⁶²

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

Using PICO[◇] sNPWT with wound fillers



Types of NPWT filler	Qualities ⁶³	Wound type ⁶³
<div>Foam</div> <div></div> <div></div> <div><p>Polyurethane foam saturated with water</p></div> <div><p>Open cell structure and hydrophobic material properties allow for efficient fluid management.</p></div>	<ul style="list-style-type: none">▪ Stimulates rapid granulation tissue formation, wound contraction▪ Removal of high volumes of exudate	<ul style="list-style-type: none">▪ Compartment syndrome▪ Acute wounds with large tissue loss▪ Post operative open abdomens and sternotomy wounds
<div>Gauze</div> <div></div> <div></div> <div><p>Gauze saturated with water</p></div> <div><p>Hydrophilic nature of fibres absorbs and retains fluid. Not optimal for management of viscous fluid.</p></div>	<ul style="list-style-type: none">▪ Easy and fast to apply on uneven, large, deep wounds▪ Reduced pain and trauma reported on removal▪ Slower initial stimulation of granulation tissue▪ Reduced post therapy contraction may reduce scarring	<ul style="list-style-type: none">▪ Cosmetic surgery▪ Skin flap preparation/stabilisation▪ Skin grafts and use over joints


Did you know?

The PICO System demonstrated superior healing outcomes, where only 13.5% of application with PICO System required a filler.*^{29,64}

*in a study of VLUs and DFUs.; Vs. tNPWT; At 12 weeks; n-161; p≤0.021

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

PICO[◇] sNPWT is compatible with foam & gauze



Clinical studies demonstrate the overall rates of healing with foam and gauze are similar⁶⁵⁻⁶⁶ however the characteristics of granulation tissue can differ

Foam: thick granulation tissue – suitable for use on wounds where scarring does not pose a problem e.g. sternotomy wounds⁶⁷

Gauze: less thick but dense granulation tissue – suitable for use where cosmetic results are of greater importance or where scar tissue may restrict movement⁶⁷

References



1. World Union of Wound Healing Societies (WUWHS) Consensus Document. Closed surgical incision management: understanding the role of NPWT. Wounds International, 2016
2. Loveluck J, Copeland T, Hill J, Hunt A, Martin R. Biomechanical modelling of the forces applied to closed incisions during single-use negative pressure wound therapy. ePlasty. 2016;16:e20.
3. Glaser DA, Farnsworth CL, Varley ES, Nunn TA, Sayad-Shah M, Breisch EA, et al. Negative pressure therapy for closed spine incisions: a pilot study. Wounds. 2012;24(11):308–16.
4. Kilpadi DV, Cunningham MR. Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system. Wound Repair Regen. 2011;19:588–596.
5. Birke-Sorensen H, et al. Evidence-based recommendations for Negative Pressure Wound Therapy: Treatment variables (pressure levels, wound filler and contact layer) - Steps towards an international consensus. Journal of plastic, reconstructive and aesthetic surgery: JPRAS 64 Suppl, S1–S16 (2011)
6. Scalise A, Calamita R, Tartaglione C, et al. Improving wound healing and preventing surgical site complications of closed surgical incisions: a possible role of Incisional Negative Pressure Wound Therapy. A systematic review of the literature. Int Wound J. 2016;13:1260–1281.
7. Shim HS, Choi JS, Kim SW. A role for postoperative negative pressure wound therapy in multitissue hand injuries. Biomed Res Int. 2018;2018
8. Malmström M, Huddleston E, Martin R. Biological effects of a disposable, canisterless negative pressure wound therapy system. ePlasty. 2014;14:e15.
9. Ma Z, Shou K, Li Z, et al. Negative pressure wound therapy promotes vessel destabilization and maturation at various stages of wound healing and thus influences wound prognosis. Exp Ther Med. 2016;11(4):1307–1317
10. Xia CY, Yu AX, Qi B, et al. Analysis of blood flow and local expression of angiogenesis associated growth factors in infected wounds treated with negative pressure wound therapy. Mol Med Rep. 2014;9(5):1749–1754
11. Saxena V, et al. Vacuum-assisted closure: microdeformations of wounds and cell proliferation. Plast Reconstr Surg. 114:1086–96.
12. Borquist O, et al. Micro and macromechanical effects on the wound bed of negative pressure wound therapy using foam and gauze. Annals of Plastic Surgery. 64: 789–93.
13. Ichioka S, et al. A technique to visualise wound bed microcirculation and the acute effect of negative pressure. Wound Rep Regen. 16(3); 460–465.
14. Karlakki S, Brem M, Giannini S, Khanduja V, Stannard J, Martin R. Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication. Bone Joint Res. 2013;2(12):276–84.
15. Galiano R, Djohan R, Shin J, Hudson D, et al. The effects of a single use canister-free Negative Pressure Wound Therapy (NPWT) System* on the prevention of postsurgical wound complications in patients undergoing bilateral breast reduction surgery. Paper presented at: The British Association of Aesthetic Plastic Surgeons (BAAP's) 30th Annual Scientific Meeting; 2014; London
16. Karlakki SL, Hamad AK, Whittall C, et al. Incisional negative pressure wound therapy dressings (iNPWTd) in routine primary hip and knee arthroplasties: A randomised controlled trial. Bone Joint Res. 2016;5(8):328–337.
17. Nherera LM, Trueman P, Karlakki SL. Cost-effectiveness analysis of single-use negative pressure wound therapy dressings (sNPWT) to reduce surgical site complications (SSC) in routine primary hip and knee replacements. Wound Repair Regen. 2017;25(3):474–482.
18. Nordmeyer M, Pauser J, Biber R, et al. Negative pressure wound therapy for seroma prevention and surgical incision treatment in spinal fracture care. Int Wound J. 2016;13(6):1176–1179.
19. Fleming CA, Kuteva M, O'Hanlon K, O'Brien G, McGreal G. Routine use of PICO dressings may reduce overall groin wound complication rates following peripheral vascular surgery. J Hosp Infect. 2018;99(1):75–80.
20. Gillespie BM, Rickard CM, Thalib L, et al. Use of Negative-Pressure Wound Dressings to Prevent Surgical Site Complications After Primary Hip Arthroplasty: A Pilot RCT. Surg Innov. 2015;22(5):488–495.
21. Galiano RD, Hudson D, Shin J, et al. Incisional Negative Pressure Wound Therapy for Prevention of Wound Healing Complications Following Reduction Mammoplasty. Plast Reconstr Surg Glob Open. 2018;6(1):e1560.
22. Pellino G, Sciaudone G, Candilio G, et al. Effects of a new pocket device for negative pressure wound therapy on surgical wounds of patients affected with Crohn's disease: a pilot trial. Surg Innov. 2014;21(2):204–212.
23. Pellino G, Sciaudone G, Candilio G, et al. Preventive NPWT over closed incisions in general surgery: does age matter? Int J Surg. 2014;12 Suppl 2:S64–S68.
24. Selvaggi F, Pellino G, Sciaudone G, et al. New Advances in Negative Pressure Wound Therapy (NPWT) for Surgical Wounds of Patients Affected with Chron's Disease. Surg Technol Int. 2014;24:83–89.
25. Lalezari S, Lee CJ, Borovikova AA, et al. Deconstructing negative pressure wound therapy. Int Wound J. 2017;14:649–657.
26. Traves KP, Studdiford JS, Pickle S, Tully AS. Edema: diagnosis and management. Am Fam Physician. 2013;88:102–110.
27. Hess CT. Checklist for factors affecting wound healing. Adv Skin Wound Care. 2011 Apr;24(4):192. doi: 10.1097/01.ASW.0000396300.04173.ec. PMID: 21422844.
28. Suh H, Lee AY, Park EJ, Hong JP. Negative pressure wound therapy on closed surgical wounds with dead space: animal study using a swine model. Ann Plast Surg. 2016;76:717–722.
29. Kirsner R, Dove C, Reyzelman A, Vayser D, Jaimes H. A Prospective, Randomised, Controlled Clinical Trial on the Efficacy of a single-use Negative Pressure Wound Therapy System, compared to Traditional Negative Pressure Wound Therapy in the Treatment of Chronic Ulcers of the Lower Extremities. Wound Repair Regen. 2019;27(5):519–529.
30. Smith+Nephew 2018. Summary of routine QA testing on MVP of PICO dressings. 2018. Internal Report. DS/18/153/R.

References



31. Hudson DA, Adams KG, Van Huyssteen A, Martin R, Huddleston EM. Simplified negative pressure wound therapy: clinical evaluation of an ultraportable, no-canister system. *Int Wound J*. 2015;12(2):195–201
32. Smith+Nephew 2019. Use of Moisture Vapour Permeability* (MVP) and Moisture Vapour Transmission rate** (MVTR) data to support product claims referring to moist wound healing. Internal Report. EO.AWM.PCSgen.001.v2.
33. Payne C, Edwards D. 2014. Application of the single use negative pressure wound therapy device (PICO) on a heterogeneous group of surgical 46-47, and traumatic wounds. ePlasty.
34. Smith+Nephew 2020. Bacterial barrier testing of the PICO \diamond Dressing. Internal Report. 2001002.
35. Rossington A. 2015. A prospective, open, non-comparative, multicentre study to evaluate the functionality and dressing performance of a new negative pressure enhanced dressing (NPED) in acute wounds. CT09/02, May
36. Morykwas MJ, Argenta LC, Shelton-Brown EI, McGuirt W. Vacuum-assisted closure: a new method for wound control and treatment: animal studies and basic foundation. *Ann Plast Surg*. 1997 Jun;38(6):553–62.
37. Young SR, Hampton S, Martin R, et al. Non-invasive assessment of NPWT using high frequency diagnostic ultrasound: oedema reduction and new tissue accumulation. *Int Wound J*. 2013;10(4):383–388.
38. Brownhill VR, Huddleston E, Bell A, et al. Pre-Clinical Assessment of Single-Use Negative Pressure Wound Therapy During In Vivo Porcine Wound Healing. *Adv Wound Care (New Rochelle)*. 2020;0(0):1–12.
39. Schwartz JA, Goss SG, Facchin F, Avdagic E, Lantis JC. Single use-negative pressure wound therapy for the treatment of chronic lower leg wounds. *J Wound Care*. 2015;24(2).
40. Dowsett C, Hampton K, Myers D, Styche T. Use of PICO™ to improve clinical and economic outcomes in hard-to-heal wounds. *Wounds International*. 2017;8(2):52–58
41. Hampton J. Providing cost-effective treatment of hard-to-heal wounds in the community through use of NPWT. *Community Wound Care*. 2015:S14–S20.
42. Kamolz LP, Andel H, Haslik W, et al. Use of subatmospheric pressure therapy to prevent burn wound progression in human: first experiences. *Burns*. 2004;30(3):253–258.
43. Karlakki SL, Hamad AK, Whittall C, Graham NM, Banerjee RD, Kuiper JH. 2016. Incisional negative pressure wound therapy dressings (iNPWTd) in routine primary hip and knee arthroplasties: A randomised controlled trial. *Bone & joint research*, 5(8), pp.328–337.
44. Canonico S, Campitiello F, Della Corte A, et al. Therapeutic possibilities of portable NPWT. Initial multidisciplinary observation with the negative pressure therapy device. *Acta Vulnologica*. 2012;10:57–66.
45. Smith+Nephew. Bacterial barrier testing (wet-wet) of PICO Dressing with a 7 day test duration against *S. marcescens*. Internal Report 1102010.
46. Smith+Nephew January 2019. Air Leak Tolerance Report: A comparison of PICO v2 (PICO 7 and PICO 14) Devices to PICO v1.6 (PICO) Devices. Internal Report. RD/19/006.
47. Sharpe A, Myers D, Searle R. Using single use negative pressure wound therapy for patients with complicated diabetic foot ulcers: an economic perspective. *Wounds UK*. 2018;14(3):80–137.
48. Myers D, Sharpe A. EP549 Service Delivery in Complex DFU Patients using Single Use NPWT - A UK Perspective. Paper presented at: EWMA; 2018; Krakow.
49. Smith+Nephew December 2018. Acoustic Testing Report: PICO v2 (PICO 7 and PICO 14) Devices. Internal Report. RD/18/131.
50. PICO \diamond 7 User Manual (English). 08493. Date of issue October 2-2017.
51. Stryja J, Staffa R, Říha D, Stryjová K, Nicielníková K. Cost-effectiveness of negative pressure wound therapy in outpatient setting. *Prolekare*. 2015.
52. Smith+Nephew July 2018. PICO \diamond 7Y Non-NPWT Wound Model Summary. Internal Report. DS.18.260.R
53. Casey C, Ambler G, Huddleston E. Consistent delivery of therapeutic negative pressure levels by a single use negative pressure wound therapy system (sNPWT) in a wound model. Poster EP079 presented at the European Wound Management Association annual meeting, June 5-7, 2019; Gothenburg, Sweden
54. Smith+Nephew 2021. PICO™ Pressure Mapping Study. Internal Report. DS/19/211/R - Part B
55. Smith+Nephew 2019. PICO \times Biomechanical Study. Internal Report. DS/19/211/R.
56. Watkins H, Brownhill R, Huddleston E. Unlocking the mode of action: why might a no-canister, ultra-portable, single use negative pressure wound therapy (sNPWT) device be clinically superior to traditional negative pressure wound therapy (tNPWT)? Poster presented at the Wounds UK Annual Conference. November 4-6, 2019; Harrogate, UK. Available at: <https://epostersonline.com/wounds2019/node>. Accessed March 2021.
57. Innocenti M, Santini M, Dreassi E, et al. Effects of cutaneous negative pressure application on perforator artery flow in healthy volunteers: a preliminary study. *J Reconstr Microsurg*. 2019;35(3):189–193.
58. Wilkinson HN, Longhorne FL, Roberts ER, Brownhill VR, Hardman MJ. Cellular benefits of single-use negative pressure wound therapy demonstrated in a novel ex vivo human skin wound model. *Wound Rep Reg*. 2021;29:298–305.
59. Smith+Nephew 2018. Summary wound model report for Opal PICO 7. Internal Report DS/18/015/R.
60. Mcmanus H, Woodmansey E. Bacterial retention within a multi-layered absorbent AIRLOCK™ Technology Single Use Negative Pressure Wound Therapy (sNPWT) dressing. Paper presented at: EWMA; 2018; Krakow, Poland.

References



61. Smith+Nephew April 2018. PICO 7Y - Waterproofness of dressings. Internal Report. DS/18/134/R.
62. Hurd T, Gilchrist B. Single use negative pressure wound therapy (sNPWT) in the community management of chronic open wounds deeper than 2cm. Paper presented at: Symposium on Advanced Wound Care/Wound Healing Society Meeting; 2020; Abu Dhabi.
63. Dowsett, et al. JCN supplement 2015, Vol 29, No 5.
64. Smith+Nephew 2021. A Prospective, Randomized, Comparative Effectiveness Study of a Single-Use, Negative Pressure Wound Therapy System (PICO™) versus a Traditional Negative Pressure Wound Therapy System (tNPWT) in the Treatment of Lower Extremity Ulcers- Dressing Application. Internal Statement. EO.AWM.PCS261.004.v1.
65. Campbell PE, et al. Retrospective clinical evaluation of gauze-based Negative Pressure Wound Therapy. Int Wound J; 5 (2): 280–286.
66. Bondokji S, et al. Clinical efficacy of a new variant of a foam- based NPWT system. J Wound Care; 20: (2) 62, 64–7
67. Malmsjö M, Borgquist O. NPWT Settings and Dressing Choices Made Easy. Wounds International (2010);1(3).
68. Molnar JA, Simpson JL, Voignier DM, Morykwas MJ, Argenta LC. Management of an Acute Thermal Injury with Subatmospheric Pressure. J Burns Wounds. 2005;4:83-92.



Advanced Wound Management

Smith & Nephew Medical Ltd

101 Hessle Road

Hull HU3 2BN, UK

T +44 (0) 1482 225181

F +44 (0) 1482 328326

°Trademark of Smith+Nephew

All Trademarks acknowledged

AWM-AWD-36940. Published Aug 2022

©2022 Smith+Nephew. GMC

[smith-nephew.com](https://www.smith-nephew.com)

For detailed product information, including indications for use, contraindications, effects, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.