


CORI[®] Surgical System with RI.KNEE: the first robotics platform indicated for use in revision total knee arthroplasty (rTKA)¹

rTKA presents unique technical challenges and decision-points for the orthopaedic surgeon, such as restoring the anatomical joint line,² addressing bone loss³ and achieving balanced flexion/extension gaps.³ Restoration of the joint line during rTKA has shown to have a positive effect on the clinical outcome^{4,5} and in one study, patients whose joint line was restored to within ±5mm of preoperative values, were 3.88 times more likely to meet or exceed the substantial clinical benefit for KOOS-JR.⁴

CORI System with RI.KNEE for rTKA enables joint line restoration and improves PROMs⁶


Joint line restoration

Enabling joint line restoration in rTKA⁶




New clinical data (retrospective case series; n=115), has demonstrated joint line restoration, to that of the contralateral native knee, in **93% of cases**, when applied across a wide range of **complex revision surgeries**

Image supplied with permission from Dr Seyler. LEGION[®] RK System with CORI Surgical System with RI.KNEE




PROMs improvement

Significantly improving PROMs, compared to pre-operative values (p<0.01)⁶




Pain decreased from 6.1 to 3.0 (p<0.01; final follow-up)




PROMIS Depression score decreased from 51.2 to 47.9 (p=0.01; final follow-up)

Post-operative outcomes


Reduced length of hospital stay and readmissions compared to conventional rTKA⁶



Mean length of hospital stay (days)



Readmission rate at **30 days**



Readmission rate at **90 days**

	Conventional rTKA		
	CORI rTKA	Study 1 ⁷	Study 2 ⁸
LoS (days)	2.3	5.6	4.6
30-day readmission	1%	--	5.7%
90-day readmission	3%	23%	13%

Figure. Comparison of CORI rTKA readmission rates to conventional methods reveals a reduction in readmissions of 4.7% at 30 days and up to a 20% reduction at 90 days

Abbreviations: KOOS-JR = Knee Injury and Osteoarthritis Outcome Score for Joint Replacement; LoS = length of stay, PROMs = patient-reported outcomes measurement; PROMIS = Patient-Reported Outcomes Measurement Information System.

References: 1. Food and Drug Administration. Available at: https://www.accessdata.fda.gov/cdrh_docs/pdf22/K220958.pdf. Accessed October 13, 2023. 2. Hofmann AA, Kurtin SM, Lyons S, Tanner AM, Bolognesi MP. Clinical and radiographic analysis of accurate restoration of the joint line in revision total knee arthroplasty. *J Arthroplasty*. 2006;21:1154–1162. 3. Graichen H. TKA revision — reasons, challenges and solutions. *J Orthop*. 2014;11:1–4. 4. Buller LT, Metzger CM, Deckard ER, Meneghini RM. The effect of joint line elevation on patient-reported outcomes after contemporary revision total knee arthroplasty. *J Arthroplasty*. 2022;37:1146–1152. 5. Sadaka C, Kabalan Z, Hoyek F, Fares GA, Lahoud J. Joint line restoration during revision total knee arthroplasty: an accurate and reliable method. *Springerplus*. 2015;4:736. 6. Cochrane NH, Kim BI, Stauffer TP, et al. Revision total knee arthroplasty with an imageless, second-generation robotic system. *J Arthroplasty*. Published online February 12, 2024. 7. Nichols CI, Vose JG. Clinical Outcomes and Costs Within 90 Days of Primary or Revision Total Joint Arthroplasty. *J Arthroplasty*. 2016;31(7):1400–1400. 8. Schairer WW, Sing DC, Vail TP, Bozic KJ. Causes and frequency of unplanned hospital readmission after total hip arthroplasty. *Clin Orthop Relat Res*. 2014;472(2):464–470.



Scan for more information

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area. For detailed product information including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.