

Anterior Cut First Surgical Technique



Introduction

The GENESIS® II Total Knee System has been designed to offer the orthopaedic surgeon solutions to address intraoperative situations. Implant function is directly related to accurate surgical technique. GENESIS II instrumentation has been developed to be an easy-to-use system that will assist the surgeon in obtaining accurate and reproducible knee alignment.

The instrumentation can be used in minimally invasive or standard exposures. While it has been the designers' objective to develop accurate, easy-to-use instrumentation, each surgeon must evaluate the appropriateness of the following technique based on his or her medical training, experience and patient evaluation.

Contributing Clinicians

Robert B. Bourne, MD, FRCSC
Chief of Orthopaedic Surgery
University Hospital
The University of Western Ontario
London, Ontario, Canada

William B. Smith, MD
Assistant Clinical Professor in Orthopaedic Surgery
Medical College of Wisconsin
Columbia Hospital
Milwaukee, Wisconsin

Steven B. Haas, MD, MPH
Associate Professor of Orthopaedic Surgery
Weill Medical College of Cornell University
Associate Chief of the Knee Service
The Hospital for Special Surgery
New York, New York

Mark A. Snyder, MD
Clinical Instructor
University of Cincinnati
Orthopaedic Surgeon
Christ Hospital
Cincinnati, Ohio

Richard S. Laskin, MD
Professor of Orthopaedic Surgery
Weill Medical College of Cornell University
Co-Chief, Knee Service
The Hospital for Special Surgery
New York, New York

Todd V. Swanson, MD
Desert Orthopaedic Center
Las Vegas, Nevada

Jan Victor, MD
Department of Orthopaedics
St. Lucas Hospital
Brugge, Belgium

Michael D. Ries, MD
Professor and Vice Chairman
University of California, San Francisco
Department of Orthopaedic Surgery
San Francisco, CA

ACF Short Technique	4
Preop Planning	6
Femoral Preparation	7
Tibial Preparation	15
Tibial Sizing	24
Posterior-Stabilized Femoral Resection	26
Resurfacing Patellar Preparation	28
Biconvex Patellar Preparation	35
Component Trialing	37
Implantation	40
Appendix A: GENESIS [®] II P-S High Flex Insert Impaction Technique for Minimally Invasive Surgery (MIS)	43
Appendix B: Articular Insert Interchangeability Chart	44
Appendix C: Anterior and Posterior Referencing	45
Appendix D: All-Polyethylene Tibial Baseplates	46
Catalog Information	47
Implant Compatibilities	53

Nota Bene:

This surgical technique is provided for educational and informational purposes only. Smith & Nephew does not provide medical advice and it is not intended to serve as such. It is the responsibility of the treating physician to determine and utilize the appropriate products and techniques according to their own clinical judgment for each of their patients. For more information on the product, including its indications for use, contraindications, and product safety information, please refer to the product's label and the Instructions for Use (IFU).

Anterior Cut First Short Technique

Femoral Preparation



Use the 9.5mm drill to open up the femoral canal and slide the valgus alignment assembly up the IM rod until it contacts the distal femur.



Place the anterior stylus tip on the lateral ridge of the anterior cortex to determine resection level.



Resect anterior cortex.



Attach distal cutting block and distal stylus to valgus alignment assembly.



Remove valgus alignment assembly and resect the distal femur.



Place the stylus tip of the femoral sizing guide on the provisional anterior cut surface. Read the size indicated by the line across the stylus shaft. If in-between sizes, choose the smaller size.



Place the correctly sized A/P cutting block on the distal femur and resect the femur.

Tibial Preparation



Extramedullary tibial alignment:

Assemble extramedullary tibial guide and place on tibia. Align guide over medial third of the tibial tubercle and parallel to the tibia.



Intramedullary tibial alignment:

Place intramedullary alignment assembly on the tibia. The alignment rod should align with the medial third of the tibial tubercle. Impact assembly.



Attach the tibial stylus to the tibial cutting block and lower the cutting block until the stylus touches the lowest point on the least affected side of the tibia. Once the resection level is determined, insert pins to secure and remove alignment assembly.



Resect the proximal tibia.



Size the tibia.

Posterior-Stabilized



Attach the PS collet to the PS housing block by tightening the gold thumbscrew, then pin to the distal femur.



Ream through the collet until the depth stop contacts the collet and then move reamer anterior and posterior until it contacts the depth stops.



Impact the housing box chisel anteriorly and posteriorly through the housing resection collet to square the corners of the housing.

Final Preparation



Prepare the patella using surgeon's preferred technique.



After trial ROM and alignment checks, select the appropriate trial fin punch and punch through the trial.



Seat the tibial implant with the tibial impactor.



Place the femoral implant on the femur and use the femoral impactor to fully seat the implant.



Place the patellar implant onto the patella and clamp onto the bone to pressurize.



Attach the articular inserter/extractor to the tibial tray (for standard inserts). Lift inserter superiorly until the anterior lip of the insert is fully seated.


Preop Planning

Determine the angle between the anatomical and the mechanical axes. This measurement will be used intraoperatively to select the appropriate valgus angle so that correct limb alignment is restored. *(Beware of misleading angles in knees with a flexion contracture or rotated lower extremities.)* The T-template provided as part of the GENESIS[®] II templates will help in this determination.

Tip: Many surgeons prefer to simply select a standard angle for the distal femoral cut (i.e., 5° 6° or 7°) based on the patient and surgical experience.

Recommended Sawblades*	
Cat. No.	Description
71512901	Stryker 2000 (Fanned)
71512905	Stryker 2000 (Straight)
71512911	Hall Powerpro (Fanned)
71512904	3M (Fanned)
Or any 0.053" or 1.35mm thickness sawblade	

*For MIS-style blocks only.



M = Mechanical Axis
A = Anatomical Axis
T = Transverse Axis
V = Vertical Axis

Femoral Preparation

Intramedullary Femoral Alignment

1. Identify the rotational reference landmarks:
 - A/P axis (as described by Whiteside)
 - Medial-lateral posterior femoral condyles
 - Epicondylar axis
2. Open the femoral canal (generally just anterior to the PCL insertion) with the 9.5mm drill (*Figure 1*).

Instrument Assembly:

- a. Attach the modular T-handle to the intramedullary rod.
- b. Select the appropriate valgus angle bushing based on preoperative measurements.
- c. Slide the bushing into the valgus alignment guide (left or right). Make sure the bushing is positioned so that “left” is facing anteriorly when operating on a left knee and “right” is facing anteriorly when operating on a right knee.
- d. Attach a quick-connect handle to the valgus alignment guide.
- e. Slide the rod through the bushing (*Figure 2*).

3. Insert the intramedullary rod into the canal. Position the valgus alignment guide until it contacts the distal femur (*Figure 3*).

Note: Do not engage the floating pins until rotation is set.




Figure 1



Figure 2




Figure 3

Femoral Preparation

Femoral Rotational Alignment Without Paddles

Rotation of the valgus alignment guide is set neutral to the posterior femoral condyles by using the landmarks described on page 7, step 1, either with or without rotational alignment paddles.

Without Paddles

1. Flex the knee to 70° - 90°.

2. Align:

A. The femoral alignment stylus (*Figure 4*) with the A/P axis. The femoral alignment template is designed such that setting it parallel to the A/P axis aligns the valgus alignment guide in neutral rotation.

Use a bovie or pen to mark the A/P axis line (which is represented by the deepest part of the trochlear groove).

The femoral alignment stylus is placed over the valgus angle bushing to guide rotational alignment (*Figure 5*). Make sure that the template is positioned so that “left” is facing out when operating on a left knee and “right” is facing out when operating on a right knee. The valgus alignment guide is placed in neutral orientation by aligning the outrigger of the template with the A/P line (*Figure 6*).

B. The posterior aspect of the valgus alignment guide is parallel to the posterior condyles.

C. The line laser-etched across the distal surface of the valgus alignment guide parallel to the epicondylar axis. The line on the valgus alignment guide is drawn such that placing it parallel to the epicondylar axis aligns the guide in neutral rotation.

3. Secure the valgus alignment assembly on the distal femur by impacting the floating spikes.

MIS Tip: The surgeon may use a pin driver or a tibial punch to gain clearance to impact the spikes.



Figure 4



Figure 5




Figure 6

With Paddles

Instrument Assembly:

Unlock the capture mechanism on the modular paddles. The arm on the paddles distracts posteriorly and rotates to either side to unlock so the anterior lip can engage the slot in the posterior aspect of the valgus alignment guide.

1. Insert the anterior lip of the paddles into the slot in the valgus alignment guide. Rotating the arm back centrally into the recess will lock the paddles onto the valgus alignment guide (Figure 7).
2. Position the paddles under the posterior condyles (Figure 8).

Tip: It may be helpful to flex the knee greater than 90°, as this will help fit the paddles under the femoral condyles.

Tip: Posterior condylar referencing may be less reliable in knees with deficient posterior condyles (e.g. severe valgus deformity). If the posterior condyles are deficient, the A/P or epicondylar axis should be used to determine alignment.

3. Secure the valgus alignment assembly on the distal femur by impacting the floating spikes.

MIS Tip: The surgeon may use a pin driver or a tibial punch to gain clearance to impact the spikes.




Figure 7



Figure 8

Femoral Preparation

Preliminary Anterior Femoral Resection

1. Release and remove the modular paddles, if applicable.

Instrument Assembly:

Place the anterior resection guide into the valgus alignment guide (Figure 9) and attach the anterior stylus to the anterior resection guide by sliding the foot into the cutting slot (Figure 10).

2. Place the anterior stylus tip on the lateral ridge of the anterior femoral cortex. Pin the anterior resection guide with a 1/8" trocar pin in any available hole and remove the anterior stylus (Figure 11).
3. Resect the anterior cortex (Figure 12).

Tip: Removing a small area of soft tissue down the bone over the distal and lateral anterior femur allows the stylus to fully seat. This will help prevent overestimating femoral sizing, which could lead to "over-stuffing" of the patellofemoral joint.




Figure 9




Figure 10



Figure 11




Figure 12

Distal Femoral Resection

Instrument Assembly:

- a. Assemble the distal femoral cutting block with the distal resection stylus by pressing the gold button and sliding the stylus until it hits a stop (*Figure 13*).
- b. The word “primary” should show through the cutting slot. This will resect the standard 9.5mm from the distal femur. For a large femur or in the case of flexion contracture, up to 7mm additional resection can be taken by sliding the cutting block proximally so that the desired resection level shows through the cutting slot.

1. Secure the distal femoral cutting block to the anterior cortex by impacting or drilling unheaded or headed pins through the holes marked “0.” Use of a third oblique pin is recommended for additional stability (*Figure 14*).




Figure 13




Figure 14

Femoral Preparation

2. Attach the slap hammer to the valgus bushing and remove the rod, distal resection stylus and valgus alignment assembly (*Figure 15*).
3. Only the distal femoral cutting block should remain on the femur (*Figure 16*).
4. Resect the distal femur (*Figure 17*), then remove the distal femoral cutting block.

Tip: To take an additional distal resection after resection, simply reposition the block through the pin holes marked +2, +4 or +6mm for the desired level of resection after removing the oblique pin.



Figure 15



Figure 17




Figure 16

Femoral Sizing

MIS Surgical Tip: To make sizing easier, you may wish to resect the tibia before further femoral preparation.

Instrument Assembly:

Attach a quick-connect handle to the femoral sizing guide.

1. Place the femoral sizing guide on the distal femur, and place the stylus tip on the provisional anterior cut (Figure 18).
2. Read the size indicated by the line across the stylus shaft (**not the retention pin**) and choose the smaller size if in-between two sizes (Figure 19).

Tip: Traditionally, surgeons using an Anterior Cut First/Anterior Referencing approach have always chosen the smaller size component between sizes. However, some surgeons choose to use the larger size, particularly when using a P-S component, as sacrificing the PCL can increase the flexion space from 4-6mm and could leave the knee loose in flexion if the smaller component is used.



Figure 18




Figure 19

Femoral Preparation

A/P Femoral Resection

1. Choose the correct size A/P cutting block and place it on the distal femur in a medialized position. M/L positioning of this block is not critical.

2. Secure the A/P cutting block first with a straight pin through either the medial, lateral or central pin hole on the distal block surface. Pin through the angled holes in the ears on the medial and lateral sides of the block as bone quality dictates to achieve stability (*Figure 20*). Remove the pin(s) from the distal face before making chamfer cuts.

Tip: The A-P cutting block should seat flush with the cut anterior and distal surfaces.

3. Complete the anterior, posterior and chamfer cuts (*Figures 21-24*). If desired, the chamfer cuts may be made through the posterior-stabilized femoral housing block or a dedicated chamfer cutting block.




Figure 20




Figure 21




Figure 22




Figure 23




Figure 24

Tibial Preparation

The system allows the surgeon to perform either extramedullary or intramedullary tibial alignment. For intramedullary tibial alignment, turn to page 19.

When using the extramedullary tibial alignment, the surgeon may use a non-spiked or spiked fixation rod. For tibial preparation using the extramedullary guide with a non-spiked fixation rod see below. For tibial preparation using the extramedullary guide with a spiked rod, turn to page 17.

Extramedullary Tibial Alignment

Instrument Assembly:

- Insert the ankle clamp into the distal end of the alignment tube and thread the locking pin into the ankle clamp (Figure 26).
- After the ankle clamp is moved into the proper position, lock into place with the gold knob.
- Choose the correct left or right tibial cutting block. Select the spiked or non-spiked fixation rod.



Figure 26

Non-spiked Fixation Rod

Instrument Assembly:

- Place the appropriate left or right tibial cutting block on top of the disc on the non-spiked fixation rod (Figure 27). Tighten the central knob to lock the block into position.
- Introduce the rod into the Extramedullary Assembly and adjust and lock the cam in the assembly.




Figure 27

Figure 28

- Place the arms of the extramedullary alignment clamp around the ankle, and adjust the distal M/L slide directly over the middle of the tibiotalar joint, which is also approximated by the second ray of the foot proximal to the malleoli (Figure 28). The cutting block on the proximal end of the assembly should be proximal to the tibial tubercle (Figure 29).




Figure 29

Tibial Preparation

2. Assess rotation of the alignment guide and slope of the cutting plane. The goal is to align the extramedullary alignment assembly rotationally so that it aligns over the medial third of the tibial tubercle and over the second toe (*Figure 30*).
3. Rotational alignment is critical due to the 3° posterior sloped cut. The slope can be adjusted according to the patient's anatomy (*Figure 31*).

Note: 3-5° of slope is built into the articular insert (depending on which insert is chosen) and 3° of slope is built into the tibial cutting block. A neutral or slightly sloped alignment should usually be chosen.

Tip: Neutral or minimally sloped alignment may be achieved by palpating the fibula followed by aligning the alignment guide parallel to the fibula. Tibial bowing and soft tissue bulk may make external tibial referencing unreliable.



Figure 30



Figure 31

Spiked Fixation Rod

Instrument Assembly:

- a. Place the spiked fixation rod through the central anterior hole in the tibial cutting guide; adjust the block and tighten the central knob to lock the block in position.
- b. Introduce the spiked fixation rod into the proximal end of the alignment assembly and adjust and lock the cam on the assembly (Figure 32).

1. Place the arms of the extramedullary alignment clamp around the ankle, and adjust the distal M/L slide directly over the middle of the tibiotalar joint which is also approximated by the second ray of the foot proximal to the malleoli (Figure 33). The cutting block on the proximal end of the assembly should be proximal to the tibial tubercle (Figure 34).
2. Impact the longer spike of the spiked fixation rod into the proximal tibia (Figure 35).



Figure 32



Figure 33



Figure 34




Figure 35

Tibial Preparation

3. Assess rotation of the alignment guide and slope of the cutting plane. The goal is to align the extramedullary alignment assembly rotationally so that it aligns over the medial third of the tibial tubercle and over the second toe (*Figure 36*).
4. Rotational alignment is critical due to the 3° posterior sloped cut. The slope can be adjusted according to the patient's anatomy (*Figure 37*). Impact the second spike to secure the assembly (*Figure 38*).

Note: 3-5° of slope is built into the articular insert (depending on which insert is chosen) and 3° of slope is built into the tibial cutting block. A neutral or slightly sloped alignment should usually be chosen.

Tip: Neutral or minimally sloped alignment may be achieved by palpating the fibula followed by aligning the alignment guide parallel to the fibula. Tibial bowing and soft tissue bulk may make external tibial referencing unreliable.



Figure 36



Figure 37




Figure 38

Intramedullary Tibial Alignment

Instrument Assembly:

- a. Insert the external rod of the Intramedullary Tibial Alignment Guide through the middle hole on the correct left or right tibial cutting block and lock the cam (*Figure 39*).
- b. Attach the T-handle to the IM rod and pass it through the cannulated alignment sleeve on the alignment assembly (*Figure 40*).

1. Make a 9.5mm pilot hole into the tibial canal (*Figure 41*) (generally 5mm medial to the midline). A preliminary resection of the tibial spine may facilitate seating of the tibial drill guide onto the proximal tibia.



Tibial Preparation

2. Slowly insert the IM rod into the tibial canal.
3. Assess rotation of the intramedullary tibial alignment guide. Rotational alignment is critical due to the 3° posterior sloped cut. The alignment rod of the intramedullary tibial alignment assembly should align with the medial third of the tibial tubercle (*Figure 42*).
4. Impact the proximal end of the cannulated alignment sleeve to drive the distal spikes into the proximal tibia to lock rotational alignment (*Figure 43*).



Figure 42




Figure 43

Tibial Resection

1. Attach the tibial stylus to the tibial cutting block by inserting the stylus foot into the cutting slot.
2. Lower the cutting block until the stylus touches the low point on the less affected side of the tibia (Figures 44). The stylus can be adjusted for a 9, 11 or 13mm tibial resection by twisting the knob on top of the stylus.
3. Pin the tibial cutting block to the tibia by inserting pins first through the central holes; then the oblique hole.

Tip: Pinning through the central holes marked 0mm with smooth pins will allow the block to be moved +2mm should additional resection be required (Figure 45).

Tip: A 9mm resection is recommended since 9mm of metal and plastic is the thinnest available component.

Tip: To do an extramedullary alignment check, place the extramedullary alignment rod through the tibial cutting block.



Figure 44




Figure 45

Tibial Preparation

4. To remove the assembly:

- a. For the intramedullary alignment assembly, use the universal extractor leaving the cutting block on the anterior tibia (*Figure 46*).
- b. For the extramedullary assembly with spiked rod, release the cam at the top of the alignment tube and use the slap hammer to remove the spiked fixation rod (*Figure 47*).
- c. The extramedullary assembly with the non-spiked rod may be left in place or removed by loosening the thumbscrew and lowering the non-spiked rod to disengage from the tibial cutting block.



Figure 46



Figure 47

5. Cut the tibia by first directing the blade in the posterior direction and then laterally (*Figure 48*).
6. Check alignment and balance with spacer block and rod (*Figures 49 & 50*). Balance ligaments in standard fashion.

Tip: Since the spacer block has one end for flexion and one for extension, ensure that the appropriate end is used.



Figure 48



Figure 49



Figure 50

Tibial Sizing

Option A – Stemless Tibial Trials

1. Attach a quick-connect handle to a stemless trial one size below the femoral component size and place on the cut tibia to assess coverage (*Figure 51*). As needed, additional sizes should be templated using the stemless trials.
2. Once the appropriate size is determined, pin the medial size of the selected stemless trial with a short headed pin.
3. Place a trial insert into the stemless tibial trial tray and perform a trial range of motion to allow the baseplate to center on the femoral trial.
(As a secondary check, the surgeon may pass the alignment rod through the quick-connect handle to assess alignment) (*Figure 52*). Pin the lateral side of the trial.

Tip: After putting the knee through a trial ROM, the surgeon should note the proper rotation of the trial tibial component on the proximal tibia and mark the tibia for future reference.

Tip: The center-line marks on the femoral and tibial trial components usually line up.

4. Using the tibial fin/stem punch, rotational alignment may be set now or at the time of trial placement.
See page 37.

Tip: In the case of sclerotic bone, first drill for the stem using the 11mm tibial drill. To avoid fracture, predrill the tibial plateau with a 1/8" drill bit.




Figure 51




Figure 52

Option B – Stemmed Tibial Trials

1. Place a tibial drill guide one size below the femoral component size on the cut tibia to assess coverage. As needed, additional sizes should be templated (Figure 53).
2. Once the tibial drill guide has been centralized on the proximal tibia, pin the drill guide in place. Retract the gold collar on the drill guide handle and insert the 11mm tibial collet.
3. With the 11mm tibial collet in place, drill with the 11mm tibial drill (Figure 54) and punch with the 11mm tibial punch (Figure 55). If a 9.5mm hole has already been made for use of the intramedullary tibial alignment assembly, you only need to utilize the 11mm tibial punch at this time.
4. Remove the tibial drill guide.
5. Place the stemmed tibial trial into the prepared hole.
6. Using the tibial fin punch, rotational alignment may be set now or at the time of trial placement. See page 40.

Tip: After putting the knee through a trial ROM, the surgeon should note the proper rotation of the trial tibial component on the proximal tibia and mark the tibia for future reference.

Tip: The center-line marks on the femoral and tibial trial components usually line up.



Figure 53



Figure 54




Figure 55

Posterior-Stabilized Femoral Resection

1. Attach the PS housing resection collet to the housing resection block by tightening the gold thumbscrew in the most anterior position (*Figure 56*).
2. The PS housing resection block must be centered on the femur, as this will determine component position.

Tip: The housing resection blocks have the same M-L dimension as the implants.

Tip: The only difference between the cruciate-retaining and the posterior-stabilized femoral components is the addition of the housing for the cam mechanism. All other box dimensions are the same. The anterior and posterior chamfer resections can be made through the posterior-stabilized housing resection block.

3. Secure with 1/8" trocar pins through the straight holes in the front of the block. If the chamfer cuts are made through this block, the angled holes in the sides of the block should be used.

Instrument Assembly:

Attach the housing reamer dome and the PS reamer sleeve to the patellar reamer shaft (*Figure 57*).




Figure 56



Figure 57

4. Ream through the housing resection collet until the automatic depth stop contacts the collet, loosen the thumbscrew and then move the reamer anterior and posterior until it contacts the automatic stop (*Figure 58*).
5. Impact the housing box chisel through the housing resection collet to square the corners of the housing. The housing box chisel should be used anteriorly and posteriorly to ensure that the full length of the box is prepared (*Figure 59*).
6. If the chamfer resections have not been made, they can now be made by cutting through the chamfer slots in the housing resection block.




Figure 58




Figure 59

Resurfacing Patellar Preparation

The surgeon can choose from a freehand cutting technique with towel clips or if desired, he or she can choose one of the following instrumented techniques.

Resection Guide Technique

1. Measure the overall thickness of the patella with the patellar calipers (*Figure 60*).
2. Subtract from this number the thickness of the GENESIS® II round resurfacing patellar component, which is 9mm.

Note: The thickness of the GENESIS II oval resurfacing patellar component varies by diameter. See the chart on page 33.

3. The guide is set at the amount of bone that should remain after cutting the patella – i.e. the difference between the original patellar thickness and 9mm. The guide is set at this level by turning the knurled knob (*Figure 61*).

For example:

A. Measure the overall thickness of the patella with the patellar calipers. For this example, the patella measures 25mm.

B. Subtract the thickness of the round resurfacing patellar component. In this example, 9mm. (25mm - 9mm = 16mm). The guide should be set at 16mm for this example.



Figure 60




Figure 61

4. Cut the patella through the dedicated saw guides (Figure 62).
5. Drill for the three pegs (Figure 63), insert the resurfacing patellar trial and remeasure. The overall thickness should be equivalent to the original thickness (Figure 64).

Reaming Technique

The reaming technique described for the biconvex patella on page 35 can be used with the resurfacing patellar implant as well. The only differences in technique between it and the biconvex are the use of the RED resurfacing depth gauge, resurfacing reamers and the resurfacing drill guides.



Figure 62




Figure 63



Figure 64

Resurfacing Patellar Preparation

Patellar Large Reamer Resurfacing Instrumentation

The objective of this technique is to resurface the articular surface of the patella with the precision of a reaming technique. The reamed patellar surface can accommodate an oval or round resurfacing patellar component.

1. Trim tissue surrounding the patella using electrocautery (bovie) (Figure 65).
2. Use a rongeur to remove osteophytes and reduce the patella to its true size (Figure 66). The bovie should also be used to release soft tissue attachments to the estimated level of resection.
3. Place the collet over the patella so that it fits snugly around the patellar diameter (Figure 67). The goal is to reduce the patella to its smallest diameter so that the smallest possible collet will fit around the entire patella. Use the patellar reamer collet as a sizing template to select the appropriately sized collet and reamer.

Tip: The collet should be resting on the soft tissue surrounding the patella. If the patella does not enter the collet evenly but instead enters at an angle, the collet may not be completely surrounding the patella, but instead resting on part of the bone. If the collet is only slightly smaller than the patella, you may trim 1-2mm of the medial and lateral edges of the patella to ensure a snug fit. If the collet is far smaller than the patella, choose the next size up and assess fit.




Figure 65



Figure 66



Figure 67

Surgeon Acknowledgement: The technique for the Patellar Large Reamer Resurfacing System was developed in conjunction with Warren Jablonsky, MD, McHenry County Orthopedics, Crystal Lake, IL.

4. Measure patellar thickness with the patellar calipers (Figure 68).

Tip: The patella should measure a minimum of 19mm before reaming to use this resurfacing technique.

Determine the design and diameter of the patellar implant to be used. A round or oval resurfacing design may be chosen. The round resurfacing patella is 9mm thick, and the depth stop for this technique prepares for 9mm resection. The oval patella's thickness is variable.

Tip: Minor adjustments may be necessary at the time of resection to accommodate the largest diameter oval patellar implants. Please see chart on page 33.

Instrument Assembly:

- Slide the correct diameter of patellar reamer collet into place on the patellar reamer guide.
- Attach the patellar reamer guide to the patella.
- Secure the patellar reamer guide on the patella by tightening the set screw.
- Attach the matching size patellar reamer dome and large patellar depth stop to the patellar reamer shaft.

5. Rotate the BLACK resurfacing patellar depth gauge around so that the hooked end or "claw" surrounds the patellar reamer shaft (Figure 69). Lower the depth stop by compressing the button until it meets the depth gauge (Figure 70). Remove the depth gauge from the assembly. Ream the patella until the depth stop engages the patellar reamer guide (Figure 71).

Tip: Excessive force on the reamer shaft may alter the depth of resection, causing overreaming.



Figure 68



Figure 69



Figure 70



Figure 71

Resurfacing Patellar Preparation

6. After reaming, the patella should have a completely flat articular surface (*Figure 72*). Measure the resected patella to ensure adequate resection (the resected patella should measure its original depth minus 9mm).
7. Drill the appropriate fixation holes for the resurfacing patellar implant using the correctly sized drill guide and resurfacing drill (*Figure 73*).
8. Place the patellar trial into the prepared patella. If desired, use the calipers to remeasure the composite thickness of bone and trial.



Figure 72



Figure 73

Oval Patellar Preparation

The oval patellar implant can be prepared for use with any resurfacing technique; however, there are a few differences in final preparation. The patella has to be implanted in the proper orientation, where the extended lateral flange will be riding on the lateral side of the femoral component.

The oval patellar implant does not have the same thickness for all sizes. This is due to the varying offset needed to obtain the correct design for the different diameters. (See the chart for sizing/ thickness options.)

1. Mark the medial facet axis of the patella superior and inferiorly with a marking pen or use the laser etch line on the sizing guide to mark the vertical ridge of the patella.
2. Measure the depth of the patella at its maximum depth centrally along the medial facet (Figures 74 & 75).

Oval Patellar Sizing Options	
Oval Resurfacing Implant	
Diameter	Thickness
29mm	8.5mm
32mm	9.0mm
35mm	9.0mm
38mm	9.5mm
41mm	10.0mm



Figure 74




Figure 75

The technique for the Oval Patella was developed in conjunction with William J. Robb III, MD, Illinois Bone and Joint Institute, Glenbrook Hospital, Evanston Northwestern Healthcare.

Resurfacing Patellar Preparation

3. Resect the patella using the preferred method.
4. Measure the diameter of the resected patella with the trial templates (*Figure 76*).
5. Centralize the thickest portion of the prosthetic patella along the line of the previously marked medial facet eminence.
6. Place the appropriate drill guide on the patellar reamer guide and clamp the guide to the patella. Drill to the measured depths (*Figure 77*).
7. Place the trial on the patella and remeasure the patella if desired (*Figures 78 & 79*).



Figure 76




Figure 77



Figure 78



Figure 79

Biconvex Patellar Preparation

Biconvex Patella

Instrument Assembly:

Determine the appropriate diameter patellar implant, and select the correctly-sized patellar reamer collet and slide it into place on the patellar reamer guide (*Figure 77*).

1. Attach the patellar reamer guide to the patella. Tighten the patellar reamer guide on the patella (*Figure 81*).
2. Use the patellar calipers to measure the thickness of the patella (*Figure 82*).



Figure 80




Figure 81




Figure 82

Biconvex Patellar Preparation

Instrument Assembly:

- a. Attach the BLUE patellar depth gauge to the reamer guide (Figure 83).
- b. Attach the matching sized patellar reamer dome and patellar depth stop to the patellar reamer shaft (Figures 84 & 85). Lower the assembly through the patellar reamer guide until the reamer dome contacts the patella.

3. Swing the patellar depth gauge around so that the "claw" surrounds the patellar reamer shaft.
4. Lower the patellar depth stop by pushing the gold button until it contacts the patellar depth gauge. The patellar depth stop will automatically lock in place (Figure 86).
5. Remove the depth gauge.
6. Ream the patella until the depth stop engages the patellar reamer guide.




Figure 83



Figure 84



Figure 85




Figure 86

Component Trialing

1. Flex the knee to 90° and insert the femoral trial using the femoral trial impactor (*Figure 87*).
2. Use the appropriate insert trial (begin with a 9mm trial) to determine stability and alignment.
3. Perform a trial range of motion. The alignment marks on the front of the femoral and tibial trials should line up (*Figure 88*). The quick-connect handle may be attached to the tibial trial and used to set the appropriate rotational alignment.

*Option: Extend the knee fully with the handle attached to the tibial trial. Pass the extramedullary rod through the handle to assess full leg alignment (*Figure 89*).*

Tip: The technique of tibial trial, then femoral trial and then trial insert works for all GENESIS® II inserts EXCEPT the dished inserts. For the deep dished, insert the trial bearing BEFORE the femoral trial.




Figure 87




Figure 88




Figure 89

Component Trialing

4. Mark correct tibial rotational alignment on the anterior tibia using a cautery knife (*Figure 90*).
5. Determine whether a porous or nonporous tibial implant will be used. Select the appropriate tibial fin punch to prepare the fins and punch through the tibial trial (*Figure 91*).

Tip: If the tibial bone is sclerotic, begin the fin slot with a burr or thin sawblade before using the fin punch to prevent tibial fracture.

6. Place the patellar trial into the prepared patella (*Figure 92*).




Figure 90



Figure 91



Figure 92

7. Perform a trial range of motion to assess patellar tracking. With cruciate-retaining knees, medial-lateral placement of the femoral trial can be adjusted to optimize patellar tracking (*Figure 93*).
8. For cruciate-retaining femorals, prepare the femoral lug holes through the femoral trial with the femoral lug punch (*Figure 94*).
Note: This step is also required for cemented posterior-stabilizing femorals using femoral lugs or Flex-Lok pegs.
9. Remove the tibial trial. Attach the end of the universal extractor to the femoral trial (*Figure 95*). Remove the femoral trial. Use a towel clip to remove the patellar trial.




Figure 93



Figure 94



Figure 95

Implantation

Tibial Implantation

1. Apply cement on the proximal tibia and/or the implant and seat the tibial implant with the tibial impactor (*Figure 96*). Remove excess cement.

Femoral Implantation

1. Mix and prepare bone cement for femoral component and distal femur. Apply to the femoral component or prepared bone, based on the surgeon's preference.

Tip: Many surgeons put cement on the bone rather than, or supplemental to, cement on the underside of the implant.

Note: If using femoral lugs or Flex-Lok pegs, screw those components into the femur prior to cement application.

2. Place the femoral implant onto the femur and use the femoral impactor to fully seat the implant (*Figure 97*).
3. Remove excess cement. Extend the knee to remove cement anteriorly without retracting the proximal soft tissue.
4. Place the tibial insert trial onto the tibial implant and extend the leg to pressurize the cement.

Tip: Place the CR tibial trial in the tibial implant tray to assist with aligning the femoral component during implantation.




Figure 96

Figure 97

Patellar Implantation

1. Assemble the patellar cement clamp to the patellar reamer guide.
2. Apply bone cement to the patella.
3. Place the patellar implant onto the patella and clamp into the bone (*Figure 98*). Remove excess cement.

Cruciate-Retaining, Dished and Posterior-Stabilized Insert Placement

1. Determine the correct articular insert thickness.
2. Clear any debris from the locking mechanism and slide the insert into the tibial baseplate engaging the locking mechanism. For the PS insert, begin insertion in flexion and extend the leg to engage the locking mechanism.
3. Attach the articular inserter/extractor to the tibial tray. Lift the inserter superiorly until the anterior lip of the articular insert is fully seated (*Figure 99*).




Figure 98



Figure 99

Implantation

P-S High Flex and C-R Deep Flex Insert Placement

1. Attach the appropriately sized bumper (either 1-2 or 3-8) to the impactor handle.
2. Position the knee in approximately 90° flexion.
3. Align the articular insert with the locking mechanism of the tibial baseplate.
4. Push the insert posteriorly until the top of the anterior rail of the baseplate is visible.
5. Place the bumper on the anterior chamfer of the insert. The mating surfaces should be very conforming (*Figures 100 & 101*).
6. Impact the handle until the insert is fully seated.

MIS Note: To use the PS High Flexion Insert in minimally invasive surgery, please see the technique described in Appendix A.




Figure 100



Figure 101

Appendix A

GENESIS[®] II P-S High Flex Insert Impaction Technique for MIS

When using the PS High Flex Insert in a minimally invasive procedure, the femoral cam mechanism is likely to prevent the insert from fully seating into the locking mechanism while the knee is in flexion. To use the PS High Flex Insert in a MIS case:

1. Flex the knee to 90° and push the insert as far back as it will go posteriorly with the knee in flexion (*Figure 102*).

Tip: Lift the distal femur to prevent scratching of the posterior condyle of the component.

2. Placing your thumb on the anterior of the insert to hold it on the baseplate (*Figure 103*), move the knee into extension.
3. Use the impactor handle with the appropriately sized bumper to fully seat the insert and engage the anterior portion of the dovetail locking mechanism (*Figure 104*).




Figure 102




Figure 103




Figure 104

Appendix B

GENESIS[®] II Articular Insert Interchangeability Chart

Cruciate-Retaining Inserts: Completely interchangeable with all size femoral components

Posterior-Stabilized (PS), Dished (DD), High Flex Posterior-Stabilized (HFPS) and Cruciate-Retaining Deep Flex (CRDF): Limited interchangeability; chart applies.

	Femoral Size								
Insert Size	1	2	3	4	5	6	7	8	9
1-2 PS, DD									
1-2 HFPS, CRDF									
3-4 PS, DD									
3-4 HFPS, CRDF									
5-6 PS, DD									
5-6 HFPS, CRDF									
7-8 PS, DD									
7-8 HFPS, CRDF									

Appendix C

Anterior and Posterior Referencing

Anterior Referencing

An anterior referencing technique is based on the anterior cortex, which serves as the primary reference point. The anterior resection is fixed while the posterior resection varies with size. Because the component will be flush against the anterior cortex, this will enable the reapproximation of the patellofemoral joint. When the sizing guide indicates the femoral implant is between two sizes, the smaller size should be selected. Choosing the smaller size results in more bone resection from the posterior condyles thereby increasing the flexion space.

Anterior Referencing	
Advantages	Disadvantages
Reapproximation of the patellofemoral joint	Knee may be loose in flexion
Reduced chance of notching the anterior cortex	

Posterior Referencing

A posterior referencing technique is based on the posterior femoral condyles which serve as the reference point. The posterior resection remains constant while the anterior resection varies with respect to the anterior cortex. Therefore, the posterior resection will equal the posterior thickness of the prosthesis, resulting in a balanced flexion-extension space. When the sizing guide indicates the femoral implant is between two sizes, the larger size should be chosen. Even though there is a slight chance in overstuffing the patellofemoral joint with a larger size, there is a reduced risk in notching the anterior cortex of the femur.

Posterior Referencing	
Advantages	Disadvantages
Balanced flexion and extension spaces	May overstuff the patellofemoral joint

Appendix D

All-Polyethylene Tibia Baseplates

All-Polyethylene Tibia baseplates are available in both Cruciate-Retaining (CR) (Figure 105a) and Posterior-Stabilized (PS) (Figure 105b).

Follow the same surgical steps as outlined in the Tibial Preparation, Sizing and Implantation sections of this document.

An All-Poly impactor is available to help seat the component (Figure 106).

After applying cement, use the All-Poly impactor to set the baseplate into place (Figure 107).




Figure 105a



Figure 105b




Figure 106




Figure 107

Catalog Information

GENESIS[®] II PS Non-Porous Femoral (CoCr)

Cat. No.	Description
71420096	GENESIS II Non-Porous PS Femoral size 1 Left
71420098	GENESIS II Non-Porous PS Femoral size 2 Left
71420100	GENESIS II Non-Porous PS Femoral size 3 Left
71420102	GENESIS II Non-Porous PS Femoral size 4 Left
71420104	GENESIS II Non-Porous PS Femoral size 5 Left
71420106	GENESIS II Non-Porous PS Femoral size 6 Left
71420108	GENESIS II Non-Porous PS Femoral size 7 Left
71420110	GENESIS II Non-Porous PS Femoral size 8 Left
71420112	GENESIS II Non-Porous PS Femoral size 1 Right
71420114	GENESIS II Non-Porous PS Femoral size 2 Right
71420116	GENESIS II Non-Porous PS Femoral size 3 Right
71420118	GENESIS II Non-Porous PS Femoral size 4 Right
71420120	GENESIS II Non-Porous PS Femoral size 5 Right
71420122	GENESIS II Non-Porous PS Femoral size 6 Right
71420124	GENESIS II Non-Porous PS Femoral size 7 Right
71420126	GENESIS II Non-Porous PS Femoral size 8 Right

GENESIS II PS Non-Porous Femoral (OXINIUM[®])

Cat. No.	Description
71421012	GENESIS II Non-Porous PS OXINIUM Femoral size 2 Left
71421013	GENESIS II Non-Porous PS OXINIUM Femoral size 3 Left
71421014	GENESIS II Non-Porous PS OXINIUM Femoral size 4 Left
71421015	GENESIS II Non-Porous PS OXINIUM Femoral size 5 Left
71421016	GENESIS II Non-Porous PS OXINIUM Femoral size 6 Left
71421017	GENESIS II Non-Porous PS OXINIUM Femoral size 7 Left
71421018	GENESIS II Non-Porous PS OXINIUM Femoral size 8 Left
71930008	GENESIS II Non-Porous PS OXINIUM Femoral size 9 Left
71421112	GENESIS II Non-Porous PS OXINIUM Femoral size 2 Right
71421113	GENESIS II Non-Porous PS OXINIUM Femoral size 3 Right
71421114	GENESIS II Non-Porous PS OXINIUM Femoral size 4 Right
71421115	GENESIS II Non-Porous PS OXINIUM Femoral size 5 Right
71421116	GENESIS II Non-Porous PS OXINIUM Femoral size 6 Right
71421117	GENESIS II Non-Porous PS OXINIUM Femoral size 7 Right
71421118	GENESIS II Non-Porous PS OXINIUM Femoral size 8 Right
71930010	GENESIS II Non-Porous PS OXINIUM Femoral size 9 Right

GENESIS II CR Non-Porous Femoral (CoCr)

Cat. No.	Description
71420000	GENESIS II Non-Porous CR Femoral Size 1 Left
71420002	GENESIS II Non-Porous CR Femoral Size 2 Left
71420004	GENESIS II Non-Porous CR Femoral Size 3 Left
71420006	GENESIS II Non-Porous CR Femoral Size 4 Left
71420008	GENESIS II Non-Porous CR Femoral Size 5 Left
71420010	GENESIS II Non-Porous CR Femoral Size 6 Left
71420012	GENESIS II Non-Porous CR Femoral Size 7 Left
71420014	GENESIS II Non-Porous CR Femoral Size 8 Left
71420016	GENESIS II Non-Porous CR Femoral Size 1 Right
71420018	GENESIS II Non-Porous CR Femoral Size 2 Right
71420020	GENESIS II Non-Porous CR Femoral Size 3 Right
71420022	GENESIS II Non-Porous CR Femoral Size 4 Right
71420024	GENESIS II Non-Porous CR Femoral Size 5 Right
71420026	GENESIS II Non-Porous CR Femoral Size 6 Right
71420028	GENESIS II Non-Porous CR Femoral Size 7 Right
71420030	GENESIS II Non-Porous CR Femoral Size 8 Right

GENESIS II CR Non-Porous Femoral (OXINIUM[®])

Cat. No.	Description
71420130	GENESIS II Non-Porous CR OXINIUM Femoral size 2 Left
71420132	GENESIS II Non-Porous CR OXINIUM Femoral size 3 Left
71420134	GENESIS II Non-Porous CR OXINIUM Femoral size 4 Left
71420136	GENESIS II Non-Porous CR OXINIUM Femoral size 5 Left
71420138	GENESIS II Non-Porous CR OXINIUM Femoral size 6 Left
71420140	GENESIS II Non-Porous CR OXINIUM Femoral size 7 Left
71420142	GENESIS II Non-Porous CR OXINIUM Femoral size 8 Left
71930007	GENESIS II Non-Porous CR OXINIUM Femoral size 9 Left
71420146	GENESIS II Non-Porous CR OXINIUM Femoral size 2 Right
71420148	GENESIS II Non-Porous CR OXINIUM Femoral size 3 Right
71420150	GENESIS II Non-Porous CR OXINIUM Femoral size 4 Right
71420152	GENESIS II Non-Porous CR OXINIUM Femoral size 5 Right
71420154	GENESIS II Non-Porous CR OXINIUM Femoral size 6 Right
71420156	GENESIS II Non-Porous CR OXINIUM Femoral size 7 Right
71420158	GENESIS II Non-Porous CR OXINIUM Femoral size 8 Right
71930009	GENESIS II Non-Porous CR OXINIUM Femoral size 9 Right

GENESIS[®] II PS Inserts

Cat. No.	Description
71420802	GENESIS II PS Inserts Size 1-2, 9mm
71420804	GENESIS II PS Inserts Size 1-2, 11mm
71420806	GENESIS II PS Inserts Size 1-2, 13mm
71420808	GENESIS II PS Inserts Size 1-2, 15mm
71420810	GENESIS II PS Inserts Size 1-2, 18mm
71420812	GENESIS II PS Inserts Size 1-2, 21mm
71420814	GENESIS II PS Inserts Size 1-2, 25mm
71420816	GENESIS II PS Inserts Size 3-4, 9mm
71420818	GENESIS II PS Inserts Size 3-4, 11mm
71420820	GENESIS II PS Inserts Size 3-4, 13mm
71420822	GENESIS II PS Inserts Size 3-4, 15mm
71420824	GENESIS II PS Inserts Size 3-4, 18mm
71420826	GENESIS II PS Inserts Size 3-4, 21mm
71420828	GENESIS II PS Inserts Size 3-4, 25mm
71420830	GENESIS II PS Inserts Size 5-6, 9mm
71420832	GENESIS II PS Inserts Size 5-6, 11mm
71420834	GENESIS II PS Inserts Size 5-6, 13mm
71420836	GENESIS II PS Inserts Size 5-6, 15mm
71420838	GENESIS II PS Inserts Size 5-6, 18mm
71420840	GENESIS II PS Inserts Size 5-6, 21mm
71420842	GENESIS II PS Inserts Size 5-6, 25mm
71420844	GENESIS II PS Inserts Size 7-8, 9mm
71420846	GENESIS II PS Inserts Size 7-8, 11mm
71420848	GENESIS II PS Inserts Size 7-8, 13mm
71420850	GENESIS II PS Inserts Size 7-8, 15mm
71420852	GENESIS II PS Inserts Size 7-8, 18mm
71420854	GENESIS II PS Inserts Size 7-8, 21mm
71420856	GENESIS II PS Inserts Size 7-8, 25mm

GENESIS II PS High Flex Inserts

Cat. No.	Description
71421501	GENESIS II PS High Flex Inserts Size 1-2, 9mm
71421502	GENESIS II PS High Flex Inserts Size 1-2, 11mm
71421503	GENESIS II PS High Flex Inserts Size 1-2, 13mm
71421504	GENESIS II PS High Flex Inserts Size 1-2, 15mm
71421505	GENESIS II PS High Flex Inserts Size 1-2, 18mm
71421506	GENESIS II PS High Flex Inserts Size 1-2, 21mm
71421507	GENESIS II PS High Flex Inserts Size 1-2, 25mm
71421508	GENESIS II PS High Flex Inserts Size 3-4, 9mm
71421509	GENESIS II PS High Flex Inserts Size 3-4, 11mm
71421510	GENESIS II PS High Flex Inserts Size 3-4, 13mm
71421511	GENESIS II PS High Flex Inserts Size 3-4, 15mm
71421512	GENESIS II PS High Flex Inserts Size 3-4, 18mm
71421513	GENESIS II PS High Flex Inserts Size 3-4, 21mm
71421514	GENESIS II PS High Flex Inserts Size 3-4, 25mm
71421515	GENESIS II PS High Flex Inserts Size 5-6, 9mm
71421516	GENESIS II PS High Flex Inserts Size 5-6, 11mm
71421517	GENESIS II PS High Flex Inserts Size 5-6, 13mm
71421518	GENESIS II PS High Flex Inserts Size 5-6, 15mm
71421519	GENESIS II PS High Flex Inserts Size 5-6, 18mm
71421520	GENESIS II PS High Flex Inserts Size 5-6, 21mm
71421521	GENESIS II PS High Flex Inserts Size 5-6, 25mm
71421522	GENESIS II PS High Flex Inserts Size 7-8, 9mm
71421523	GENESIS II PS High Flex Inserts Size 7-8, 11mm
71421524	GENESIS II PS High Flex Inserts Size 7-8, 13mm
71421525	GENESIS II PS High Flex Inserts Size 7-8, 15mm
71421526	GENESIS II PS High Flex Inserts Size 7-8, 18mm
71421527	GENESIS II PS High Flex Inserts Size 7-8, 21mm
71421528	GENESIS II PS High Flex Inserts Size 7-8, 25mm

GENESIS® II CR Inserts

Cat. No.	Description
71420480	GENESIS II CR Articular Insert Size 1-2, 9mm
71420482	GENESIS II CR Articular Insert Size 1-2, 11mm
71420484	GENESIS II CR Articular Insert Size 1-2, 13mm
71420486	GENESIS II CR Articular Insert Size 1-2, 15mm
71420488	GENESIS II CR Articular Insert Size 1-2, 18mm
71420490	GENESIS II CR Articular Insert Size 3-4, 9mm
71420492	GENESIS II CR Articular Insert Size 3-4, 11mm
71420494	GENESIS II CR Articular Insert Size 3-4, 13mm
71420496	GENESIS II CR Articular Insert Size 3-4, 15mm
71420498	GENESIS II CR Articular Insert Size 3-4, 18mm
71420500	GENESIS II CR Articular Insert Size 5-6, 9mm
71420502	GENESIS II CR Articular Insert Size 5-6, 11mm
71420504	GENESIS II CR Articular Insert Size -6, 13mm
71420506	GENESIS II CR Articular Insert Size -6, 15mm
71420508	GENESIS II CR Articular Insert Size 5-6, 18mm
71420510	GENESIS II CR Articular Insert Size 7-8, 9mm
71420512	GENESIS II CR Articular Insert Size 7-8, 11mm
71420514	GENESIS II CR Articular Insert Size 7-8, 13mm
71420516	GENESIS II CR Articular Insert Size 7-8, 15mm
71420518	GENESIS II CR Articular Insert Size 7-8, 18mm

GENESIS II CR Deep Flex Inserts

Cat. No.	Description
71421531	GENESIS II CR Deep Flexion Insert Size 1-2, 9mm
71421532	GENESIS II CR Deep Flexion Insert Size 1-2, 11mm
71421533	GENESIS II CR Deep Flexion Insert Size 1-2, 13mm
71421534	GENESIS II CR Deep Flexion Insert Size 1-2, 15mm
71421535	GENESIS II CR Deep Flexion Insert Size 1-2, 18mm
71421536	GENESIS II CR Deep Flexion Insert Size 3-4, 9mm
71421537	GENESIS II CR Deep Flexion Insert Size 3-4, 11mm
71421538	GENESIS II CR Deep Flexion Insert Size 3-4, 13mm
71421539	GENESIS II CR Deep Flexion Insert Size 3-4, 15mm
71421541	GENESIS II CR Deep Flexion Insert Size 3-4, 18mm
71421542	GENESIS II CR Deep Flexion Insert Size 5-6, 9mm
71421543	GENESIS II CR Deep Flexion Insert Size 5-6, 11mm
71421544	GENESIS II CR Deep Flexion Insert Size -6, 13mm
71421545	GENESIS II CR Deep Flexion Insert Size -6, 15mm
71421546	GENESIS II CR Deep Flexion Insert Size 5-6, 18mm
71421547	GENESIS II CR Deep Flexion Insert Size 7-8, 9mm
71421548	GENESIS II CR Deep Flexion Insert Size 7-8, 11mm
71421549	GENESIS II CR Deep Flexion Insert Size 7-8, 13mm
71421551	GENESIS II CR Deep Flexion Insert Size 7-8, 15mm
71421552	GENESIS II CR Deep Flexion Insert Size 7-8, 18mm

GENESIS II CR Deep Dish Inserts

Cat. No.	Description
71420754	GENESIS II Dished Insert Size 1-2, 9mm
71420756	GENESIS II Dished Insert Size 1-2, 11mm
71420758	GENESIS II Dished Insert Size 1-2, 13mm
71420760	GENESIS II Dished Insert Size 1-2, 15mm
71420762	GENESIS II Dished Insert Size 1-2, 18mm
71420764	GENESIS II Dished Insert Size 1-2, 21mm
71420766	GENESIS II Dished Insert Size 3-4, 9mm
71420768	GENESIS II Dished Insert Size 3-4, 11mm
71420770	GENESIS II Dished Insert Size 3-4, 13mm
71420772	GENESIS II Dished Insert Size 3-4, 15mm
71420774	GENESIS II Dished Insert Size 3-4, 18mm
71420776	GENESIS II Dished Insert Size 3-4, 21mm
71420778	GENESIS II Dished Insert Size 5-6, 9mm
71420780	GENESIS II Dished Insert Size 5-6, 11mm
71420782	GENESIS II Dished Insert Size 5-6, 13mm
71420784	GENESIS II Dished Insert Size 5-6, 15mm
71420786	GENESIS II Dished Insert Size 5-6, 18mm
71420788	GENESIS II Dished Insert Size 5-6, 21mm
71420790	GENESIS II Dished Insert Size 7-8, 9mm
71420792	GENESIS II Dished Insert Size 7-8, 11mm
71420794	GENESIS II Dished Insert Size 7-8, 13mm
71420796	GENESIS II Dished Insert Size 7-8, 15mm
71420798	GENESIS II Dished Insert Size 7-8, 18mm
71420800	GENESIS II Dished Insert Size 7-8, 21mm

GENESIS[®] II PS All-Poly Tibial Baseplate

Cat. No.	Description
71420352	GENESIS II PS All-Poly Tibia Size 1, 9mm, Left
71420354	GENESIS II PS All-Poly Tibia Size 1, 11mm, Left
71420356	GENESIS II PS All-Poly Tibia Size 1, 13mm, Left
71420358	GENESIS II PS All-Poly Tibia Size 1, 15mm, Left
71420360	GENESIS II PS All-Poly Tibia Size 2, 9mm, Left
71420362	GENESIS II PS All-Poly Tibia Size 2, 11mm, Left
71420364	GENESIS II PS All-Poly Tibia Size 2, 13mm, Left
71420366	GENESIS II PS All-Poly Tibia Size 2, 15mm, Left
71420368	GENESIS II PS All-Poly Tibia Size 3, 9mm, Left
71420370	GENESIS II PS All-Poly Tibia Size 3, 11mm, Left
71420372	GENESIS II PS All-Poly Tibia Size 3, 13mm, Left
71420374	GENESIS II PS All-Poly Tibia Size 3, 15mm, Left
71420376	GENESIS II PS All-Poly Tibia Size 4, 9mm, Left
71420378	GENESIS II PS All-Poly Tibia Size 4, 11mm, Left
71420380	GENESIS II PS All-Poly Tibia Size 4, 13mm, Left
71420382	GENESIS II PS All-Poly Tibia Size 4, 15mm, Left
71420384	GENESIS II PS All-Poly Tibia Size 5, 9mm Left
71420386	GENESIS II PS All-Poly Tibia Size 5, 11mm, Left
71420388	GENESIS II PS All-Poly Tibia Size 5, 13mm, Left
71420390	GENESIS II PS All-Poly Tibia Size 5, 15mm, Left
71420392	GENESIS II PS All-Poly Tibia Size 6, 9mm, Left
71420394	GENESIS II PS All-Poly Tibia Size 6, 11mm, Left
71420396	GENESIS II PS All-Poly Tibia Size 6, 13mm, Left
71420398	GENESIS II PS All-Poly Tibia Size 6, 15mm, Left
71420400	GENESIS II PS All-Poly Tibia Size 7, 9mm, Left
71420402	GENESIS II PS All-Poly Tibia Size 7, 11mm, Left
71420404	GENESIS II PS All-Poly Tibia Size 7, 13mm, Left
71420406	GENESIS II PS All-Poly Tibia Size 7, 15mm, Left
71420408	GENESIS II PS All-Poly Tibia Size 8, 9mm, Left
71420410	GENESIS II PS All-Poly Tibia Size 8, 11mm, Left
71420412	GENESIS II PS All-Poly Tibia Size 8, 13mm, Left
71420414	GENESIS II PS All-Poly Tibia Size 8, 15mm, Left
71420416	GENESIS II PS All-Poly Tibia Size 1, 9mm, Right
71420418	GENESIS II PS All-Poly Tibia Size 1, 11mm, Right
71420420	GENESIS II PS All-Poly Tibia Size 1, 13mm, Right
71420422	GENESIS II PS All-Poly Tibia Size 1, 15mm, Right
71420424	GENESIS II PS All-Poly Tibia Size 2, 9mm, Right
71420426	GENESIS II PS All-Poly Tibia Size 2, 11mm, Right
71420428	GENESIS II PS All-Poly Tibia Size 2, 13mm, Right
71420430	GENESIS II PS All-Poly Tibia Size 2, 15mm, Right
71420432	GENESIS II PS All-Poly Tibia Size 3, 9mm, Right
71420434	GENESIS II PS All-Poly Tibia Size 3, 11mm, Right
71420436	GENESIS II PS All-Poly Tibia Size 3, 13mm, Right
71420438	GENESIS II PS All-Poly Tibia Size 3, 15mm, Right

Cat. No.	Description
71420440	GENESIS II PS All-Poly Tibia Size 4, 9mm, Right
71420442	GENESIS II PS All-Poly Tibia Size 4, 11mm, Right
71420444	GENESIS II PS All-Poly Tibia Size 4, 13mm, Right
71420446	GENESIS II PS All-Poly Tibia Size 4, 15mm, Right
71420448	GENESIS II PS All-Poly Tibia Size 5, 9mm, Right
71420450	GENESIS II PS All-Poly Tibia Size 5, 11mm, Right
71420452	GENESIS II PS All-Poly Tibia Size 5, 13mm, Right
71420454	GENESIS II PS All-Poly Tibia Size 5, 15mm, Right
71420456	GENESIS II PS All-Poly Tibia Size 6, 9mm, Right
71420458	GENESIS II PS All-Poly Tibia Size 6, 11mm, Right
71420460	GENESIS II PS All-Poly Tibia Size 6, 13mm, Right
71420462	GENESIS II PS All-Poly Tibia Size 6, 15mm, Right
71420464	GENESIS II PS All-Poly Tibia Size 7, 9mm, Right
71420466	GENESIS II PS All-Poly Tibia Size 7, 11mm, Right
71420468	GENESIS II PS All-Poly Tibia Size 7, 13mm, Right
71420470	GENESIS II PS All-Poly Tibia Size 7, 15mm, Right
71420472	GENESIS II PS All-Poly Tibia Size 8, 9mm, Right
71420474	GENESIS II PS All-Poly Tibia Size 8, 11mm, Right
71420476	GENESIS II PS All-Poly Tibia Size 8, 13mm, Right
71420478	GENESIS II PS All-Poly Tibia Size 8, 15mm, Right

GENESIS II Tibial Baseplate (Cemented)

Cat. No.	Description
71420160	GENESIS II Cemented Tibial Baseplate Size 1 Left
71420162	GENESIS II Cemented Tibial Baseplate Size 2 Left
71420164	GENESIS II Cemented Tibial Baseplate Size 3 Left
71420166	GENESIS II Cemented Tibial Baseplate Size 4 Left
71420168	GENESIS II Cemented Tibial Baseplate Size 5 Left
71420170	GENESIS II Cemented Tibial Baseplate Size 6 Left
71420172	GENESIS II Cemented Tibial Baseplate Size 7 Left
71420174	GENESIS II Cemented Tibial Baseplate Size 8 Left
71931923	GENESIS II Cemented Tibial Baseplate Size 9 Left
71420176	GENESIS II Cemented Tibial Baseplate Size 1 Right
71420180	GENESIS II Cemented Tibial Baseplate Size 2 Right
71420182	GENESIS II Cemented Tibial Baseplate Size 3 Right
71420184	GENESIS II Cemented Tibial Baseplate Size 4 Right
71420186	GENESIS II Cemented Tibial Baseplate Size 5 Right
71420188	GENESIS II Cemented Tibial Baseplate Size 6 Right
71420190	GENESIS II Cemented Tibial Baseplate Size 7 Right
71420191	GENESIS II Cemented Tibial Baseplate Size 8 Right
71931716	GENESIS II Cemented Tibial Baseplate Size 9 Right

GENESIS® II CR All-Poly Tibial Baseplate

Cat. No.	Description
71420224	GENESIS II CR All-Poly Tibia Size 1, 9mm, Left
71420226	GENESIS II CR All-Poly Tibia Size 1, 11mm, Left
71420228	GENESIS II CR All-Poly Tibia Size 1, 13mm, Left
71420230	GENESIS II CR All-Poly Tibia Size 1, 15mm, Left
71420232	GENESIS II CR All-Poly Tibia Size 2, 9mm, Left
71420234	GENESIS II CR All-Poly Tibia Size 2, 11mm, Left
71420236	GENESIS II CR All-Poly Tibia Size 2, 13mm, Left
71420238	GENESIS II CR All-Poly Tibia Size 2, 15mm, Left
71420240	GENESIS II CR All-Poly Tibia Size 3, 9mm, Left
71420242	GENESIS II CR All-Poly Tibia Size 3 11mm, Left
71420244	GENESIS II CR All-Poly Tibia Size 3, 13mm, Left
71420246	GENESIS II CR All-Poly Tibia Size 3, 15mm, Left
71420248	GENESIS II CR All-Poly Tibia Size 4, 9mm, Left
71420250	GENESIS II CR All-Poly Tibia Size 4, 11mm, Left
71420252	GENESIS II CR All-Poly Tibia Size 4, 13mm, Left
71420254	GENESIS II CR All-Poly Tibia Size 4, 15mm, Left
71420256	GENESIS II CR All-Poly Tibia Size 5, 9mm, Left
71420258	GENESIS II CR All-Poly Tibia Size 5, 11mm, Left
71420260	GENESIS II CR All-Poly Tibia Size 5, 13mm, Left
71420262	GENESIS II CR All-Poly Tibia Size 5, 15mm, Left
71420264	GENESIS II CR All-Poly Tibia Size 6, 9mm, Left
71420266	GENESIS II CR All-Poly Tibia Size 6, 11mm, Left
71420268	GENESIS II CR All-Poly Tibia Size 6, 13mm, Left
71420270	GENESIS II CR All-Poly Tibia Size 6, 15mm, Left
71420272	GENESIS II CR All-Poly Tibia Size 7, 9mm, Left
71420274	GENESIS II CR All-Poly Tibia Size 7, 11mm, Left
71420276	GENESIS II CR All-Poly Tibia Size 7, 13mm, Left
71420278	GENESIS II CR All-Poly Tibia Size 7, 15mm, Left
71420280	GENESIS II CR All-Poly Tibia Size 8, 9mm, Left
71420282	GENESIS II CR All-Poly Tibia Size 8, 11mm, Left
71420284	GENESIS II CR All-Poly Tibia Size 8, 13mm, Left
71420286	GENESIS II CR All-Poly Tibia Size 8, 15mm, Left
71420288	GENESIS II CR All-Poly Tibia Size 1, 9mm, Right
71420290	GENESIS II CR All-Poly Tibia Size 1, 11mm, Right
71420292	GENESIS II CR All-Poly Tibia Size 1, 13mm, Right
71420294	GENESIS II CR All-Poly Tibia Size 1, 15mm, Right
71420296	GENESIS II CR All-Poly Tibia Size 2, 9mm, Right
71420298	GENESIS II CR All-Poly Tibia Size 2, 11mm, Right
71420300	GENESIS II CR All-Poly Tibia Size 2, 13mm, Right
71420302	GENESIS II CR All-Poly Tibia Size 2, 15mm, Right
71420304	GENESIS II CR All-Poly Tibia Size 3, 9mm, Right
71420306	GENESIS II CR All-Poly Tibia Size 3, 11mm, Right
71420308	GENESIS II CR All-Poly Tibia Size 3, 13mm, Right
71420310	GENESIS II CR All-Poly Tibia Size 3, 15mm, Right

Cat. No.	Description
71420312	GENESIS II CR All-Poly Tibia Size 4, 9mm, Right
71420314	GENESIS II CR All-Poly Tibia Size 4, 11mm, Right
71420316	GENESIS II CR All-Poly Tibia Size 4, 13mm, Right
71420318	GENESIS II CR All-Poly Tibia Size 4, 15mm, Right
71420320	GENESIS II CR All-Poly Tibia Size 5, 9mm, Right
71420322	GENESIS II CR All-Poly Tibia Size 5, 11mm, Right
71420324	GENESIS II CR All-Poly Tibia Size 5, 13mm, Right
71420326	GENESIS II CR All-Poly Tibia Size 5, 15mm, Right
71420328	GENESIS II CR All-Poly Tibia Size 6, 9mm, Right
71420330	GENESIS II CR All-Poly Tibia Size 6, 11mm, Right
71420332	GENESIS II CR All-Poly Tibia Size 6, 13mm, Right
71420334	GENESIS II CR All-Poly Tibia Size 6, 15mm, Right
71420336	GENESIS II CR All-Poly Tibia Size 7, 9mm, Right
71420338	GENESIS II CR All-Poly Tibia Size 7, 11mm, Right
71420340	GENESIS II CR All-Poly Tibia Size 7, 13mm, Right
71420342	GENESIS II CR All-Poly Tibia Size 7, 15mm, Right
71420344	GENESIS II CR All-Poly Tibia Size 8, 9mm, Right
71420346	GENESIS II CR All-Poly Tibia Size 8, 11mm, Right
71420348	GENESIS II CR All-Poly Tibia Size 8, 13mm, Right
71420350	GENESIS II CR All-Poly Tibia Size 8, 15mm, Right

GENESIS® II Round Resurfacing Patella

Cat. No.	Description
71420580	GENESIS II Round Resurfacing Patella 26mm
71420574	GENESIS II Round Resurfacing Patella 29mm
71420576	GENESIS II Round Resurfacing Patella 32mm
71420578	GENESIS II Round Resurfacing Patella 35mm
71926225	GENESIS II Round Resurfacing Patella 38mm
71926226	GENESIS II Round Resurfacing Patella 41mm

GENESIS II Oval Resurfacing Patella

Cat. No.	Description
71421029	GENESIS II Oval Resurfacing Patella 29mm
71421032	GENESIS II Oval Resurfacing Patella 32mm
71421035	GENESIS II Oval Resurfacing Patella 35mm
71421038	GENESIS II Oval Resurfacing Patella 38mm
71421041	GENESIS II Oval Resurfacing Patella 41mm

GENESIS II Biconvex Resurfacing Patella

Cat. No.	Description
71420566	GENESIS II Biconvex Resurfacing Patella 23mm
71420568	GENESIS II Biconvex Resurfacing Patella 26mm
71420570	GENESIS II Biconvex Resurfacing Patella 29mm
71420572	GENESIS II Biconvex Resurfacing Patella 32mm

GENESIS II PROFIX® Stem

Cat. No.	Description
76547001	GENESIS II PROFIX Mobile Bearing Stem 14mm x 33mm

GENESIS II Femoral Flex-Lok Peg

Cat. No.	Description
71420063	GENESIS II Femoral Flex-Lok Peg

GENESIS II Primary Femoral Lug

Cat. No.	Description
71420999	GENESIS II Femoral Lug

GENESIS II Tibial Offset Couplers

Cat. No.	Description
71422002	GENESIS II Tibial Offset Coupler with Sleeve 2mm
71422004	GENESIS II Tibial Offset Coupler with Sleeve 4mm

GENESIS II Long Stem

Cat. No.	Description
71420628	GENESIS II Long Stem 10mm x 100mm
71420630	GENESIS II Long Stem 12mm x 100mm
71420632	GENESIS II Long Stem 14mm x 100mm
71420634	GENESIS II Long Stem 16mm x 100mm
71420636	GENESIS II Long Stem 18mm x 100mm
71420638	GENESIS II Long Stem 20mm x 100mm
71420640	GENESIS II Long Stem 22mm x 100mm
71420642	GENESIS II Long Stem 24mm x 100mm
71420647	GENESIS II Long Stem 10mm x 150mm
71420648	GENESIS II Long Stem 14mm x 150mm
71420649	GENESIS II Long Stem 12mm x 150mm
71420650	GENESIS II Long Stem 16mm x 150mm
71421310	Revision Press-fit Stem 10mm x 100mm
71421312	Revision Press-fit Stem 12mm x 100mm
71421314	Revision Press-fit Stem 14mm x 100mm
71421316	Revision Press-fit Stem 16mm x 100mm
71421318	Revision Press-fit Stem 10mm x 150mm
71421320	Revision Press-fit Stem 12mm x 150mm
71421322	Revision Press-fit Stem 14mm x 150mm
71421324	Revision Press-fit Stem 16mm x 150mm

GENESIS[®] II PS Implant Construct

Recommended Combination

Femoral	Insert	Tibial Baseplate	Patella
GENESIS II PS Non-Porous Femoral (CoCr)	GENESIS II PS Insert	GENESIS II Baseplate (Cemented)	GENESIS II Round Resurfacing Patella
GENESIS II PS Non-Porous Femoral (OXINIUM [®])	GENESIS II PS High Flex Insert	GENESIS II PS All-Poly Tibial Baseplate	GENESIS II Oval Resurfacing Patella
			GENESIS II Biconvex Resurfacing Patella

Compatibility with Additional Components

	GENESIS II PS Non-Porous Femoral (CoCr)	GENESIS II Tibial Baseplate (Cemented)
Optional Component	GENESIS II Femoral Flex-Lok Peg GENESIS II Primary Femoral Lug	GENESIS II Hemi-Stepped Tibial Wedge GENESIS II/PROFIX[®] Metaphyseal Tibial Stem (76547001 Only) GENESIS II Long Stem GENESIS II Tibial Offset Coupler

GENESIS[®] II Compatibility

GENESIS II PS Non-Porous Femoral (CoCr) GENESIS II PS Non-Porous Femoral (OXINIUM [®])		GENESIS II PS Insert GENESIS II PS High Flex Insert GENESIS II PS Constrained Insert*	
Patella	Inserts/Baseplate	Tibial Baseplate	Femoral
GENESIS II Round Resurfacing Patella	GENESIS II PS Insert	GENESIS II Tibial Baseplate (Cemented)	GENESIS II PS Non-Porous Femoral (CoCr)
GENESIS II Oval Resurfacing Patella	GENESIS II PS High Flex Insert	LEGION Revision Tibial Baseplate (Cemented)	GENESIS II PS Non-Porous Femoral (OXINIUM)
GENESIS II Biconvex Resurfacing Patella	GENESIS II PS All Poly Tibial Baseplate	LEGION Porous HA Tibial Baseplate with Holes	LEGION PS Non-Porous Femoral (CoCr)
	LEGION [®] PS High Flex Insert	LEGION Porous HA Tibial Baseplate without Holes	LEGION PS Narrow Non-Porous Femoral (CoCr)
	ANTHEM [®] PS High Flex Insert	ANTHEM Tibial Baseplate**	LEGION PS Femoral (OXINIUM)
			LEGION PS Narrow Femoral (OXINIUM)
			LEGION RK Constrained Femoral (CoCr)
			LEGION RK Constrained Femoral (OXINIUM)
			ANTHEM PS Femoral (CoCr)**
			ANTHEM PS Narrow Femoral (CoCr)**

* Components must be used with the LEGION Femorals (CoCr and OXINIUM) and tibial stems.

** Component used with GENESIS II PS High Flex Insert.

Compatibility Table

GENESIS [®] II Component	Compatible Component	Size
GENESIS II PS Non-Porous Femoral (CoCr)	GENESIS II PS Insert	1-8, 9-25 mm
	GENESIS II PS High Flex Insert	1-8, 9-25 mm
	LEGION [®] PS High Flex Insert	1-8, 9-21 mm
	ANTHEM [®] PS High Flex Insert	1-8, 9-18 mm
GENESIS II PS Non-Porous Femoral (OXINIUM [®])	GENESIS II PS All-Poly Tibial Baseplate	1-8 RT/LT, 9-15mm
	GENESIS II Resurfacing Patella	26-41 mm
	GENESIS II Oval Resurfacing Patella	29-41 mm
	GENESIS II Biconvex Resurfacing Patella	23-32 mm
GENESIS II PS Non-Porous Femoral (CoCr)	GENESIS II Femoral Flex-Lok Peg	N/A
	GENESIS II Primary Femoral Lug	N/A

GENESIS II PS Insert	GENESIS II PS Non-Porous Femoral (CoCr)	1-8 RT/LT
	GENESIS II PS Non-Porous Femoral (OXINIUM)	2-9 RT/LT
	LEGION RK Constrained Femoral (CoCr)	2-8 RT/LT
	LEGION RK Constrained Femoral (OXINIUM)	2-8 RT/LT
	LEGION PS Non-Porous Femoral (CoCr)	2-8 RT/LT
GENESIS II PS High Flex Insert	LEGION PS Narrow Non-Porous Femoral (CoCr)	3-6 RT/LT
	LEGION PS Femoral (OXINIUM)	2-8 RT/LT
	LEGION PS Narrow Femoral (OXINIUM)	3-6 RT/LT
	ANTHEM PS Femoral (CoCr)**	3-8 RT/LT
GENESIS II PS Constrained Insert (PE)*	ANTHEM PS Narrow Femoral (CoCr)**	1-6 RT/LT
	GENESIS II Tibial Baseplate (Cemented)	1-9 RT/LT
	LEGION Revision Tibial Baseplate (Cemented)	1-8 RT/LT
	LEGION Porous HA Tibial Baseplate with Holes	2-8 RT/LT
	LEGION Porous HA Tibial Baseplate without Holes	2-8 RT/LT
	ANTHEM Tibial Baseplate**	1-8 RT/LT

* Components must be used with the LEGION Femorals (CoCr and OXINIUM) and tibial stems.

** Component used with GENESIS II PS High Flex Insert.

GENESIS [®] II Component	Compatible Component	Size
GENESIS II Tibial Baseplate (Cemented)	GENESIS II PS Insert	1-8, 9-25 mm
	GENESIS II PS High Flex Insert	1-8, 9-25 mm
	GENESIS II PS Constrained Insert*	1-8, 11-30 mm
	LEGION [®] PS Constrained Insert*	1-8, 9 mm
	LEGION PS High Flex Insert	1-8, 9-21 mm
	ANTHEM [®] PS High Flex Insert	1-8, 9-18 mm
	GENESIS II Hemi-Stepped Tibial Wedge	1-8, 10 mm and 15 mm
	GENESIS II/PROFIX [®] Metaphyseal Tibial Stem	14 x 33 mm
	GENESIS II Long Stem	10-24, 100 mm and 150 mm
	GENESIS II Tibial Offset Coupler	2 mm and 4 mm
GENESIS II PS All-Poly Tibial Baseplate	GENESIS II PS Non-Porous Femoral (CoCr)	1-8 RT/LT
	GENESIS II PS Non-Porous Femoral (OXINIUM)	2-9 RT/LT
	LEGION PS Non-Porous Femoral (CoCr)	2-8RT/LT
	LEGION PS Narrow Non-Porous Femoral (CoCr)	3-6RT/LT
	LEGION PS Femoral (OXINIUM)	2-8RT/LT
	LEGION PS Narrow Femoral (OXINIUM)	3-6RT/LT

* Components must be used with the LEGION Femorals (CoCr and OXINIUM) and tibial stems.

GENESIS[®] II CR Implant Construct

Recommended Combination

Femoral	Insert	Tibial Baseplate	Patella
GENESIS II CR Non-Porous Femoral (CoCr) GENESIS II CR Non-Porous Femoral (OXINIUM[®])	GENESIS II CR Insert GENESIS II CR Deep Dish Insert GENESIS II CR High Flex (Deep Flex) Insert	GENESIS II Baseplate (Cemented) GENESIS II CR All-Poly Tibial Baseplate	GENESIS II Round Resurfacing Patella GENESIS II Oval Resurfacing Patella GENESIS II Biconvex Resurfacing Patella

Compatibility with Additional Components

	GENESIS II CR Non-Porous Femoral (CoCr)	GENESIS II Tibial Baseplate (Cemented)
Optional Component	GENESIS II Femoral Flex-Lok Peg	GENESIS II Hemi-Stepped Tibial Wedge GENESIS II/PROFIX[®] Metaphyseal Tibial Stem (76547001 Only) GENESIS II Long Stem GENESIS II Tibial Offset Coupler

GENESIS II° Compatibility

GENESIS II CR Non-Porous Femoral (CoCr) GENESIS II CR Femoral (OXINIUM°)		GENESIS II CR Insert GENESIS II CR Deep Dish Insert GENESIS II CR High Flex (Deep Flex) Insert	
Patella	Inserts/Baseplate	Tibial Baseplate	Femoral
GENESIS II Round Resurfacing Patella	GENESIS II CR Insert	GENESIS II Tibial Baseplate (Cemented)	GENESIS II CR Non-Porous Femoral (CoCr)
GENESIS II Oval Resurfacing Patella	GENESIS II CR Deep Dish Insert	LEGION Revision Tibial Baseplate (Cemented)	GENESIS II CR Femoral (OXINIUM°)
GENESIS II Biconvex Resurfacing Patella	GENESIS II CR All Poly Tibial Baseplate	LEGION Porous HA Tibial Baseplate with Holes	LEGION CR Non-Porous Femoral (CoCr)
	GENESIS II CR High Flex (Deep Flex) Insert	LEGION Porous HA Tibial Baseplate without Holes	LEGION CR Narrow Non-Porous Femoral (CoCr)
	LEGION° CR Deep Dish Insert	LEGION Porous HA Tibial Baseplate without Holes	LEGION CR Porous Femoral (CoCr)
	LEGION CR High Flex Insert	ANTHEM Tibial Baseplate†	LEGION CR Porous + HA Femoral (CoCr)
	ANTHEM° CR High Flex Insert		LEGION CR Femoral (OXINIUM°)
			LEGION CR Narrow Femoral (OXINIUM°)
			ANTHEM CR Femoral (CoCr)†
			ANTHEM CR Narrow Femoral (CoCr)†

† Component used with GENESIS II CR High Flex (Deep Flex) Insert.

Compatibility Table

GENESIS [®] II Component	Compatible Component	Size
GENESIS II CR Non-Porous Femoral (CoCr) GENESIS II CR Non-Porous Femoral (OXINIUM [®])	GENESIS II CR Insert	1-8, 9-18 mm
	GENESIS II CR High Flex (Deep Flex) Insert	1-8, 9-18 mm
	GENESIS II CR Deep Dish Insert	1-8, 9-21 mm
	LEGION [®] CR Deep Dish Insert	1-8, 9-21 mm
	LEGION CR High Flex Insert	1-8, 9-18 mm
	ANTHEM [®] CR High Flex Insert	1-8, 9-18 mm
	GENESIS II CR All-Poly Tibial Baseplate	1-8 RT/LT, 9-15mm
	GENESIS II Resurfacing Patella	26-41 mm
	GENESIS II Oval Resurfacing Patella	29-41 mm
	GENESIS II Biconvex Resurfacing Patella	23-32 mm
GENESIS II CR Non-Porous Femoral (CoCr)	GENESIS II Femoral Flex-Lok Peg	N/A

GENESIS II CR Insert GENESIS II CR High Flex Insert GENESIS II CR Deep Dish Insert	GENESIS II CR Non-Porous Femoral (CoCr)	1-8 RT/LT
	GENESIS II CR Femoral (OXINIUM [®])	2-9 RT/LT
	LEGION CR Non-Porous Femoral (CoCr)	2-8 RT/LT
	LEGION CR Narrow Non-Porous Femoral (CoCr)	3-6 RT/LT
	LEGION CR Porous Femoral (CoCr)	2-8 RT/LT
	LEGION CR Porous + HA Femoral (CoCr)	2-8 RT/LT
	LEGION CR Femoral (OXINIUM)	2-8 RT/LT
	LEGION CR Narrow Femoral (OXINIUM)	3-6 RT/LT
	ANTHEM CR Femoral (CoCr) [†]	3-8 RT/LT
	ANTHEM CR Narrow Femoral (CoCr) [†]	1-6 RT/LT
	GENESIS II Tibial Baseplate (Cemented)	1-9 RT/LT
	LEGION Revision Tibial Baseplate (Cemented)	1-8 RT/LT
	LEGION Porous HA Tibial Baseplate with Holes	2-8 RT/LT
	LEGION Porous HA Tibial Baseplate without Holes	2-8 RT/LT
	ANTHEM Tibial Baseplate [†]	1-8 RT/LT

[†] Component used with GENESIS II CR High Flex (Deep Flex) Insert.

GENESIS [®] II Component	Compatible Component	Size
GENESIS II Tibial Baseplate (Cemented)	GENESIS II CR Insert	1-8, 9-18 mm
	GENESIS II CR High Flex (Deep Flex) Insert	1-8, 9-18 mm
	GENESIS II CR Deep Dish Insert	1-8, 9-21 mm
	LEGION [®] CR Deep Dish Insert	1-8, 9-21 mm
	LEGION CR High Flex Insert	1-8, 9-18 mm
	ANTHEM [®] CR High Flex Insert	1-8, 9-18 mm
	GENESIS II Hemi-Stepped Tibial Wedge	1-8, 10 mm and 15 mm
	GENESIS II/PROFIX [®] Mobile Bearing Tibial Stem (76547001 Only)	14 x 33 mm
	GENESIS II Long Stem	10-24, 100 and 150 mm
	GENESIS II Tibial Offset Coupler	2 mm and 4 mm
GENESIS II CR All-Poly Tibial Baseplate	GENESIS II CR Non-Porous Femoral (CoCr)	1-8 RT/LT
	GENESIS II CR Femoral (OXINIUM [®])	2-9 RT/LT
	LEGION CR Non-Porous Femoral (CoCr)	2-8 RT/LT
	LEGION CR Narrow Non-Porous Femoral (CoCr)	3-6 RT/LT
	LEGION CR Porous Femoral (CoCr)	2-8 RT/LT
	LEGION CR Porous + HA Femoral (CoCr)	2-8 RT/LT
	LEGION CR Femoral (OXINIUM)	2-8 RT/LT
	LEGION CR Narrow Femoral (OXINIUM)	3-6 RT/LT

[illegible]



Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith & Nephew representative or distributor if you have questions about the availability of Smith & Nephew products in your area.

Smith & Nephew, Inc.

1450 Brooks Road
Memphis, TN 38116
USA

www.smith-nephew.com

Telephone: 1-901-396-2121
Information: 1-800-821-5700
Orders/Inquiries:
1-800-238-7538

®Trademark of Smith & Nephew.
All trademarks acknowledged.

©2022 Smith & Nephew, Inc.
All rights reserved.
29596 V2 71281337 REVA 10/22