

# IV3000 Application

IV3000 Frame Delivery Dressing 6cm x 7cm

Suitable for dressing peripheral catheters


## FRAME DELIVERY

### Step 1




Clean and thoroughly dry the application area in accordance with clinical protocols. The sterile strips supplied on the backing paper can be used on the catheter wings or body for improved security, to secure tubing or to stabilise catheter lumens.

### Step 2



Remove the central window marked 1 at scored area

### Step 3



Remove the printed backing paper 2 marked from the dressing


### Step 4



Holding the frame, centre the dressing over the catheter insertion site.

Note you have full visibility of the insertion site for a more accurate application.

### Step 5



Place the dressing on the catheter hub and adhere the dressing to the patient's skin. Mould the film around the catheter. Remove the paper frame using the tab.

### Step 6



Smooth down the edges of the dressing, ensuring the dressing is secure and the patient is comfortable.

Finally after recording relevant information on the date label, place it on the patient's skin or the dressing (without obscuring the view of the catheter site).




## Removal

Take a corner of the film and stretch it parallel to the patient's skin while stabilising the patient's skin and catheter with the other hand.

Release the film and repeat the process.

## Ordering codes



Product code	Dressing size	Quantity per carton	Recommended indication
59410082	6cm x 7cm	100	Peripheral
59410882	10cm x 12cm	50	Central

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

Smith & Nephew Pty Ltd  
Australia  
T +61 2 9857 3999 F +61 2 9857 3900  
www.smith-nephew.com/australia  
Customer Service  
T 13 13 60 F 1800 671 000

Smith & Nephew Ltd  
New Zealand  
T +64 9 820 2840 F +64 9 820 2841  
www.smith-nephew.com/new-zealand  
Customer Service  
T 0800 657 799 (Surgical) T 0800 807 663 (Wound)