SmithNephew

Use of NAVIO^o Surgical System for unicompartmental knee arthroplasty (UKA) resulted in a significantly reduced total revision rate compared with conventional UKA at short-term follow up

Plus points

Significantly lower revision rate with NAVIO

No specific complications

Significantly increased

Overview

- Single centre, retrospective study performed between January 2013 and December 2018 comparing the use of NAVIO UKA and conventional UKA
 - 200 NAVIO UKAs (mean age, 66.7 years)
 - 191 conventional UKAs (mean age, 67.1 years)
 - Mean follow-up was 22.5 months for NAVIO UKA and 30.2 months for conventional UKA (p<0.001)

Results

- NAVIO UKA had a significantly reduced total revision rate compared to conventional UKA at last follow-up (4 vs 11%, p=0.014; Figure)
 - Revision due to malalignment was significantly lower with NAVIO UKA compared to conventional UKA (0 vs 5.2%, p=0.002)
- No specific complications associated with use of NAVIO, in particular, no issues due to the use of navigation pins
- Total reoperation rate (without implant removal) was reduced with NAVIO UKA compared to conventional UKA at last follow-up (6.5 vs 9.4%)
- At the last follow-up, KSS functional score was significantly higher with NAVIO UKA compared to conventional UKA (92.8 vs 88.4, p=0.01)
- No significant difference in duration of surgery (NAVIO UKA, 81 min; conventional UKA, 76 min)

- Data were collected preoperatively and at 2, 6, 12 months and at last follow-up
 - Revisions, intraoperative and postoperative complications, functional and radiological results were collected



Figure. Total revision rate (%) of NAVIO UKA and conventional UKA at last follow-up

www.smith-nephew.com/education

Conclusions

Citation

*Mergenthaler G, Batailler C, Lording T, et al. Is robotic-assisted unicompartmental knee arthroplasty a safe procedure? A case control study. ESSKA. 2020. [Epub ahead of print]

Available at: European Society of Sports Traumatology, Knee Surgery, Arthroscopy

Smith & Nephew, Inc., 1450 Brooks Road, Memphis, TN 38116, USA. 27577-en V1 11/20. Published November 2020. ©2020 Smith+Nephew. All rights reserved. °Trademark of Smith+Nephew. All Trademarks acknowledged.