

+ Evidence in focus

Compendium of peer-reviewed clinical evidence

JOURNEY^{II}
Total Knee Arthroplasty

June 2022

Smith+Nephew



JOURNEY[◇] II TKA Compendium of Evidence

Key evidence

All studies

Key outcomes

Function

Recovery

Patient
satisfaction

Survivorship

PDF navigation tips

Interactivity has been included throughout the compendium to aid navigation



will take you to this slide



will take you to a summary of each outcome



will take you to a grid summary of the 51 peer-reviewed studies highlighted in this compendium. Key evidence is represented in **blue** and supporting evidence in **turquoise**.

Full summary

will take you to the Evidence in Focus full summary of the relevant study


Published paper

will take you to the official online journal website on which the original paper is published

Arrows at the bottom of each page will take you to previous or subsequent pages, as indicated

What are the issues with conventional TKA?

A “Forgotten Joint” is the desired goal for both patients and surgeons following TKA. However, the majority of patients do not obtain a normal feeling knee post-TKA with high numbers of patients experiencing reduced function and dissatisfaction:



80% patients feel that their joint is artificial¹



>50% patients report a degree of functional limitation^{1,2}




20% patients are dissatisfied³

Studies show that following conventional TKA, patients demonstrate abnormalities in their gait and functional kinematics of the knee:⁴⁻⁸

Velocity



Stride length



Flexion during walking



Quadricep muscle strength following conventional TKA fails to recover to the same level as those of healthy subjects, resulting in reduced physical functioning of the knee^{9,10}




Conventional TKAs have been successful in their aim to relieve pain and provide long term survivorship.^{11,12} However, conventional TKAs do not live up to all patients functional expectations,² leaving a high percentage feeling dissatisfied.³

Kinematics and patient satisfaction — how do they relate?

Fluoroscopic analysis of weight-bearing motion has led to a better understanding of TKA kinematics. During flexion, TKA kinematic patterns can vary considerably from the normal knee.¹³

Paradoxical motion (Figure), where there is a sudden anterior translation of the femur relative to the tibia, has been observed in many studies and is related to reduced knee flexion and quadriceps efficiency.¹³



An implant designed to better replicate the shape and position of the normal knee could improve knee function and subsequently patient satisfaction.

For the first time, a relationship between kinematic patterns of TKA and patient satisfaction has been demonstrated.^{14,15}

Kinematic patterns for one weight-bearing and one non-weight-bearing motion were analysed for JOURNEY II BCS and two posterior-stabilised TKAs. Whilst no differences were observed during non-weight bearing activity, during weight-bearing motion, poor patient-reported outcomes were associated with:¹⁴

Pronounced paradoxical anterior motion (medial side)

Less stable medial compartment in midflexion

Less posterior translation in deep flexion (lateral side)

“Reproduction of optimal kinematic patterns during TKA could be instrumental in improving patient satisfaction.”¹⁴



Why choose JOURNEY II TKA for your patients?

Shape



Position

Motion

Unlike conventional TKA design, the shape and position of the JOURNEY II TKA have been designed to replicate the normal knee


Shape¹⁶⁻¹⁹

Replication of anatomic asymmetric femoral and tibial profiles



Position^{16,19,20}


Restoration of native anterior/posterior (A/P) starting position and the anatomic 3° varus joint line




Motion^{16-18,20-23}

Replication of native femoral rollback and axial rotation


Normal knee



JOURNEY II TKA knee



Why choose JOURNEY II TKA for your patients?



*Compared to non-JOURNEY II knees. †As observed in JOURNEY II BCS and JOURNEY II CR.



Key outcome: Improved function

43 studies

Postoperative ROM is one of the most important factors influencing patient satisfaction following TKA, as a full ROM is critical for a patient to perform daily activities.⁴¹

JOURNEY II TKA has demonstrated better kinematics,^{17,24,25} ROM²⁶⁻²⁹ and gait^{24,30} compared to conventional TKAs.

Improved knee flexion^{24,27-29}



Post-TKA:

3 months

Significantly improved knee flexion during walking ($p < 0.01$) with JOURNEY II CR compared to Attune™ CR* TKA.²⁴

A significantly increased mean ROM with JOURNEY II BCS compared to conventional PS TKAs:

+6° ($p = 0.04$)²⁷

+6° ($p = 0.002$)²⁸

+23° ($p < 0.001$)²⁹



1 year

Improved gait^{24,30} and stair climbing⁴²

Post-TKA:

3 months


Significantly greater increase in external rotation whilst walking with JOURNEY II CR compared to Attune™ CR* TKA ($p < 0.01$).²⁴

6 months

Significant improvements in gait with JOURNEY II BCS compared to LEGION^o PS TKA ($p = 0.03$).³⁰

1 year

Significantly more patients able to climb stairs with JOURNEY II BCS compared to LEGION PS TKA ($p < 0.05$).⁴²



Better kinematic outcomes^{16,17}

Compared to the normal knee after 1 year post-TKA:



A similar axial rotation in early and late flexion with JOURNEY II BCS¹⁶



A similar AP translation and axial rotation in mid flexion with JOURNEY II XR¹⁷

*Manufactured by DePuy Synthes Orthopaedics, Warsaw, IN, USA.





Key outcome: Smoother recovery

16 studies

Pain relief, function and return to activities are some of the most common patient expectations post-TKA.²

JOURNEY II TKA* has been shown to result in speedier recovery^{29,31,32} and improved muscle activation and muscle strength in the early recovery period,^{24,30} compared to conventional implants.

Improved quadricep activation and strength^{24,30}

Improved muscle activation and muscle strength in the early recovery period with JOURNEY II BCS and JOURNEY II CR, compared to LEGION^o PS TKA and Attune™† CR, respectively.^{24,30}

Improving quadriceps function is important for limiting post-TKA functional deficits⁹




Speedier recovery^{31,32}

Compared with patients receiving conventional TKA, JOURNEY II TKA patients are:


More likely to be discharged home (p<0.001)³¹

Discharged from hospital sooner (p<0.0001)³²



Higher function outcome measures^{26,29}

Compared to conventional TKA patients, JOURNEY II TKA patients have demonstrated higher scores in functional outcome measures at 6 and 12 weeks post-TKA.^{26,29}



Similar levels of patient satisfaction and patient-reported outcomes between JOURNEY II BCS and THA patients at 3 months post-op³³



*As observed in JOURNEY II BCS and JOURNEY II CR. †Manufactured by DePuy Synthes Orthopaedics, Warsaw, IN, USA.





Key outcome: Higher patient satisfaction

 17 studies

The primary determinant of patient satisfaction is the fulfilment of patient expectations, of which pain relief, improved knee function and return to sports are the most common.²

JOURNEY II TKA has demonstrated improvements in patient-reported outcome measures post-TKA,^{26,29,33-37} compared to pre-TKA.

Improved WOMAC scores²⁶

Significant reductions in WOMAC scores with JOURNEY II CR compared to P.F.C. Sigma™ CR* at 6 (p=0.018) and 12 months (p=0.008) post-TKA.²⁶



WOMAC scores⁴³, KSS⁷⁵, and KOOS⁷⁶ directly correlate to patient satisfaction

Better KSS^{25,26,29,36} and KOOS³⁵

Significantly better KSS outcomes with:

JOURNEY II CR compared to P.F.C. Sigma™ CR* at 3, 6, and 24 months post-TKA (p<0.05) and at 12 months post-TKA (p<0.001)²⁶

JOURNEY II BCS compared to Persona™† PS TKA (p<0.05)²⁵ and to LEGION^o PS TKA (p<0.001)²⁹ at 12 months post-TKA

JOURNEY II BCS at 24 months post-TKA compared to pre-TKA (p<0.01)³⁶



Quality of life³³

Patients receiving JOURNEY II BCS have reported similar levels of satisfaction as those receiving THA at 3 months and 1 year post-op, with no significant differences in patient quality of life.³³

Return to work and sports³³

One study reported a
2 month

median time to return to work, activities of daily living and sporting activities³³



Significantly higher KOOS with JOURNEY II XR^o compared to pre-TKA scores at 6 weeks, 3 months, 6 months and 12-months post-TKA (p<0.05).³⁵

Improving patient satisfaction is important for maintaining quality of care and patient loyalty which is linked to reimbursement to payors and healthcare providers^{44,45}





Key outcome: High survivorship

6 studies

Revision TKA is a technically challenging procedure and is associated with a high risk of complications.⁴⁶

Studies on survivorship of JOURNEY II BCS, JOURNEY II CR and JOURNEY II XR^o have shown favourable survivorship outcomes, with high early^{35,38} and mid-term^{39,40} survivorship.

High early survivorship

Early results from prospective studies on both JOURNEY II CR (n=174)³⁸ and JOURNEY II XR (n=165)³⁵ have demonstrated a low risk of revision:

<1.2% Risk of revision at 1 year^{35,38}


Revision TKA poses significant economic costs to the healthcare industry.⁴⁶

\$75,000
Average estimated US cost of revision TKA⁴⁷


High mid-term survivorship

JOURNEY II BCS has demonstrated high mid-term survivorship*, comparable to AOANJRR cemented PS TKA control.^{39,40}

Retrospective, international study of **2,059 JOURNEY II BCS TKAs**³⁹



Prospective, multicentre, US study of **209 JOURNEY II BCS TKAs**⁴⁰



JOURNEY II BCS resulted in:
<1% major revisions, lower than registry control^{39,40}
and a lower revision rate in <55 year olds compared to registry control (3 vs 7%)³⁹

*Kaplan-Meier (K-M) approach was used to obtain estimates of implant survival.





Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study²⁴

Di Benedetto P, Vidi D, Colombo A, Buttironi MM, Cainero V, Causero A. *Acta Biomed.* 2019;90:91–97

Independent pilot study of:

12 JOURNEY II CR
Mean age: 70.3 years

14 Attune™ CR*
Mean age: 71.8 years

Follow up: **3 months**

Assessed:

Kinematics were assessed at flexion-extension internal-external rotation and abduction-adduction

Muscle strength and activation using electromyography


KSS

KOOS



Results

Compared to Attune™ CR* at 3 months post-TKA, JOURNEY II CR patients demonstrated:



Better muscle activation timing, comparable to the physiological knee

37.5 **35.4** A numerically higher KSS satisfaction score (maximum score: 40)

84.9 **76.2** A numerically higher KOOS quality of life score (maximum score: 100)
JOURNEY II CR Attune CR

Conclusion

JOURNEY II CR patients experienced significantly better flexion, external rotation and muscle strength during walking at 3 months post-TKA compared to Attune™ CR* TKA. In addition, JOURNEY II CR patients demonstrated physiological activation timing of most muscles.

*Manufactured by DePuy Synthes Orthopaedics, Warsaw, IN, USA.



In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy¹⁶

Grieco TF, Sharma A, Dessinger GM, Cates HE, Komistek RD. *J Arthroplasty*. 2018;33:565–571

Retrospective, single-surgeon study of:

40 JOURNEY II BCS
Mean age: 69.8 years

10 normal knees
Mean age: 57.4 years

Mean follow up: **14.3 months**

Assessed:

Kinematics assessed during a weight-bearing deep knee bend at 30° increments from full extension to 120° of knee flexion to determine whether the dual cam-post mechanism is able to replicate the cruciate ligament (ACL and PCL) function



Results



0–30°

JOURNEY II BCS subjects exhibited similar patterns of femoral rollback and axial rotation compared with normal knee subjects



30–60°

JOURNEY II BCS subjects experienced minimal anterior-posterior motions and axial rotation, whereas normal knees continued to roll back and externally rotate



60–90°

JOURNEY II BCS resumed posterior motion with axial rotation increasing in a normal-like fashion after 90°

Conclusion

JOURNEY II BCS exhibited normal-like kinematic patterns and moved as designed under in vivo observation. Similarities in early and late kinematic patterns between the two groups suggest the dual cam-post design and asymmetric articular geometries of the JOURNEY II BCS adequately replicate ACL and PCL function of the normal knee.



A comparison of rollback ratio between bicruciate substituting total knee arthroplasty and Oxford unicompartmental knee arthroplasty²¹

Iriuchishima T and Ryu K. *J Knee Surg.* 2018;31:568–572

Retrospective, single-surgeon study of:

64 JOURNEY II BCS Mean age: 71.3 years	50 Oxford™ Partial Knee* Mean age: 73.8 years	50 normal knees N/A
---	--	------------------------



Follow up: **6–9 months**

Assessed:


When subjects had recovered their range of knee flexion (6–9 months post-op), lateral radiographs in active flexion were taken to measure rollback ratio and flexion angle




Results



JOURNEY II BCS



Oxford™ UKA*



Normal knee

No significant difference in rollback ratio and no significant difference in knee flexion angle among the three groups

A significant correlation between rollback ratio and knee flexion angle among the three groups ($p=0.002$)

Conclusion

JOURNEY II BCS showed no significant difference in rollback ratio when compared with UKA or asymptomatic control knees. The results suggest that JOURNEY II BCS reproduces native anterior and posterior cruciate function and native knee rollback.

*Manufactured by Zimmer Biomet, Warsaw, Indiana, USA.



The bicruciate substituting knee design and initial experience²⁹

Nodzo SR, Carroll KM, Mayman DJ. *Tech Orthop.* 2018;33:37–41

Retrospective, single-surgeon, study of:

100 JOURNEY II BCS 100 LEGION^o Total Knee System
Mean age: 51.0 years

Mean follow up: **1 year**

Assessed at pre-TKA and 6 weeks and 1 year post-TKA:


ROM

KSS



Results

JOURNEY II BCS patients, compared to PS TKA patients, demonstrated:



Significantly greater mean ROM at 1 year post-TKA (119° vs 96°; $p < 0.0001$)



Numerically higher mean KSS score at 6 weeks post-TKA (88 vs 73)



Significantly higher mean KSS scores at 1 year post-TKA (89 vs 81; $p < 0.001$)

Conclusion

JOURNEY II BCS led to significant improvements in ROM and patient-reported outcomes at 1 year post-TKA, compared with standard PS TKA. The authors noted that the results suggest that the more anatomic design of the implant, which is intended to replicate a more normal knee position and kinematic patterns, may be responsible for the improved flexion and patient satisfaction in JOURNEY II BCS patients.



In vivo knee kinematics: how important are the roles of femoral geometry and the cruciate ligaments?¹⁷

Smith LA, Nachtrab J, LaCour M, Cates H, Freeman MG, Komistek RD. *J Arthroplasty*. 2021;36:1445–1454

Kinematic study of:

40 JOURNEY II CR
Mean age: 68.8 years

10 JOURNEY II XR[®]
Mean age: 62.3 years

10 normal knees
Mean age: 57.4 years

Mean follow up: **24.8 months (JOURNEY II CR)**
and **16 months (JOURNEY II XR)**

Assessed:

Kinematics were measured at full extension and at 30° increments to full flexion, during a weight-bearing deep knee bend



Results



Early flexion
0–30°

JOURNEY II XR demonstrated a similar magnitude of posterior femoral rollback (PFR) to normal knees



Mid flexion
30–60°

JOURNEY II XR demonstrated similar anterior-posterior translation or axial rotation compared to normal knees



Deeper flexion
60–90°

No significant difference with JOURNEY II XR or CR in anterior-posterior translation compared to normal knees



Late flexion
90°+

Conclusion

In early and mid flexion, JOURNEY II XR patients experienced more normal-like kinematic patterns compared to JOURNEY II CR patients, demonstrating the role of ACL retention and subtle changes in femoral geometry in kinematic outcomes. However, the authors stated that JOURNEY II CR subjects achieved increased amounts of lateral condyle PFR and axial rotation compared to previously studied CR TKAs.



A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY II bi-cruciate stabilizing knee system and total hip arthroplasty patients³³

Snyder MA, Sympson A, Gregg J, Levit A. *Orthop Trauma Prosth.* 2018;3:5–10

Retrospective review of clinically matched data from a regional US total joint registry (Cincinnati, Ohio):

48 JOURNEY II BCS
Mean age: 58.3 years

48 THA
Mean age: 55.9 years

Mean follow up: **1 year**

Assessed at 3 months and 1 year post-op:

Patient satisfaction

UCLA activity scores

EQ-5D-5L



Results

JOURNEY II BCS patients, compared to THA patients, reported:



No significant difference in overall satisfaction at 3 months ($p=0.398$) or 1 year ($p=0.590$)



Significantly greater median UCLA activity scores at 3 months (8 vs 7; $p=0.028$) and 1 year (8 vs 7; $p<0.001$; maximum score: 10)



Significantly greater median EQ-5D-5L scores at 3 months (90 vs 80; $p<0.001$; maximum score: 100); no significant difference at 1 year

Conclusion

JOURNEY II BCS TKA patients reported similar levels of satisfaction and patient reported outcomes as THA patients at 3 months and 1 year post-op. The kinematic design advancements of JOURNEY II BCS demonstrate improvements in patient satisfaction following TKA compared to past TKA procedures, via comparison with THAs.



Anatomical bi-cruciate retaining TKA improves gait ability earlier than bi-cruciate stabilized TKA based on triaxial accelerometry data: a prospective cohort study⁴⁸

Amemiya K, Kaneko T, Omata M, Igarashi T, Takada K, Ikegami H, Musha Y. *AP-SMART*. 2021;25:35–41

Independent, single-surgeon, prospective study of:

10 JOURNEY II XR^o
Mean age: 75.0 years

15 JOURNEY II BCS
Mean age: 74.4 years

Follow up: **1.5–3 months**

Assessed at pre-TKA and at 6 weeks and 3 months post-TKA:

Walking time, number of steps, walking cycle time, coefficients of variability (CV) of double leg support time, gait velocity and stride length during a 10m walk test



Results


JOURNEY II XR patients, compared to JOURNEY II BCS patients, experienced:



Significantly improved gait ability at 6 weeks post-TKA ($p < 0.05$), with a significant reduction in walking time (9.2 vs 13.7s; $p = 0.01$) and in number of steps (18.1 vs 23.6; $p = 0.02$)



Significantly greater improvement in gait stability from pre-TKA to 6 weeks post-TKA (-27.2 vs 85.8; $p = 0.04$)



Significantly greater velocity (66.2 vs 45.8m/min; $p = 0.01$) and stride length (56.2 vs 43.2cm; $p = 0.01$)

Conclusion

JOURNEY II XR patients experienced significantly improved gait ability compared to JOURNEY II BCS at 6 weeks post-TKA, with a significant improvement in gait stability compared to pre-TKA.



Midterm performance of a guided-motion bicruciate-stabilized total knee system: results from the international study of over 2000 consecutive primary total knee arthroplasties³⁹

Harris AI, Christen B, Malcorps JJ, O'Grady CP, Kopjar B, Sensiba PR, Vandenneucker H, Huang BK, Cates HE, Hur J, Marra DA. *J Arthroplasty*. 2019;34:S201–S208

Retrospective, multicentre, international case series study of:

2,059 JOURNEY II BCS

Mean age, 64.3 years

Mean follow up: **3.87 years**

Assessed at last follow-up:

Survivorship compared to all other PS cemented TKA in the AOANJRR



Results

At 5 years, JOURNEY II BCS, compared to the AOANJRR control, resulted in:



Similar high survivorship rate (96.4 vs 95.9%)



Lower percentage of major revisions, accounting for <1% of JOURNEY II BCS TKAs and defined as all revisions involving tibial and/or femoral component removal (29.9 vs 41.6%)



Lower revision rate in patients <55 year olds (3 vs 7%)

Conclusion

The JOURNEY II BCS knee system performs favourably when compared to PS TKA control from the AOANJRR, particularly in patients <55 years of age. At up to 6.1 years post TKA, less than 1% of all JOURNEY II BCS TKAs required major revision.



Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty³⁰

Hyodo K, Kanamori A, Kadone H, Takahashi T, Kajiwara M, Yamazaki M. *Arthroplast Today*. 2020;6:338–342

A retrospective analysis of:

12 JOURNEY II BCS TKA
Mean age: 69.4 years

12 LEGION^o TKA
Mean age: 70.0 years

Mean length of study: **6 months**

Assessed at 6 months post-TKA:

Gait and motion capture analysis with a force plate and electromyogram of the lower limb muscles



Results

JOURNEY II BCS TKA, compared to conventional TKA, demonstrated:




Significantly longer step length (p=0.03)



Significantly faster gait speed (p=0.03)

0.55Nm/kg
0.34Nm/kg



Significantly greater maximum knee extension moment (p=0.04)

■ JOURNEY II TKA ■ Conventional TKA



Significantly greater maximum knee flexion angle during initial stance phase (p=0.04)

Conclusion

JOURNEY II BCS patients demonstrate faster gait speed, longer step length and greater knee extension moment, compared to conventional PS TKA. This suggests that the quadriceps muscles are more effectively activated and that anterior stability function of the anterior cruciate ligament can be achieved with JOURNEY II BCS TKA compared to conventional PS TKA.



In vivo kinematics of bicruciate-retaining total knee arthroplasty with anatomical articular surface under high-flexion conditions⁴⁹

Kono K, Inui H, Tomita T, Yamazaki T, Taketomi S, Tanaka S. *J Knee Surg.* 2021;34:452–459

A single-center case series of:

17 JOURNEY II XR[®]
Mean age: 72.0 years

Mean length of study: **7.6 months**

Assessed:

In vivo kinematics during high-flexion activities

Fluoroscopy during squatting and cross-legged sitting motions, starting with legs fully extended




Results

JOURNEY II XR kinematics were similar to the normal knee:

During cross-legged sitting: slight femoral internal rotation from early- flexion to mid-flexion, followed by femoral external rotation beyond mid- flexion

Significantly larger external rotation from 80 to 110° of flexion during squatting compared to cross-legged sitting ($p < 0.05$)



Conclusion

JOURNEY II XR patients exhibited kinematics which were consistent with those reported in the literature for normal knees and which differed depending on high-flexion activity.



Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics²⁶

Lutes W and Fitch D. Presented at: 39th SICOT Orthopaedic World Congress; 2018; Montreal, Canada

Retrospective, single-surgeon study of:		Assessed at 3, 6, 12 and 24 months post-TKA:	
52 JOURNEY II CR Mean age: 67.3 years	60 P.F.C. Sigma™ CR* Mean age: 70.2 years	KSS	✓
Mean follow up: 2 years		WOMAC scores	✓
		ROM (only assessed up to 12 months)	✓

Results:

JOURNEY II CR patients, compared to P.F.C. Sigma™ CR* patients, reported:



Significantly greater KSS scores at 3 (69.5 vs 63.0), 6 (84.4 vs 70.1), 12 (93.0 vs 86.1) and 24 (96.4 vs 91.7) months post-TKA ($p < 0.05$)



Significantly lower WOMAC scores at 6 (17.8 vs 24.6) and 12 (12.4 vs 18.5) months post-TKA ($p < 0.05$)



Significantly greater change in ROM from baseline at 3 (-4.4 vs -10.1), 6 (5.8 vs -1.8) and 12 (11.4 vs 4.0) months post-TKA ($p < 0.05$)

Conclusion

JOURNEY II CR patients reported significant improvements in functional outcomes compared to P.F.C. Sigma™ CR* TKA.



Hospital-related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patient³¹

Mayman DJ, Patel AR, Carroll KM. Poster presented at: ISPOR Symposium; 2018; Baltimore, Maryland, USA

A retrospective cohort study with real world evidence of:

1,692 JOURNEY II BCS TKA
Mean age: 64.1 years

1,692 other TKA devices
Mean age: 63.9 years

Assessed:

Hospital-related clinical and economic outcomes (reported as 2016 US dollars)

1:1 propensity score matching to control for patients and provider characteristics



Results:

JOURNEY II BCS, compared to other TKA knees, was associated with:



Patients receiving JOURNEY II BCS were also:



Conclusion

JOURNEY II BCS led to a significantly lower total hospital cost and significantly shorter hospital stay, and patients were less likely to be readmitted within 30 days. Patients receiving this device were also more likely to be discharged to home when compared with patients undergoing primary TKA with other posterior-stabilized TKA systems.



Evaluation of anteroposterior accelerometric change after bi-cruciate stabilized total knee arthroplasty and posterior stabilized total knee arthroplasty²⁵

Tomite T, Saito H, Kijima H, Ishikawa N, Hatakeyama Y, Tazawa H, Miyakoshi N, Shimada Y. *Knee*. 2021;32:121–130

Independent, single-surgeon, prospective study of:

30 JOURNEY II BCS
Mean age: 77.9 years

30 PERSONA™ PS*
Mean age: 77.3 years

Mean follow up: **1 year**

Assessed at pre-TKA and 12 months post-TKA:

ROM

KSS

FTA


Lateral X-rays of the standing extended knee

Accelerometer data



Results


JOURNEY II BCS patients, compared to PS TKA patients, demonstrated:



Significantly lower posterior offset ratio (2.1 vs 17.9%, respectively), with the AP positioning of the femur and tibia close to that of the normal knee



Lower anteroposterior (AP) acceleration on the femoral side of the knee



Conclusion

JOURNEY II BCS resulted in better functional kinematics, closer positioning to that of the normal knee on lateral X-ray, and lower anteroposterior acceleration on the femoral side compared to PS TKA at 12 months post-TKA.

All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Amemiya K, et al. 2021 ⁴⁸	Anatomical bi-cruciate retaining TKA improves gait ability earlier than bi-cruciate stabilized TKA based on triaxial accelerometry data: A prospective cohort study	AP-SMART		XR, BCS	✓	✓			
Biały M, et al. 2021 ⁵¹	A comparison of the JOURNEY II ^o bi-cruciate stabilized total knee system and Genesis II cruciate-retaining implant	Physiotherapy Review		BCS	✓	✓	✓		
Boese K, et al. 2019 ³⁵	Early clinical and patient-reported results of a bi-cruciate retaining total knee implant: six-month results of a prospective multicentre study of 149 primary TKAs	EORS Congress		XR	✓	✓	✓	✓	
Christen B, et al. 2018 ⁵²	Second-generation bi-cruciate stabilized total knee system has a lower reoperation and revision rate than its predecessor	Arch Orthop Trauma Surg		BCS		✓			
Coomer S, et al. 2021 ⁵³	Determining patella function in non-implanted knees having functional cruciate ligaments and subjects having a bi-cruciate stabilized total knee arthroplasty	Orthopaedic Proceedings		BCS	✓				

All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function 	Recovery 	Patient Satisfaction 	Survivorship 	Health Economics 
Di Benedetto P, et al. 2019 ²⁴	Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study	Acta Biomed		CR	✓	✓	✓		
Grieco TF, et al. 2018 ¹⁶	In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy	J Arthroplasty		BCS	✓				
Harris AI, et al. 2018 ³⁷	Short-term safety and effectiveness of a second-generation motion-guided total knee system	Arthro Today		BCS	✓	✓	✓		
Harris AI, et al. 2019 ³⁹	Midterm performance of a guided-motion bicruciate-stabilized total knee system: results from the international study of over 2000 consecutive primary total knee arthroplasties	J Arthroplasty		BCS				✓	
Harris AI, et al. 2019 ⁴⁰	Guided motion total knee arthroplasty system: five-year outcomes of the prospective multicenter US study	EFFORT Congress		BCS				✓	

All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function 	Recovery 	Patient Satisfaction 	Survivorship 	Health Economics 
Harris AI, et al. 2019 ⁵⁴	Guided motion total knee arthroplasty (TKA) system in younger patients has a lower revision rate than registry controls: results from the international multicenter study with up to 6 Years follow-up	EKS Congress		BCS				✓	
Harris AI, et al. 2019 ⁵⁵	Guided motion total knee arthroplasty (TKA) in patients with BMI of 40kg/m ² or more: results from the international multicentre study of 2,059 primary TKAs with up to 6 years follow-up	EKS Congress		BCS				✓	
Heir S, et al. 2019 ³⁶	Clinical and functional outcomes of a second-generation guided motion total knee arthroplasty system: two-year results of a prospective multicentre study	EKS Congress		BCS	✓	✓	✓		
Hino K, et al. 2018 ⁵⁶	Bi-cruciate substituting total knee arthroplasty provides varus–valgus stability throughout the midflexion range	Knee		BCS	✓				
Hommel H, et al. 2017 ⁵⁷	Good early results obtained with a guided-motion implant for total knee arthroplasty: A consecutive case series	Open Orthop J		BCS	✓				

All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Hyodo K, et al. 2020 ³⁰	Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty	Arthroplasty Today		BCS	✓	✓			
Inui H, et al. 2018 ⁵⁸	The relationship between anteroposterior stability and medial-lateral stability of the bi-cruciate stabilized total knee arthroplasty	Knee		BCS	✓				
Inui H, et al. 2019 ⁵⁹	The relationship between soft-tissue balance and intraoperative kinematics of guided motion total knee arthroplasty	J Knee Surg		BCS	✓				
Iriuchishima T, et al. 2018 ²¹	A comparison of rollback ratio between bicruciate substituting total knee arthroplasty and Oxford unicompartmental knee arthroplasty	J Knee Surg		BCS	✓				
Iriuchishima T, et al. 2019 ⁶⁰	Bicruciate substituting total knee arthroplasty improves stair climbing ability when compared with cruciate-retain or posterior stabilizing total knee arthroplasty	Indian J Orthop		BCS	✓				

All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Ishibashi T, et al. 2020 ⁴²	Kinematics of bicruciate and posterior stabilized total knee arthroplasty during deep knee flexion and stair climbing	J Orthop Res		BCS	✓				
Ishida K, et al. 2017 ⁶¹	Comparison of intra-operative navigation-based kinematics between bi-cruciate-stabilised total knee arthroplasty (TKA) and conventional posterior-stabilised TKA	Orthop Proceedings		BCS	✓				
Itou J, et al. 2021 ⁶²	Anterior prominence of the femoral condyle varies among prosthesis designs and surgical techniques in total knee arthroplasty	BMC Musculoskeletal Disorders		BCS	✓				
Kage T, et al. 2021 ¹⁵	The association between in vivo knee kinematics and patient-reported outcomes during squatting in bicruciate-stabilized total knee arthroplasty	J Knee Surg		BCS	✓		✓		
Kaneko T, et al. 2017 ¹⁹	Bi-cruciate substituting total knee arthroplasty improved medio-lateral instability in midflexion range	J Orthop		BCS	✓		✓		

All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Kaneko T, et al. 2018 ⁶³	The influence of compressive forces across the patello-femoral joint on patient reported outcome after bi-cruciate stabilized total knee arthroplasty	Bone Joint J		BCS			✓		
Kaneko T, et al. 2020 ⁶⁴	The influence of tibiofemoral joint forces on patient-reported outcome measurements after bicruciate stabilized total knee arthroplasty	J Orthop Surg		BCS	✓		✓		
Kiyohara M, et al. 2021 ⁶⁵	Comparison of in vivo knee kinematics before and after bicruciate-stabilized total knee arthroplasty during squatting	BMC Musculoskeletal Disorders		BCS	✓	✓			
Kono K, et al. 2019 ⁶⁶	Bicruciate-stabilised total knee arthroplasty provides good functional stability during high-flexion weight-bearing activities	Knee Surg Sports Traumatol Arthosc		BCS	✓				
Kono K, et al. 2021 ⁴⁹	In vivo kinematics of bicruciate-retaining total knee arthroplasty with anatomical articular surface under high-flexion conditions	J Knee Surg		XR	✓				

All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Kono K, et al. 2021 ⁶⁷	In vivo kinematics and cruciate ligament forces in bicruciate-retaining total knee arthroplasty	Scientific Reports		XR	✓	✓			
Kopjar B, et al. 2019 ³⁸	Clinical and functional outcomes of JOURNEY ^o II CR total knee system. Interim results of an ongoing, prospective, multicenter study	ISTA Congress		CR	✓		✓	✓	
Kosse NM, et al. 2018 ²⁷	Minor adaptations in implant design bicruciate-substituted total knee system improve maximal flexion	EFFORT Congress		BCS	✓				
Lutes W, et al. 2018 ²⁶	Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics	SICOT OW Congress		CR	✓	✓	✓		
Mayman DJ, et al. 2018 ³¹	Hospital related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patients	ISPOR Symposium		BCS		✓			✓
Murakami K, et al. 2018 ¹⁸	In vivo kinematics of gait in posterior-stabilized and bicruciate-stabilized total knee arthroplasties using image-matching techniques	Int Orthop		BCS	✓				

All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Murakami K, et al. 2018 ²²	Knee kinematics in bi-cruciate stabilized total knee arthroplasty during squatting and stair-climbing activities	J Orthop		BCS	✓				
Murakami K, et al. 2018 ⁶⁸	Preoperative tibial mechanical axis orientation and articular surface design influence on the coronal joint line orientation relative to the ground during gait after total knee arthroplasties	Knee Surg Sports Traumatol Arthrosc		BCS	✓				
Moewis P, et al. 2020 ⁶⁹	Retention of posterior cruciate ligament alone may not achieve physiological knee joint kinematics after total knee arthroplasty: a retrospective study	J Bone Joint Surg Am		BCS, CR	✓		✓		
Nodzo SR, et al. 2018 ²⁹	The bicruciate substituting knee design and initial experience	Tech Orthop		BCS	✓	✓	✓		
Oikonomidis L, et al. 2020 ⁷⁰	The Journey bicruciate knee replacement: design modifications yield better early functional results and reduce complications	J Knee Surg		BCS	✓		✓		

All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function 	Recovery 	Patient Satisfaction 	Survivorship 	Health Economics 
Patel AR, et al. 2019 ³²	Hospital related clinical and economical outcomes of two premium knee systems in total knee arthroplasty (TKA) patients	ISPOR Symposium		BCS		✓			✓
Salzman M, et al. 2017 ⁷¹	Does postoperative mechanical axis alignment have an effect on clinical outcome of primary total knee arthroplasty? A retrospective cohort study	Open Orthop J		BCS	✓		✓		
Smith AL, et al. 2021 ¹⁷	In vivo knee kinematics: how important are the roles of femoral geometry and the cruciate ligaments?	J Arthroplasty		CR, XR	✓				
Snyder MA, et al. 2018 ³³	A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY ^o II bi-cruciate stabilizing knee system and total hip arthroplasty patients	Orthop Trauma Prosth		BCS	✓	✓	✓		
Takubo A, et al. 2017 ²⁸	Comparison of muscle recovery following bicruciate substituting versus posterior stabilized total knee arthroplasty in an Asian population	J Knee Surg		BCS	✓				

All evidence

Select the study icon to see the study overview (if applicable).

Authors/ Year	Study Title	Journal/ Source	Useful Links	JOURNEY II	Function	Recovery	Patient Satisfaction	Survivorship	Health Economics
									
Tomite T, et al. 2016 ⁷²	Gait analysis of conventional total knee arthroplasty and bicruciate stabilized total knee arthroplasty using a triaxial accelerometer	Case Report Orthop		BCS	✓				
Tomite T, et al. 2021 ²⁵	Evaluation of anteroposterior accelerometric change after bi-cruciate stabilized total knee arthroplasty and posterior stabilized total knee arthroplasty	Knee		BCS	✓				
Tria Jr AJ, et al. 2021 ⁷³	Bicruciate retaining total knee arthroplasty: Short-term clinical outcomes using a novel prosthesis	J Long Term Eff Med Implants		XR	✓				
West JA, et al. 2019 ³⁴	Clinical outcomes and patient satisfaction after total knee arthroplasty: a follow-up of the first 50 cases by a single surgeon	J Int Med Res		BCS	✓	✓	✓		
Zambianchi F, et al. 2018 ⁷⁴	Changes in total knee arthroplasty design affect in vivo kinematics in a redesigned total knee system: A fluoroscopy study	Clin Biomech		BCS	✓				

Smith & Nephew, Inc
1450 Brooks Road
Memphis,
TN 38116
USA

www.smith-nephew.com

◇ Trademark of Smith+Nephew
All Trademarks acknowledged
© 2022 Smith+Nephew
17058 V4 0622. Published June 2022.

Developed by Evidence Communications,
Global Clinical & Medical Affairs

Smith+Nephew

References

- Collins M, Lavigne M, Girard J, Vendittoli PA. Joint perception after hip or knee replacement surgery. *Orthop Traumatol Surg Res.* 2012;98:275–280.
- Noble PC, Gordon MJ, Weiss JM, Reddix RN, Conditt MA, Mathis KB. Does total knee replacement restore normal knee function? *Clin Orthop Relat Res.* 2005;431:157–165.
- Scott CEH, Howie CR, MacDonald D, Biant LC. Predicting dissatisfaction following total knee replacement. *J Bone Joint Surg Am.* 2010;92-B(9):1253–1258.
- Lee A, Park J, Lee S. Gait analysis of elderly women after total knee arthroplasty. *J Phys Ther Sci.* 2015;27:591–595.
- McClelland JA, Webster KE, Feller JA, Menz HB. Knee kinematics during walking at different speeds in people who have undergone total knee replacement. *The Knee.* 2010;18:151–5.
- Dorr LD, Ochsner JL, Gronley J, Perry J. Functional Comparison of Posterior Cruciate-Retained versus Cruciate-Sacrificed Total Knee Arthroplasty. *Clin Orthop Relat Res.* 1988;236:36–43.
- Kramers-de Quervain IA, Stussi E, Muller R, Drobny T, Munzinger U, Gschwend N. Quantitative Gait Analysis After Bilateral Total Knee Arthroplasty With Two Different Systems Within Each Subject. *J Arthroplasty.* 1997;12(2):168–179.
- Saari T, Tranberg R, Züchner R, Uvehammer J, Karrholm J. Changed gait pattern in patients with total knee arthroplasty but minimal influence of tibial insert design. *Acta Orthop.* 2005;76(2):253–260.
- Thomas AC, Stevens-Lapsley JE. Importance of attenuating quadriceps activation deficits after total knee arthroplasty. *Exerc Sport Sci Rev.* 2012;40:95–101.
- Lester DK, Shantharam R, Zhang K. Dynamic electromyography after cruciate-retaining total knee arthroplasty revealed a threefold quadriceps demand compared with the contralateral normal knee. *J Arthroplasty.* 2013;28(4):557–562.
- Birdsall PD, Hayes JH, Cleary R, Pinder IM, Moran CG, Sher JL. Health outcome after total knee replacement in the very elderly. *J Bone Joint Surg Br.* 1999;81:660–662.
- Evans JT, Evans JP, Walker RB, Blom AW, Sayers A, Whitehouse MR. How long does a hip replacement last? A systematic review and meta-analysis of case series and national registry reports with more than 15 years of follow-up. *Lancet.* 2019;393:647–654.
- Angerame MR, Holst DC, Jennings JM, Komistek RD, Dennis DA. Total knee arthroplasty kinematics. *J Arthroplasty.* 2019;34:2502–2510.
- Van Onsem S, Verstraete M, Van Eenoo W, Van Der Straten C, Victor J. Are TKA kinematics during closed kinetic chain exercises associated with patient-reported outcomes? A preliminary analysis. *Clin Orthop Relat Res.* 2020;478:255–263.
- Kage T, Inui, H, Tomita T, et al. The association between in vivo knee kinematics and patient-reported outcomes during squatting in bicruciate-stabilized total knee arthroplasty. *J Knee Surgery.* 2021. doi:10.1055/s-0041-1723763.
- Grieco TF, Sharma A, Dessinger GM, Cates HE, Komistek RD. In vivo kinematic comparison of a bicruciate stabilized total knee arthroplasty and the normal knee using fluoroscopy. *J Arthroplasty.* 2018;33:565–571.
- Smith LA, Nachtrab J, LaCour M, Cates H, Freeman MG, Komistek RD. In vivo knee kinematics: how important are the roles of femoral geometry and the cruciate ligaments? *J Arthroplasty.* 2020. doi: <https://doi.org/10.1016/j.arth.2020.10.020>.
- Murakami K, Hamai S, Okazaki K, et al. In vivo kinematics of gait in posterior-stabilized and bicruciate-stabilized total knee arthroplasties using image matching techniques. *Int Orthop.* 2018;42:2573–2581.
- Kaneko T, Kono N, Mochizuki Y, Hada M, Toyoda S, Musha Y. Bi-cruciate substituting total knee arthroplasty improved medio-lateral instability in mid-flexion range. *J Orthop.* 2017;14:201–206.
- Catani F, Ensini A, Belvedere C, et al. In vivo kinematics and kinetics of a bi-cruciate substituting total knee arthroplasty: a combined fluoroscopic and gait analysis study. *J Orthop Res.* 2009;27(12):1569–1575.
- Iriuchishima T, Ryu K. A Comparison of rollback ratio between bicruciate substituting total knee arthroplasty and Oxford unicompartmental knee arthroplasty. *J Knee Surg.* 2018;31:568–572.
- Murakami K, Hamai S, Okazaki K, et al. Knee kinematics in bicruciate stabilized total knee arthroplasty during squatting and stair-climbing activities. *J Orthop.* 2018;15:650–654.
- Carpenter RD, Brilhault J, Majumdar S, Ries MD. Magnetic resonance imaging of in vivo patellofemoral kinematics after total knee arthroplasty. *Knee.* 2009;16(5):332–336.
- Di Benedetto P, Vidi D, Colombo, Buttironi MM, Cainero V, Causero A. Pre-operative and post-operative kinematic analysis in total knee arthroplasty. A pilot study. *Acta Biomed.* 2019;90:91–97.
- Tomite, T, Saito, H, Kijima H, et al. Evaluation of anteroposterior accelerometric change after bi-cruciate stabilized total knee arthroplasty and posterior stabilized total knee arthroplasty. *Knee.* 2021;32:121–130.
- Lutes W, Fitch D. Comparison of functional outcomes following total knee arthroplasty with a conventional implant design or one designed to mimic natural knee kinematics. Presented at: 39th SICOT Orthopaedic World Congress; October 10–13, 2018; Montreal, Canada.
- Kosse NM, Heesterbeek PJC, Defoort KC, Wymenga AB, van Hellemond G. Minor adaptations in implant design bicruciate-substituted total knee system improve maximal flexion. Poster presented at: 2nd World Arthroplasty Congress; 19–21 April, 2018; Rome, Italy.
- Takubo A, Ryu K, Iriuchishima T, Tokuhashi Y. Comparison of muscle recovery following bicruciate substituting versus posterior stabilized total knee arthroplasty in an Asian population. *J Knee Surg.* 2017;30:725–729.
- Nodzo SR, Carroll KM, Mayman DJ. The bicruciate substituting knee design and initial experience. *Tech Orthop.* 2018;33:37–41.
- Hyodo K, Kanamori A, Kadone H, Takahashi T, Kajiwara M, Yamazaki M. Gait analysis comparing kinematic, kinetic, and muscle activation data of modern and conventional total knee arthroplasty. *Arthroplast Today.* 2020;6:338–342.
- Mayman DJ, Patel AR, Carroll KM. Hospital related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patients. Poster presented at: ISPOR Symposium; May 19–23, 2018; Baltimore, Maryland, USA.
- Patel AR, Delhogue G. Hospital related clinical and economical outcomes of two premium knee systems in total knee arthroplasty (TKA) patients. Poster presented at: ISPOR Symposium; May 18–22, 2019; New Orleans, Louisiana, USA.
- Snyder MA, Simpson A, Gregg J, Levit A. A comparison of patient reported outcomes between total knee arthroplasty patients receiving the JOURNEY II bi-24. cruciate stabilizing knee system and total hip arthroplasty patients. *Orthop Trauma Prosth.* 2018; 2018; 3:5–10.
- West JA, Scudday T, Anderson S, Amin NH. Clinical outcomes and patient satisfaction after total knee arthroplasty; a follow-up of the first 50 cases by a single surgeon. *J Int Med Res.* 2019;47:1667–1676.
- Boese K, MacDonald J, Huang W, et al. Early clinical and patient-reported results of a bi-cruciate retaining total knee implant: six-month results of a prospective multicentre study of 149 primary TKAs. Poster presented at: European Orthopaedic Research Society (EORS); October 2–4, 2019; Maastricht, the Netherlands.
- Heir S, Catani F, van Hellemond G, et al. Clinical and functional outcomes of a second generation guided motion total knee arthroplasty system: two-year results of a prospective multicentre study. Abstract number P37 presented at: European Knee Society; May 2–3, 2019; Valencia, Spain.

References (cont.)

37. Harris AI, Luo TD, Lang JE, Kopjar B. Short-term safety and effectiveness of a second generation motion-guided total knee system. *Arthroplast Today*. 2018;4:240–243.
38. Kopjar B, Archan P, MacDonald J, et al. Clinical and functional outcomes of JOURNEY II CR total knee system. Interim results of an ongoing, prospective, multicenter study. Abstract presented at: International Society for Technology in Arthroplasty Congress (ISTA); October 2–5, 2019; Toronto, Canada.
39. Harris AI, Christen B, Malcorps JJ, et al. Mid-term performance of a guided motion bicruciate stabilized total knee system; results from the international study of over 2,000 consecutive primary total knee arthroplasties. *J Arthroplasty*. 2019;34:S201–S208.
40. Harris A, Luo TD, Lang JE, et al. Guided motion total knee arthroplasty system: five-year outcomes of the prospective multicentre US study. Presentation number 2761 presented at: EFORT; June 5–7, 2019; Lisbon, Portugal.
41. Matsuda S, Kawahara S, Okazaki K, Tashiro Y, Iwamoto Y. Postoperative alignment and ROM affect patient satisfaction after TKA. *Clin Orthop Relat Res*. 2013;471:127–133.
42. Ishibashi T, Tomita T, Yamazaki T, Tsuji S, Yoshikawa H, Sugamoto K. Kinematics of bicruciate and posterior stabilized total knee arthroplasty during deep knee flexion and stair climbing. *J Orthop Res*. 2020; doi: 10.1002/jor.24773.
43. Walker LC, Clement ND, Bardgett M, et al. The WOMAC score can be reliably used to classify patient satisfaction after total knee arthroplasty. *Knee Surg Sports Traumatol Arthrosc*. 2018;26:3333–3341.
44. Lyu H, Wick E, Housman M, Freischlag JA, Makary MA. Patient satisfaction as a possible indicator of quality surgical care. *JAMA Surg*. 2013;148:362–367.
45. Otani K, Waterman B, Faulkner K, et al. Patient satisfaction: focusing on “excellent”. *J Healthc Manag*. 2009;54:93–102.
46. Hamilton DF, Howie CR, Burnett R, Simpson AHRW, Patton JT. Dealing with the predicted increase in demand for revision total knee arthroplasty: challenges, risks and opportunities. *Bone Joint J*. 2015;97-b:723–728.
47. Delanois RE, Mistry JB, Gwam CU, Mohamed NS, Choksi US, Mont MA. Current epidemiology of revision total knee arthroplasty in the United States. *J Arthroplasty*. 2017;32:2663–2668.
48. Amemiya K, Kaneko T, Omata M, et al. Anatomical bi-cruciate retaining TKA improves gait ability earlier than bi-cruciate stabilized TKA based on triaxial accelerometry data: a prospective cohort study. *AP-SMART*. 2021;25:35–41.
49. Kono K, Inui H, Tomita T, Yamazaki T, Taketomi S, Tanaka S. In vivo kinematics of bicruciate-retaining total knee arthroplasty with anatomical articular surface under high-flexion conditions. *J Knee Surg*. 2021; 34 4:452–459.
50. Husain A, Lee GC. Establishing realistic patient expectations following total knee arthroplasty. *J Am Acad Orthop Surg*. 2015;23:707–713.
51. Biały M, Corte RD, Dec J, Pierzchała A, Gnat R. Evaluation of functional status in patients after total knee arthroplasty: a comparison of the JOURNEY II bi-cruciate stabilized total knee system and Genesis II cruciate-retaining implant. *Physiotherapy review*. 2021;25;3:24–34. doi:10.5114/phr.2021.109028.
52. Christen B, Kopjar B. Second-generation bi-cruciate stabilized total knee system has a lower reoperation and revision rate than its predecessor. *Arch Orthop Trauma Surg*. 2018;138:1591–1599.
53. Coomer S, LaCour M, Khasian M, Cates H, Komistek R. Determining patella function in non-implanted knees having functional cruciate ligaments and subjects having a bicruciate-stabilised total knee arthroplasty. *Orthopaedic Proceedings*. 2021;103:4–4.
54. Harris A, O’Grady C, Sensiba PR, et al. Guided motion total knee arthroplasty (TKA) system in younger patients has a lower revision rate than registry controls: results from the international multicenter study with up to 6 years follow-up. Abstract number P35 presented at: European Knee Society; May 2–3, 2019; Valencia, Spain.
55. Harris A, O’Grady C, Sensiba PR, et al. Guided motion total knee arthroplasty (TKA) in patients with BMI of 40kg/m² or more: results from the international multicentre study of 2,059 primary TKAs with up to 6 years follow-up. Abstract number P36 presented at: European Knee Society; May 2–3, 2019; Valencia, Spain.
56. Hino K, Kutsuna T, Watamori K, et al. Bi-cruciate substituting total knee arthroplasty provides varus–valgus stability throughout the midflexion range. *Knee*. 2018;25:897–902.
57. Hommel H, Wilke K. Good early results obtained with a guided-motion implant for total knee arthroplasty: A consecutive case series. *Open Orthop J*. 2017;11:51–56.
58. Inui H, Taketomi S, Yamagami R, Kawaguchi K, Nakazato K, Tanaka S. The relationship between anteroposterior stability and medial-lateral stability of the bi-cruciate stabilized total knee arthroplasty. *Knee*. 2018;25:1247–1253.
59. Inui H, Taketomi S, Yamagami R, Shirakawa N, Kawaguchi K, Tanaka S. The relationship between soft-tissue balance and intraoperative kinematics of guided motion total knee arthroplasty. *J Knee Surg*. 2019;32:91–96.
60. Iriuchishima T, Ryu K. Bicruciate substituting total knee arthroplasty improves stair climbing ability when compared with cruciate-retain or posterior stabilizing total knee arthroplasty. *Indian J Orthop*. 2019. DOI:10.4103/ortho.IJOrtho_392_18.
61. Ishida K, Shibanuma N, Toda A, et al. Comparison of intra-operative navigation-based kinematics between bi-cruciate-stabilised total knee arthroplasty (TKA) and conventional posterior stabilised TKA. *Orthopaedic proceedings*. 2017;99B:Suppl 4.
62. Itou J, Kuwashima U, Itoh M, Okazaki K. Anterior prominence of the femoral condyle varies among prosthesis designs and surgical techniques in total knee arthroplasty. *BMC musculoskeletal disorders*. 2021;22:1–8.
63. Kaneko T, Kono N, Mochizuki Y, et al. The influence of compressive forces across the patellofemoral joint on patient reported outcome after bi-cruciate stabilized total knee arthroplasty. *Bone Joint J*. 2018; 100-B:1585–1591.
64. Kaneko T, Kono N, Mochizuki Y, et al. The influence of tibiofemoral joint forces on patient-reported outcome measurements after bicruciate stabilized total knee arthroplasty. *J Orthop Surg*. 2020;28:1–2.
65. Kiyohara M, Hamai S, Gondo H, et al. Comparison of in vivo knee kinematics before and after bicruciate-stabilized total knee arthroplasty during squatting. *BMC Musculoskeletal Disorders*. 2021;22:1–10.
66. Kono K, Inui H, Tomita T, et al. Bicruciate-stabilised total knee arthroplasty provides good functional stability during high-flexion weight-bearing activities. *Knee Surg Sports Traumatol Arthrosc*. 2019;27:2096–2103.
67. Kono K, Inui H, Tomita T, et al. In vivo kinematics and cruciate ligament forces in bicruciate-retaining total knee arthroplasty. *Scientific Reports*. 2021;11:1–10.
68. Murakami K, Hamai S, Okazaki K, et al. Preoperative tibial mechanical axis orientation and articular surface design influence on the coronal joint line orientation relative to the ground during gait after total knee arthroplasties. *Knee Surg Sports Traumatol Arthrosc*. 2018;26:3368–3376.
69. Moewis P, Duda GN, Trepczynski A, et al. Retention of posterior cruciate ligament alone may not achieve physiological knee joint kinematics after total knee arthroplasty: a retrospective study. *J Bone Joint Surg Am*. 2020; doi: 10.2106/JBJS.20.00024.
70. Oikonomidis L, Santini AJA, Davidson JS, Banks JV, Phillipson A, Pope J. The JOURNEY bicruciate knee replacement: design modifications yield better early functional results and reduce complications. *J Knee Surg*. 2020; doi: 10.1055/s-0040-1718599.
71. Salzmann M, Fennema P, Becker R, Hommel H. Does postoperative mechanical axis alignment have an effect on clinical outcome of primary total knee arthroplasty? A retrospective cohort study. *Open Orthop J*. 2017;11:1330–1336.
72. Tomite T, Saito H, Aizawa T, Kijima H, Miyakoshi N, Shimada Y. Gait analysis of conventional total knee arthroplasty and bicruciate stabilized total knee arthroplasty using a triaxial accelerometer. *Case Report Orthop*. 2016;6875821.
73. Tria Jr AJ, Bateman DK, Preston JS, Diamond KB. Bicruciate retaining total knee arthroplasty: short-term clinical outcomes using a novel prosthesis. *J long-term eff med implants*. 2021;31:3.
74. Zambianchi F, Fiacchi F, Lombardi V et al. Changes in total knee arthroplasty design affect in vivo kinematics in a redesigned total knee system: A fluoroscopy study. *Clin Biomech*. 2018;54:92–102.
75. Goh GS, Bin Abd Razak HR, Tay DKJ, Lo N, Yeo S. Early post-operative oxford knee score and knee society score predict patient satisfaction 2 years after total knee arthroplasty. *Arch Orthop Trauma Surg*. 2021;141:129–137. doi: 10.1007/s00402-020-03612-2.
76. Wendelspiess S, Kaelin R, Vogel N, Rychen T, Arnold MP. No difference in patient-reported satisfaction after 12 months between customised individually made and off-the-shelf total knee arthroplasty. *Knee Surg Sports Traumatol Arthrosc*. 2022. doi: 10.1007/s00167-022-06900-z.

Abbreviations

ACL	Anterior cruciate ligament	PS	Posterior stabilised
AOANJRR	Australian Orthopaedic Association National Joint Replacement Registry	ROM	Range of motion
BCS	Bicruciate stabilised	THA	Total hip arthroplasty
CR	Cruciate retaining	TKA	Total knee arthroplasty
FJS-12	Forgotten joint score-12	TKS	Total knee system
ITB	Iliotibial band	UCLA	University of California Los Angeles
KOOS	Knee injury and osteoarthritis outcome score	UKA	Unicompartmental knee arthroplasty
KSS	Knee society score	WOMAC	The Western Ontario and McMaster Universities Arthritis Index
LOS	Length of stay	XR	Bicruciate retaining
PFR	Posterior femoral rollback		

JOURNEY II TKA: further information

Learn more:



Access additional
JOURNEY II TKA
resources and training at
+ Education Unlimited

Explore similar resources:



Robotics compendium
of clinical evidence



Contact your sales representative for further information