




Inflatable tamponade designed to create a fast, single treatment solution with a low-profile and self-lubricating properties

 **smith&nephew**  
**RAPID RHINO<sup>®</sup>**  
Epistaxis Products

# RAPID RHINO<sup>◇</sup> Epistaxis Products


Offering a full range of epistaxis solutions

Ultra low profile and self-lubricating hydrocolloid fabric provides for unsurpassed<sup>1</sup> ease of insertion and removal.




## Features

- Easy to insert and remove
- Quickly conforms to nasal anatomy
- Provides gentle and even compression to areas of epistaxis
- Designed to prevent adherence to tissue or blood clots upon removal
- Available in a wide range of sizes



RAPID RHINO Nasal Pac  
(Uninflated 5.5cm anterior shown)




RAPID RHINO with  
Carboxymethylcellulose (CMC)


## Comparison

A comparison of RAPID RHINO<sup>®</sup> epistaxis products to Rhino Rocket<sup>2</sup> for the treatment of epistaxis in the emergency department<sup>1</sup>: RAPID RHINO epistaxis products provide greater patient comfort at insertion and at removal.


### Ease of insertion



### Ease of removal




### Mean pain scores




■ RAPID RHINO  
■ Rhino Rocket


## RAPID RHINO product usage directions




1 Soak in sterile water for a FULL 30 seconds.




2 Insert along superior aspect of the hard palate until the blue indicator is past the nares.




3 Using a 20ml syringe, inflate the Rapid Rhino device with AIR only. Monitor the pilot cuff for direct tactile feedback; Stop inflation when the pilot cuff becomes rounded and feels firm when squeezed.



4 Inflate the cuff to provide a gentle, low-pressure tamponade delivering the CMC fabric directly to the bleed site.



5 Reassess after 15-20 minutes; reinflate to ensure proper pressure (if necessary) and tape to patient's cheek away from the upper lip.



6 Removal should occur 24-72 hours after treatment.

## Unilateral models



RR 450 nasal tamponade – 4.5cm inflatable epistaxis tamponade device



RR 530 non-inflatable tamponade RAPID-PAC® epistaxis tamponade device is 5.5cm offers a non-inflatable option for epistaxis



RR 550 inflatable tamponade – 5.5cm inflatable epistaxis tamponade device



RR 551 inflatable tamponade with airway – 5.5cm inflatable epistaxis tamponade device that contains an airway



RR 750 inflatable tamponade – 7.5cm inflatable epistaxis tamponade device



RR 751 inflatable tamponade with airway – 7.5cm inflatable epistaxis tamponade device that contains an airway



RR 900 – 9cm inflatable tamponade device with two balloons to control intranasal hemorrhage originating from the sphenopalatine artery

## Bilateral models



RR 552 5.5cm bilateral epistaxis device; RR 555 5.5cm bilateral epistaxis device with an airway



RR 752 7.5cm bilateral epistaxis device; RR 755 7.5cm bilateral epistaxis device with an airway

### References

- 1 Singer AJ, Blanda M, Cronin K, LoGiudice-Khwaja M, Gulla J, Bradshaw J, Katz. A Comparison of nasal tampons for the treatment of epistaxis in the emergency department: a randomized controlled trial. *Ann Emerg Med.* 2005; 45:134–139
- 2 Rhino Rocket is a registered trademark of Shippert Medical Technologies Corp. and is used for reference purposes only.

## Ordering information

### RAPID RHINO® Inflatable Tamponade Devices

Reference #	Description
RR 450	4.5cm Inflatable Unilateral without airway
RR 530	5.5cm Non-Inflatable Rapid-Pac®
RR 550	5.5cm Inflatable Unilateral without airway
RR 551	5.5cm Inflatable Unilateral with airway
RR 552	5.5cm Inflatable Bilateral without airway
RR 555	5.5cm bilateral epistaxis device with airway
RR 750	7.5cm Inflatable Unilateral without airway
RR 751	7.5cm Inflatable Unilateral with airway
RR 752	7.5cm Inflatable Bilateral without airway
RR 755	7.5cm Inflatable Bilateral with airway
RR 900	9cm Inflatable Dual Balloon

Supporting healthcare professionals for over 150 years



ArthroCare Corporation  
7000 West William Cannon Drive  
Austin, TX 78735 USA

[www.smith-nephew.com](http://www.smith-nephew.com)  
Order Entry: +1 800 343 5717  
Order Entry Fax: +1 888 994 2782

®Trademark of Smith & Nephew.  
©2016 Smith & Nephew. All rights reserved.  
P/N RR117 Rev. M 08/16