

+ Evidence in focus

Smith+Nephew

Real world implementation of the infection management (IM) pathway:
clinical cases involving management of wounds with biofilm and local infections



What is the infection management pathway?


A simple, evidence-based, easy-to-use tool for specialists and non-specialists to help improve patient outcomes by:¹

- **Simplifying differentiation and diagnosis** of biofilms and local infections
- **Guiding early and appropriate management** and selection of antimicrobials, supporting antimicrobial stewardship
- **Promoting effective and consistent care** among specialists and non-specialists
- **Prompting referral** and facilitating communication among the multidisciplinary team



Helps to identify signs and symptoms of infection¹


Uses a comprehensive **A B C D E approach** to wound assessment as recommended in the T.I.M.E. clinical decision support tool (CDST) and provides an illustrated guide to the key clinical signs and symptoms of infection^{1,2}



Biofilm and covert (subtle) infections

Overt (classic) and spreading or systemic infections

Guides appropriate management based on presence of biofilm versus local wound infection¹



Please refer to full infection management pathway for supporting references and footnotes

Guides use of antimicrobial dressings and multiple treatment approaches¹



Use of antimicrobial dressings is recommended for a minimum of two weeks

After two weeks, re-evaluate:


- *Discontinue* if signs and symptoms of infection have resolved
- *Continue* with antimicrobial if wound is progressing but there are still signs and symptoms
- *Consider* an alternative antimicrobial and refer to an appropriate specialist if no improvement



Step down, step up approach

Use of multiple treatment approaches can be scaled down or up depending on how well the infection is responding to treatment and how well wound healing is progressing

Applying the infection management pathway: specific cases



Note:

Some products used in these case studies may not be available in all countries




A chronic non-healing ulcer^{5,6}

PRESENTATION

- **74-year-old male** referred from ICU with a stalled wound (4 months)
- **Relevant medical history:**
Coronary artery disease, diabetes mellitus, hypertension, end-stage renal disease, stroke, anaemia
- **Prior treatments:**
Systemic antibiotics for osteomyelitis, local wound dressings of Aquacel™ Ag+ (ConvaTec, Canada) and Mepilex™ (Mölnlycke Health Care, Sweden)
- **Wound dimensions (length, width, depth):**
7.0×6.0×0.4cm



Biofilm: Chronic non-healing ulcer^{5,6}



REFER TO INFECTION MANAGEMENT PATHWAY

- Delayed healing despite optimal wound/patient management
- Copious amount of serous drainage, concurrent inflammation
- Hypergranulated
- Bleeding, friable granulation
- Malodorous
- 50% granulating, 50% slough




**Suspected
biofilm**

TREATMENT

- **Repeated aggressive debridement**
- **IODOSORB[◇] Ointment/Gel**
- **ALLEVYN[◇] Gentle Border Dressing**
- **Daily dressing changes** due to high volume of exudate
- **Off loading** with air mattress and positioning



Biofilm: Chronic non-healing ulcer^{5,6}



DAY 7 – REFER TO INFECTION MANAGEMENT PATHWAY

- Cleaner wound bed
- Healthy pink colour
- Non-friable
- Flattening periwound margins
- Highly exudating
- 70% granulating, 30% slough



**Suspected
biofilm**


DAY 7 – CHANGE TREATMENT

- Wound dimensions
(length, width, depth): **7.0×4.0×0.2cm**

**Change to IODOSORB[◇] Powder
for increased absorbency**




Biofilm: Chronic non-healing ulcer^{5,6}



DAY 13

100% granulation tissue
Wound dimensions (length, width, depth):
6.0×4.0×0.2cm
**Continue IODOSORB[◇] Powder plus
ALLEVYN[◇] Gentle Border Dressing**



DAY 17

95% granulation tissue,
5% epithelialisation tissue
Wound dimensions (length, width, depth):
6.0×4.0×0.2cm



**Positive
signs of
progression
to healing**





Pressure ulcer (trochanter area)^{5,6}


PRESENTATION

- **75-year-old male with a stalled wound** (>6 months)
- **Relevant medical history:**
Parkinson's disease, osteomyelitis, previous treatment for colon cancer, anaemia
- **Prior treatments:**
Systemic antibiotics for osteomyelitis, PHMB ribbon and foam dressing
- **Wound dimensions (length, width, depth):**
4.5×4.5×2.0cm

PHMB, polyhexamethylene biguanide



Local infection: Pressure ulcer (trochanter)^{5,6}



REFER TO INFECTION MANAGEMENT PATHWAY

- Erythema
- Warmth
- Oedema/swelling
- Purulent discharge, increasingly malodorous
- Pain
- 60% granulating, 40% slough




Local wound infection

TREATMENT

- **ACTICOAT[◇] FLEX 3 Dressing**
- **ALLEVYN[◇] Gentle Border Dressing**
- **Antibiotics** (osteomyelitis) to help prevent sepsis
- **Off loading** with air mattress and positioning



Local infection: Pressure ulcer (trochanter)^{5,6}



**Positive
signs of
progression
to healing**




DAY 7 – RE-ASSESSMENT

- Deeper wound, sloughy tissue reaccumulating
- More (hyper)granulation
- Heavy exudate
- Malodorous
- 70% granulating, 30% slough
- Wound dimensions (length, width, depth):
5.0×4.5×3.0cm



**Continue
treatment**
(add
analgesia)

Local infection: Pressure ulcer (trochanter)^{5,6}



Positive signs of progression to healing

DAY 11

Highly exudating,
malodorous, pain

90% granulation,
10% slough

Wound dimensions
(length, width, depth):
4.8×4.5×3.0cm



DAY 26

Highly exudating, odour
improved, less painful

Wound dimensions
(length, width, depth):
3.5×3.5×2.0cm



DAY 47

Wound dimensions
(length, width, depth):
3.0×3.0×1.5cm





Haematoma (lower leg)^{5,6}

PRESENTATION

- 61-year-old female, traumatic wound, lateral gaiter area (left leg) 5 days post injury. Referred to acute facility with suspected compartment syndrome and Fournier gangrene
- **Relevant medical history:**
Diabetes mellitus, lymphoedema
- **Wound dimensions (length, width, depth):**
7.0×6.0×0.4cm




11 DAYS AFTER PRESENTATION



- **Wound dimensions after debridement (length, width):**
20×10cm

Complex (mixed) infection: Haematoma (lower leg)^{5,6}



Treatment approach **before** formal referral to the infection management pathway

TREATMENT (DAY 6)

Daily enzymatic debridement with antimicrobial dressing

Wound dimensions
(length, width):
25x17cm



DAY 12

Continued daily enzymatic debridement with antimicrobial dressing

Wound dimensions
(length, width):
24x15cm




DAY 13

Continued daily enzymatic debridement with antimicrobial dressing

Wound dimensions
(length, width): 24x15cm



Complex (mixed) infection: Haematoma (lower leg)^{5,6}



REFER INFECTION MANAGEMENT PATHWAY


- Erythema
- Warmth
- Oedema/swelling
- Purulent discharge, increasing malodour, friable tissue
- Pain
- Delayed wound healing



Local
wound
infection

DAY 15

- **ACTICOAT[®] FLEX 3 Dressing**
- **ALLEVYN[®] LIFE Foam Dressing**
- 25% granulating, 75% slough
- Wound dimensions (length, width): 24×15cm




DAY 21

- 30% granulation, 40% slough, 30% necrotic tissue
- Wound dimensions (length, width): 24×15cm



Improvement
in local
infection,
wound healing
stalled

Complex (mixed) infection: Haematoma (lower leg)^{5,6}



REFER INFECTION MANAGEMENT PATHWAY

- Delayed wound healing
- Serous drainage with concurrent inflammation
- Friable hypergranulation
- Wound breakdown and enlargement
- Increasing pain and malodour



**Suspect
biofilm**

DAY 21


- **Aggressive debridement**
- **IODOSORB[®] Powder, ALLEVYN LIFE Foam Dressing, compression**
- Wound dimensions (length, width): 24x15cm



DAY 38 (ASSESSED AT DAYS 28 & 35)

- Wound dimensions (length, width): 21x13cm
(Day 28: 23.5x14cm;
Day 35: 22x13.5cm)

Day 35




Day 38



**Positive
signs of
progression
to healing**



Complex (mixed) infection: Haematoma (lower leg)^{5,6}



Step down, step up
(transitions between multiple
therapies and standard of care)⁴

DAY 43 – RE-ASSESSMENT

- Signs and symptoms of biofilm resolved
- **STEP UP** treatment with **PICO[®] Single Use Negative Pressure Wound Therapy System** (2 weeks)
- 80% granulation, 20% slough
- Wound dimensions (length, width): 20x12cm



DAY 64

- **STEP DOWN** treatment to **standard wound care, foam dressings**
- 95% granulation tissue, 5% yellow fibrin
- Wound dimensions (length, width): 18x10cm





Guided differentiation
between local infection
and biofilm



Simplified
dressing selection



**Facilitated
communication**
between clinicians



Improved understanding of the
different management approaches to
local infection and **biofilm**



Eased
decision
making



1. Dowsett C, Bellingeri A, Carville K, Garten A, Woo K. A route to more effective infection management: the infection management pathway. *Wounds Int.* 2020;11(3):20–27.
2. Moore Z, Dowsett C, Smith G, et al. TIME CDST: an updated tool to address the current challenges in wound care. *J Wound Care.* 2019;28(3):154-161.
3. Ayello EA, Carville K, Fletcher J, et al. Appropriate use of silver dressings in wounds. An expert working group consensus. *Wounds International.* 2012. Available from: <https://www.woundsinternational.com/> Last accessed 18 January 2021.
4. Schultz G, Bjarnsholt T, James GA, et al. Consensus guidelines for the identification and treatment of biofilms in chronic nonhealing wounds. *Wound Repair Regen.* 2017;25(5):744-757.
5. Woo K. Implementing the new infection management pathway to optimise outcomes: real-world case series. *Wounds International.* 2020;11(4):50–57.
6. Woo K. Implementing the infection management pathway: case series review. AWM-AWC-27378. September 2020. Presentation at the S+N Webinar: Practical hints and tips for understanding and managing infected chronic wounds. 21 October 2020. Available at: <https://www.smith-nephew.com/education>.

Developed by Evidence Communications,
Global Clinical & Medical Affairs
www.smith-nephew.com/education

Smith+Nephew does not provide medical advice. It is the responsibility of the treating physician to determine and utilise the appropriate products and techniques according to their own clinical judgment for each of their patients.

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

Smith+Nephew

♦Trademark of Smith+Nephew
All Trademarks acknowledged
©March 2021 Smith+Nephew
AWM-AWC-28965

Advanced Wound Management
Smith & Nephew Medical Ltd
101 Hessle Road
Hull HU3 2BN, UK

T +44 (0) 1482 225181
F +44 (0) 1482 328326
smith-nephew.com
GMC1314