

Smith+Nephew

ENGAGE[◇]
Cementless Partial Knee System

Clinical and Economic
Benefits Summary



Executive summary

Only partial knee available in the US that is cleared for cementless use

Only system that creates compression at the tibial interface using a patented 4th component, the Tibial Anchor. This limits micromotion and interface stability to promote osteointegration.

Eliminating cement from the procedure has potential clinical and economic benefits.

Removing cementing steps results in up to 15 minutes less time in the OR.

Cost saving due to not having to purchase cement or accessories.

Eliminates the risk of tibial loosening due to cement failure.

Reduces the incidence of third-body wear due to loose cement debris in the joint.



Clinical summary

Eliminating cement from the unicompartmental knee replacement (UKA) procedure presents potential clinical benefits

- It is estimated that up to 1 in 5 traditional UKAs have excess or loose cement floating in the joint¹
- Loose cement particles present in UKA leads to increases in polyethylene wear rates by 10X²
- Polyethylene wear debris causes osteolysis. Aseptic loosening resulting from osteolysis is the most common failure mode of joint arthroplasty³⁻⁶

Reduction in revision rates are possible with cementless partial knee replacements

- Cementless UKA have a lower risk of post-operative loosening compared to those using cement fixation^{7,8}
- The risk of revision related to aseptic loosening is cut in half for cementless partial knee replacements compared to cemented UKA constructs at 10 years⁷

Possible reduction in surgical time (up to 15 min) by eliminating cementing steps results in potential benefits for the patient⁹⁻¹²

- Reduced chance of infection¹³
- Fewer tourniquet related complications¹⁴

Cementless UKAs have less chance of loosening than cementless total knee replacements (TKA)

- Implant aseptic loosening failures occurred nearly 50% less in cementless UKA than TKA⁶



Product performance summary

Measured improvement in product performance is demonstrated for the ENGAGE[◇] Cementless Partial Knee System over traditional keel-based tibial trays (Figure 1)




Figure 1. Comparison of product performance of an Engage tibial construct with tibial anchor to a traditional tibial tray utilizing keel fixation.

The only partial knee system that produces implant compression at the tibial tray-to-bone interface.^{15,20}

Economic analysis

The positive financial impacts are significant; and result in direct and indirect cost savings per surgery to hospital or outpatient surgery center.

Direct savings

Significant reduction in OR time

- Clinical studies for cementless UKA demonstrated up to 15 minutes saving per case⁹⁻¹²
- Literature estimates the average cost for OR time including staff and instrumentation at approximately \$100 per min²¹
- Possible Direct Saving = \$100 / min x 15min = > \$1,500 per surgery

Eliminate cost of bone cement

- Average price for bone cement with antibiotics = \$264 per packet²²
- Possible Direct Savings = \$264 per surgery

Indirect savings

UKA may result in lower unplanned patient re-admissions reducing CMS-related penalties

- 30 day re-admission rates were cut in half for UKA (2.2%) when compared to TKA (4%)²³
- >67% reduction in re-admission risk for UKA as compared to TKA²⁴
- Knee Society Functional Score has been shown to be higher in cementless UKA vs. cemented UKA at 5 years (p=0.003)⁹
- Cementless fixation avoids technical errors related to cementing, like inadequate cementation technique, presence of loose fragments, or excess cement causing impingement²⁵
- Registry data indicates the risk of revision related to aseptic loosening is cut in half for cementless UKA compared to cemented UKA constructs after 10 years⁷

Item	Potential savings per case (\$US)
Reduced OR Time	\$1500
Eliminate Cement	\$264
Total Savings	\$1764





Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area.

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, Tennessee 38116
USA

www.smith-nephew.com

®Trademark of Smith+Nephew
All Trademarks acknowledged
©2022 Smith+Nephew
36637 V1 09/22

References

1. Hauptmann SM, Weber P, Glaser P, et al. Free Bone cement fragments after minimally invasive unicompartmental knee arthroplasty: an underappreciated problem. *Knee Surg Sports Traumatol Arthrosc.* 2008;16(8):770-775.
2. Paulus AC, Franke PM, Kraxenberger M, et al. PMMA Third-Body Wear after Unicompartmental Knee Arthroplasty Decouples the UHMWPE Wear Particle Generation in Vitro. *Biomed Res Int.* 2015.
3. HG, Semlitsch M. Reactions of the articular capsule to wear products of artificial joint prostheses. *J Biomed Mater Res.* 1977;11(2):157-164.
4. Archibeck MJ, Jacobs JJ, Roebuck KA, Glant TT. The Basic Science of Periprosthetic Osteolysis. *Instr Course Lect.* 2000;82-A(10):1478-1489.
5. Hallab NJ, Jacobs JJ. Biologic Effects of Implant Debris. *Bull NYU Hosp Jt Dis.* 2009;67(2):182-188.
6. Van der List JP, Sheng DL, Kleeblad LJ, Chawla H, Pearle AD. Outcome of cementless unicompartmental and total knee arthroplasty: A systematic review. *Knee.* 2017;24(3):497-507.
7. Mohammad H, Matharu G, Judge A, Murray D. Comparison of the 10-year outcomes of cemented and cementless unicompartmental knee replacements: data from the National Joint Registry for England, Wales, Northern Ireland and the Isle of Man. *Acta Orthop.* 2020;91(1):76-81.
8. Murray D, Mohammad H, Matharu G, Mellon S, Judge A. A Comparison of The Outcomes of Cemented and Cementless Oxford Unicompartmental Knee Arthroplasty: A Propensity-Matched Cohort Study of 10 836 Knees. *Orthop Proceedings.* 2018.
9. Pandit H, Liddle AD, Kendrick BJL, et al. Improved Fixation in Cementless Unicompartmental Knee Replacement. *J Bone Joint Surg.* 2013;95(15):1365-1372.
10. Stempin R, Kaczmarek W, Stempin K, Dutka J. Midterm Results of Cementless and Cemented Unicompartmental Knee Arthroplasty with Mobile Meniscal Bearing: A Prospective Cohort Study. *Open Orthop J.* 2017;11:1173-1178.
11. Kerens B, Schotanus MGM, Boonen B, et al. Cementless versus cemented Oxford unicompartmental knee arthroplasty: early results of a non-designer user group. *Knee Surg Sports Traumatol Arthrosc.* 2017;25(3):703-709.
12. Akan B, Karagüven D, Güçlü B, et al. Cemented versus uncemented Oxford unicompartmental knee arthroplasty: is there a difference?. *Adv orthop.* 2013;2013.
13. Jämsen E, Varonen M, Huhtala H, et al. Incidence of Prosthetic Joint Infections After Primary Knee Arthroplasty. *J Arthroplasty.* 2010;25(1):87-92.
14. Olivecrona C, J Lasse, Benson La, Benson Li, Blomfeldt R. Tourniquet time affects postoperative complications after knee arthroplasty. *International Orthopaedics (SICOT).* 2013;37:827-832.
15. Engage 2019. Anchor Compression Test Report. 101-09912-004-01.
16. Engage 2018. Anchor Fixation Report. 101-09912-001-01.
17. ENGAGE 2019. Tibia Tray Holding Power on ENGAGE's Unicompartmental Knee System, Pullout Testing. 10-3182.
18. Engage 2018. Static Tibial Cadaver Strength Report. 101-09912-009-01.
19. Engage 2019. Engage Unicompartmental Knee System Biomechanical Testing. 1906527.000_0138.
20. Smith+Nephew 2022. 35898 V1 71282199 REVO ENGAGE Cementless Partial Knee Surgical Technique 0822.
21. Ting N, Moric M, Della Valle C, Levine B. Use of Knotless Suture for Closure of Total Hip and Knee Arthroplasties. A Prospective, Randomized Clinical Trial. *J Arthroplasty.* 2012;27(10):1783-1788.
22. Orthopedic Network News. 2019 Hip and Knee Implant Review. 2019;30(3).
23. Drager J, Hart A, Khalik JA, et al. Shorter Hospital Stay and Lower 30-Day Readmission After Unicompartmental Knee Arthroplasty Compared to Total Knee Arthroplasty. *J Arthroplasty.* 2016;31(2):356-361.
24. Hansen EN, Ong KL, Lau E, Kurtz SM, Lonner JH. Unicompartmental Knee Arthroplasty Has Fewer Complications but Higher Revision Rates Than Total Knee Arthroplasty in a Study of Large United States Databases. *J Arthroplasty.* 2019;34(8):1617-1625.
25. Campi S, Pandit HG, Dodd CAF, Murray DW. Cementless fixation in medial unicompartmental knee arthroplasty: a systematic review. *Knee Surg Sports Traumatol Arthrosc.* 2017;25(3):736-745.