OXINIUM-Technologie in der Knieendoprothetik **Smith+Nephew**

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OXINIUM[◇]
Oxidiertes Zirkonium

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Literaturangaben

1. Klug A, Gramlich Y, Rudert M, et al. The projected volume of primary and revision total knee arthroplasty will place an immense burden on future health care systems over the next 30 years. *Knee Surgery, Sports Traumatology, Arthroscopy*. 2020;15:1-12. **2.** Sloan M, Premkumar A, Sheth NP. Projected Volume of Primary Total Joint Arthroplasty in the U.S., 2014 to 2030. *J Bone Joint Surg Am*. 2018;100:1455-1460. **3.** Ackerman IN, Bohensky MA, Zomer E, et al. The projected burden of primary total knee and hip replacement for osteoarthritis in Australia to the year 2030. *BMC Musculoskeletal Disorders*. 2019;23(1):90. **4.** Davis ET, Pagkalos J, Kopjar B. Bearing surface and survival of cementless and hybrid total hip arthroplasty in the National Joint Registry of England, Wales, Northern Ireland and the Isle of Man. *Journal of Bone Joint Surgery*. 2020;5(2):pe0075. **5.** Peters RM, Van Steenberghe LN, Stevens M, et al. The effect of bearing type on the outcome of total hip arthroplasty. *Acta Orthopaedica*. 2018; 89(2):163-169. **6.** Atrey A, Ancarani C, Fitch D, Bordini B. Impact of bearing couple on long-term component survivorship for primary cementless total hip replacement in a large arthroplasty registry. Poster presented at: Canadian Orthopedic Association; June 20–23, 2018; Victoria, British Columbia, Canada. **7.** Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR) Hip, Knee & Shoulder Arthroplasty: 2020 Annual Report. **8.** Innocenti M, Matassi F, Carulli C, Nistri L, Civinini C. Oxidized zirconium femoral component for TKA: A follow-up note of a previous report at a minimum of 10 years. *The Knee*. 2014;21:858–861. **9.** Fernandez-Fernandez R, Cruz-Pardos A, Garcia-Rey E. Revision Total Hip Arthroplasty: Epidemiology and Causes. In: Rodríguez-Merchán E. *Revision Total Joint Arthroplasty*. Springer, 2020. **10.** Lewis PL, Robertsson O, Graves SE, et al. Variation and trends in reasons for knee replacement revision: a multi-registry study of revision burden. *Acta Orthopaedica*. 2020. DOI: 10.1080/17453674.2020.1853340. **11.** Hunter G, Dickinson J, Herb B, et al. Creation of oxidized zirconium orthopaedic implants. *Journal of ASTM International*. 2005;2:1-14. **12.** Long M, Riester L, Hunter G. Nano-hardness Measurements of Oxidized Zr-2.5Nb and Various Orthopaedic Materials. Abstract presented at: 24th Annual Meeting of the Society for Biomaterials. April 22-26, 1998, San Diego, California. **13.** Parikh A, Hill P, Hines G, Pawar V. Wear of conventional and highly crosslinked polyethylene liners during simulated fast walking/jogging. Poster presented at: 55th Annual Meeting of the Orthopaedic Research Society, 2009. Poster no. 2340. **14.** Parikh A, Hill P, Pawar V, Sprague J. Long-term Simulator Wear Performance of an Advanced Bearing Technology for THA. Poster presented at: 2013 Annual Meeting of the Orthopaedic Research Society. Poster no. 1028. **15.** Papannagari R, Hines G, Sprague J, Morrison M. Long-term wear performance of an advanced bearing technology for TKA. Poster presented at: 2011 Annual Meeting of the Orthopaedic Research Society. Poster no. 1141. **16.** Aldinger P, Williams T, Woodard E. Accelerated Fretting Corrosion Testing of Zirconia Toughened Alumina Composite Ceramic and a New Composition of Ceramicised Metal Femoral Heads. Poster presented at: 2017 Annual Meeting of the Orthopaedic Research Society. Poster no. 1037. **17.** Smith+Nephew 2016. OR-16-127. **18.** 2005 ASM International Engineering Materials Achievement Award. **19.** Dalal A, Pawar V, McAllister K, Weaver C, Hallab NJ. Orthopaedic implant cobalt-alloy particles produce greater toxicity and inflammatory cytokines than titanium alloy and zirconium alloy-based particles in vitro, in human osteoblasts, fibroblasts, and macrophages. *J Biomed Mater Res Part A*. 2012;100A:2147-2158. **20.** Pawar V, Jones B, Sprague J, Salehi A, Hunter G. Acidic Fretting Tests of Oxidized Zr-2.5Nb, CoCr and SS Femoral Heads. Paper presented at: Materials & Processes for Medical Devices Conference, ASMI, St. Paul, Minnesota, 2004. **21.** Mayman DJ, Patel AR, Carroll KM. Hospital related clinical and economic outcomes of a bicruciate knee system in total knee arthroplasty patients. Poster presented at: ISPOR Symposium; May 19-23, 2018; Baltimore, Maryland, USA. **22.** Patrick C, Delhougne G, Patel A. Retrospective Analysis of Oxidized Zirconium Bearing Surface in Hip Replacement 90-Day Episode Claims. ISPOR EU Poster, 2019. **23.** Civinini R, Carulli C, Matassi F, et al. The Survival of Total Knee Arthroplasty: Current Data from Registries on Tribology. *HSS Journal*. 2017;13:28-31. **24.** Civinini R, Matassi F, Carulli C, et al. Clinical Results of Oxidized Zirconium Femoral Component in TKA. A Review of Long-Term Survival. *HSS Journal*. 2017;13:32-34. **25.** Patel AR, Delhougne G. Hospital related clinical and economic outcomes of two premium knee system in total knee arthroplasty patients. Poster presented at: ISPOR Symposium; May 22, 2019; New Orleans, LA, USA. **26.** ASTM International Standard Specification for Wrought Zirconium-2.5 Niobium Alloy for Surgical Implant Applications (UNS R60901) Designation: F 2384 – 10. **27.** National Joint Registry for England, Wales and Northern Ireland: LEGION revision OXINIUM (with revision tibial) implant summary report. 18 May 2022. Available at: <https://bit.ly/3w8TEI9>. **28.** Hunter G, Long M. Abrasive wear of oxidized Zr-2.5 Nb, CoCrMo, and Ti-6 Al-4 V against bone cement. Abstract presented at: Sixth World Biomaterials Congress, 2000. **29.** ASTM International Standard Specification for Cobalt-28 Chromium-6 Molybdenum Alloy Castings and Casting Alloy for Surgical Implants (UNS R30075): Designation: F 75 – 12. **30.** Glover AW, Santini AJA, Davidson JS, Pope JA. Mid- to long-term survivorship of oxidized zirconium total knee replacements performed in patients under 50 years of age. *The Knee*. 2018;25:617-622. **31.** Ahmed I, Salmon LJ, Waller A, et al. Total knee arthroplasty with an oxidized zirconium femoral component: Ten Year Survivorship Analysis. *Bone Joint J*. 2016;98-B:58–64. **32.** The Orthopaedic Data Evaluation Panel (ODEP). www.odep.org.uk. Accessed April 12th, 2021. **33.** Lee KY, Slavinsky JP, Ries MD, Blumenkrantz G, Majumdar S. Magnetic Resonance Imaging of In Vivo Kinematics After Total Knee Arthroplasty. *Journal of Resonance Imaging*. 2005;21:172–178. **34.** Hallab N, Merritt K, Jacobs JJ. Metal Sensitivity in Patients with Orthopaedic Implants. *JBJS*. 2001;83(3):428-436.

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©2021 Smith & Nephew, Inc.
Smith & Nephew GmbH
Smith & Nephew Schweiz AG