 **smith&nephew**

**LEGION<sup>◇</sup>**

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## Bibliografia

1. Rand J: Augmentation of a total knee arthroplasty with a modular metal wedge. *Journal of Bone and Joint Surgery*. 77-A-2: 266-268, February, 1995.
2. Haas SD; Insall, JN; Montgomery III W; Windsor R E: Revision total knee arthroplasty with use of modular components with stems inserted without cement. *Journal of Bone and Joint Surgery*, 77-A-11: 1700-1707, December, 1995.
3. Bugbee WD, Ammeen DJ, Engh GA: Does implant selection affect outcome of revision knee arthroplasty? *The Journal of Arthroplasty*. Volume 16, Issue 5, p. 581-585, August, 2001.
4. Jones RE: Management of complex revision problems with a modular total knee system. *Orthopedics*. 19-9: 802-804, September, 1996.
5. Sculco TP: Augmentation options in total knee replacement. *Orthopedics*. 19-9: 800-801, September, 1996.
6. Smith & Nephew Research report. OR-16-008.
7. Data on file Smith & Nephew.
8. J.D. Bobyn, S.A. Hacking, S.P. Chan, K.K. Toh, J.J. Krygier and M. Tanzer, "Characterization of a new porous tantalum biomaterial for reconstructive
9. Koval, Kummer: Internal Smith & Nephew documents and AAOS 2004.
10. Gonzalez, et al.: Radius der curvatur des femur in der sagittalebene. *Gegebenaurs morphol. Jahrb – Leipzig*. 135 (1989) 5, S. 761-764.
11. Stephenson, Seedhom: Modeling femoral curvature in the sagittal plane: A cadaveric study. *Proc Instn Mech Engrs*. Volume 215(2), Part H, p. 221-228, 2001.
12. Harper, Carson: Curvature of the femur and the proximal entry point for an intramedullary rod. *Clinical Orthopaedics and Related Research*. Number 220, p. 155-161, July, 1987.
13. Noble, et al.: The effect of aging on the shape of the proximal femur. *Clinical Orthopaedics and Related Research*. Number 316, p. 31-44.
14. Zuber, et al.: Form und dimension der markhohle menschlicher femora in hinblick suf die passing von marknagelimplantaten. *Unfallchirurg*. Volume 91(7), p. 314-319.
15. Onoue, et al.: Treatment of the femoral shaft fracture with a curved heat-treated COP clover-leaf nail. *International Orthopedics (SICOT)*, Volume 3, p. 203-210, 1979.
16. Trudell MB, et al.: Anterior femoral curvature revisited: Race assessment from the femur. *Journal of Forensic Sciences*. Volume 44 (4) p. 700-707, July, 1999.
17. Walensky NA: A study of anterior femoral curvature in man. *Anat Rec*. Volume 151, p. 559-570, April, 1965.
18. Long M, Riestler L, and Hunter GB: Nano-hardness measurements of oxidized Zr-2.5Nb and various orthopaedic materials. *Trans. Society for Biomaterials*, Minneapolis, MN, p. 528, 1998.
19. Poggie RA, Wert JJ, Mishra AK, and Davidson JA: Friction and wear characteristics of UHMWPE in reciprocating sliding contact with Co-Cr, Ti-6Al-4V, and zirconia implant bearing surfaces. In *Wear and Friction of Elastomers, ASTM STP 1145*, Philadelphia, PA, 1992.
20. Hunter GB, and Long M: Abrasive wear of oxidized Zr-2.5Nb, CoCrMo, and Ti-6Al-4V against bone cement. 6th World Biomaterials Congress. *Trans. Society for Biomaterials*, Minneapolis, MN, p. 835, 2000.
21. Spector M, Ries MD, Bourne RB, Sauer W, Long M, and Hunter GB: UHMWPE wear performance of oxidized zirconium total knee femoral components. Scientific exhibit, AAOS, 2001.
22. Ries MD, Salehi A, Widding K, Hunter GB: Polyethylene wear performance of oxidized zirconium and cobalt chrome knee components under abraded conditions. *The Journal of Bone & Joint Surgery*. Volume 83-A, Number 2, p. 129.
23. Hallab N, Meritt K, and Jacobs JJ: Metal sensitivity in patients with orthopaedic implants. *The Journal of Bone & Joint Surgery*. Volume 83-A, Number 3, p. 428, March 1, 2001.
24. ASTM International Standard Specification for Wrought Zirconium-2.5Niobium Alloy for Surgical Implant Applications (UNS R60901) Designation: F 2384 – 05 and Standard Specification for Cobalt-28 Chromium-6 Molybdenum Alloy Castings and Casting Alloy for Surgical Implants (UNS R30075): Designation: F 75 – 07
25. Westrich GH, Haas, SB, Insall JN, Franchie A: Resection specimen analysis of proximal tibial anatomy based upon one hundred total knee replacement specimens. *Journal of Arthroplasty*. 10:1, 1995.
26. Yoshioka Y, Siu D, Scudamore RA, Cooke TDV: Tibial anatomy and functional axes. *Journal of Orthopaedic Research*. 7, 1:132-137, 1989.
27. Stulberg BN, Dombrowski RM, Froimson M, Easley K: Computer tomography analysis of proximal tibial coverage. *Clinical Orthopaedics*. 311:148-156, 1995.
28. Smith JR, Hofmann AA: Morphology of the proximal tibia in the arthritic knee. Scientific exhibit, AAOS, 1992.
29. Engh, et al.: Tibial baseplate wear: A major source of debris with contemporary modular knee implants. Scientific exhibit, AAOS, 2000.
30. Wevers HW, Simurda M, Griffin M, arrel J: Improved fit y asymmetric tibial prosthesis for total knee arthroplasty. *Medical Engineering Physics*. 16:297-300, 1994.
31. Hitt K, Shurman JR, Greene K, McCarthy J, Moskal J, Hoeman T, Mont MA: Anthropometric measurements of the human knee: Correlation to the sizing of current knee arthroplasty systems. *The Journal of Bone and Joint Surgery (A)*. Volume 85, p. 115-122, 2003.
32. Bargren JH, Day WH, Freeman MAR, Swanson SAV: Mechanical properties of four nonhinged knee prostheses. Proceedings of the 24th Annual Meeting of ORS. Dallas, TX, February, 1978.
33. Rielly D, Walker, PS, Ben-Dov M, Ewald FC: Effects of tibial components on load transfer in the upper tibia. *Clinical Orthopaedics and Related Research*. 165:273-282, 1982.
34. Walker PS, Thatcher J, Ewald FC, Milden J: Variables affecting the fi ation of tibial components. *Engineering in Medicine*. 11, 2:83-87, 1982.
35. Figgie H, Davy D, Heiple K, Hart R: Load-bearing capacity of the tibial component of the total condylar knee prosthesis. *Clinical Orthopaedics and Related Research*. 183:288-297, 1984.
36. Courtland G, Belniak R, Plowman M, Hopfer S, Knight J, Sunderman W: Intraarticular carcinogenesis bioassays of CoCrMo and Ti-6Al-4V alloys in rats. *The Journal of Arthroplasty*. Volume 10 (1). p. 75-80, February 1995.
37. Condit MA, Ismaili SK, Alexander JW, Noble PC: Backside wear of modular ultra-high molecular weight polyethylene tibial inserts. *The Journal of Bone and Joint Surgery (A)*. Volume 86, p. 1031-1037, 2004.
38. Cooper M, Holbrook J: Fatigue and s tatic testing of the GENESIS II posterior stabilized tibial insert locking mechanism. Smith & Nephew Richards Inc., Technical Report, OR-95-50, April, 1995.
39. Daigle K, Thomas R: Medial-lateral fatigue testing of the GENESIS II constrained tibial insert. Smith & Nephew Richards Inc., Technical Report, OR-96-45, May, 1996.
40. Daigle K, Thomas R: Posterior/anterior fatigue and static testing of the GENESIS II constrained tibial insert locking mechanism. Smith & Nephew Richards Inc., Technical Report, OR-96-48, May, 1996.
41. Condit MA, Thompson MT, Usrey MM, Ismaili SK, Noble PC: Backside wear of polyethylene tibial inserts: Mechanism and magnitude of material loss. *The Journal of Bone and Joint Surgery (A)*. Volume 87, p. 326-331, 2005.
42. Wasielewski RC, et al.: Tibial insert undersurface as a contributing source of polyethylene wear debris. *Clinical Orthopaedics and Related Research*. Number 345, p. 53-59, 1994.

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©2021 Smith & Nephew.
Smith & Nephew, Inc.
Smith & Nephew S.r.l.