

TRIGEN[®]
META-NAIL[®]

Semi-extended Instrument Set



Table of contents

Introduction	1
TRIGEN° META-NAIL° Tibial Nail Specifications	2
Surgical Technique	3
Patient positioning	3
Prepare the Entry Tube.....	3
Incision and entry point: Lateral subluxation of the patella	4
Entry portal acquisition.....	5
Entry portal	5
Alternative technique: Entry portal.....	6
Reduce the fracture.....	7
Determine the implant length	8
Unreamed technique	9
Reamed technique.....	9
Nail assembly.....	10
Nail insertion	11
Check nail depth	12
Lock the screw	13
Proximal locking	14
Proximal locking: Compression	15
Nail extraction: optional	16
An alternative method for extraction	17
Catalog information	18

Nota Bene

The following technique is for informational and educational purposes only. It is not intended to serve as medical advice. It is the responsibility of treating physicians to determine and utilize the appropriate products and techniques, according to their own clinical judgment, for each of their patients. For more information on the product, including its Indications for Use, contraindications, cleaning, sterilization and product safety information, please refer to the product's label and the Instructions for Use (IFU) for the product.


Introduction

The TRIGEN® META-NAIL® Tibial Nail System is designed specifically to repair fractures of the proximal and distal thirds of the tibia, including the shaft, stable and unstable fractures, nonunions, malunions, and for the prophylactic nailing of impending pathological fractures.

The META-NAIL System can be used with or without the SURESHOT® Distal Targeting System. If using the SURESHOT Distal Targeting System, be sure to read and understand the TRIGEN SURESHOT Distal Targeting System User Manual. Only trained operators are allowed to use the TRIGEN SURESHOT Distal Targeting System.



TRIGEN[◇] META-NAIL[◇] Tibial Nail Specifications



M/L view
Figure 1

Specifications	TRIGEN META-NAIL Tibia (8.5mm)	TRIGEN META-NAIL Tibia
Material	Ti6Al4V	Ti6Al4V
Diameter	8.5mm	10, 11.5 & 13mm
Lengths	16-50cm*	16-50cm*
Nail Color	Grey	Gold
Cross Section	Round	Round
Proximal Diameter (driving end)	12mm	12mm (10, 11.5 diameter) 13mm (13 diameter)
Distal Diameter (non-driving end)	8.5mm	10, 11.5 and 13mm
Smallest Through Diameter	4.8mm	5.0mm
Wall Thickness	1.9mm	2.3mm (10) 3.0mm (11.5) 2.3mm (13)
Guide Bolt Thread	5/16-24 UNF	5/16-24 UNF
Screw Diameter	4.5mm	5.0mm
Screw Color	Grey	Gold
Major Diameter	4.5mm	5.0mm
Minor Diameter (core)	4.0mm	4.3mm
Screw Lengths	25-65mm	25-110mm
Hex Size	4.7mm	4.7mm
Alternative Hexdrivers	RT Femoral & Recon 7.0mm Cannulated Screw PERI-LOC [®] Locking Screw	RT Femoral & Recon 7.0mm Cannulated Screw PERI-LOC 4.7mm Hexdriver PROFIX [®] 4.7mm Hexdriver
Alternative Modes	No	No
Proximal Locking (Driving End)		
Static Lock Locations/Orientations	17mm/45° Screw Locked w/META-NAIL Cap 23mm/45° Threaded w/bushing 30mm/25° Threaded 40mm/25° Threaded	
Static Locking Hole Dimensions	Threaded 4.3mm minor diameter Threaded 4.7mm major diameter	Threaded 4.5mm minor diameter Threaded 5.3mm major diameter
Compression/Dynamic Slot Location	10mm	10mm
Compression/Dynamic Slot Diameter/Length	4.7mm/7mm	5.3mm/7mm
Degree of Proximal Bend (Herzog)	10°	10°
Proximal Bend Location	27mm	27mm
Distal Locking (Non-Driving End)		
Static Lock Locations/Orientations	25mm/M/L 15mm/A/P 5mm/M/L (Threaded)	25mm/M/L 15mm/A/P 5mm/M/L (Threaded)
Static Locking Hole Dimensions	4.7mm**	5.3mm**
Degree of Distal Bend	2°	2°
Distal Bend Location	60mm	60mm

Note: These views are not to scale and should be used as a pictorial representation only.

*Set does not include all sizes.

**Most distal hole threaded.

Surgical Technique

Patient positioning

Position the patient supine on a radiolucent table with the unaffected limb extended away from the affected limb.

For the semi-extended technique, the operative limb should be in 10–20° of flexion.

Caution: If using the TRIGEN® SURESHOT® Distal Targeting System, verify that there are no metal objects in the immediate targeting area. Metal interference will cause the SURESHOT System to be inaccurate.




Figure 2

Prepare the Entry Tube

1. Inspect the Entry Tube to ensure that it is not damaged, bent or chipped. Any flaws in the tube can lead to damage of the surrounding tendons and tissues.
2. Assemble the Entry Tube (71654526) to the Entry Portal Handle (71674092) by pulling back the black trigger of the Entry Portal Handle and inserting the Entry Tube into the handle ring.
3. Turn the tube until it clicks into a locked position.
4. Place the Guide Pin Sleeve (71654527) through the tube until it locks into both the end of the sleeve and Entry Portal Handle ring.

Optional: Suction can be applied to the Entry Portal Handle.




Figure 3

Incision and entry point: Lateral subluxation of the patella

1. Make a midline skin incision approximately 1.5 inches in length from the upper pole of the patella to the middle of the patella (Figure 4). Make a second deep incision medial to the patella, cutting the superior two-thirds of the medial retinaculum, but leaving the cuff intact.
2. Extend the incision 1–2cm into the quadriceps tendon (Figure 5). This incision is through the medial one-third of the quadriceps tendon as outlined in the image.
3. Sublux the patella laterally and pass the fully assembled entry tube through the skin incision until it is sitting on the tibial plateau (Figure 6).



Figure 4




Figure 5



Figure 6

Entry portal acquisition

1. Use the femoral condyles as guides for proper tube placement, making sure that the Guide Pin Sleeve sits at the appropriate entry point (Figure 7).
2. Locate the entry point just medial to the lateral tibial eminence in the (anteroposterior) A/P view, and in line with the anterior cortex and intramedullary canal in the lateral view (Figure 8).
3. Attach the 3.2mm Brad Point Guide Pin (71674130 or 71631436) to the drill via the Mini Connector (71631186) and insert the guide pin into the Guide Pin Sleeve (Figure 9).



Figure 7

Note: If the guide pin insertion is suboptimal, insert the Semi-extended Honeycomb (71654533) and rotate the Honeycomb within the Entry Tube to the desired location. Insert another 3.2mm Brad Point Guide Pin. Avoid over-insertion of the guide pin as this can establish a false trajectory and subsequently cause fracture malalignment.




Figure 8




Figure 9

Entry portal

1. After guide pin placement, remove the Guide Pin Sleeve and carefully slide the Entry Tube down and onto the tibial plateau (Figure 10).
2. Attach the 12.5mm Entry Reamer (71631116) to a power drill. Advance the reamer over the guide pin and through the Entry Tube into the tibia to a depth of 4–6cm. Maintain alignment to avoid penetration of the posterior cortex.
3. Check position of the entry reamer via radiographic imaging and then remove the Entry Reamer and Guide Pin.



Figure 10

Alternative technique: Entry portal

1. With the Entry Tube (Figure 11) and Entry Portal Handle (Figure 12) still in the joint space, attach the T-handle (71674076 or 71674576) (Figure 13) to the Cannulated Awl (Figure 14).

Note: Introduce the 3.2mm T-handle Trocar (71674074) (Figure 15) into the back of the assembly prior to insertion in order to prevent awl slippage and the accumulation of cortical bone within the cannulation.

2. Insert the awl into the proximal tibia to a depth of 4–6cm.



Entry Tube

Cat. No. 71654526

Figure 11



Entry Portal Handle

Cat. No. 71674092

Figure 12



T-handle

Cat. No. 71674076 or 71674576

Figure 13



Cannulated Awl

Cat. No. 71674000

Figure 14



3.2mm T-handle Trocar

Cat. No. 71674074

Figure 15

Reduce the Fracture

1. Open the Gripper (71674080).
2. Insert the smooth back end of the 3.0mm Ball Tip Guide (71631626) into the front of the Gripper and gently close the trigger grip (Figure 16).
3. Connect the Reducer and Reducer Connector (71674077) so that the words "Slot Orientation" are in line with the opening at the Reducer's tip (Figure 17).




Figure 16



Figure 17

4. Connect the assembly to the T-handle.
5. Advance the Reducer through the Entry Tube into the intramedullary canal and use the curved tip to direct the 3.0mm Ball Tip Guide Rod past the fracture into the region of the distal epiphyseal scar (Figure 18).



Figure 18

Note: The guide rod should be center-center in the A/P and lateral views.

Reduce removal

1. Once the guide rod reaches the desired depth, detach the Gripper and remove the Reducer from the tibial canal.
2. To maintain guide rod position within the canal, slide the Obturator (71674078) into the back of the T-handle during extraction (Figure 19).




Figure 19

Determine the implant length

1. Confirm that the Ruler opens easily. If it does not, adjust the thumb-wheel connection at the end to ensure free movement.
2. After the Reducer has been removed, reconfirm the guide rod placement within the distal tibia.
3. Slide the Ruler (71674079) over the guide rod until the metal tip contacts the proximal tibia (Figure 20).
4. Confirm the guide rod position in the window at the opposite end of the Ruler as shown in order to ensure accurate implant measurement (Figure 21).
5. Push down on the top of the Ruler until it contacts the 3.0mm Ball Tip Guide Rod.

Note: Ensure that the ruler is in contact with the Ball Tip Guide Rod.

6. Read the implant length from the exposed calibrations at the end of the Ruler.
7. Confirm fracture reduction to ensure that the implant length is not underestimated. Reference the fibula for accurate fracture distraction or compression.




Figure 20



Figure 21

Unreamed technique

1. Use radiographic templating to determine nail size (Figure 22).

Note: The appropriate-diameter implant provides translational fill within the isthmus of the intramedullary canal.

2. To help avoid implant incarceration during insertion, select a nail approximately 1.0–1.5mm less than the narrowest canal measurement on the lateral radiograph.

Note: The 8.2mm diameter of the Reducer provides an initial determination of the canal width in small-diameter tibias.




Figure 22: TRIGEN® META-NAIL® Tibial Nail Radiographic Template, Cat. No. 71180810

Reamed technique

1. Use radiographic templating to determine nail size.
2. Use the 9.0mm Front Cutting Reamer Head (71118231) and Flexible Reamer Shaft (71118200) to ream the intramedullary canal sequentially in half-millimeter increments to a size 1.0–1.5mm larger than the selected nail size (Figure 23).
3. To ensure Guide Rod placement during reaming, insert the Obturator (71674078) into the back of the Reamer unit during retraction.
4. Continue to confirm Guide Rod placement in the distal tibia throughout reaming.

Note: Periodically move the reamer back and forth in the canal to clear debris from the cutting flutes.



Figure 23

Nail assembly

1. Use the Semi-extended Guide Bolt (71654525) to attach the META-NAIL[®] Semi-extended Drill Guide (71654524) to the nail.
 2. Tighten with the Guide Bolt Wrench (71631140) and T-handle. The nail is correctly aligned when:
 - The vertical black line on the posterior side of the insertion barrel aligns with the line on the posterior side of the nail.
 - The “A” on the anterior side of the nail aligns with the “A” on the anterior side of the insertion barrel (Figure 24).
 - The apex of the nail’s proximal Herzog Bend faces posteriorly and the drill guide is oriented anteriorly (Figure 25).

Note: The bevel on the front of the nail marks the connection to the drill guide and can be seen in the lateral view as a means to determine proximal insertion depth.

3. Attach the META-NAIL Anterior Drop (71654501) to the drill guide.

To verify the accuracy of the target, insert a gold 9.0mm Drill Sleeve (71631152) and silver 4.0mm Drill Sleeve (71674083) into the drop (Figure 26), and then pass a 4.0mm Long Pilot Drill (71631110) through the assembly.

Note: An incorrectly attached nail will not target.

Note: If using the SURESHOT[®] Distal Targeting System, refer to the TRIGEN[®] SURESHOT Distal Targeting System’s User Manual for the field accuracy check instructions.



Figure 24




Figure 25




Figure 26

Nail insertion

Remove the Entry Tube and Handle

1. Remove the Anterior Drop and attach the Cannulated Impactor-Medium (71675081) to the drill guide.
2. Orient the drill guide assembly in the A/P position.
3. Use light blows with the Slotted Hammer (71674082) to advance the nail over the guide rod.

Note: If excessive force is required to insert the nail, additional reaming of the intramedullary canal may be required.

4. Verify fracture reduction as the nail crosses the fracture site by paying close attention to rotation, length, alignment, distraction and/or shortening of the affected limb.
5. Check the final nail position in both the A/P and lateral views to confirm correct alignment.

Note: To obtain provisional fixation of the proximal tibial fragments, reattach the Anterior Drop and insert the gold 9.0mm Drill Sleeve (71631152) and silver 4.0mm Drill Sleeve (71674083) into one of the two Blocking Screw holes in the Drop, then pass a 4.0mm Long Pilot Drill (71631110)* through the drill sleeves.




Figure 27

*The 4.0mm Long Pilot Drill (71631110) is interchangeable with 4.0mm AO Long Drill (71631121)

Check nail depth

Check the proximal nail depth

In the lateral view, confirm nail position by observing the nail/drill guide junction (Figure 28).

Note: If compression or dynamic locking is desired, countersink the nail approximately 10mm in order to avoid implant prominence.

Check the distal nail depth

1. In the A/P and lateral views, confirm that the nail has been inserted to the desired depth (Figure 29).

Note: Optimal insertion depth, which will allow room for screw insertion below the fracture, is essential. Distal third tibial fractures require at least two locking screws to maintain stability.

2. Remove the Guide Rod once the nail is fully seated.
3. Reattach the Anterior Drop.
4. Following nail insertion, confirm that the nail and drill guide are securely connected. Hammering can loosen the Guide Bolt.



Figure 28



Figure 29

Lock the screw

Determine the length of the locking screw

There are three options available to determine the locking screw lengths:

- Gold 9.0mm Drill Sleeve, silver 4.0mm Drill Sleeve and 4.0mm Long Pilot Drill* (Figure 30A).
- Screw Depth Gauge (71631189) (Figure 30B).
- Screw Length Sleeve (71674085) and 4.0mm Short Drill (71631117)** (Figure 30C).

Note: For information regarding Drill Depth Measurement Software (DDM) using the TRIGEN[®] SURESHOT[®] Targeting System, refer to the TRIGEN SURESHOT Targeting System User Manual. The DDM software can be used for distal locking only.




Figure 30

Insert the Locking Screws

Proximal locking options include three statically locked threaded holes and one slot that allows for both fracture compression and/or dynamization. These are targeted through the orange and blue color-coded holes on the Anterior Drop.

Distal locking options include three statically locked holes: two mediolateral (M/L) and one A/P. The most distal M/L hole is threaded for additional stability.

Gold 5.0mm locking screws are compatible with 10mm, 11.5mm and 13mm diameter nails. Grey 4.5mm locking screws are compatible with 8.5mm diameter nails.

Note: Do not use the 4.0mm Short Step Drill (71641123) when drilling for a grey 4.5mm locking screw. Its diameter transitions from 4.0mm to 4.7mm and will drill too large a hole in the near cortex, which may compromise locking screw purchase.

Note: If using the SURESHOT Distal Targeting System, be sure to remove the probe prior to proximal locking.




Figure 31: META-NAIL[®] Tibial Nail with all four proximal screws inserted

*The 4.0mm Long Pilot Drill (71631110) is interchangeable with 4.0mm AO Long Drill (71631121)

**The 4.0mm Short Drill (71631117) is interchangeable with 4.0mm AO Short Drill (71631123)

Proximal locking

Static locking

1. Make a small incision at the site of screw entry.
2. Insert the gold 9.0mm Drill Sleeve and silver 4.0mm Drill Sleeve through the static slot (STAT) on the Anterior Drop (Figure 32) down to bone.
3. Drill both cortices with the 4.0mm Long Pilot Drill*.
4. Measure for screw length using either the calibrations on the 4.0mm Long Pilot Drill* or by removing the 4.0mm Drill Sleeve and using the Screw Depth Gauge.
5. Attach the appropriate-length screw to the end of the Medium Hexdriver (71631066).
6. Use power to insert the screw through the gold 9.0mm Drill Sleeve until the laser-etched ring on the Hexdriver reaches the back of the Drill Sleeve.
7. Attach the T-handle to the Hexdriver and tighten the screw by hand.

Dynamic locking

1. With the nail countersunk approximately 10mm, make a small incision at the site of screw entry and insert the gold 9.0mm Drill Sleeve and silver 4.0mm Drill Sleeve through the dynamic slot (DYN) on the Anterior Drop down to bone (Figures 33 and 34).
2. Drill both cortices with the 4.0mm Long Pilot Drill*.
3. Use the techniques described in the “Lock the Screw” section to measure and insert the screw.

Note: Use the “TIBIA” and “ALL NAILS” holes marked on the META-NAIL[®] Anterior Drop.

Note: If the screw is left in the dynamic mode, a nail cap cannot be used as it will push against the locking screw.



Figure 32



Figure 33



Figure 34

*The 4.0mm Long Pilot Drill (71631110) is interchangeable with 4.0mm AO Long Drill (71631121)

Proximal locking: Compression

There are two methods that can be used to compress the fracture:

Compression Driver method

1. Countersink the nail a minimum of 10mm (Figure 35).
2. Lock the nail distally.
3. Reduce the fracture as much as possible to maximize the advantage of the compression driver.
4. Insert a screw through the proximal side of dynamic slot as described in the “Dynamic locking” section.
5. Remove the gold 9.0mm Drill Sleeve and Medium Hexdriver.
6. Attach the Universal Compression Driver (71654528) to the T-handle and thread it through the guide bolt into the top of the nail until it contacts the most proximal 5.0mm locking screw (Figure 36).
7. Turn the Compression Driver clockwise to advance the driver and drive the locking screw distally, which will compress the fracture up to 7mm (Figure 37).

Note: If there are screws in the other proximal holes, the fracture cannot be compressed.

Note: Use fluoroscopy to check dynamic screw position when compressing the fracture.

8. Once the fracture gap is adequately compressed, lock the nail with up to three additional proximal static screws while the Anterior Drop is still attached to the drill guide.

Note: If the Compression Driver is progressed after the fracture is fully compressed, the locking screw will begin to bend. In extreme cases where excess force is applied, the screw may break upon removal.



Figure 35




Figure 36




Figure 37

Nail Cap Set Screw method

1. Lock the nail distally, then fully insert the dynamic locking screw as previously described.
2. Remove the Drill Guide/Anterior Drop assembly.
3. Insert the Nail Cap Set Screw (71656000) into the top of the nail and advance with the Medium Hexdriver/T-handle assembly until the fracture is compressed and the Nail Cap Set Screw has fully engaged the locking screw (Figures 38a and 38b).

Note: If there are screws in the other proximal holes, the fracture cannot be compressed.

Note: Use fluoroscopy to check dynamic screw position when compressing the fracture.

Note: For information on the use of Blocking Screws, refer to the TRIGEN® META-NAIL® Tibial Nail Surgical Technique (71181610).




Figure 38a



Figure 38b

Nail extraction: *optional*

Standard technique

1. If implanted, remove the Nail Cap or Nail Cap Set Screw by using the Medium Hexdriver/T-Handle assembly.
2. Remove all the distal locking screws, and all but one of the proximal locking screws.
3. Thread the Cannulated Impactor-Medium (71675081) or Cannulated Impactor-Long (71631185)* into the back of the Disposable Nail Extractor (71631320)**.
4. Thread the assembly into the top of the nail.
5. Remove the remaining proximal locking screw.
6. Use the Slotted Hammer with a back-slapping motion to extract the nail (Figure 39).




Figure 39

*The Cannulated Impactor-Long is located in the original TRIGEN® Instrument Set (71631326)

**The Disposable Nail Extractor (71631320) is interchangeable with the Large Nail Extractor located in the original TRIGEN Instrument Set (71631326) and the HFN Instrument Set (71700001)

Percutaneous technique

In the absence of a Nail Cap or Nail Cap Set Screw:

1. Use the Medium Hexdriver/T-handle assembly to remove all distal locking screws and all but one of the proximal locking screws.
2. Under fluoroscopy, insert a 3.2mm Brad Point Tip Guide Pin (71674130 or 71631436) into the top of the nail either using power or by hand (Figure 40).
3. Make a 2cm incision around the pin.
4. Advance the 12.5mm Entry Reamer over the pin and into the top of the nail to remove any bony ingrowth (Figure 41).

Note: The flared portion of the Entry Reamer enters the top of the nail.

5. Thread the Cannulated Impactor-Medium or Cannulated Impactor-Long* (71631185) into the back of the Disposable Nail Extractor** (71631320).
6. Thread the assembly into the top of the nail.
7. Remove the remaining proximal locking screw.
8. Extract the nail with a back-slapping motion.




Figure 40




Figure 41

An alternative method for extraction

Guide rod jamming technique

1. Advance the end of a 3.0mm Ball Tip Guide Rod through the end of the nail.
2. Insert the smooth end of a 2.0mm Ball Tip Guide Rod (71751146) in the same manner.
3. With both guide rods in place, attach the Gripper to the end of the 3.0mm Ball Tip Guide Rod.
4. Pull the Gripper back so that it wedges the ball tip against the 2.0mm Guide Rod.
5. Backslap against the Gripper with the Slotted Hammer to extract the nail.

Guide rods	
Cat. No.	Description
71631626	3.0mm x 1000mm Ball Tip Guide Rod
71751146	2.0mm x 600mm Ball Tip Guide Rod

Additional removal items	
Cat. No.	Description
115074	Large Extractor Hook*
115073	Small Extractor Hook*

*Available sterile packed. For nail removal only; do not use for nail insertion.

**Located in RUSSELL-TAYLOR Extraction Kit (Set #7508), available through Loaners.


Catalog information



TRIGEN[®] Base Instrument Set


Set No. 71674012

Cat. Item	Description	Qty	Cat. Item	Description	Qty
71129401	Small Outer Case	1	71674078	Obturator	1
71129402	Lid for Outer Case	1	71674079	Ruler	1
71674021	TRIGEN Base Tray	1	71674080	Gripper	1
71631066	Medium Hexdriver	1	71674081	Impactor	1
71631068	Short Hexdriver	1	71674082	Slotted Hammer	1
71631116	12.5mm Entry Reamer	1	71674083	4.0mm Drill Sleeve	2
71631140	Guide Bolt Wrench	1	71674084	Screwdriver Release Handle	1
71631152	9.0mm Drill Sleeve	2	71674085	Screw Length Sleeve	1
71631161	Multipurpose Driver	1	71674092	Entry Portal Handle	1
71631186	Mini Connector	1	71671212	TRIGEN Reamer Set, Optional	1
71631189	Screw Depth Gauge	1	7111-8200	SCULPTOR [®] Flexible Reamer, Optional	1
71674000	Cannulated Awl	1	71631130	Flexible Reamer Extender, Optional	1
71674060	Entry Portal Tube	1	71641123	TRIGEN 4.0mm Diaphyseal Drill	1
71674074	3.2mm T-handle Trocar	1	71631121	4.0mm Long AO Pilot Drill, 333mm, Disposable	2
71674075	Honeycomb	1	71671123	4.0mm Short AO Pilot Drill, 161mm, Disposable	1
71674576	T-handle	1	71674130	3.2mm x 343mm Brad Point Tip Guide Pin, Disposable	3
71674077	Reducer	1	71631070	Long Hexdriver, Optional	1
71674077	Reducer Connector	1	71751153	AO Mini Connector, Optional	1
			71631187	Trinkle to Mini Connector, Optional	1




TRIGEN[®] META-NAIL[®]
Instrument Set
Set No. 71654002

Cat. Item	Description	Qty
71654501	META-NAIL Anterior Drop	1
71654502	META-NAIL Drill Guide	1
71654503	META-NAIL Extension Drill Guide	1
71654505	Extension Guide Bolt (23mm)	2
71654506	Guide Bolt Long (51mm)	2
71654520	Long Screw Length Sleeve	1
71654554	Cannulated Impactor-Short	1
71631025	Large Nail Extractor	1



Blocking Screw
Instrument Set (Optional)
Set No. 71654001

Cat. Item	Description	Qty
71654515	Blocking Screw Device	1
71654509	Tibial Blocking Screw Attachment	1
71654522	11.0mm T-handle Awl	1
71654511	8.5mm/10mm Blocking Screw Cartridge	1
71654513	11.5mm/13mm Blocking Screw Cartridge	1
71654514	Offset Blocking Screw Cartridge	1
71654523	Blocking Screw Alignment Pin	2
71634508	Retrograde Femoral Blocking Screw Attachment	1



TRIGEN[®] META-NAIL[®]
Semi-extended Instrument Set
Set No. 71654004

META-NAIL Disposables
Set No. 71654003

Cat. Item	Description	Qty
71654524	Semi-extended Drill Guide	1
71654525	Semi-extended Guide Bolt	2
71654527	Semi-extended Guide Pin Sleeve	1
71654533	Semi-extended Honeycomb	1
71654526	Semi-extended Entry Tube	1
71654554	Cannulated Impactor - Short	1
71674092	Entry Portal Handle	1

Cat. Item	Description
71631110	4.0mm Long Pilot Drill*
71631117	4.0mm Short Drill**
71631626	3.0mm x 1000mm Ball Tip Guide Rod
71674130 or 71631436	3.2mm Brad Point Guide Pin
71654528	Universal Compression Driver
71631320	Disposable Nail Extractor***

*The 4.0mm Long Pilot Drill (71631110) is interchangeable with 4.0mm AO Long Drill (71631121)

**The 4.0mm Short Drill (71631117) is interchangeable with 4.0mm AO Short Drill (71631123)

***The Disposable Nail Extractor (71631320) is interchangeable with the Large Nail Extractor (71631278) located in the original TRIGEN Instrument Set (71631326) and the HFN Instrument Set (71700001)

Implants

TRIGEN® Low Profile Screws 4.5mm and 5.0mm

Set No. 71642005

Set No. 71642000

Cat. Item	Length
71645420	4.5mm x 20mm
71645422	4.5mm x 22.5mm
71645425*	4.5mm x 25mm
71645427*	4.5mm x 27.5mm
71645430*	4.5mm x 30mm
71645432*	4.5mm x 32.5mm
71645435*	4.5mm x 35mm
71645437*	4.5mm x 37.5mm
71645440*	4.5mm x 40mm
71645442*	4.5mm x 42.5mm
71645445*	4.5mm x 45mm
71645447*	4.5mm x 47.5mm
71645450*	4.5mm x 50mm
71645452*	4.5mm x 52.5mm
71645455*	4.5mm x 55mm
71645457*	4.5mm x 57.5mm
71645460*	4.5mm x 60mm
71645462*	4.5mm x 62.5mm
71645465*	4.5mm x 65mm

Cat. Item	Length
71645020	5.0mm x 20mm
71645022	5.0mm x 22.5mm
71645025*	5.0mm x 25mm
71645027*	5.0mm x 27.5mm
71645030*	5.0mm x 30mm
71645032*	5.0mm x 32.5mm
71645035*	5.0mm x 35mm
71645037*	5.0mm x 37.5mm
71645040*	5.0mm x 40mm
71645042*	5.0mm x 42.5mm
71645045*	5.0mm x 45mm
71645047*	5.0mm x 47.5mm
71645050*	5.0mm x 50mm
71645052*	5.0mm x 52.5mm
71645055*	5.0mm x 55mm
71645057*	5.0mm x 57.5mm
71645060*	5.0mm x 60mm
71645062*	5.0mm x 62.5mm
71645065*	5.0mm x 65mm
71645067*	5.0mm x 67.5mm
71645070*	5.0mm x 70mm
71645072*	5.0mm x 72.5mm
71645075*	5.0mm x 75mm



4.5mm



5.0mm

TRIGEN META-NAIL® 8.5mm Tibial

Set No. 71653002

Cat. Item	Length	Cat. Item	Length
71655024*	24cm	71655035*	35cm
71655026*	26cm	71655036*	36cm
71655028*	28cm	71655037*	37cm
71655029	29cm	71655038*	38cm
71655030*	30cm	71655039	39cm
71655031	31cm	71655040	40cm
71655032*	32cm	71655041	41cm
71655033*	33cm	71655042	42cm
71655034*	34cm		



*Contained in the standard implant set

TRIGEN° META-NAIL° 10mm Tibial
Set No. 71653000

Cat. Item	Length
71655124	24cm
71655126	26cm
71655128*	28cm
71655129	29cm
71655130*	30cm
71655131	31cm
71655132*	32cm
71655133	33cm
71655134*	34cm
71655135*	35cm
71655136*	36cm
71655137*	37cm
71655138*	38cm
71655139*	39cm
71655140*	40cm
71655141	41cm
71655142	42cm



TRIGEN META-NAIL 11.5mm Tibial
Set No. 71653001

Cat. Item	Length
71655224	24cm
71655226	26cm
71655228	28cm
71655229	29cm
71655230*	30cm
71655231	31cm
71655232*	32cm
71655233	33cm
71655234*	34cm
71655235*	35cm
71655236*	36cm
71655237*	37cm
71655238*	38cm
71655239*	39cm
71655240*	40cm
71655241	41cm
71655242	42cm



TRIGEN META-NAIL 13mm Tibial

Cat. Item	Length
71655318	18cm
71655320	20cm
71655322	22cm
71655324	24cm
71655326	26cm
71655328	28cm
71655329	29cm
71655330	30cm
71655331	31cm
71655332	32cm
71655333	33cm
71655334	34cm
71655335	35cm
71655336	36cm
71655337	37cm
71655338	38cm
71655339	39cm
71655340	40cm
71655341	41cm
71655342	42cm
71655343	43cm
71655344	44cm
71655346	46cm
71655348	48cm
71655350	50cm



Nail Cap Set Screw
Cat. Item 71656000



TRIGEN Nail Caps

Cat. Item	Length
71634000	0mm
71634005	5mm
71634010	10mm
71634015	15mm
71634020	20mm

*Contained in the standard implant set

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith & Nephew representative or distributor if you have questions about the availability of Smith & Nephew products in your area.

Supporting healthcare professionals for over 150 years

Caution: *U.S. Federal law restricts these devices to sale by or on the order of a physician.*

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, TN 38116
USA

www.smith-nephew.com
Telephone: 1-901-396-2121
Information: 1-800-821-5700
Orders/Inquiries: 1-800-238-7538

®Trademark of Smith & Nephew.
All trademarks acknowledged.

©2022 Smith & Nephew.
All rights reserved.
01361 V8 71181612 REVG 11/22