


REGENETEN[®] Bioinductive Implant

Introduction

The REGENETEN Implant has a large and growing evidence base in rotator cuff repair.¹⁻¹³ The REGENETEN Implant supports the body's natural healing response to facilitate new tendon-like tissue growth and change the course of rotator cuff tear progression.^{1,2,12-16}

The REGENETEN Implant is resorbed and shows tendon-like tissue at 6 months^{*1,13,14}



With a substantial body of literature evaluating its use in rotator cuff repair (>650 patients across >10 published studies),^{†1-13} the REGENETEN Implant has been associated with:



Lower re-tear rates than reported for standard repair techniques^{4-6,10,17}



Rapid post-operative recovery^{9,13,18}



Significant and clinically meaningful post-operative improvements in clinical outcomes, compared to pre-operative levels^{3,9,19}



Low complication rates^{2,3,5-10,13}



Effective treatment of rotator cuff tears can be challenging. Conventional treatment of PT tears commonly involves removing healthy tissue²⁰ and results in lengthy rehabilitation,²¹ whilst repair of large FT tears has been associated with re-tear rates >50%.²²

29 REGENETEN Implant publications



- 2 Randomised Controlled Trials (RCTs)^{5,13}
- 15 Human Case Series^{‡2,6-8,10-12,15,19,23,24-28}
- 2 Human Histological Case Studies^{1,4}
- 5 Systematic Literature Reviews (SLRs)²⁹⁻³³
- 2 Cost Analyses^{34,35}
- 3 Registries^{3,9,18}

Clinical studies length of follow-up



>800
total patients studied

*On human biopsy (n=1) and in vivo sampling. †As of July 2024. ‡Representing 11 independent cohorts.

Overview of REGENETEN[◇] Bioinductive Implant clinical study publications in rotator cuff repair

Study	Design	Tear	n= (of patients)	Follow-up	Key findings
<p>Ruiz Ibán MA, et al. Augmentation of a transosseous equivalent repair in posterosuperior nonacute rotator cuff tears with a bioinductive collagen implant decreases the retear rate at one year. A randomised controlled trial.</p> <p><i>Arthroscopy.</i> 2024;40(6):1760–1773.</p> <p>Available at (open access): Arthroscopy: The Journal of Arthroscopic & Related Surgery</p>	Randomised controlled trial	FT	124 (enrolled), 122 (at follow-up) (REGENETEN Implant group n=60; control group n=62)	1 year	<p>At 12 months follow-up, compared with repair alone, repair of medium and large full-thickness tears augmented with the REGENETEN Implant demonstrated:</p> <ul style="list-style-type: none"> • Significantly lower re-tear rate (8.3 vs 25.8%; p=0.01) • Significantly better tendon integrity (defined as Sugaya grades 1–3; 91.7 vs 74.2%; p=0.03) • A three times lower risk of re-tear (RR=0.32; 95% CI: 0.13–0.83) • No differences in clinical outcome measures between groups • No difference in the number of serious or minor complications
<p>Camacho Chacón JA, et al. An isolated bioinductive repair vs sutured repair for full-thickness rotator cuff tears: 2-year results of a double blinded, randomized controlled trial.</p> <p><i>J Shoulder Elbow Surg.</i> 2024 May 09. 2024;33(9):1894–1904.</p> <p>Available at: Journal of Shoulder and Elbow Surgery</p>	Randomised controlled trial	FT	60 (REGENETEN Implant group n=30; control group n=30)	2 years	<p>Patients with full-thickness rotator cuff tears, with an intact stable rotator cable, who received an IBR using the REGENETEN Implant at 2 years post-operatively demonstrated:</p> <ul style="list-style-type: none"> • Superior collagen organisation, without inflammation, seen on biopsy (p<0.0001) <ul style="list-style-type: none"> – Control group showed poorly organised collagen in 80% of patients • Significantly thicker tendon at 6 months, compared to the control group, maintained throughout the follow-up period (p<0.0001) • Healed tendon at 12 months (all patients had 100% tendon gap fill-in) • Significantly better ASES and CMS scores at all timepoints, compared to the control group (p<0.0001) <ul style="list-style-type: none"> – Significantly lower VAS pain score at 6 months, compared to the control group (p<0.0001) – No difference in VAS pain score between groups at 24 months • Significantly faster return to work and fewer days in a sling, compared to the control group (p<0.0001) • Significantly fewer physical therapy sessions, compared to the control group (p<0.0001)

Study	Design	Tear	n= (of patients)	Follow-up	Key findings
<p>Bushnell BD, et al. Rotator cuff repair with a bioinductive bovine collagen implant has a low incidence of post-operative stiffness: review of 406 shoulders. Poster presented at the Annual Academy of Orthopaedic Surgeons (AAOS) Annual Congress; March 7–11 2023; Las Vegas, Nevada, USA.</p>	Case series	PT + FT	406 cases (368 patients, some patients had a bilateral rotator cuff repair)	32-month average follow-up (range: 2–102 months)	<ul style="list-style-type: none"> A low incidence of post-operative stiffness (3.2%; 13/406) with a minimum of 2 months follow-up: <ul style="list-style-type: none"> – 1.5% (4/269) in full-thickness tears – 6.6% (9/137) in partial-thickness tears A low incidence of cases requiring secondary intervention for post-operative stiffness (2.7%; 11/406) with a minimum of 2 months follow-up: <ul style="list-style-type: none"> – Arthroscopic lysis of adhesions (n=10) – Isolated manipulation (n=1) Most of the patients who experienced stiffness (11/13) had one or more risk factors for stiffness including: <ul style="list-style-type: none"> – Diabetes – Smoking/tobacco use – Thyroid disease – Workers' compensation status
<p>McIntyre LF, et al. Resorbable bioinductive collagen implant is cost effective in the treatment of rotator cuff tears. <i>Arthrosc Sports Med Rehabil.</i> 2023;5(2):e367–374. Available at: Arthroscopy, Sports Medicine and Rehabilitation</p>	Economic analysis	FT	100	N/A	<ul style="list-style-type: none"> In the base case analysis, the REGENETEN Implant + SoC resulted in improved healing/re-tear rates, at an increased cost, compared with SoC alone When analysed by tear size, the cost-effectiveness of the REGENETEN Implant + SoC relative to SoC alone improved with increasing tear size In the scenario analysis including the impact on return to work, the REGENETEN Implant + SoC was cost-saving compared with SoC alone, saving \$469,017 per 100 treated patients with full-thickness rotator cuff tears One-way sensitivity analysis demonstrated increased cost-effectiveness, compared with the base case analysis, in patients at a higher risk of re-tear due to the following risk factors: age >60 years, hypertension, alcohol consumption and obesity
<p>Camacho-Chacon JA, et al. Bioinductive collagen implants facilitate tendon regeneration in rotator cuff tears. <i>J Exp Orthop.</i> 2022;9(1):53. Available at: Journal of Experimental Orthopaedics</p>	Case series	PT + FT	30 (enrolled), 29 (at follow-up)	1 year	<ul style="list-style-type: none"> No evidence of inflammatory, scarring or ischaemic changes in histological examination of all samples at 6 months post-operatively <ul style="list-style-type: none"> – No evidence of any remaining REGENETEN Implant MRI demonstrated complete healing in 27 patients after 6 months; mean tendon thickness increased significantly by 1.84mm (from 4.18mm pre-operatively to 6.02mm; p=0.001) <ul style="list-style-type: none"> – All patients showed filling of the defect at 6 months with 90% completely filled and 10% with greater than 50% of the defect filled – No evidence of tear propagation at 1 year – Neotendon was indistinguishable from native tendon on MRI in all samples at 6 months and 1 year Patients experienced statistically significant improvements vs pre-operative values in VAS pain score (p=0.003), ASES (p=0.001) and CMS (p=0.001) at 6 months post-operatively, which were sustained at 1 year

Study	Design	Tear	n= (of patients)	Follow-up	Key findings
Bushnell BD, et al. Two-year outcomes with a bioinductive collagen implant used in augmentation of arthroscopic repair of full-thickness rotator cuff tears: Final results of a prospective multi-center study. <i>JSES Int.</i> 2022. 2022;31(12):2532–2541. Available at: Journal of Shoulder and Elbow	Case series	FT	115	2 years	<ul style="list-style-type: none"> Between baseline and 2 year follow-up, the average mean tendon thickness increased by 12.5% for medium tears and 17.1% for large tears Overall re-tear rate of medium tears; 4.6% at 3 months (3/66) and 10.6% at 1 year and 2 years (7/66) Overall re-tear rate of large tears; 20.4% at 3 months (10/49), 24.5% at 1 year (12/49) and 28.6% at 2 years (14/49) Significantly lower re-tear rate in patients who had repair with double-row (12/91, 13.2%) technique at 2 years, compared to single-row technique (p=0.0061) <ul style="list-style-type: none"> Re-tear rates with single-row (9/24, 37.5%) technique at 2 years More than 90% of patients had significant post-operative improvements in both medium and large tears in ASES Shoulder and CMS scores that exceeded respective MCIDs (p<0.001) Most patients (97.1%) surveyed were satisfied with the procedure; 100% of patients surveyed would recommend the procedure to a friend Nine re-operations were required: seven for persistent symptoms, one for superficial infection and one for inflammation and osteopenia
Bushnell BD, et al. Treatment of partial-thickness rotator cuff tear repairs with a resorbable bioinductive bovine collagen implant: 1-year results from a prospective multi-center registry. <i>Orthop J Sports Med.</i> 2021;9(8): 23259671211027850. Available at (open access): Orthopaedic Journal of Sports Medicine	Registry	PT	272 (enrolled), 227 (at follow-up)	1 year	<ul style="list-style-type: none"> All PROMs were significantly improved at 3, 6 and 1 year from pre-operative values, except VR-12 MCS, which was significant at 1 year (p<0.05) <ul style="list-style-type: none"> MCIDs for ASES, SANE and WORC scores were met or exceeded by >90% of patients at 1 year In higher grade tears (≥grade 2) and compared with the augmented takedown and repair group (n=26), isolated bioinductive repair (REGENETEN Implant without surgical repair; n=201) resulted in: <ul style="list-style-type: none"> Significantly better ASES, SANE and WORC scores at 2 and 6 weeks (p<0.05) No significant differences in 1 year outcomes, except VR-12 PCS, which was significantly improved with isolated bioinductive repair (p=0.0213) Significantly less sling time (19.1 vs 34.3 days; p<0.0001) and faster return to non-overhead sports (72.2 vs 128.9 days; p=0.0192); no other significant differences in recovery

Study	Design	Tear	n= (of patients)	Follow-up	Key findings
<p>McIntyre LF, et al. Full-thickness rotator cuff tears can be safely treated with a resorbable bioinductive bovine collagen implant: one-year results of a prospective, multicenter registry. <i>Arthrosc Sports Med Rehabil.</i> 2021;3(5):e1473-e1479.</p> <p>Available at (open access): Arthroscopy, Sports Medicine and Rehabilitation</p>	Registry	FT	210 (enrolled), 192 (at follow-up)	1 year	<ul style="list-style-type: none"> At 6 months and 1 year, ASES, SANE, VR-12 PCS and WORC were significantly improved from pre-operative values (p<0.001) <ul style="list-style-type: none"> MCIDs were met or exceeded by 90.5% of patients for ASES, 84.3% of patients for SANE and 87.2% of patients for WORC were observed VR-12 MCS was significantly improved at 6 months (p=0.002); no significant difference to pre-operative values at 1 year An ad-hoc analysis demonstrated similar results at 1 year regardless of tear size Mean duration of post-operative recovery (days): sling time, 36.3 (n=188); return to driving, 24.0 (n=135); return to work, 48.4 (n=128); return to non-overhead sports, 105.4 (n=71), return to overhead sports, 131.7 (n=42) 18 patients required re-operation;* re-tear (n=11), infection[†] (n=3), shoulder stiffness/adhesive capsulitis (n=3), bursitis (n=1) or implant displacement after a fall (n=1)
<p>Schlegel TF, et al. Isolated bioinductive repair of partial thickness rotator cuff tears using a resorbable bovine collagen implant: Two-year radiologic and clinical outcomes from a prospective multicentre study. <i>J Shoulder Elbow Surg.</i> 2021;30(8):1938–1948.</p> <p>Available at: Journal of Shoulder and Elbow Surgery</p>	Case series	PT	33 (enrolled), 31 (at follow-up)	2 years	<ul style="list-style-type: none"> At 2 years, 87% of tears had reduced in size by >50% from pre-operative measurements No compliant patient progressed to a full-thickness tear; one non-compliant patient progressed ~1 month after surgery Significant increase in mean tendon thickness from pre-operative values to 2 years: 1.2mm in intermediate-grade tears; 1.8mm in high-grade tears (p≤0.012) Improvements in CMS met or exceeded MCIDs for all patients with intermediate-grade tears and 94% of patients with high-grade tears No significant difference in clinical outcomes, tear healing or tendon thickness based on tear location
<p>Dai A, et al. Collagen-based bioinductive implant for treatment of partial thickness rotator cuff tears. <i>Bull Hosp Jt Dis.</i> 2020;78(3):195–201.</p> <p>Available at (open access): Bulletin of the Hospital for Joint Diseases</p>	Case series	PT	30 (eligible), 24 (at follow-up)	Mean: 19.1 months	<ul style="list-style-type: none"> Significant increase in mean ASES score by 22.5 points (45.6 pre-operatively to 68.1 at 19.1 months post-operatively; p=0.001) Significant decrease in mean VAS pain from 8.3 pre-operatively to 3.8 at 19.1 months (p<0.001) Improvements in ASES score and VAS pain from pre-operative values were greatest in articular-sided and intrasubstance tears Mean patient satisfaction score was 7.5/10 Mean tendon thickness significantly increased by 0.8mm (5.7mm preoperatively to 6.5mm at 9.9 months post-operatively; p=0.007) No implant-related complications; one patient had a traumatic re-tear 4 months post-operatively

*Re-operations can be attributed to multiple reasons.

[†]Includes two recurrent infections, which were resolved with treatment.

Study	Design	Tear	n= (of patients)	Follow-up	Key findings
Thon SG, et al. REGENETEN bio-inductive collagen scaffold for rotator cuff tears: indications, technique, clinical outcomes, and review of current literature. <i>Ann Jt.</i> 2020;5:41. Available at (open access): Annals of Joint	Review	PT + FT	251	Mean: 14.9 months (10.8–32.0 months)	<p>Review article including Bokor 2015, Bokor 2016, Schlegel 2018, Thon 2019, and McIntyre 2019, reporting:</p> <ul style="list-style-type: none"> • 92–94% patient satisfaction in partial thickness tears (2 studies) • 94–100% healing rate in partial thickness tears (2 studies) • 89–91% patient satisfaction in full thickness tears (2 studies) • 96–100% healing rate in full thickness tears (2 studies) • 3.9% reoperation rate (10/251; 5 studies) • 5.9% (15/251) failure rate (lack of patient satisfaction, lack of tendon healing, or need for reoperation) based on 5 studies • 9.9% (25/251) complication rate (5 studies)
Micheloni GM, et al. Bio-inductive implant for rotator cuff repair: our experience and technical notes. <i>Acta Biomed.</i> 2020;91:14-S. Available at (open access): Acta Biomedica	Case series	PT + FT	4	6 weeks	<ul style="list-style-type: none"> • In 4 patients (1 PT tear receiving isolated bioinductive repair, 3 FT tears receiving a REGENETEN Implant to augment surgical repair), no complications were found at 6 weeks follow-up • Increase in procedure duration by 10 minutes

Study	Design	Tear	n= (of patients)	Follow-up	Key findings
<p>McIntyre LF, et al. Patient-reported outcomes after use of a bioabsorbable collagen implant to treat partial and full-thickness rotator cuff tears. <i>Arthroscopy.</i> 2019;35(8):2262–2271.</p> <p>Available at: Arthroscopy</p>	Registry	PT + FT	203 (eligible), 173 (at follow-up)	1 year	<p>For patients with partial-thickness tears (n=90):</p> <ul style="list-style-type: none"> • MCIDs were achieved in VAS pain from 2 weeks (5.3 at baseline to 3.3 at 2 weeks and 1.1 at 1 year; p<0.001) and ASES score from 6 weeks (47.0 at baseline to 60.6 at 6 weeks and 85.6 at 1 year; p<0.001) • Mean VAS pain and ASES score at 3 months were significantly improved compared to published results from patients undergoing transtendon or takedown repair without a REGENETEN Bioinductive Implant (p<0.05) • Significant improvements were observed from 6 weeks in SANE score (42.5 at baseline to 59.4 at 6 weeks and 86.6 at 1 year; p<0.001), VR-12 PCS (35.8 at baseline to 39.1 at 6 weeks and 49.7 at 1 year; p<0.002) and WORC Index (38.2 at baseline to 53.5 at 6 weeks and 84.4 at 1 year; p<0.001) • Time in sling and time to return to driving and sport were found to be shorter than published results for standard surgical techniques <p>For patients with full-thickness tears (n=83):</p> <ul style="list-style-type: none"> • MCIDs were achieved in VAS pain from 2 weeks (5.2 at baseline to 3.7 at 2 weeks and 1.2 at 1 year; p<0.001) and ASES score from 3 months (45.5 at baseline to 68.4 at 3 months and 83.8 at 1 year; p<0.001) • Significant improvements from 6 weeks in WORC index (35.0 at baseline to 41.1 at 6 weeks and 80.1 at 1 year; p<0.001), and from 3 months in SANE score (39.2 at baseline to 63.3 at 3 months and 80.7 at 1 year; p<0.001) and VR-12 PCS (34.5 at baseline to 40.8 at 3 months and 45.7 at 1 year; p<0.001) • Shorter time to return to driving and sport compared to published results for standard surgical techniques
<p>Thon SG, et al. Evaluation of healing rates and safety with a bioinductive collagen patch for large and massive rotator cuff tears; 2-year safety and clinical outcomes. <i>Am J Sports Med.</i> 2019;47(8):1901–1908.</p> <p>Available at: American Journal of Sports Medicine</p>	Case series	FT	23	2 years	<ul style="list-style-type: none"> • No implant-related adverse events • Tendon healing in 22/23 patients (96%) on ultrasound at 2 years • Treatment success in 21/23 patients (91%) at 2 years; one healing failure and one failure due to progression of glenohumeral osteoarthritis • Mean tendon thickness on ultrasound increased from 6.29mm at 3 months to 7.72mm at 1 year, decreasing to 7.28mm at 2 years

Study	Design	Tear	n= (of patients)	Follow-up	Key findings
<p>Bokor DJ, et al. Healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a highly-porous collagen implant: a 5-year clinical and MRI follow-up. <i>Muscles Ligaments Tendons J.</i> 2019;9(3):338–347.</p> <p>Available at (open access): Muscles, Ligaments and Tendons Journal</p>	Case series	PT	11	5 years	<ul style="list-style-type: none"> 8/11 patients (73%) had no decline in tendon integrity between 2 and 5 years Mean tendon thickness significantly decreased between 2 and 5 years (5.9 vs 5.2mm; p=0.0012), but remained significantly greater than pre-operative values (4.3mm; p<0.0001) Significant improvements from baseline in ASES total score and ASES pain score were sustained to 5 years (p≤0.01) and were not significantly different to 2 year values
<p>Schlegel TF, et al. Radiologic and clinical evaluation of a bioabsorbable collagen implant to treat partial-thickness tears: a prospective multicenter study. <i>J Shoulder Elbow Surg.</i> 2018;27(2):242–251.</p> <p>Available at: Journal of Shoulder and Elbow Surgery</p>	Case series	PT	33	1 year	<ul style="list-style-type: none"> Reduction in tear size of ≥1 grade from baseline in 31/33 patients (94%) at 1 year <ul style="list-style-type: none"> 8/33 patients (24%) had no visible defect Significant increase in mean tendon thickness on 1 year MRI versus pre-operative values (p<0.0001) Improvements in ASES pain, ASES shoulder index and CMS were significant at 1 year compared with pre-operative values (p<0.0001) <ul style="list-style-type: none"> Improvements in each score were greater than twice the respective MCIDs 30/33 patients (94%) were satisfied with the results of their procedure at 1 year
<p>Arnoczky SP, et al. Histologic evaluation of biopsy specimens obtained after rotator cuff repair augmented with a highly porous collagen implant. <i>Arthroscopy.</i> 2017;33(2):278–283.</p> <p>Available at: Arthroscopy</p>	Case series	FT	7	Up to 9 months	<ul style="list-style-type: none"> At the earliest time period (5 weeks), the biopsy showed host cell ingrowth and early collagen formation At 3 months, there was increased collagen formation, maturation and organization on the surface of the implant By 6 months, the implant was no longer visible, with new tendon-like tissue and oriented collagen indicative of functional loading No evidence of foreign body or inflammatory reactions at any time point

Study	Design	Tear	n= (of patients)	Follow-up	Key findings
<p>Bokor DJ, et al. Evidence of healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a collagen implant: a 2-year MRI follow-up. <i>Muscles Ligaments Tendons J.</i> 2016;6(1):16–25.</p> <p>Available at (open access): Muscles, Ligaments and Tendons Journal</p>	Case series	PT	13	2 years	<ul style="list-style-type: none"> At 3 months after surgery there was a significant increase in mean tendon thickness of 2.2mm at 3 months versus pre-operative values ($p<0.0001$) At 2 years, new tissue was indistinguishable from underlying tissue in 12/13 patients (92%) Tendon thickness at 2 years was significantly greater ($p<0.0001$) than pre-operative values At 1 year, all assessable patients had a reduction in tear size of ≥ 1 grade, with complete tear disappearance in 7 of 10 patients with a measurable pre-operative tear size Significant improvement in clinical scores over the 2 year follow-up period from pre-operative values: CMS, Constant-Murley pain score, ASES total score and ASES pain score ($p\leq 0.01$) Outcomes were satisfactory for 12/13 patients (92%) at 2 years
<p>Bokor DJ, et al. Preliminary investigation of a biological augmentation of rotator cuff repairs using a collagen implant: a 2-year MRI follow-up. <i>Muscles Ligaments Tendons J.</i> 2015;5(3):144–150.</p> <p>Available at (open access): Muscles, Ligaments and Tendons Journal</p>	Case series	FT	9	2 years	<ul style="list-style-type: none"> No MRI evidence of re-tear or gap formation, with the integrity of all repaired tendons intact at 2 years Significant increase in mean tendon thickness versus published values at 3, 6, 12 and 24 months ($p<0.01$), with an average of 2mm new tissue over the bursal surface New tissue rapidly matured, improved in quality and was indistinguishable from the native tendon by 1 year From 1 to 2 years, tendon thickness slightly decreased, likely reflecting continued functional remodeling Significant improvement in clinical scores at 2 years versus pre-operative measures: CMS and Constant-Murley pain score (both $p<0.001$); ASES score and ASES pain score (both $p<0.001$) Outcomes were satisfactory for 8/9 patients (89%) at 2 years

References

1. Arnoczky SP, Bishai SK, Schofield B, et al. Histologic evaluation of biopsy specimens obtained after rotator cuff repair augmented with a highly porous collagen implant. *Arthroscopy*. 2017;33(2):278–283. <https://doi.org/10.1016/j.arthro.2016.06.047>
2. Bokor DJ, Sonnabend D, Deady L, et al. Evidence of healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a collagen implant: a 2-year MRI follow-up. *Muscles Ligaments Tendons J*. 2016;6(1):16–25. <https://pubmed.ncbi.nlm.nih.gov/27331028/>
3. Bushnell BD, Bishai SK, Krupp RJ, et al. Treatment of partial-thickness rotator cuff tears with a resorbable bioinductive bovine collagen implant: 1-year results from a prospective multicenter registry. *Orthop J Sports Med*. 2021;9(8):23259671211027850. <https://doi.org/10.1177/23259671211027850>
4. Camacho-Chacon JA, Cuenca-Espierrez J, Roda-Rojo V, et al. Bioinductive collagen implants facilitate tendon regeneration in rotator cuff tears. *J Exp Orthop*. 2022;9(1):53. <https://doi.org/10.1186/s40634-022-00495-7>
5. Ruiz Ibán MA, Navlet MG, Marco SM, et al. Augmentation of a transosseous-equivalent repair in posterosuperior non-acute rotator cuff tears with a bioinductive collagen implant decreases the re-tear rate at one year. A randomised controlled trial. *Arthroscopy*. 2024;40(6):1760–1773. <https://doi.org/10.1016/j.arthro.2023.12.014>
6. Thon SG, O'Malley L 2nd, O'Brien MJ, Savoie FH 3rd. Evaluation of healing rates and safety with a bioinductive collagen patch for large and massive rotator cuff tears: 2-year safety and clinical outcomes. *Am J Sports Med*. 2019;47(8):1901–1908. <https://doi.org/10.1177/0363546519850795>
7. Bokor DJ, Sonnabend D, Deady L, et al. Preliminary investigation of a biological augmentation of rotator cuff repairs using a collagen implant: a 2-year MRI follow-up. *Muscles Ligaments Tendons J*. 2015;5(3):144–150. <https://pubmed.ncbi.nlm.nih.gov/26605186/>
8. Dai A, Campbell A, Bloom D, Baron S, Begly J, Meislin R. Collagen-based bioinductive implant for treatment of partial thickness rotator cuff tears. *Bull Hosp Jt Dis*. 2020;78(3):195–201. <https://pubmed.ncbi.nlm.nih.gov/32857027/>
9. McIntyre LF, McMillan S, Trenhaile SW, Bishai SK, Bushnell BD. Full-thickness rotator cuff tears can be safely treated with a resorbable bioinductive bovine collagen implant: one-year results of a prospective, multicenter registry. *Arthrosc Sports Med Rehabil*. 2021;3(5):e1473–e1479. <https://doi.org/10.1016/j.asmr.2021.07.009>
10. Bushnell BD, Connor PM, Harris HW, Ho CP, Trenhaile SW, Abrams JS. Two-year outcomes with a bioinductive collagen implant used in augmentation of arthroscopic repair of full-thickness rotator cuff tears: final results of a prospective multi-center study. *J Shoulder Elbow Surg*. 2022. 2022;31(12):2532–2541. <https://doi.org/10.1016/j.jse.2022.05.025>
11. Micheloni GM, Salmasso G, Zecchinato G, Giaretta S, Barison E, Momoli A. Bio-inductive implant for rotator cuff repair: our experience and technical notes. *Acta Biomed*. 2020;91(14-5):e2020004. <https://doi.org/10.23750/abm.v91i14-s.10824>
12. Schlegel TF, Abrams JS, Bushnell BD, Brock JL, Ho CP. Radiologic and clinical evaluation of a bioabsorbable collagen implant to treat partial-thickness tears: a prospective multicenter study. *J Shoulder Elbow Surg*. 2018;27(2):242–251. <https://doi.org/10.1016/j.jse.2017.08.023>
13. Camacho Chacón JA, Roda Rojo V, Martínez AM, et al. An isolated bioinductive repair vs sutured repair for full-thickness rotator cuff tears: 2-year results of a double blinded, randomized controlled trial. *J Shoulder Elbow Surg*. 2024;33(9):1894–1904. <https://doi.org/10.1016/j.jse.2024.03.043>
14. Van Kampen C, Arnoczky S, Parks P, et al. Tissue-engineered augmentation of a rotator cuff tendon using a reconstituted collagen scaffold: a histological evaluation in sheep. *Muscles Ligaments Tendons J*. 2013;3(3):229–235. <https://pubmed.ncbi.nlm.nih.gov/24367785/>
15. Bokor DJ, Sonnabend DH, Deady L, et al. Healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a highly-porous collagen implant: a 5-year clinical and MRI follow-up. *Muscles Ligaments Tendons J*. 2019;9(3):338–347. <http://dx.doi.org/10.32098/mltj.03.2019.07>
16. McElvany MD, McGoldrick E, Gee AO, Neradilek MB, Matsen FA, 3rd. Rotator cuff repair: published evidence on factors associated with repair integrity and clinical outcome. *Am J Sports Med*. 2015;43(2):491–500. <https://doi.org/10.1177/0363546514529644>
17. Hein J, Reilly JM, Chae J, Maerz T, Anderson K. Retear rates after arthroscopic single-row, double-row, and suture bridge rotator cuff repair at a minimum of 1 year of imaging follow-up: a systematic review. *Arthroscopy*. 2015;31(11):2274–2281. <https://doi.org/10.1016/j.arthro.2015.06.004>
18. McIntyre LF, Bishai SK, Brown PB 3rd, Bushnell BD, Trenhaile SW. Patient-reported outcomes after use of a bioabsorbable collagen implant to treat partial and full thickness rotator cuff tears. *Arthroscopy*. 2019;35(8):2262–2271. <https://doi.org/10.1016/j.arthro.2019.02.019>
19. Schlegel TF, Abrams JS, Angelo RL, Getelman MH, Ho CP, Bushnell BD. Isolated bioinductive repair of partial-thickness rotator cuff tears using a resorbable bovine collagen implant: two-year radiologic and clinical outcomes from a prospective multicenter study. *J Shoulder Elbow Surg*. 2021;30(8):1938–1948. <https://doi.org/10.1016/j.jse.2020.10.022>
20. Lo IKY, Burkhart SS. Transendon arthroscopic repair of partial-thickness, articular surface tears of the rotator cuff. *Arthroscopy*. 2004;20(2):214–220. <https://doi.org/10.1016/j.arthro.2003.11.042>
21. Millett PJ, Wilcox RB 3rd, O'Holleran JD, Warner JJR. Rehabilitation of the rotator cuff: an evaluation-based approach. *J Am Acad Orthop Surg*. 2006;14(11):599–609. <https://doi.org/10.5435/00124635-200610000-00002>
22. Rashed MS, Cooper C, Cook J, et al. Increasing age and tear size reduce rotator cuff repair healing rate at 1 year. *Acta Orthop*. 2017;88(6):606–611. <https://doi.org/10.1080/17453674.2017.1370844>
23. Bushnell BD, Piller CP, Hicks JS, Jarvis B, Jarvis RC, Baudier RS. Rotator cuff repair with a bioinductive bovine collagen implant has a low incidence of post-operative stiffness: review of 406 shoulders. Poster presented at: the Annual Academy of Orthopaedics Surgeons (AAOS) Annual Congress; March 7–11 2023; Las Vegas, Nevada, USA.
24. Yatseta Yatseta O, Escobar PVF. All-inside bioinductive collagen implant augmentation for the repair of myotendinous rotator cuff tears: surgical technique and report of preliminary results. *Acta Scientific Orthopaedics*. 2022;5(11):96–101. <https://actascientific.com/ASOR/ASOR-05-0607.php>
25. Zhang T, Ajayi A, Hajjar M, Fleckenstein CM, Nolan J, Hasan SS. Arthroscopic repair of retracted large and massive rotator cuff tears with and without augmentation with a bio-Inductive collagen implant reveals substantial and comparable clinical improvement. *Arthroscopy*. 2024;40(5):1434–1442. <https://doi.org/10.1016/j.arthro.2023.10.024>
26. Yeazell S, Lutz A, Bohon H, et al. Increased stiffness and reoperation rate in partial rotator cuff repairs treated with a bovine patch: a propensity-matched trial. *J Shoulder Elbow Surg*. 2022;31(6):S131–S135. <https://doi.org/10.1016/j.jse.2022.02.003>
27. Ting RS, Loh YC, Rosenthal R, et al. Revision rotator cuff repair with versus without an arthroscopically inserted onlay bioinductive implant in workers' compensation patients. *Orthop J Sports Med*. 2023;11(6):23259671231175883. <https://doi.org/10.1177/23259671231175883>
28. Bushnell BD, Connor PM, Harris HW, Ho CP, Trenhaile SW, Abrams JS. Retear rates and clinical outcomes at 1 year after repair of full-thickness rotator cuff tears augmented with a bioinductive collagen implant: a prospective multicenter study. *JSES Int*. 2020;5(2):228–237. <https://doi.org/10.1016/j.jseint.2020.10.020>
29. Priyavadhana S, Ashraf M, Sakthivelan V, Purudappa PP, Mounasamy V, Sambandam S. Non-acute rotator cuff tear: repair augmented with reconstituted absorbable collagen scaffold (RACS)- systematic review. *Arch Bone Jt Surg*. 2021;9(2):135–140. <https://doi.org/10.22038/abjs.2020.52764.2612>
30. Mandalia K, Mousad A, Welborn B, et al. Scaffold- and graft-based biological augmentation of rotator cuff repair: an updated systematic review and meta-analysis of preclinical and clinical studies for 2010–2022. *J Shoulder Elbow Surg*. 2023;32(9):1784–1800. <https://doi.org/10.1016/j.jse.2023.03.031>
31. Warren JR, Domingo-Johnson ER, Sorensen AA, Cheng AL, Latz KH, Cil A. Bioinductive patch as an augmentation for rotator cuff repair, a systematic review and meta-analysis. *J Shoulder Elbow Surg*. 2024;S1058-2746(24)00448–8. <https://doi.org/10.1016/j.jse.2024.05.002>
32. Longo UG, Marino M, Lalli A, Bandini B, Giannarelli D, Ruiz Ibán MÁ. Arthroscopic management of isolated partial-thickness rotator cuff tears. *Knee Surg Sports Traumatol Arthrosc*. 2024;32(9):2358–2375. <https://doi.org/10.1002/ksa.12326>
33. Thon SG, Belk JW, Bravman JT, McCarty EC, Savoie FH. Regenerated bio-inductive collagen scaffold for rotator cuff tears: indications, technique, clinical outcomes, and review of current literature. *Ann Joint*. 2020;5:41. <https://aoj.amegroups.org/article/view/5816/html>
34. McIntyre LF, Nherera LM, Schlegel TF. Resorbable bioinductive collagen implant is cost effective in the treatment of rotator cuff tears. *Arthrosc Sports Med Rehabil*. 2023;5(2):e367–e374. <https://doi.org/10.1016/j.asmr.2023.01.002>
35. Rognoni C, Nherera LM, Garofalo R, et al. Economic evaluation of a bioinductive implant for the repair of rotator cuff tears compared with standard surgery in Italy. *Adv Ther*. 2023;40(12):5271–5284. <https://doi.org/10.1007/s12325-023-02686-9>

Abbreviations

ASES = American Shoulder and Elbow Surgeons; CI = confidence interval; CMS = Constant-Murley score; FT = full-thickness; IBR = isolated bioinductive repair; MCIIDs = minimal clinically important differences; NA = not applicable; PROMs = patient-reported outcome measures; PT = partial-thickness; QALY = quality-adjusted life years; RR = relative risk; SANE = single assessment numeric evaluation; SoC = standard of care; VAS = visual analogue scale; VR-12 MCS = Veterans RAND 12 Item Health Survey Mental Component Score; VR-12 PCS = Veterans RAND 12 Item Health Survey Physical Component Score; WORC = Western Ontario Rotator Cuff.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Smith+Nephew representative or distributor if you have questions about the availability of Smith+Nephew products in your area. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.