

Always innovating, never imitating

In over 20 years of clinical use, TAYLOR SPATIAL FRAME has been used to treat more than 132,000 patients with deformity or traumatic injury in over 50 countries around the world.



95% success rate in acute trauma¹



91% of patients with non-union or malunion achieved union¹



86% of patients did not require further surgery for complications¹

Smith+Nephew

TAYLOR SPATIAL
FRAME[◇]

External Fixator

Adults

TAYLOR SPATIAL FRAME[◇] External Fixator

1996 Orthopaedic Surgeon J. Charles Taylor collaborated with Smith+Nephew to develop the TAYLOR SPATIAL FRAME (TSF[◇]) External Fixator.

Dr. Taylor took mathematical algorithms already employed by the aerospace and automotive industries, and married them to Professor Ilizarov's principles of Distraction Osteogenesis to produce the first-of-its-kind hexapod for limb reconstruction.

The TSF construct is two rings attached to bone and connected by six telescoping struts. A prescription for strut adjustment is generated by web-based software allowing correction of deformity in six axes simultaneously, and achieving reduction to within 1mm and 1°.

More than

200

clinical publications reference
TAYLOR SPATIAL FRAME¹

Every year we connect

500

surgeons around the world
with master faculty at
our industry-leading
instructional courses

Smith+Nephew has led the
education of orthopaedic
surgeons in circular fixation
since the first trip to Kurgan in

1988

We help you push the boundaries in limb restoration and allow your patients to rediscover the joy of **Life Unlimited.**

TAYLOR SPATIAL FRAME[◊] External Fixator is the **most widely used hexapod in the world.**



95% success rate in acute trauma¹



91% of patients achieved non-union or malunion consolidation¹



86% of patients were managed without need for further surgery¹



In over 20 years of clinical use, TSF has been used to treat more than 132,000 patients with deformity or traumatic injury in over 50 countries around the world.



95% success rate in acute trauma¹

Challenge

Complex trauma such as open fractures and periarticular fractures are often not amenable to treatment with traditional internal fixation methods.^{18,20-22}

TAYLOR SPATIAL FRAME[◊] Solution

TAYLOR SPATIAL FRAME can be used to stabilise complex fractures without causing further injury to compromised tissue. Fine Wire fixation facilitates restoration of the joint surface. The stability of TSF[◊] allows for near immediate weight-bearing.^{21-23,30}

Upon completion of treatment, no hardware is left behind.


TSF works in harmony with the body's ability to re-build bone.

Studies that assessed the use of TSF for acute trauma (fractures) in adults¹

Study	Number of patients	Patients with successfully realised treatment goals for consolidation (%)
Ahearn 2014	21 adult patients with fractures	21/21 (100%)
Menakaya 2014	37 patients (mean age 45 yr) with fractures	37/37 (100%)
O'Neill 2016	15 adult patients with fractures	13/15 (86.7%)
Rampurada 2008	26 adult patients with fractures	23/26 (88.5%)
Sala 2017	20 adult patients with fractures	17/20 (85%)

Smith+Nephew 2019. Systematic literature review with meta-analysis of TSF clinical effectiveness. Internal report. EO/TRAUMA/TSF/001/v6.

Open Pilon Fracture – TSF applied with bony deformity to allow soft tissue closure. Gradual TSF correction performed as soft tissue heals.



Images courtesy of J. Tracy Watson MD



91% of patients with non-union or malunion achieved union¹



Challenge

Many factors may interfere with normal fracture repair, resulting in non-union with persistent motion between fracture fragments. Fractures of the tibia have the highest incidence of nonunion¹⁷

TAYLOR SPATIAL FRAME[◇] Solution

Fractures that have failed to heal can be resected and stabilised with TAYLOR SPATIAL FRAME using the Ilizarov method.


Normal limb length and alignment can be achieved¹

Studies looking at outcomes for non-unions and malunions treated with TSF¹

Study	Number of patients	Patients with successfully realised treatment goals for union (%)
Arvesen 2017	37 adult patients	35/37 (94.6%)
Khunda 2016	40 adult and adolescent patients (mean age 39.5 yr)	39/40 (97.5%)
Napora 2017	75 patients (mean age 45.7 yr)	70/75 (93.3%)
Rozbruch 2008	38 adult and adolescent patients (mean age 43 yr)	27/38 (71.1%)

Smith+Nephew 2019. Systematic literature review with meta-analysis of TSF clinical effectiveness. Internal report. EO/TRAUMA/TSF/001/v6.

Thirty-one-year-old male with non-union following treatment with unilateral external fixator for crush injury.



Images courtesy of S. Robert Rozbruch MD



86% of patients did not require further surgery to manage complications¹

Challenge

External Fixation is often associated with more complications than traditional internal fixation methods


TAYLOR SPATIAL FRAME[◇] Solution

- 70% of patients in four studies reported no complications
- Most complications are managed conservatively and remotely
- Complications such as pin site infection are usually managed conservatively without compromising the success of the treatment

Proportion of patients who experienced no complications, or required a non-operative/operative intervention to resolve complications. Based on four Adult studies across all indications, that included patient-specific complication results, with more than 10 patients.¹

Study Author	Number of patients with operative interventions	Number of patients with non-operative interventions	Number of patients without complications
Ahearn 2014	0	8	13
Arvesen 2017	4	1	32
Rampurada 2008	2	10	14
Rozbruch 2008	11	0	27

Smith+Nephew 2019. Systematic literature review with meta-analysis of TSF clinical effectiveness. Internal report. EO/TRAUMA/TSF/001/v6.



Systematic literature review with meta-analyses of TSF[◇] clinical effectiveness EO/TRAUMA/TSF/001/v6

Nanocrystalline Silver has been shown to act as an effective antimicrobial barrier.⁹⁻¹²

ACTICOAT[◇] EXFIX contains a nanocrystalline silver layer that has been shown to be effective against over 150 pathogens*¹³⁻¹⁶


ACTICOAT rapidly kills bacteria in as little as 30 minutes*¹³⁻¹⁶

*As demonstrated *in vitro*



In over 20 years of clinical use, TSF[◇] has been used to treat more than 132,000 patients with deformity or traumatic injury in over 50 countries around the world.

Forty-five-year-old female with congenital fibrous dysplasia and left leg length discrepancy of 78mm




Treated with TAYLOR SPATIAL FRAME. Union achieved after 4.5 months.
Images courtesy of S. Robert Rozbruch MD

“...About a year ago, Taylor’s computer-interface patent expired. As expected, competitors drooled with anticipation of tapping into the deformity correction market. A number of new circular fixators – each with an integrated computer program – appeared on fabric-covered tables at the world’s orthopedic meetings. The products’ champions lauded their imagined advantages over the TAYLOR SPATIAL FRAME™: a tweak here, a nip there, and yes, the strut colored markers are prettier too!


Rest assured, dear reader: many moons will pass before the proponents of the aforementioned gadgets will accumulate the wisdom contained herein...”

Stuart A. Green, MD Clinical Professor,
Orthopedic Surgery University
of California, Irvine

Foreword to “the Art of Limb Alignment:
TAYLOR SPATIAL FRAME” RIAO,
Sinai Hospital of Baltimore



SPATIALFRAME.com




How TAYLOR SPATIAL FRAME[◇] External Fixator Works

With streamlined instrumentation and innovative hardware, the TAYLOR SPATIAL FRAME External Fixator offers the maximum benefits of a circular fixator without the complexity of traditional Ilizarov methods.^{2,3,5}





Angled Pin Connector

Angled Pin Connectors allow for the placement of a steering pin, which provides stability in oblique fractures.^{5,6}



Rancho with Post

Smith+Nephew offers the most comprehensive range of Rancho cubes and posts since their design by Dr. Stuart A. Green in 1991. Ranchos are designed to be modular, for use with 4mm, 5mm, and 6mm Half Pins. Threaded Rancho Posts are designed to simplify the connection of Pins to Rings.




Half Pin

Hydroxyapatite Coated Half Pins have a tapered minor – constant major diameter for improved bi-cortical purchase.



Circular frame

The 7-hole Tab offers more options for fixation and flexibility for Strut attachment.



Wires

Wires are designed to be minimally invasive and to allow for stable fixation in small fragments.

Drill tip wires reduce heat generation and chance of thermal necrosis.^{7*}

**tested at 700rpm in vitro*

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, Tennessee 38116
USA

www.smith-nephew.com

®Trademark of Smith+Nephew
All Trademarks acknowledged
©2020 Smith & Nephew, Inc.
17127 V3 07/20

References

1. Smith & Nephew 2019. Systematic literature review with meta-analysis of TSF clinical effectiveness. Internal report. EO/TRAUMA/TSF/001/v6. Complete Bibliography listed below. 2. Feldman DS, Shin SS, Madan S, Koval KJ. Correction of Tibial Malunion and Nonunion With Six-Axis Analysis Deformity Correction Using the Taylor Spatial Frame. *J Orthop Trauma*. 2003;17(8):549-554. 3. Moroni A, Cadossi M, Romagnoli M, Faldini C, Giannini S. A Biomechanical and Histological Analysis of Standard Versus Hydroxyapatite-Coated Pins for External Fixation. *J Biomed Mater Res Part B: Appl Biomater* 86B: 417-421, 2008 4. HA Coated Half Pins Design Rationale. Smith & Nephew 7108-0616 08/04. 5. Taylor C. Dynamic interfragmentary compression in oblique fractures stabilized with half pin external fixation: The Steerage Pin. Paper presented at: Annual Meeting of the American Academy of Orthopaedic Surgeons 1994; New Orleans, Louisiana. 6. Lenarz C, Bledsoe G, Watson JT. Circular External Fixation Frames with Divergent Half Pins. *Clin Orthop Relat Res* (2008) 466:2933-2939 7. Livingstone J, Hartsell Z. Evaluation of heat generated with drill tip k-wires. Smith & Nephew 7118-1446 07/09. 8. www.smarttrak.com Market data 9. Holder IA, Durkee P, Supp AP, Boyce ST. Assessment of a silver-coated barrier dressing for potential use with skin grafts on excised burns. *Burns* 29 (2003) 445-448 10. Thomas S, McCubbin P. A comparison of the antimicrobial effects of four silver-containing dressings on three organisms. *Journal of Wound Care* Vol 12. No 3. March 2003 11. Strohal R, Schelling M, Takacs M, Jurecka W, Gruber U, Offner F. Nanocrystalline Silver Dressings as an Efficient Anti-MRSA Barrier: A New Solution to an Increasing Problem. *Journal of Hospital Infection* (2005) 60, 226-230 12. Edwards-Jones V. Antimicrobial and barrier effects of silver against methicillin-resistant *Staphylococcus aureus*. *J of Wound Care* Vol 15. No 7. July 2006. 13. Wright JB, Lam K, Burrell RE. Wound Management in an era of increasing bacterial antibiotic resistance: a role for topical silver treatment. *Amer J of Infection Control* 1998; 26:6 572-577 14. Wright JB, Lam K, Hansen D, Burrell RE. Efficacy of Topical silver against fungal burn wound pathogens. *Amer J of Infection Control* 1999; 27: 344-350 15. Wright JB, Hansen D, Burrell RE. The comparative efficacy of two antimicrobial barrier dressings: In-vitro Examination of Two Controlled Release of Silver Dressings. *WOUNDS: A Compendium of Clinical Research and Practice* Vol 10 No 6 Nov/Dec 1998. 16. Yin HQ, Langford R, Burrell RE. Comparative Evaluation of the Antimicrobial Activity of ACTICOAT Antimicrobial Barrier Dressing. *Journal of Burn Care & Rehabilitation* May/June 1999 195-200 17. Catagni MA. Treatment of Fractures, Nonunions, and bone Loss of the Tibia with the Ilizarov Method. Smith & Nephew 2007.

The following Clinical Publications were included in the Systematic literature review with meta-analysis of TSF clinical effectiveness. Internal report. EO/TRAUMA/TSF/001/v6

1. Ahearn N, Oppy A, Halliday R, et al. The outcome following fixation of bicondylar tibial plateau fractures. *Bone Joint J*. 2014;96-B(7):956-962. 2. Lahoti O, Findlay I, Abhishetty N. Intentional deformation and closure of soft tissue defect in open tibial fractures with Taylor spatial frame: 'sat nav' technique. *Bone Joint J*. 2013;95-B((Supp 20)):10. 3. Menakaya C, Rigby A, Hadland Y, Barron E, Sharma H. Fracture healing following high energy tibial trauma: Ilizarov versus Taylor Spatial Frame. *Ann R Coll Surg Engl*. 2014;96(2):106-110. 4. O'Neill B, Fox C, Molloy A, O'hEireamhoin S, Moore D. The use of circular external fixators in the management of lower limb trauma in Dublin: a single surgeon's 20-year experience. *Ir J Med Sci*. 2016;185(1):133-138. 5. Rampurda A, Madan S, Tadross T. Treatment of complex tibial plateau and distal tibial fractures with Taylor spatial frame: experience in a district general hospital. *Eur J Orthop Surg Traumatol*. 2008;18:521. 6. Sala F, Thabet A, Capitani P, Bove F, Abdelgawad A, Lovisetti G. Open supracondylar-intercondylar fractures of the femur treatment with taylor spatial frame. *J Orthop Trauma*. 2017;31(10):546-553. 7. Sharma H, Nunn T. Conversion of open tibial IIIb to IIIa fractures using intentional temporary deformation and the Taylor Spatial Frame. *Strategies Trauma Limb Reconstr*. 2013;8(2):133-140. 8. Arvesen J, Watson J, Israel H. Effectiveness of treatment for distal tibial nonunions with associated complex deformities using a hexapod external fixator. *J Orthop Trauma*. 2017;31(2):e43-e48. 9. Khunda A, Al-Maiyah M, Eardley W, Montgomery R. The management of tibial fracture non-union using the Taylor Spatial Frame. *J Orthop*. 2016;13(4):360-363. 10. Napora J, Weinberg D, Eagle B, Kaufman B, Sontich J. Hexapod frame stacked transport for tibial infected nonunions with bone loss: analysis of use of adjunctive stability. *J Orthop Trauma*. 2017;31(7):393-399. 11. Rozbruch S, Pugsley J, Fragomen A, Ilizarov S. Repair of tibial nonunions and bone defects with the Taylor Spatial Frame. *J Orthop Trauma*. 2008;22(2):88-95. 12. Schoenleber S, Hutson J. Treatment of hypertrophic distal tibia nonunion and early malunion with callus distraction. *Foot Ankle Int*. 2015;36(4):400-407. 13. Ashfaq K, Fragomen A, Nguyen J, Rozbruch S. Correction of proximal tibia varus with external fixation. *J Knee Surg*. 2012;25(5):375-384. 14. Baumgartner H, Weber A. Segmental bone transport and deformity correction in patients with osteomyelitis using the Taylor spatial frame (TSF). *Archives of Orthopaedic and Trauma Surgery*. 2017;137(11):1602. 15. Docquier P, Rodriguez D, Mousny M. Three-dimensional correction of complex leg deformities using a software assisted external fixator. *Acta Orthop Belg*. 2008;74(6):816-822. 16. Horn D, Fragomen A, Rozbruch S. Supramalleolar osteotomy using circular external fixation with six-axis deformity correction of the distal tibia. *Foot Ankle Int*. 2011;32(10):986-993. 17. Manggala Y, Angthong C, Primadhi A, Kungwan S. The deformity correction and fixator-assisted treatment using Ilizarov versus Taylor spatial frame in the foot and ankle. *Orthop Rev (Pavia)*. 2018;9(4):114-117. 18. Nakase T, Kitano M, Kawai H, et al. Distraction osteogenesis for correction of three-dimensional deformities with shortening of lower limbs by Taylor Spatial Frame. *Arch Orthop Trauma Surg*. 2009;129(9):1197-1201. 19. Robinson P, Papanna M, Somanchi B, Khan S. High tibial osteotomy in medial compartment osteoarthritis and varus deformity using the Taylor spatial frame: early results. *Strat Traum Limb Recon*. 2011;6(3):137-145. 20. Rozbruch S, Segal K, Ilizarov S, Fragomen A, Ilizarov G. Does the Taylor Spatial Frame accurately correct tibial deformities? *Clin Orthop Relat Res*. 2010;468(5):1352-1361. 21. Sökücü S, Karakoyun Ö, Arikani Y, Küçükçaya M, Kabukcuoğlu Y. Efficacy of the Taylor spatial frame in the treatment of deformities around the knee. *Acta Orthop Traumatol Turc*. 2013;47(2):86-90. 22. Tawari G, Maheshwari R, Madan S. Extra-articular deformity correction using Taylor spatial frame prior to total knee arthroplasty. *Strategies Trauma Limb Reconstr*. 2018;Mar 20;doi: 10.1007/s11751-11018-10310-11755 23. Thiryayi W, Naqui Z, Khan S. Use of the Taylor spatial frame in compression arthrodesis of the ankle: a study of 10 cases. *J Foot Ankle Surg*. 2010;49(2):182-187.

*Cover image of Johnnie Yellock used with permission of Dr Joseph Hsu, Atrium Health Carolinas Medical Center.



Johnnie Yellock