

# REGENETEN<sup>†</sup> Bioinductive Implant

## Introduction

The REGENETEN Implant has a large and growing evidence base in rotator cuff repair.<sup>1-12</sup> By promoting the formation of new tendon-like tissue,<sup>1,2,6,12,13</sup> the REGENETEN Implant changes the course of rotator cuff tear progression.<sup>2,12,14,15</sup>

The REGENETEN Implant is resorbed and replaced by tendon-like tissue over 6 months<sup>1,13,†</sup>

5 weeks      6 months


With a substantial body of literature evaluating its use in rotator cuff repair (>650 patients across >10 published studies),<sup>1-12\*</sup> the REGENETEN Implant has been associated with:

<p>Lower re-tear rates for FT tears versus standard repair techniques<sup>5-7,10</sup></p>	<p>Rapid post-operative recovery<sup>9,15,16</sup></p>
<p>Significant and clinically meaningful post-operative improvements in clinical outcomes<sup>3,9,17</sup></p>	<p>Low complication rates<sup>2,3,5-11</sup></p>

Effective treatment of rotator cuff tears can be challenging. Conventional treatment of PT tears commonly involves removing healthy tissue<sup>18</sup> and results in lengthy rehabilitation,<sup>19</sup> whilst repair of large FT tears has been associated with re-tear rates >50%.<sup>16</sup>

**18**  
Studies

- 2 Randomised Controlled Trials (RCTs)
- 9 Human Case Series
- 2 Human Histological Case Studies
- 1 Systematic Literature Review (SLR)
- 1 Cost Analysis
- 3 Registries



**>650**  
total patients studied

\*As of September 2022.  
† On human biopsy (n=1) and in-vivo sampling.

## Overview of REGENETEN<sup>◇</sup> Bioinductive Implant clinical study publications in rotator cuff repair

Study	Design	Tear	n=	Follow-up	Key findings
<p>Ruiz Ibán MA, et al. <b>The Effect on healing rate of the addition of a bioinductive implant to a rotator cuff repair.</b></p> <p>Preliminary report presented at: The European Society for Surgery of the Shoulder and Elbow (SECEC) Annual Congress; September 7–9, 2022; Dublin, Ireland.</p> <p>Please note: These are interim results, final results will be published by the end of 2023.</p>	Randomised controlled trial	FT	Enrolled: 60 (n=30 REGENETEN group, n=30 control group)  Interim follow-up: 57 (29/28)	1 year	<p>Interim results from this study, show that rotator cuff repair with REGENETEN (REGENETEN group) compared to standard rotator cuff repair alone (control group) demonstrated:</p> <ul style="list-style-type: none"> <li>Significantly lower re-tear rates in the REGENETEN group vs control group at 1 year follow-up (3.5% vs 25.0%, p=0.022)</li> <li>No differences in post-operative complications between groups</li> <li>The failure rate at the musculotendinous junction was significantly lower in the REGENETEN group (3.5% vs 22.0%, p=0.044)</li> <li>Post-operative fatty infiltration was lower in the REGENETEN group compared to the control group (10.0% vs 25.0%)</li> </ul>
<p>Ferreira Barros A, et al. <b>Use of bio inductive bovine collagen patch augmentation for full thickness cuff tears - 12-month follow up results of an ongoing prospective randomised trial.</b></p> <p>Preliminary report presented at: The European Society for Surgery of the Shoulder and Elbow (SECEC) Annual Congress; September 7–9, 2022; Dublin, Ireland.</p>	Randomised controlled trial	FT	Enrolled: 56 (n=27 REGENETEN group, n=29 control group)	1 year	<ul style="list-style-type: none"> <li>Lower re-tear rates in the REGENETEN group vs control group at 1 year follow-up (7% vs 13%)</li> <li>Patients in the REGENETEN group experienced greater improvement in function and pain scores from pre-operative levels to 1 year post-operatively: <ul style="list-style-type: none"> <li>CMS increased on average by 40 points compared to the control group which increased by 30 (REGENETEN group average scores; 49 pre-operatively to 89 1 year post-operatively vs control group; 52 pre-operatively to 82 1 year post-operatively)</li> <li>Greater reduction in VAS pain score was seen in the REGENETEN group vs control (6.8 vs 6.2) (REGENETEN group average scores; 7.5 pre-operatively to 0.7 1 year post-operatively vs control group; 7.2 at pre-operatively to 1.0 1 year post-operatively)</li> </ul> </li> <li>Three cases of shoulder stiffness/adhesive capsulitis (two in the REGENETEN group, one in the control group) were reported</li> </ul>

Study	Design	Tear	n=	Follow-up	Key findings
Camacho-Chacon JA, et al. <b>Bioinductive collagen implants facilitate tendon regeneration in rotator cuff tears.</b> <i>J Exp Orthop.</i> 2022;9(1):53. Available at: <a href="#">Journal of Experimental Orthopaedics</a>	Case series	PT + FT	30 (enrolled), 29 (at follow-up)	1 year	<ul style="list-style-type: none"> <li>No evidence of inflammatory, scarring or ischaemic changes in histological examination of all samples at 6 months post-operatively <ul style="list-style-type: none"> <li>No evidence of any remaining REGENETEN Implant</li> </ul> </li> <li>MRI demonstrated complete healing in 27 patients after 6 months; mean tendon thickness increased significantly by 1.84mm (from 4.18mm pre-operatively to 6.02mm; p=0.001) <ul style="list-style-type: none"> <li>All patients showed filling of the defect at 6 months with 90% completely filled and 10% with greater than 50% of the defect filled</li> <li>No evidence of tear propagation at 1 year</li> <li>Neotendon was indistinguishable from native tendon on MRI in all samples at 6 months and 1 year</li> </ul> </li> <li>Patients experienced statistically significant improvements vs pre-operative values in VAS pain score (p=0.003), ASES (p=0.001) and CMS (p=0.001) at 6 months post-operatively, which were sustained at 1 year</li> </ul>
Bushnell BD, et al. <b>Two-year outcomes with a bioinductive collagen implant used in augmentation of arthroscopic repair of full-thickness rotator cuff tears: Final results of a prospective multi-center study.</b> <i>JSES Int.</i> 2022. [Epub ahead of print]. Available at: <a href="#">Journal of Shoulder and Elbow</a>	Case series	FT	115	2 years	<ul style="list-style-type: none"> <li>Between baseline and 2 year follow-up, the average mean tendon thickness increased by 12.5% for medium tears and 17.1% for large tears</li> <li>Overall re-tear rate of medium tears; 4.6% at 3 months (3/66) and 10.6% at 1 year and 2 years (7/66)</li> <li>Overall re-tear rate of large tears; 20.4% at 3 months (10/49), 24.5% at 1 year (12/49) and 28.6% at 2 years (14/49)</li> <li>Significantly lower re-tear rate in patients who had repair with double-row (12/91, 13.2%) technique at 2 years, compared to single-row technique (p=0.0061) <ul style="list-style-type: none"> <li>Re-tear rates with single-row (9/24, 37.5%) technique at 2 years</li> </ul> </li> <li>More than 90% of patients had significant post-operative improvements in both medium and large tears in ASES Shoulder and CMS scores that exceeded respective MCIDs (p&lt;0.001)</li> <li>Most patients (97.1%) surveyed were satisfied with the procedure; 100% of patients surveyed would recommend the procedure to a friend</li> <li>Nine re-operations were required: seven for persistent symptoms, one for superficial infection and one for inflammation and osteopenia</li> </ul>

Study	Design	Tear	n=	Follow-up	Key findings
<p>Nherera, L. <b>Early economic analysis of the REGENETEN Bioinductive Implant in the treatment of full-thickness rotator cuff tears.</b></p> <p>Poster presented at: ISPOR Conference, May 15–18, 2022; Washington, D.C., USA.</p> <p>Available at: <a href="#">Value in Health</a></p>	Economic analysis	FT	NA	NA	<ul style="list-style-type: none"> <li>In the base case analysis, the REGENETEN Implant + SoC was dominant and resulted in reduced total costs and improved clinical outcomes, compared to SoC alone</li> <li>Sensitivity analyses showed that the REGENETEN Implant + SoC remained the dominant intervention from a societal perspective in FT tears of all sizes, and cost-effectiveness improved with tear size: <ul style="list-style-type: none"> <li>Medium: \$31,451 per healed tear, \$152,469 per QALY</li> <li>Large: \$18,555 per healed tear, \$89,953 per QALY</li> <li>Massive: \$4,275 per healed tear, \$20,725 per QALY</li> </ul> </li> <li>A one-way sensitivity analysis demonstrated increased cost-effectiveness, in comparison to the base case analysis, in patients at a higher risk of re-tear due to the following risk factors: age &gt;60 years, hypertension, alcohol consumption and obesity</li> </ul>
<p>Bushnell BD, et al. <b>Treatment of partial-thickness rotator cuff tear repairs with a resorbable bioinductive bovine collagen implant: 1-year results from a prospective multi-center registry.</b> <i>Orthop J Sports Med.</i> 2021;9(8): 23259671211027850.</p> <p>Available at (open access): <a href="#">Orthopaedic Journal of Sports Medicine</a></p>	Registry	PT	272 (enrolled), 227 (at follow-up)	1 year	<ul style="list-style-type: none"> <li>All PROMs were significantly improved at 3, 6 and 1 year from pre-operative values, except VR-12 MCS, which was significant at 1 year (p&lt;0.05) <ul style="list-style-type: none"> <li>MCIDs for ASES, SANE and WORC scores were met or exceeded by &gt;90% of patients at 1 year</li> </ul> </li> <li>In higher grade tears (≥grade 2) and compared with the augmented takedown and repair group (n=26), isolated bioinductive repair (REGENETEN Implant without surgical repair; n=201) resulted in: <ul style="list-style-type: none"> <li>Significantly better ASES, SANE and WORC scores at 2 and 6 weeks (p&lt;0.05)</li> <li>No significant differences in 1 year outcomes, except VR-12 PCS, which was significantly improved with isolated bioinductive repair (p=0.0213)</li> <li>Significantly less sling time (19.1 vs 34.3 days; p&lt;0.0001) and faster return to non-overhead sports (72.2 vs 128.9 days; p=0.0192); no other significant differences in recovery</li> </ul> </li> </ul>
<p>McIntyre LF, et al. <b>Full-thickness rotator cuff tears can be safely treated with a resorbable bioinductive bovine collagen implant: one-year results of a prospective, multicenter registry.</b> <i>Arthrosc Sports Med Rehabil.</i> 2021;3(5):e1473-e1479.</p> <p>Available at (open access): <a href="#">Arthroscopy, Sports Medicine and Rehabilitation</a></p>	Registry	FT	210 (enrolled), 192 (at follow-up)	1 year	<ul style="list-style-type: none"> <li>At 6 months and 1 year, ASES, SANE, VR-12 PCS and WORC were significantly improved from pre-operative values (p&lt;0.001) <ul style="list-style-type: none"> <li>MCIDs were met or exceeded by 90.5% of patients for ASES, 84.3% of patients for SANE and 87.2% of patients for WORC were observed</li> <li>VR-12 MCS was significantly improved at 6 months (p=0.002); no significant difference to pre-operative values at 1 year</li> </ul> </li> <li>An ad-hoc analysis demonstrated similar results at 1 year regardless of tear size</li> <li>Mean duration of post-operative recovery (days): sling time, 36.3 (n=188); return to driving, 24.0 (n=135); return to work, 48.4 (n=128); return to non-overhead sports, 105.4 (n=71), return to overhead sports, 131.7 (n=42)</li> <li>18 patients required re-operation†; re-tear (n=11), infection* (n=3), shoulder stiffness/adhesive capsulitis (n=3), bursitis (n=1) or implant displacement after a fall (n=1)</li> </ul>

\*Includes two recurrent infections, which were resolved with treatment.

† Re-operations can be attributed to multiple reasons.

Study	Design	Tear	n=	Follow-up	Key findings
Schlegel TF, et al. <b>Isolated bioinductive repair of partial thickness rotator cuff tears using a resorbable bovine collagen implant: Two-year radiologic and clinical outcomes from a prospective multicentre study.</b> <i>J Shoulder Elbow Surg.</i> 2021;30(8):1938-1948. Available at: <a href="#">Journal of Shoulder and Elbow Surgery</a>	Case series	PT	33 (enrolled), 31 (at follow-up)	2 years	<ul style="list-style-type: none"> <li>At 2 years, 87% of tears had reduced in size by &gt;50% from pre-operative measurements</li> <li>No compliant patient progressed to a full-thickness tear; one non-compliant patient progressed ~1 month after surgery</li> <li>Significant increase in mean tendon thickness from pre-operative values to 2 years: 1.2mm in intermediate-grade tears; 1.8mm in high-grade tears (p&lt;0.012)</li> <li>Improvements in CMS met or exceeded MCIDs for all patients with intermediate-grade tears and 94% of patients with high-grade tears</li> <li>No significant difference in clinical outcomes, tear healing or tendon thickness based on tear location</li> </ul>
Dai A, et al. <b>Collagen-based bioinductive implant for treatment of partial thickness rotator cuff tears.</b> <i>Bull Hosp Jt Dis.</i> 2020;78(3):195-201. Available at (open access): <a href="#">Bulletin of the Hospital for Joint Diseases</a>	Case series	PT	30 (eligible), 24 (at follow-up)	1 year	<ul style="list-style-type: none"> <li>Significant increase in mean ASES score by 22.5 points (45.6 pre-operatively to 68.1 19.1 months post-operatively; p=0.001)</li> <li>Significant decrease in mean VAS pain from 8.3 pre-operatively to 3.8 at 19.1 months (p&lt;0.001)</li> <li>Improvements in ASES score and VAS pain from pre-operative values were greatest in articular-sided and intrasubstance tears</li> <li>Mean patient satisfaction score was 7.5/10</li> <li>Mean tendon thickness significantly increased by 0.8mm (5.7mm preoperatively to 6.5mm 9.9 months post-operatively; p=0.007)</li> <li>No implant-related complications; one patient had a traumatic re-tear 4 months post-operatively</li> </ul>
Thon SG, et al. <b>REGENETEN bio-inductive collagen scaffold for rotator cuff tears: indications, technique, clinical outcomes, and review of current literature.</b> <i>Ann Jt.</i> 2020;5:41. Available at (open access): <a href="#">Annals of Joint</a>	Review	PT + FT	251	Mean: 14.9 months (10.8–32.0 months)	<p>Review article including Bokor 2015, Bokor 2016, Schlegel 2018, Thon 2019, and McIntyre 2019, reporting:</p> <ul style="list-style-type: none"> <li>92–94% patient satisfaction in partial thickness tears (2 studies)</li> <li>94–100% healing rate in partial thickness tears (2 studies)</li> <li>89–91% patient satisfaction in full thickness tears (2 studies)</li> <li>96–100% healing rate in full thickness tears (2 studies)</li> <li>3.9% reoperation rate (10/251; 5 studies)</li> <li>5.9% (15/251) failure rate (lack of patient satisfaction, lack of tendon healing, or need for reoperation) based on 5 studies</li> <li>9.9% (25/251) complication rate (5 studies)</li> </ul>

Study	Design	Tear	n=	Follow-up	Key findings
<p>Micheloni GM, et al. <b>Bio-inductive implant for rotator cuff repair: our experience and technical notes.</b> <i>Acta Biomed.</i> 2020;91:14-5.</p> <p>Available at (open access): <a href="#">Acta Biomedica</a></p>	Case series	PT + FT	4	6 weeks	<ul style="list-style-type: none"> <li>In 4 patients (1 PT tear receiving isolated bioinductive repair, 3 FT tears receiving a REGENETEN Implant to augment surgical repair), no complications were found at 6 weeks follow-up</li> <li>Increase in procedure duration by 10 minutes</li> </ul>
<p>McIntyre LF, et al. <b>Patient-reported outcomes after use of a bioabsorbable collagen implant to treat partial and full-thickness rotator cuff tears.</b> <i>Arthroscopy.</i> 2019;35(8):2262-2271.</p> <p>Available at: <a href="#">Arthroscopy</a></p>	Registry	PT + FT	203 (eligible), 173 (at follow-up)	1 year	<p>For patients with partial-thickness tears (n=90):</p> <ul style="list-style-type: none"> <li>MCIDs were achieved in VAS pain from 2 weeks (5.3 at baseline to 3.3 at 2 weeks and 1.1 at 1 year; p&lt;0.001) and ASES score from 6 weeks (47.0 at baseline to 60.6 at 6 weeks and 85.6 at 1 year; p&lt;0.001)</li> <li>Mean VAS pain and ASES score at 3 months were significantly improved compared to published results from patients undergoing transtendon or takedown repair without a REGENETEN Bioinductive Implant (p&lt;0.05)</li> <li>Significant improvements were observed from 6 weeks in SANE score (42.5 at baseline to 59.4 at 6 weeks and 86.6 at 1 year; p&lt;0.001), VR-12 PCS (35.8 at baseline to 39.1 at 6 weeks and 49.7 at 1 year; p&lt;0.002) and WORC Index (38.2 at baseline to 53.5 at 6 weeks and 84.4 at 1 year; p&lt;0.001)</li> <li>Time in sling and time to return to driving and sport were found to be shorter than published results for standard surgical techniques</li> </ul> <p>For patients with full-thickness tears (n=83):</p> <ul style="list-style-type: none"> <li>MCIDs were achieved in VAS pain from 2 weeks (5.2 at baseline to 3.7 at 2 weeks and 1.2 at 1 year; p&lt;0.001) and ASES score from 3 months (45.5 at baseline to 68.4 at 3 months and 83.8 at 1 year; p&lt;0.001)</li> <li>Significant improvements from 6 weeks in WORC index (35.0 at baseline to 41.1 at 6 weeks and 80.1 at 1 year; p&lt;0.001), and from 3 months in SANE score (39.2 at baseline to 63.3 at 3 months and 80.7 at 1 year; p&lt;0.001) and VR-12 PCS (34.5 at baseline to 40.8 at 3 months and 45.7 at 1 year; p&lt;0.001)</li> <li>Shorter time to return to driving and sport compared to published results for standard surgical techniques</li> </ul>

Study	Design	Tear	n=	Follow-up	Key findings
<p>Thon SG, et al. <b>Evaluation of healing rates and safety with a bioinductive collagen patch for large and massive rotator cuff tears; 2-year safety and clinical outcomes.</b> <i>Am J Sports Med.</i> 2019;47(8):1901-1908.</p> <p>Available at: <a href="#">American Journal of Sports Medicine</a></p>	Case series	FT	23	2 years	<ul style="list-style-type: none"> <li>No implant-related adverse events</li> <li>Tendon healing in 22/23 patients (96%) on ultrasound at 2 years</li> <li>Treatment success in 21/23 patients (91%) at 2 years; one healing failure and one failure due to progression of glenohumeral osteoarthritis</li> <li>Mean tendon thickness on ultrasound increased from 6.29mm at 3 months to 7.72mm at 1 year, decreasing to 7.28mm at 2 years</li> </ul>
<p>Bokor DJ, et al. <b>Healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a highly-porous collagen implant: a 5-year clinical and MRI follow-up.</b> <i>Muscles Ligaments Tendons J.</i> 2019;9(3):338-347.</p> <p>Available at (open access): <a href="#">Muscles, Ligaments and Tendons Journal</a></p>	Case series	PT	11	5 years	<ul style="list-style-type: none"> <li>8/11 patients (73%) had no decline in tendon integrity between 2 and 5 years</li> <li>Mean tendon thickness significantly decreased between 2 and 5 years (5.9 vs 5.2mm; p=0.0012), but remained significantly greater than pre-operative values (4.3mm; p&lt;0.0001)</li> <li>Significant improvements from baseline in ASES total score and ASES pain score were sustained to 5 years (p≤0.01) and were not significantly different to 2 year values</li> </ul>
<p>Schlegel TF, et al. <b>Radiologic and clinical evaluation of a bioabsorbable collagen implant to treat partial-thickness tears: a prospective multicenter study.</b> <i>J Shoulder Elbow Surg.</i> 2018;27(2):242-251.</p> <p>Available at: <a href="#">Journal of Shoulder and Elbow Surgery</a></p>	Case series	PT	33	1 year	<ul style="list-style-type: none"> <li>Reduction in tear size of ≥1 grade from baseline in 31/33 patients (94%) at 1 year <ul style="list-style-type: none"> <li>8/33 patients (24%) had no visible defect</li> </ul> </li> <li>Significant increase in mean tendon thickness on 1 year MRI versus pre-operative values (p&lt;0.0001)</li> <li>Improvements in ASES pain, ASES shoulder index and CMS were significant at 1 year compared with pre-operative values (p&lt;0.0001) <ul style="list-style-type: none"> <li>Improvements in each score were greater than twice the respective MCIDs</li> </ul> </li> <li>30/33 patients (94%) were satisfied with the results of their procedure at 1 year</li> </ul>

Study	Design	Tear	n=	Follow-up	Key findings
<p>Arnoczky SP, et al. <b>Histologic evaluation of biopsy specimens obtained after rotator cuff repair augmented with a highly porous collagen implant.</b> <i>Arthroscopy</i>. 2017;33(2):278-283.</p> <p>Available at: <a href="#">Arthroscopy</a></p>	Case series	FT	7	Up to 9 months	<ul style="list-style-type: none"> <li>At the earliest time period (5 weeks), the biopsy showed host cell ingrowth and early collagen formation</li> <li>At 3 months, there was increased collagen formation, maturation and organization on the surface of the implant</li> <li>By 6 months, the implant was no longer visible, with new tendon-like tissue and oriented collagen indicative of functional loading</li> <li>No evidence of foreign body or inflammatory reactions at any time point</li> </ul>
<p>Bokor DJ, et al. <b>Evidence of healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a collagen implant: a 2-year MRI follow-up.</b> <i>Muscles Ligaments Tendons J</i>. 2016;6(1):16-25.</p> <p>Available at (open access): <a href="#">Muscles, Ligaments and Tendons Journal</a></p>	Case series	PT	13	2 years	<ul style="list-style-type: none"> <li>At 3 months after surgery there was a significant increase in mean tendon thickness of 2.2mm at 3 months versus pre-operative values (<math>p &lt; 0.0001</math>)</li> <li>At 2 years, new tissue was indistinguishable from underlying tissue in 12/13 patients (92%)</li> <li>Tendon thickness at 2 years was significantly greater (<math>p &lt; 0.0001</math>) than pre-operative values</li> <li>At 1 year, all assessable patients had a reduction in tear size of <math>\geq 1</math> grade, with complete tear disappearance in 7 of 10 patients with a measurable pre-operative tear size</li> <li>Significant improvement in clinical scores over the 2 year follow-up period from pre-operative values: CMS, Constant-Murley pain score, ASES total score and ASES pain score (<math>p \leq 0.01</math>)</li> <li>Outcomes were satisfactory for 12/13 patients (92%) at 2 years</li> </ul>
<p>Bokor DJ, et al. <b>Preliminary investigation of a biological augmentation of rotator cuff repairs using a collagen implant: a 2-year MRI follow-up.</b> <i>Muscles Ligaments Tendons J</i>. 2015;5(3):144-150.</p> <p>Available at (open access): <a href="#">Muscles, Ligaments and Tendons Journal</a></p>	Case series	FT	9	2 years	<ul style="list-style-type: none"> <li>No MRI evidence of re-tear or gap formation, with the integrity of all repaired tendons intact at 2 years</li> <li>Significant increase in mean tendon thickness versus published values at 3, 6, 12 and 24 months (<math>p &lt; 0.01</math>), with an average of 2mm new tissue over the bursal surface</li> <li>New tissue rapidly matured, improved in quality and was indistinguishable from the native tendon by 1 year</li> <li>From 1 to 2 years, tendon thickness slightly decreased, likely reflecting continued functional remodeling</li> <li>Significant improvement in clinical scores at 2 years versus pre-operative measures: CMS and Constant-Murley pain score (both <math>p &lt; 0.001</math>); ASES score and ASES pain score (both <math>p &lt; 0.001</math>)</li> <li>Outcomes were satisfactory for 8/9 patients (89%) at 2 years</li> </ul>

## References

1. Arnoczky SP, Bishai SK, Schofield B, et al. Histologic evaluation of biopsy specimens obtained after rotator cuff repair augmented with a highly porous collagen implant. *Arthroscopy*. 2017;33:278–283.
2. Bokor DJ, Sonnabend D, Deady L, et al. Evidence of healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a collagen implant: a 2-year MRI follow-up. *Muscles Ligaments Tendons J*. 2016;6(1):16–26.
3. Bushnell BD, Bishai SK, Krupp RJ, et al. Treatment of partial-thickness rotator cuff tears with a resorbable bioinductive bovine collagen implant: 1-year results from a prospective multicenter registry. *Orthop J Sports Med*. 2021;9:23259671211027850.
4. Camacho-Chacon JA, Cuenca-Espierrez J, Roda-Rojo V, et al. Bioinductive collagen implants facilitate tendon regeneration in rotator cuff tears. *J Exp Orthop*. 2022;9:53.
5. Ruiz Iban M. The effect on healing rate of the addition of a bioinductive implant to a rotator cuff repair. Preliminary report presented at: European Society for Surgery of the Shoulder and Elbow (SECEC) Annual Congress; September 7–9, 2022; Dublin, Ireland.
6. Thon SG, O'Malley L 2nd, O'Brien MJ, Savoie FH 3rd. Evaluation of healing rates and safety with a bioinductive collagen patch for large and massive rotator cuff tears: 2-year safety and clinical outcomes. *Am J Sports Med*. 2019;47:1901–1908.
7. Bokor DJ, Sonnabend D, Deady L, et al. Preliminary investigation of a biological augmentation of rotator cuff repairs using a collagen implant: a 2-year MRI follow-up. *Muscles Ligaments Tendons J*. 2015;5:144–150.
8. Dai A, Campbell A, Bloom D, Baron S, Begly J, Meislin R. Collagen-based bioinductive implant for treatment of partial thickness rotator cuff tears. *Bull Hosp Jt Dis*. 2020;78:195–201.
9. McIntyre LF, McMillan S, Trenhaile SW, Bishai SK, Bushnell BD. Full-thickness rotator cuff tears can be safely treated with a resorbable bioinductive bovine collagen implant: one-year results of a prospective, multicenter registry. *Arthrosc Sports Med Rehabil*. 2021;3:e1473–e1479.
10. Bushnell BD, Connor PM, Harris HW, Ho CP, Trenhaile SW, Abrams JS. Retear rates and clinical outcomes at 1 year after repair of full-thickness rotator cuff tears augmented with a bioinductive collagen implant: a prospective multicenter study. *JSES Int*. 2021;5:228–237.
11. Micheloni GM, Salmaso G, Zeccchinato G, Giaretta S, Barison E, Momoli A. Bio-inductive implant for rotator cuff repair: our experience and technical notes. *Acta Biomed*. 2020;91:e2020004.
12. Schlegel TF, Abrams JS, Bushnell BD, Brock JL, Ho CP. Radiologic and clinical evaluation of a bioabsorbable collagen implant to treat partial-thickness tears: a prospective multicenter study. *J Shoulder Elbow Surg*. 2018;27:242–251.
13. Van Kampen C, Arnoczky S, Parks P, et al. Tissue-engineered augmentation of a rotator cuff tendon using a reconstituted collagen scaffold: a histological evaluation in sheep. *Muscles Ligaments Tendons J*. 2013;3:229–235.
14. Bokor DJ, Sonnabend DH, Deady L, et al. Healing of partial-thickness rotator cuff tears following arthroscopic augmentation with a highly-porous collagen implant: a 5-year clinical and MRI follow-up. *Muscles Ligaments Tendons J*. 2019;9:338–347.
15. McIntyre LF, Bishai SK, Brown PB 3rd, Bushnell BD, Trenhaile SW. Patient-reported outcomes after use of a bioabsorbable collagen implant to treat partial and full thickness rotator cuff tears. *Arthroscopy*. 2019;35:2262–2271.
16. Rashid MS, Cooper C, Cook J, et al. Increasing age and tear size reduce rotator cuff repair healing rate at 1 year. *Acta Orthop*. 2017;88:606–611.
17. Schlegel TF, Abrams JS, Angelo RL, Getelman MH, Ho CP, Bushnell BD. Isolated bioinductive repair of partial-thickness rotator cuff tears using a resorbable bovine collagen implant: two-year radiologic and clinical outcomes from a prospective multicenter study. *J Shoulder Elbow Surg*. 2021;30:1938–1948.
18. Lo IK, Burkhart SS. Transtendon arthroscopic repair of partial-thickness, articular surface tears of the rotator cuff. *Arthroscopy*. 2004;20:214–220.
19. Millett PJ, Wilcox RB, O'Holleran JD, Warner JJP. Rehabilitation of the rotator cuff: an evaluation-based approach. *J Am Acad Orthop Surg*. 2006;14:599–609.

### Abbreviations

ASES = American Shoulder and Elbow Surgeons, CMS = Constant-Murley score, FT = full-thickness, MCIDs = minimal clinically important differences, MRI = magnetic resonance imaging, NA = not applicable, PROMs = patient-reported outcome measures, PT = partial-thickness, QALY = quality-adjusted life years, SANE = single assessment numeric evaluation, SoC = standard of care, VAS = visual analogue scale, VR-12 MCS = Veterans RAND 12 Item Health Survey Mental Component Score, VR-12 PCS = Veterans RAND 12 Item Health Survey Physical Component Score, WORC = Western Ontario Rotator Cuff