



+ Evidence in focus

# POLAR3<sup>◇</sup> clinical evidence summary

**POLAR3**  
Total Hip Solution

February 2024


**Smith+Nephew**



# Product summary: POLAR3<sup>◇</sup>

5 studies reporting on POLAR3


POLAR3 delivers excellent short- to mid-term survivorship



POLAR3 is associated with strong clinical, functional and radiological outcomes

- Lower revision risk due to aseptic loosening of either stem or cup<sup>4</sup>
- Low postoperative complication rate when used in bikini incision DAA THA<sup>3</sup>
- Strong medium-term functional and radiographic outcomes in younger patients<sup>5</sup>

POLAR3 results in excellent PROMs and better than class average patient satisfaction



\*This study included 865 cases of bikini incision DAA THA using POLARSTEM/CPCS<sup>◇</sup> stem with R3 and OXINIUM<sup>◇</sup>.


Abbreviations: CPCS = collarless polished cemented stem; DAA = direct anterior approach; EQ-5D = EuroQol-5 Dimensions; EQ-VAS = EuroQol Visual Analogue Scale; HHS = Harris Hip Score; NICE = National Institute for Health and Care Excellence; NJR = National Joint Registry; OHS = Oxford Hip Score; PROMs = patient reported outcome measures; THA = total hip arthroplasty.




Product summary: **POLARSTEM**

9 studies reporting on POLARSTEM

POLARSTEM delivers excellent short- to mid-term survivorship




POLARSTEM demonstrates excellent clinical outcomes



POLARSTEM results in improved PROMs and patient satisfaction

POLARSTEM results in significantly improved average 6-month OHS, EQ-5D and EQ-VAS scores compared with other uncemented stems ( $p < 0.001$ ),<sup>9</sup> and improvements from pre- to postoperative HHS and WOMAC scores<sup>6,8</sup>



\*A pre-entry benchmark introduced for implants being introduced through beyond compliance.

Abbreviations: EQ-5D = EuroQol-5 Dimensions; EQ-VAS = EuroQol Visual Analogue Scale; HHS = Harris Hip Score; NJR = National Joint Registry; OHS = Oxford Hip Score; PROMs = patient reported outcome measures; THA = total hip arthroplasty; WOMAC = Western Ontario and McMaster Universities Osteoarthritis Index.




# Product summary: OXINIUM<sup>+</sup>/XLPE

6 studies reporting on OXINIUM/XLPE

OXINIUM/XLPE demonstrates superior mid- to long-term survivorship compared with other bearing combinations

Highest survivorship of all bearing combinations analysed across multiple registries (AOANJRR, NJR, LRO)<sup>\*12-14</sup>




OXINIUM/XLPE results in superior clinical and health economic outcomes compared with non-OXINIUM bearings

- Lower mean hospital LoS<sup>15,16</sup>
- Lower 30-day readmission rate<sup>15,16</sup>
- Greater number of patients discharged to home health care<sup>15,16</sup>
- Lower likelihood of blood transfusion<sup>15</sup>

OXINIUM/XLPE is associated with:

- A **38%** lower revision risk from 1 year compared with Metal/XLPE<sup>12</sup>
- A **19%** lower revision risk compared with Metal/PE<sup>14</sup>
- A **60%** lower risk of revision relative to the NICE benchmark at 10 years<sup>2,13</sup>



OXINIUM/XLPE is associated with a significant **\$602** reduction in 90-day episode of care costs when compared with non-OXINIUM bearings (p=0.003)<sup>16</sup>

\*Period of registry analysis for AOANJRR: 1999–2022; NJR: 2004–2016; LRO: 2007–2016. †OXINIUM/(XL)PE includes grouped data for OXINIUM bearings on XLPE or PE.


Abbreviations: AOANJRR = Australian Orthopaedic Association National Joint Replacement Registry; Metal/PE = metal on polyethylene; Metal/XLPE = metal on crosslinked polyethylene; NICE = National Institute for Health and Care Excellence; NJR = National Joint Registry; OXINIUM/XLPE = ceramicised metal on XLPE; LRO = Landelijke Registratie Orthopedische Implantaten (Dutch Arthroplasty Register).



# Product summary: R3<sup>◇</sup>

8 studies reporting on R3

R3 delivers excellent short- to mid-term survivorship



R3 demonstrates excellent clinical outcomes

- Low rates of intraoperative complications<sup>11</sup>
- Low rate of postoperative complications<sup>18</sup>
- Strong radiographic outcomes<sup>21</sup>
- Low liner dissociation rate when compared with Pinnacle<sup>17</sup>

**ODEP 13A\***  
R3 uncemented cup<sup>10</sup>

**1.4% lower cumulative 14-year revision rate compared with all other uncemented cups<sup>20</sup>**

R3 results in improved PROMs and patient satisfaction

R3 results in significantly improved average 6-month OHS, EQ-5D and EQ-VAS scores compared with other uncemented cups ( $p < 0.001$ ),<sup>20</sup> and improvements from pre- to postoperative HHS and WOMAC scores<sup>5,19</sup>

**Patients who reported feeling much better/a little better post-THA<sup>20</sup>**

96.4%

96.0%

p=n.s.

**Patients who rated their post-THA satisfaction as 'excellent/very good'<sup>20</sup>**

81.1%

78.9%

p<0.001

■ R3    ■ Other uncemented cups

\*A pre-entry benchmark introduced for implants being introduced through beyond compliance.

Abbreviations: EQ-5D = EuroQol-5 Dimensions; EQ-VAS = EuroQol Visual Analogue Scale; HHS = Harris Hip Score; NJR = National Joint Registry; n.s. = not significant; OHS = Oxford Hip Score; PROMs = patient reported outcome measures; THA = total hip arthroplasty; WOMAC = Western Ontario and McMaster Universities Osteoarthritis Index.



# POLAR Total Hip Solution delivers high survivorship at 11 years, higher than class average patient satisfaction and improvements in PROMS<sup>4</sup>



## Overview

- Bespoke implant report produced by the UK NJR summarising usage and outcomes associated with POLAR3, the combination of POLARSTEM<sup>o</sup>, OXINIUM<sup>o</sup> femoral head, highly bearing and R3<sup>o</sup> cup
- The analysis is based on data collected by the NJR and on PROMS data collected by NHS Digital
- Reported POLAR3 usage between July 2007 and September 2023 for:
  - 25,619 THR
  - 23,087 total patients
  - 497 implanting surgeons at 104 centres

## Considerations

The data used for this analysis was obtained from the NJR Supplier Feedback System. The HQIP and/or the NJR take no responsibility for the accuracy, currency, reliability and correctness of any data used or referred to in this report, nor for the accuracy, currency, reliability and correctness of links or references to other information sources and disclaims all warranties in relation to such data, links and references to the maximum extent permitted by legislation.


## Results



**ODEP<sup>10</sup>**  
 10A\* POLARSTEM  
 Cementless Stem  
 13A\* R3 Cementless Cup

Patients who received POLAR3...

Were **significantly more likely to be satisfied** with their THR and to consider their procedure a success (p<0.001), after 6 months, compared to those who received all other THR with cementless stems



Achieved **significantly better average 6-month PROM scores** compared to patients who received other types of cementless stems (p<0.001)

PROMs	POLAR3	All other THR with NJR cementless stems	p value
Oxford Hip Score	23.1 (22.9–23.3)	21.7 (21.6–21.7)	p<0.001
EQ-5D	0.472 (0.466–0.479)	0.443 (0.442–0.444)	
EQ-VAS	14.7 (14.3–15.2)	12.3 (12.3–12.4)	

Adjusted health gain scores (95% confidence interval) at 6 months after THR. Adjusted scores correspond to the NHS Digital version 3 case-mix-adjustment model. Using case-mix adjusted scores allows for a more accurate comparison between groups by taking into account variations in patient characteristics.

**Conclusions**  
 POLAR3 delivers 11-year survivorship of 97.8%. It also delivers significantly higher patient satisfaction, success outcomes and improvements in PROMs compared to the class average for cementless stems in all other THR patients.

\*A pre-entry benchmark introduced for implants being introduced through beyond compliance.

Abbreviations: EQ-5D= EuroQol-5 Dimensions; EQ-VAS= EuroQol Visual Analogue Scale; HQIP=Healthcare Quality Improvement Partnership; NJR=national joint registry; ODEP= Orthopaedic Data Evaluation Panel; PROMS=patient reported outcome measures; THR=total hip replacements; XLPE=crosslinked polyethylene.

# Comparison of survivorship of cementless constructs in primary total hip arthroplasty THA

Available at: [NJR 20th Annual Report 2023](#)  

## Overview

- The NJR reported on primary THA between 1 April 2003 and 31 December 2022, with a follow-up at maximum of 19 years
- Cementless THA accounted for 37% of all primary THA (n=1,448,541)
- All cementless THA prostheses with at least 10 years follow-up, ≥10,000 procedures and >250 implants at risk included in this analysis

## Acknowledgments

We thank the patients and staff of all the hospitals in England, Wales and Northern Ireland who have contributed data to the National Joint Registry. We are grateful to the HQIP, the NJR Steering Committee and staff at the NJR Centre for facilitating this work. The views expressed represent those of the authors and do not necessarily reflect those of the NJR Steering Committee or the HQIP who do not vouch for how the information is presented.

## Results



**98.07%**  
survivorship with  
POLARSTEM®/R3<sup>®</sup>

**57%**  
↓  
Lower revision risk  
with POLARSTEM/R3  
compared to class average  
for cementless hip  
combinations at 10 years<sup>1</sup>

**ODEP<sup>10</sup>**  
10A\* POLARSTEM  
Cementless Stem  
13A\* R3 Cementless Cup




Figure. Cumulative revision of THA prostheses with at least 10 years follow-up, ≥10,000 procedures and >250 implants at risk

## Conclusions

POLARSTEM with R3 was associated with a 1.93% revision rate at 10 years. This revision rate corresponds to a 57% reduction in revision risk compared with the class average revision rate for cementless constructs in the NJR.

\*A pre-entry benchmark introduced for implants being introduced through beyond compliance.

Abbreviations: HQIP=Healthcare Quality Improvement Partnership; NJR=national joint registry; ODEP= Orthopaedic Data Evaluation Panel; THA=total hip arthroplasty.

# Comparison of survivorship of OXINIUM<sup>®</sup>/XLPE with other bearing combinations in primary total hip arthroplasty<sup>12</sup>

Available at: [AOANJRR Annual Report, 2023](#)  



## Overview

- The AOANJRR reported on 10 bearing surfaces, 8 of which have been used in >5,000 procedures
- Bearing surfaces include three types of femoral head (metal, ceramic, and ceramicised metal [OXINIUM]) and four types of acetabular articular surface (XLPE, non XLPE, ceramic, and metal)

## Considerations

Comparing the rates of revision for these bearings, Ceramicised Metal/XLPE\* has the lowest rate of revision at 15 years. As in previous years, the Registry urges caution in the interpretation of this result. This bearing is a single company product, used with a small number of femoral stem and acetabular component combinations. This may have a confounding effect on the outcome, making it unclear if the lower rate of revision is an effect of the bearing surface or reflects the limited combinations of femoral and acetabular prostheses. Tables and graphs have been reproduced in exact and complete form.

## Results




Number at Risk	0 Yr	1 Yr	3 Yrs	5 Yrs	10 Yrs	15 Yrs	20 Yrs
Ceramic/Ceramic	107126	101283	90732	78607	41948	14629	2623
Ceramic/Non XLPE	9962	9108	7463	6050	3541	2150	784
Ceramic/XLPE	143653	122064	86968	56755	15458	3377	286
Metal/Metal >32mm	14424	14063	13216	11981	9306	4165	100
Metal/Metal ≤32mm	5143	5021	4840	4654	3954	2651	787
Metal/Non XLPE	35807	34413	32041	29306	20949	11722	3369
Metal/XLPE	196515	181105	153647	124529	56134	15691	1385
Ceramicised Metal/XLPE	34564	30720	24139	18672	8012	2230	0



OXINIUM/XLPE has the **lowest revision rate** of all modern bearing combinations

### HR adjusted for age and gender

Ceramic/Ceramic vs Metal/XLPE	Entire Period: HR=0.99 (0.96, 1.03), p=0.770
Ceramic/Non XLPE vs Metal/XLPE	0 - 2Yr: HR=1.16 (1.01, 1.32), p=0.029 2Yr - 3.5Yr: HR=1.48 (1.15, 1.91), p<0.002 3.5Yr - 5Yr: HR=0.85 (0.58, 1.24), p<0.388 5Yr - 8Yr: HR=0.150 (1.18, 1.89), p<0.001 8Yr+: HR=2.67 (2.37, 3.01), p<0.001
Ceramic/XLPE vs Metal/XLPE	0 - 2Yr: HR=1.01 (0.97, 1.06), p=0.562 2Yr+: HR=0.77 (0.72, 0.82), p<0.001
Metal/Metal ≤32mm vs Metal/XLPE	0 - 2Wk: HR=1.27 (0.96, 1.68), p=0.092 2Wk - 1Mth: HR=0.46 (0.31, 0.67), p<0.001 1Mth - 9Mth: HR=0.96 (0.79, 1.15), p=0.636 9Mth - 1.5Yr: HR=2.81 (2.38, 3.33), p<0.001 1.5Yr - 2Yr: HR=4.40 (3.66, 5.29), p<0.001 2Yr - 3Yr: HR=6.46 (5.72, 7.29), p<0.001 3Yr - 8Yr: HR=9.45 (8.90, 10.02), p<0.001 8Yr - 10Yr: HR=5.92 (5.28, 6.63), p<0.001 10Yr - 12Yr: HR=4.86 (4.29, 5.51), p<0.001 12Yr+: HR=3.35 (2.97, 3.77), p<0.001
Metal/Non XLPE vs Metal/XLPE	Entire Period: HR=1.42 (1.30, 1.56), p<0.001
Ceramicised Metal/XLPE vs Metal/XLPE	0 - 6Mth: HR=1.16 (1.06, 1.27), p<0.001 6Mth - 1Yr: HR=1.04 (0.85, 1.28), p=0.694 1Yr+: HR=0.62 (0.56, 0.69), p<0.001

30% 

From 1 year, OXINIUM/XLPE demonstrates significantly lower risk of revision compared to Metal/XLPE (p<0.001)

## Conclusions

OXINIUM/XLPE has the highest survivorship of all bearing combinations up to 19 years.

\*The term 'Ceramicised Metal/XLPE' is equivalent to 'OXINIUM/XLPE', and the term 'Non XLPE' is equivalent to 'CPE'.

Abbreviations: AOANJRR=Australian Orthopaedic Association National Joint Replacement Registry; CPE=conventional polyethylene; Metal/XLPE=metal on crosslinked polyethylene; XLPE=crosslinked polyethylene.



### References

1. National Joint Registry for England, Wales and Northern Ireland. 20th Annual Report, 2023. Available at: <https://reports.njrcentre.org.uk>. Accessed December 8, 2023.
2. National Institute for Health and Care Excellence (NICE). Total hip replacement and resurfacing arthroplasty for end-stage arthritis of the hip. Available at: <https://pathways.nice.org.uk/pathways/joint-replacement/hip-replacement#content=viewnode%3Anodes-surgical-procedures>. Accessed November 11, 2021.
3. Alva A, Nizam I, Gogos S. Minimizing complications in bikini incision direct anterior approach total hip arthroplasty: a single surgeon series of 865 cases. *J Exp Orthop*. 2018;8:1.
4. National Joint Registry for England, Wales and Northern Ireland: Polarstem Cementless (Oxinium/XLPE/R3 cup) implant summary report. October 20, 2023. Copy available from Smith & Nephew on request.
5. Wade R, Shah KA. Functional and radiological outcome of uncemented total hip arthroplasty in young adults –5 year follow-up. *J Orthop*. 2019;18:237–239.
6. Willburger RE, Heukamp M, Lindenlauv TE, Peterlein CD, Schuttler KF. Excellent midterm survival and functional outcomes of a fully hydroxyapatite-coated cementless stem: first results of a prospective multicenter study. *Arthroplast Today*. 2020;6:201–205.
7. Assaf A, Manara JR, Teoh KH, Evans AR. Mid-term clinical results of the cementless R3 cup and Polarstem total hip arthroplasty. *Eur J Orthop Surg Traumatol*. 2019;29:827–833.
8. Cypres A, Fiquet A, Girardin P, et al. Long-term outcomes of a dual-mobility cup and cementless triple-taper femoral stem combination in total hip replacement: a multicenter retrospective analysis. *J Orthop Surg Res*. 2019;14:376.
9. National Joint Registry for England, Wales and Northern Ireland: Polarstem Cementless implant summary report. 8 November 2023. Copy available from Smith & Nephew on request.
10. Orthopaedic Data Evaluation Panel (ODEP). Available at <http://www.odep.org.uk>. Accessed December 8, 2023.
11. Ahmad A, Mirza Y, Evans A, Teoh K. A comparative study between uncemented and hybrid total hip arthroplasty in octogenarians. *J Arthroplasty*. 2018;33:3719–3723.
12. Australian Orthopaedic Association National Joint Replacement Registry Hip, Knee & Shoulder Arthroplasty: 2023 Annual Report. Available at: <https://aoanjrr.sahmri.com/annual-reports-2023>. Accessed October 24, 2023.
13. Davis ET, Pagkalos J, Kopjar B. Effect of bearing surface on survival of cementless and hybrid total hip arthroplasty. *JBJS OA*. 2020;5:e0075.
14. Peters RM, Van Steenberghe LN, Stevens M, Rijk PC, Bulstra SK, Zijlstra WP. The effect of bearing type on the outcome of total hip arthroplasty. *Acta Orthop*. 2018;89:163–169.
15. Duncan S, Patel A, Delhogue G, Patrick C. Can the choice of cementless implants and bearings during total hip arthroplasty have an impact on the overall costs within a bundled payment model? *J Hip Surg*. 2020;4:66–76.
16. Patrick C, Delhogue G, Patel AR. Retrospective analysis of oxidized zirconium bearing surface in hip replacement 90-day episode claims. Poster presented at: International Society for Pharmacoeconomics and Outcomes Research (ISPOR); May 18–22, 2019; New Orleans, Los Angeles, USA.
17. Gwynne-Jones D, Memon A. Acetabular liner dissociation: a comparative study of two contemporary uncemented acetabular components. *Arthroplast Today*. 2020;6:354–359.
18. Yeroushalmi D, Singh V, Maher N, Gabor JA, Zuckerman JD, Schwarzkopf R. Excellent mid-term outcomes with a hemispheric titanium porous-coated acetabular component for total hip arthroplasty: 7–10 year follow-up. *Hip Int*. 2021 Aug 19 doi:10.1177/11207000211040181.
19. Wilson AI, Turgeon TR, Gascoyne RW, Della Valle CJ, McCalden RW. Midterm results of a contemporary, porous-coated acetabular system in patients undergoing primary total hip replacement for degenerative hip disease: a prospective, multicenter study. *J Arthroplasty*. 2020;35:1862–1867.
20. National Joint Registry for England, Wales and Northern Ireland: R3 Cementless Cup implant summary report. 8 November 2023. Copy available from Smith & Nephew on request.
21. Teoh KH, Whitham RDJ, Golding DM, Wong JF, Lee PYF, Evans AR. R3 cup does not have a high failure rate in conventional bearings: a minimum of 5-year follow-up. *J Arthroplasty*. 2018;33:460–463.

Implant Summary Reports (ISRs) are produced by the National Joint Registry of England, Wales and Northern Ireland and summarises usage and outcomes associated with the specific combinations of products. These analyses are based on data collected by the NJR and PROMs data collected by NHS Digital. Disclaimer relating to the use of ISRs: The data used for this analysis was obtained from the NJR Supplier Feedback System. The HQIP and/or the NJR take no responsibility for the accuracy, currency, reliability and correctness of any data used or referred to in this report, nor for the accuracy, currency, reliability and correctness of links or references to other information sources and disclaims all warranties in relation to such data, links and references to the maximum extent permitted by legislation.

We thank the patients and staff of all the hospitals in England, Wales and Northern Ireland who have contributed data to the National Joint Registry. We are grateful to the HQIP, the NJR Steering Committee and staff at the NJR Centre for facilitating this work. The views expressed represent those of the authors and do not necessarily reflect those of the NJR Steering Committee or the HQIP, who do not vouch for how the information is presented

Smith & Nephew, Inc  
1450 Brooks Road  
Memphis,  
TN 38116,  
USA

Smith & Nephew  
Orthopaedics AG, Baar,  
Switzerland is the  
legal manufacturer  
for POLARSTEM<sup>®</sup>

<sup>®</sup>Trademark of Smith+Nephew  
All Trademarks acknowledged  
© 2024 Smith+Nephew  
38567 V2 0224.  
Published February 2024.

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

Developed by Evidence Communications,  
Global Clinical & Medical Affairs  
[www.smith-nephew.com](http://www.smith-nephew.com)

**Smith+Nephew**