

Application guide

IODOSORB Cadexomer Iodine Ointment promotes wound bed preparation by removing barriers to healing from the wound.*¹⁻⁸ With a simple and easy application and removal,^{1,6} IODOSORB Ointment is clinically proven to treat infection in chronic wounds.^{1,8-11}

Applying IODOSORB Ointment



1. Wash hands thoroughly before and after use
- Using gloved hands, clean the wound and the surrounding area with either a gentle stream of water or saline
- DO NOT DRY** the wound surface



- Apply sufficient IODOSORB Ointment to a sterile dressing to form a 3mm deep layer, which conforms to the shape of the wound
- Alternatively, apply IODOSORB Ointment to the wound and smooth into place with a gloved finger, to the shape of the wound and depth of 3mm



- Cover the wound completely with an appropriate secondary dressing
- Apply compression bandaging when appropriate
- After application, the remaining contents of the tube should be discarded due to loss of sterility

Changing IODOSORB Ointment



- IODOSORB Ointment should be changed when it has become saturated with wound fluid and all the iodine has been released
- This is indicated by loss of colour, usually two to three times a week. If the wound is discharging heavily, daily changes may be needed



- If necessary, soak the dressing for a few minutes, and then remove
- Gently remove the IODOSORB Ointment using a stream of water or saline



- Gently blot any excess fluid, leaving the wound surface slightly moist, before reapplying IODOSORB Ointment

*e.g. slough, pus, debris, and biofilm.

IODOSORB[◇] Cadexomer Iodine Ointment

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use.

Advanced Wound Management
Smith & Nephew Medical Ltd
Croxley Park, Building 5,
Lakeside, Hatters Lane, Watford,
Hertfordshire WD18 8YE, UK

www.smith-nephew.com

T +44 (0) 1923 477100
F +44 (0) 1923 477101

◇Trademark of Smith+Nephew
All Trademarks acknowledged
©April 2023 Smith+Nephew
AWM-AWC-39438 | GMC1705



IODOSORB Cadexomer Iodine Ointment

Code	Size	Carton
66001297	20g	2
66001298	10g	4
66001299	40g	1

References: 1. Skog E, Arnesjo B, Troeng T, et al. *British Journal of Dermatology*. 1983;109:77-83. 2. Troeng T, Skog E, Arnesjo B, et al. Stuttgart: Schattauer Verlag; 1983. 3. Holloway GA, Johansen KH, Barnes RW, Pierce GE. *The Western Journal of Medicine*. 1989;151(1):35-38. 4. Malone M, Johani K, Jensen SO, et al. *J Antimicrob Chemother*. 2017;72(7):2093-2101. 5. Hansson C, Persson L-M, Stenquist B, et al. *International Journal of Dermatology*. 1998;37(390-396). 6. Moberg S, Hoffman L, Grennert M-L, Holst A. *Journal of American Geriatrics Society*. 1983;31(6):462-465. 7. Harcup JW, Saul PA. *The British Journal of Clinical Practise*. 1986;40(9):360-364. 8. Lantis JC. Paper presented at: WUWHs; 2016; Florence. 9. Hillstrom L. *Acta Chir Scand Suppl*. 1988;544:53-56. 10. Schwartz JA, Lantis JC 2nd, Gendics C, et al. *Int Wound J*. 2013;10(2):193-199. 11. Malone M, Schwarzer S, Radzieta M, et al. *Int Wound J*. 2019;16(6):1477-1486.