

+ Less pain, less post-op bleeding^{1-4*}

An evidence-based look at **COBLATION[®]**
Technology Intracapsular Tonsillectomy (CIT)

Smith+Nephew



*When compared to a traditional tonsillectomy

Total tonsillectomies are considered the status quo⁵, but there is another way

New evidence has revealed that COBLATION Intracapsular Tonsillectomy (CIT) may offer:

- Pain free **4.2 days sooner**^{6,7*}
- Less time to analgesia free **4.1 days sooner**^{6*}
- Speedier return to normal activity **2.8 days sooner**^{6*}
- Low primary post tonsillectomy hemorrhage (PTH) **Under 1%**¹⁻⁴
- Happy Parents **99% of parents recommend CIT**⁷



0%
readmission
for pain¹⁰


.....
In a recent large study of CIT patients, there were no delayed discharges or readmissions for pain

*When compared to a traditional tonsillectomy

+ There is a new normal

COBLATION[®] Intracapsular Tonsillectomy (CIT)

Since the underlying arteries are not exposed, the risk of PTH is substantially reduced, leading to a safer procedure for patients, with less burden on the healthcare system as a whole.



+ Less pain means less medication consumption

There is a decreased need for post-operative analgesic medication after CIT, shortening the time patients are on pain medication by 4.1 days on average.⁸




4.1 days
less medication
consumption⁹

+ Less pain for your patients...and for you

- CIT patients are on average pain-free **4.2 days** sooner than total tonsillectomy patients.⁸
- CIT results in a less painful recovery throughout the healing process.⁸

Less painful recoveries lead to happier patients, fewer calls for pain control and more time for you to focus on your practice^{2,9}



Less pain
at
1 and 7
days⁸

*When compared to a traditional tonsillectomy

+ Less time hurting, more time playing

On average, CIT patients return to a normal diet and normal activity about 3 days sooner⁸



Normal dinners
3.5 days sooner⁸



Back to work
2.8 days earlier⁸



Back to school
2.8 days earlier⁸



+ Low Post Tonsillectomy Hemorrhage (PTH)

Accepted PTH rates of 3% from total tonsillectomies are no longer the norm.² CIT results in a PTH of under 1% helping to reduce the morbidity of the procedure¹⁻³

A recent Meta-analysis found:

Odds of a secondary hemorrhage are reduced by



79%
with intracapsular tonsillectomy¹⁰

+ 99% of parents would recommend a CIT⁹



Just like patients have a choice in ENTs, you have a choice in what post-op complications you accept

0% PTH.⁷

In a recent study, **1,045** patients underwent Tonsillectomies. Of those who had a **Total Tonsillectomy**, **6.4%** experienced PTH. Those who had **Intracapsular Tonsillectomy**, **0%** had PTH⁷.

⁸When compared to a traditional tonsillectomy

+ A COBLATION[®] Technology Intracapsular Tonsillectomy

Why COBLATION Technology for Intracapsular Tonsillectomy?

Unlike a total tonsillectomy, which removes all tonsil tissue including the capsular plane, therefore exposing the larger blood vessels located in the extracapsular space, an intracapsular tonsillectomy aims to remove all tonsil tissue by ablating the tonsil tissue to the tonsillar capsule.⁹ This provides a protective shield from large nerves and blood vessels, typically found in the muscle below the tonsil bed.

COBLATION technology involves the creation and application of a high-energy field called “glow discharge plasma.” This plasma ablates tissue through a chemical process as highly energized particles in the plasma break down molecules in the tissue.

- Operates at lower temperatures than other radio frequency (RF) based technologies¹¹
- The 100 – 200 microns plasma field allows for precise removal of soft tissue¹¹
- Integrated suction and hemostasis



Smith+Nephew offers a full portfolio of COBLATION wands for a COBLATION Intracapsular Tonsillectomy (CIT) procedure.

For further product information, please consult your Smith and Nephew ENT sales representative.



Learn more at [smith-nephew.com](https://www.smith-nephew.com)

Sports Medicine
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810

www.smith-nephew.com
T +978 749 1000
US Customer Service:
+1 800 343 5717

°Trademark of Smith+Nephew.
All trademarks acknowledged.
©2021 Smith+Nephew. All rights
reserved. Printed in USA.28643 V1 11/21

References

1. Kim JS et al. Can intracapsular tonsillectomy be an alternative to classical tonsillectomy? A meta-analysis. *Otolaryngology - Head and Neck Surgery*. 2017;157(2):178-89.
2. Francis DO et al. Postoperative bleeding and associated utilization following tonsillectomy in children: A systematic review and meta-analysis. *Otolaryngology - Head and Neck Surgery*. 2017;156(3):442-55.
3. Amin N et al. Coblation intracapsular tonsillectomy in children: A prospective study of 1257 consecutive cases with long-term follow-up. *Clinical Otolaryngology*. 46(6):1184-1192, November 2021.
4. Khemani S et al. Coblation® intracapsular tonsillectomy in children with recurrent tonsillitis: Initial experience. HYPERLINK "<https://www.sciencedirect.com/science/journal/01655876>" *International Journal of Pediatric Otorhinolaryngology*. HYPERLINK "<https://www.sciencedirect.com/science/journal/01655876/135/supp/C>" Volume 135, August 2020, 110113.
5. Hoey AW et al. Coblation intracapsular tonsillectomy (tonsillotomy) in children: A prospective study of 500 consecutive cases with long-term follow-up. *Clinical Otolaryngology*. 2017;42(6):1211-7.
6. Stalder KR et al. Coblation in Otolaryngology. SPIE 4949, *Lasers in Surgery: Advanced Characterization, Therapeutics, and Systems XIII*, (12 September 2003); doi: 10.1117/12.488349.
7. Mitchell RB et al. Clinical practice guideline: Tonsillectomy in children (Update). *Otolaryngology - Head and Neck Surgery*. 2019;160(1 Suppl):S1-S42.
8. Post-Tonsillectomy Pain Management for Children: Education for Caregivers. ENThealth. [Online] [Cited: January 18, 2021.] https://www.enthealth.org/be_ent_smart/post-tonsillectomy-pain-management-for-children-education-for-caregivers/.
9. Partial tonsillectomy using COBLATION versus alternative tonsillectomy techniques: A systematic literature review with meta-analysis.
10. Huoh et al. Current status and future trends: Pediatric intracapsular tonsillectomy in the United States. *The Laryngoscope*. 2020;131(Supple 2):S1-S9. EA/ENT/COBLATION/002/V2.
11. Zhang LY et al. Tonsillectomy or tonsillotomy? A systematic review for paediatric sleep-disordered breathing. *International Journal of Pediatric Otorhinolaryngology*. 2017;103:41-50.