Publication summary: Smith LA, et al. J Arthroplasty (2020)\*

# **Smith**Nephew

Femoral geometry and the cruciate ligaments may play a significant role in achieving more normal-like kinematic patterns post-total knee arthroplasty (TKA)

### 🕂 Plus points

experienced more normal-like kinematic patterns compared to The authors stated JOURNEY II CR knees demonstrated an improvement in lateral femoral rollback and axial rotation compared to previous studies on CR knees

• All TKA subjects were selected as a result of having a well-

Each patient was asked to perform a weight-bearing deep knee

functioning TKA, with a Knee Society Score (KSS) ≥80

Kinematics were measured at full extension and at 30°

bend (full extension to full flexion; Figure)

# **Overview**

- First study to assess the impact of the ACL, PCL and femoral condylar geometry on kinematic patterns following TKA and comparing this to the normal knee
  - 40 JOURNEY II CR knees (average follow up, 24.8 months; mean age, 68.8 years)
  - 10 JOURNEY II XR knees (average follow up, 16 months; mean age, 62.3 years)
  - 10 normal knees (mean age, 57.4 years)

### Results

Compared to the normal knee:

Early flexion (0-30°; ACL function)

- JOURNEY II XR demonstrated a similar magnitude of posterior femoral rollback (PFR)
- JOURNEY II CR demonstrated PFR, but significantly less than the normal knee (lateral, p=0.004; medial, p=0.002)
- Both JOURNEY II XR and CR exhibited external femoral rotation, but to a lower magnitude

Mid flexion (30-60°; ACL/PCL translation)

- JOURNEY II XR demonstrated no statistical difference in anterior-posterior translation or axial rotation
- JOURNEY II CR demonstrated significantly less anteriorposterior translation (lateral, p<0.004; medial, p<0.03) and axial rotation p=0.01)
  - The authors stated that JOURNEY II CR did not experience the significant amount of paradoxical sliding seen in other PCR TKAs

Deeper flexion (60-90°; PCL dominant)

increments to full flexion

- No significant difference with JOURNEY II XR or CR in anteriorposterior translation
- JOURNEY II XR and CR experienced similar axial rotation patterns where the femur externally rotated relative to the tibia

Late flexion (90°+)

JOURNEY II XR and CR demonstrated no significant differences in anterior-posterior translation or axial rotation



Figure. Stages of weight-bearing knee bend which were analysed using fluoroscopy

www.smith-nephew.com/education

# Conclusions

#### Citation

\*Smith LA, Nachtrab J, LaCour M, Cates H, Freeman MG, Komistek RD. In Vivo Knee Kinematics: How Important Are the Roles of Femoral Geometry and the Cruciate Ligaments? J Arthroplasty. 2020: doi: https://doi.org/10.1016/j.arth.2020.10.020. Available at: Journal of Arthroplasty

Smith & Nephew, Inc., 1450 Brooks Road, Memphis, TN 38116, USA. 27996-en V1 1220. Published December 2020. ©2020 Smith+Nephew. All rights reserved. °Trademark of Smith+Nephew. All Trademarks acknowledged.