

LAPAROSCOPES

Surgical Optics

Visualization

Scope Portfolio

Smith+Nephew's new 4KO^o Arthroscopes and Laparoscopes are designed for use with either LENS HD or LENS 4K Surgical Imaging System to increase brightness^{1,2*}, optimize color transmission^{1,2,3**} and improve overall image quality.¹⁻¹⁰







+ Visualization

- Exceptional 4K UHD image quality^{7,8}
- Reduced instances of image clipping on surgical display when using 4KO VAS Scopes compared to HD VAS Scopes⁹
- Crisp Clear Image
- Over 1 billion colors create vivid images and color reproduction^{7,8}
- 4KO & HD Scopes provide improved image quality for both LENS HD and LENS 4K

+ Multispecialty


- Scopes for variety of procedures including general surgery and bariatrics
- Choice of Direct View and VAS coupling
- Direct View Scopes compatible with competitive visualization systems

+ Durability

- Stainless steel construction for durability and sterilization
- Sapphire distal lenses and rear ocular windows for scratch resistance
- Withstands autoclave, EtO and STERRAD™ or STERIS™
- 40+ years of product service experience
- 98% of service requests addressed within 1 day on average

New design features

- Round tip designed to minimize the risk of iatrogenic injury
- Newly designed VAS focus ring designed to improve ease of adjustment during surgical use





^{*}Applies to 4mm DV and VAS Scopes **Applies to 4mm DV and VAS scopes as well as 2.7mm VAS Scopes

♣ 4KO[⋄] Scopes

Optimized Imaging and Design

4KO video arthroscopes have

21%

minimum increase in resolution as compared to S+N high definition video arthroscopes⁴ 4KO direct view arthroscopes have

30%

minimum increase in resolution as compared to S+N high definition direct view arthroscopes⁴ **4KO** arthroscopes have

73%

minimum increase in line pairs per field than the specification of S+N high definition scopes⁴



Before
LENS HD image with S+N
4mm 30° HD scope



After
LENS HD image with S+N
4mm 30° **4KO** scope



Ordering information

Small Joint		
Reference	Description	
7208133	HD VAS, 1.9mm x 30° x 62mm, non-A/C	
72205143	HD VAS, 2.7mm x 30° x 120mm	
72205144	HD VAS, 2.7mm x 30° x 67mm	
72205145	HD VAS, 2.7mm x 70° x 70mm	
4184	HD DV, 1.9mm x 30° x 67mm, non-A/C	
7205681	HD DV, 2.7mm x 30° x 119mm	
7205682	HD DV, 2.7mm x 30° x 67mm	
72200733	HD DV, 2.7mm x 70° x 119mm	
3948	HD DV, 2.7mm x 70° x 70mm, non-A/C	
4144	HD DV, 4.0mm x 30° x 106mm, non-A/C	
Ear + Nose +	Throat	
72205072	HD DV, 4.0mm x 0° x 185mm	
72205073	HD DV, 4.0mm x 30° x 185mm	
72205074	HD DV, 4.0mm x 30° x 185mm, reverse post	
72205075	HD DV, 4.0mm x 45° x 185mm	
72205076	HD DV, 4.0mm x 70° x 185mm	
72205077	HD DV, 4.0mm x 70° x 185mm, reverse post	
Shoulder + H	lip + Knee	
72205149	4KO VAS 4.0mm x 30° x 160mm	
72205150	4KO VAS 4.0mm x 70° x 160mm, round tip	
72202958	HD VAS 4.0mm x 70° x 157mm	
72202963	HD VAS 4.0mm x 30° x 160mm, reverse post, non-A/C	
72202964	HD VAS 4.0mm x 70° x 160mm, reverse post, non-A/C	
72205153	4KO DV 4.0mm x 0° x 160mm	
72205154	4KO DV 4.0mm x 30° x 160mm	
72205155	4KO DV 4.0mm x 70° x 160mm, round tip	

Shoulder + Hip + Knee		
Reference	Description	
72205156	4KO DV 4.0mm x 45° x 160mm	
72201845	HD DV 4.0mm 30° x 157mm, reverse post	
72202088	HD DV 4.0mm x 70° x 160mm	
7208132	HD DV 4.0mm x 30° x 175mm	
72201597	HD DV HIP Scope 4.0mm x 30° x 278mm	
72201598	HD DV HIP Scope 4.0mm x 70° x 280mm	
Laparoscope	es	
72205157	4KO Laparoscope, 5.5mm x 0° x 300mm	
72205158	4KO Laparoscope, 5.5mm x 30° x 300mm	
72205159	4KO Laparoscope, 5.5mm x 45° x 300mm	
72205163	4KO Laparoscope, 10mm x 0° x 330mm	
72205164	4KO Laparoscope, 10mm x 30° x 330mm	
72205165	4KO Laparoscope, 10mm x 45° x 330mm	
Bariatric Lap	paroscopes	
72205160	4KO Bariatric Laparoscope, 5.5mm x 0° x 424mm	
72205161	4KO Bariatric Laparoscope, 5.5mm x 30° x 424mm	
72205162	4KO Bariatric Laparoscope, 5.5mm x 45° x 425mm	
72205166	4KO Bariatric Laparoscope, 10mm x 0° x 455mm	
72205167	4KO Bariatric Laparoscope, 10mm x 30° x 455mm	
72205168	4KO Bariatric Laparoscope, 10mm x 45° x 455mm	

Learn more at smith-nephew.com

Sports Medicine Smith & Nephew, Inc. 150 Minuteman Road Andover, MA 01810 www.smith-nephew.com T +978 749 1000 US Customer Service: +1 800 343 5717 °Trademark of Smith+Nephew. All trademarks acknowledged. ©2020 Smith+Nephew. All rights reserved. Printed in USA. 24976 V1 11/20

References