

＋ ハイクオリティ絆創膏

ハイドロ ジェントルエイド[◇]



皮膚の保護、浅い創傷に




高吸収



クッション性



低刺激粘着
貼り直し可能



痛みやわらぐふんわりパッド

創傷にあたるパッド部分がポリウレタンフォームです。クッション性に優れたふんわりパッドが肌によくなじみ、保護します。



ポリウレタンフォームが滲出液を吸収、ゲル化しません。

最長7日間¹⁾使用

※傷を観察し、滲出液の量に応じて適宜交換をしてください。



剥がす時の痛みが少ない

低刺激のシリコンゲル粘着剤を採用^{2,3,4)}。剥がす時に感じる痛みが少なく^{1,5)}、しわになっても貼り直しができます。




※痛みの感じ方には個人差があります。

粘着が落ちにくいので、しわになったりズレた時に貼り直しができます。

イメージ図

過剰な滲出液を素早く蒸散




細菌の侵入を防ぎます

滲出液を吸収・保持




防水仕様^{6,7)}

イメージ図



低刺激、低アレルギー性のシリコンゲル粘着剤



従来のアクリル系粘着剤

使用例 擦り傷・切り傷・刺し傷・掻き傷・靴ずれなどの浅い傷および周囲皮膚の保護



靴ずれなど



擦り傷や掻き傷など



スキン-テアなど



擦り傷など



骨突出部位など


ハイドロ ジェントルエイド[◇]

一般医療機器 救急絆創膏
届出番号:
13B1X10222WC0016

製品番号	サイズ	パッド部サイズ	枚/箱	標準価格
66800833	5×5 cm	3×3 cm	10	¥1,780
66800834	7.5×7.5 cm	5×5 cm	10	¥2,180
66800835	10×10 cm	7.5×7.5 cm	10	¥2,690
66800840	15×15 cm	12.5×12.5 cm	10	¥6,040
66800836	5.5×12 cm	10	¥2,610	
66800837	8×15 cm	5.5×12.5 cm	10	¥4,000
66801012	10×20 cm	7.2×17.2 cm	10	¥4,860

使用目的又は効果 傷及び周囲皮膚の被覆及び保護 禁忌・禁止* 再使用禁止

*詳細は添付文書をご参照ください。



References

1. Hurd, T et al. A multi-centre-in-market evaluation of ALLEVYN Gentle Border. Wounds UK 2009;5(3):32-44. 2. Stephen Haynes, J et al. Clinical evaluation of ALLEVYN Gentle Border Lite within one health trust. BJN;2011;20:S36,S38-42. 3. Denyer JE. Wound management for children with epidermolysis bullosa. Dermatol Clin. 2010 Apr;28(2):257-64. 4. Walton G. Product Safety Statement 196"ALLEVYN Gentle Border Lite", Smith & Nephew Ltd, Hull 2016. 5. K. Vowden et al. An open prospective randomized, multi-centre clinical evaluation of a hydrocellularsilicone foam dressing in the management of exuding chronic and acute wounds. Presented at the EWMA, Brussels 2011. 6. L. Tompkins. ALLEVYN Gentle Border Lite Dressing Physical Properties. Smith & Nephew Medical Ltd, Hull DS/10/100/R1, Aug 2010. 7. H. Dharma. An open, prospective, randomized, within volunteer, comparison of ALLEVYN Gentle Border and Mepilex Border. Smith & Nephew Medical Ltd, Hull OR-DOF/011, 24 Aug 2010.

製造販売元

スミス・アンド・ネフュー株式会社 ウンドマネジメント事業部

〒105-5114 東京都港区浜松町二丁目4番1号 TEL.03-5403-8830

札幌営業所 TEL.011-736-9596 FAX.011-736-9755
 仙台営業所 TEL.022-276-6726 FAX.022-276-6652
 東京営業所 TEL.03-5403-8915 FAX.03-5403-8916
 名古屋営業所 TEL.052-221-1407 FAX.052-221-1420
 大阪営業所 TEL.06-6399-3303 FAX.06-6399-3993
 広島営業所 TEL.082-502-7322 FAX.082-502-7325
 福岡営業所 TEL.092-452-0142 FAX.092-452-0146

<https://www.smith-nephew.com/ja-jp>

◇Trademark of Smith+Nephew
 © 2019-2023 Smith+Nephew

販売店

HS-5
 202311-3