

# Safety Data Sheet

## SoloSite™ Conformable Wound Dressing

### 1. Identification

#### 1.1. Product identifier

**Product Identity** SoloSite™ Conformable Wound Dressing  
**Alternate Names** SoloSite™ Conformable Wound Dressing

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

**Intended use** Topical Dressing  
**Application Method** See Technical Data Sheet.

#### 1.3. Details of the supplier of the safety data sheet

**Company Name** Smith & Nephew  
970 Lake Carillon Drive, Suite 110  
St. Petersburg, FL 33716

#### Emergency

**Customer Service: Smith & Nephew** 1-800-876-1261

### 2. Hazard(s) identification

#### 2.1. Classification of the substance or mixture

Skin Sens. 1;H317 May cause an allergic skin reaction.

#### 2.2. Label elements

Using the Toxicity Data listed in section 11 and 12 the product is labeled as follows.



### Warning

H317 May cause an allergic skin reaction.

#### [Prevention]:

P261 Avoid breathing dust / fume / gas / mist / vapors / spray.  
P272 Contaminated work clothing should not be allowed out of the workplace.  
P280 Wear protective gloves / eye protection / face protection.

#### [Response]:

P302+352 IF ON SKIN: Wash with plenty of soap and water.  
P313 Get medical advice / attention.  
P321 Specific treatment (see information on this label).  
P333+313 If skin irritation or a rash occurs: Get medical advice / attention.

# Safety Data Sheet

## SoloSite™ Conformable Wound Dressing

P363 Wash contaminated clothing before reuse.

### [Storage]:

No GHS storage statements

### [Disposal]:

P501 Dispose of contents / container in accordance with local / national regulations.

## 3. Composition/information on ingredients

This product contains the following substances that present a hazard within the meaning of the relevant State and Federal Hazardous Substances regulations.

Ingredient/Chemical Designations	Weight %	GHS Classification	Notes
Glycerol CAS Number: 0000056-81-5	10 - 25	Not Classified	[1][2]
Cellulose, carboxymethyl ether, sodium salt CAS Number: 0009004-32-4	1.0 - 10	Not Classified	[1]
1-[1,3-bis(hydroxymethyl)-2,5-dioxoimidazolidin-4-yl]-1,3-bis(hydroxymethyl)urea CAS Number: 0078491-02-8	0.10 - 1.0	Skin Sens. 1;H317 Eye Irrit. 2;H319	[1]

In accordance with paragraph (i) of §1910.1200, the specific chemical identity and/or exact percentage (concentration) of composition has been withheld as a trade secret.

[1] Substance classified with a health or environmental hazard.

[2] Substance with a workplace exposure limit.

[3] PBT-substance or vPvB-substance.

\*The full texts of the phrases are shown in Section 16.

## 4. First aid measures

### 4.1. Description of first aid measures

#### General

In all cases of doubt, or when symptoms persist, seek medical attention.  
Never give anything by mouth to an unconscious person.

#### Inhalation

Remove to fresh air, keep patient warm and at rest. If breathing is irregular or stopped, give artificial respiration. If unconscious place in the recovery position and obtain immediate medical attention. Give nothing by mouth.

#### Eyes

Irrigate copiously with clean water for at least 15 minutes, holding the eyelids apart and seek medical attention.

#### Skin

No first aid should be needed.

#### Ingestion

If conscious, give victim water. Do not induce vomiting. Contact a poison control center for instructions

### 4.2. Most important symptoms and effects, both acute and delayed

#### Overview

##### Acute Effects

**Eye:** Eye irritation.

**Ingestion:** Gastrointestinal irritation.

See section 2 for further details.

#### Skin

May cause an allergic skin reaction.

# Safety Data Sheet

## SoloSite™ Conformable Wound Dressing

### 5. Fire-fighting measures

#### 5.1. Extinguishing media

Carbon dioxide, dry chemical, foam. Use a water spray to cool exposed containers.

#### 5.2. Special hazards arising from the substance or mixture

Hazardous decomposition: Oxides of carbon, phenol

Avoid breathing dust / fume / gas / mist / vapors / spray.

#### 5.3. Advice for fire-fighters

Self-Contained Breathing Apparatus and protective clothing should be worn.

Do not release runoff from fire control methods to sewers or waterways.

**ERG Guide No.** ----

### 6. Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

Put on appropriate personal protective equipment (see section 8).

#### 6.2. Environmental precautions

Do not allow spills to enter drains or waterways.

Use good personal hygiene practices. Wash hands before eating, drinking, smoking or using toilet. Promptly remove soiled clothing and wash thoroughly before reuse.

#### 6.3. Methods and material for containment and cleaning up

**Small Spills:** Wipe up spill. Rinse with water.

**Large spills:** Add salt to break up polymer, then mop or wipe up spill and place in proper container for disposal.

### 7. Handling and storage

#### 7.1. Precautions for safe handling

See section 2 for further details. - [Prevention]:

#### 7.2. Conditions for safe storage, including any incompatibilities

Handle containers carefully to prevent damage and spillage.

Incompatible materials: Strong oxidizing agents, strong bases

Store in a dry place at < 77°F (25°C)

See section 2 for further details. - [Storage]:

#### 7.3. Specific end use(s)

No data available.

# Safety Data Sheet

## SoloSite™ Conformable Wound Dressing

### 8. Exposure controls and personal protection

#### 8.1. Control parameters

##### Exposure

CAS No.	Ingredient	Source	Value
0000056-81-5	Glycerol	OSHA	TWA 15 mg/m <sup>3</sup> (total) TWA 5 mg/m <sup>3</sup> (resp)
		ACGIH	TWA: 3 mg/m <sup>3</sup> (respirable) 10 mg/m <sup>3</sup> (mist)
		NIOSH	No established RELs
		Supplier	No Established Limit
0009004-32-4	Cellulose, carboxymethyl ether, sodium salt	OSHA	No Established Limit
		ACGIH	No Established Limit
		NIOSH	No Established Limit
		Supplier	No Established Limit
0078491-02-8	1-[1,3-bis(hydroxymethyl)-2,5-dioximidazolidin-4-yl]-1,3-bis(hydroxymethyl)urea	OSHA	No Established Limit
		ACGIH	No Established Limit
		NIOSH	No Established Limit
		Supplier	No Established Limit

#### Carcinogen Data

CAS No.	Ingredient	Source	Value
0000056-81-5	Glycerol	OSHA	Select Carcinogen: No
		NTP	Know n: No; Suspected: No
		IARC	Group 1: No; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;
0009004-32-4	Cellulose, carboxymethyl ether, sodium salt	OSHA	Select Carcinogen: No
		NTP	Know n: No; Suspected: No
		IARC	Group 1: No; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;
0078491-02-8	1-[1,3-bis(hydroxymethyl)-2,5-dioximidazolidin-4-yl]-1,3-bis(hydroxymethyl)urea	OSHA	Select Carcinogen: No
		NTP	Know n: No; Suspected: No
		IARC	Group 1: No; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;

#### 8.2. Exposure controls

##### Respiratory

None required for normal use.

##### Eyes

None required for normal use.

##### Skin

None required for normal use.

##### Engineering Controls

Provide adequate ventilation. Where reasonably practicable this should be achieved by the use of local exhaust ventilation and good general extraction. If these are not sufficient to maintain concentrations of particulates and any vapor below occupational exposure limits suitable respiratory protection must be worn.

##### Other Work Practices

Use good personal hygiene practices. Wash hands before eating, drinking, smoking or using toilet. Promptly remove soiled clothing and wash thoroughly before reuse.

See section 2 for further details. - [Prevention]:

# Safety Data Sheet

## SoloSite™ Conformable Wound Dressing

### 9. Physical and chemical properties

<b>Appearance</b>	Gauze impregnated with translucent Gel
<b>Odor</b>	Odorless
<b>Odor threshold</b>	Not Measured
<b>pH</b>	5.3 - 7.0 @ 25°C
<b>Melting point / freezing point</b>	Not determined
<b>Initial boiling point and boiling range</b>	Not applicable
<b>Flash Point</b>	> 200°F / > 93°C TCC
<b>Evaporation rate (Ether = 1)</b>	< 1 (nBuAc = 1)
<b>Flammability (solid, gas)</b>	Not Applicable
<b>Upper/lower flammability or explosive limits</b>	<b>Lower Explosive Limit:</b> Not known <b>Upper Explosive Limit:</b> Not known
<b>Vapor pressure (Pa)</b>	> 10 mmHg
<b>Vapor Density</b>	Not applicable
<b>Specific Gravity</b>	Not determined
<b>Solubility in Water</b>	Insoluble
<b>Partition coefficient n-octanol/water (Log Kow)</b>	Not Measured
<b>Auto-ignition temperature</b>	Not determined
<b>Decomposition temperature</b>	Not available
<b>Viscosity (cSt)</b>	Not determined
<b>VOC Content</b>	Not available

#### 9.2. Other information

No other relevant information.

### 10. Stability and reactivity

#### 10.1. Reactivity

Hazardous Polymerization will not occur.

#### 10.2. Chemical stability

Stable under normal circumstances.

#### 10.3. Possibility of hazardous reactions

No data available.

#### 10.4. Conditions to avoid

Keep away from oxidizing agents, excessive heat, and ignition sources.

#### 10.5. Incompatible materials

Strong oxidizing agents, strong bases

#### 10.6. Hazardous decomposition products

Oxides of carbon, phenol

# Safety Data Sheet

## SoloSite™ Conformable Wound Dressing

### 11. Toxicological information

#### Acute toxicity

Ingredient	Oral LD50, mg/kg	Skin LD50, mg/kg	Inhalation Vapor LC50, mg/L/4hr	Inhalation Dust/Mist LC50, mg/L/4hr	Inhalation Gas LC50, ppm
Glycerol - (56-81-5)	12,600.00, Rat - Category: NA	10,000.00, Rabbit - Category: NA	No data available	No data available	No data available
Cellulose, carboxymethyl ether, sodium salt - (9004-32-4)	No data available	No data available	No data available	No data available	No data available
1-[1,3-bis(hydroxymethyl)-2,5-dioximidazolidin-4-yl]-1,3-bis(hydroxymethyl)urea - (78491-02-8)	No data available	No data available	No data available	No data available	No data available

Note: When no route specific LD50 data is available for an acute toxin, the converted acute toxicity point estimate was used in the calculation of the product's ATE (Acute Toxicity Estimate).

Classification	Category	Hazard Description
Acute toxicity (oral)	---	Not Applicable
Acute toxicity (dermal)	---	Not Applicable
Acute toxicity (inhalation)	---	Not Applicable
Skin corrosion/irritation	---	Not Applicable
Serious eye damage/irritation	---	Not Applicable
Respiratory sensitization	---	Not Applicable
Skin sensitization	1	May cause an allergic skin reaction.
Germ cell mutagenicity	---	Not Applicable
Carcinogenicity	---	Not Applicable
Reproductive toxicity	---	Not Applicable
STOT-single exposure	---	Not Applicable
STOT-repeated exposure	---	Not Applicable
Aspiration hazard	---	Not Applicable

### 12. Ecological information

#### 12.1. Toxicity

No additional information provided for this product. See Section 3 for chemical specific data.

#### Aquatic Ecotoxicity

Ingredient	96 hr LC50 fish, mg/l	48 hr EC50 crustacea, mg/l	Er C50 algae, mg/l
Glycerol - (56-81-5)	67,500.00, Oncorhynchus mykiss	10,000.00, Daphnia magna	Not Available
Cellulose, carboxymethyl ether, sodium salt - (9004-32-4)	Not Available	Not Available	Not Available
1-[1,3-bis(hydroxymethyl)-2,5-dioximidazolidin-4-yl]-1,3-	Not Available	Not Available	Not Available

# Safety Data Sheet

## SoloSite™ Conformable Wound Dressing

bis(hydroxymethyl)urea - (78491-02-8)

### 12.2. Persistence and degradability

There is no data available on the preparation itself.

### 12.3. Bioaccumulative potential

Not Measured

### 12.4. Mobility in soil

No data available.

### 12.5. Results of PBT and vPvB assessment

This product contains no PBT/vPvB chemicals.

### 12.6. Other adverse effects

No data available.

## 13. Disposal considerations

### 13.1. Waste treatment methods

Observe all federal, state and local regulations when disposing of this substance.

## 14. Transport information

	<b>DOT (Domestic Surface Transportation)</b>	<b>IMO / IMDG (Ocean Transportation)</b>	<b>ICAO/IATA</b>
<b>14.1. UN number</b>	Not Applicable	Not Regulated	Not Regulated
<b>14.2. UN proper shipping name</b>	Not Regulated	Not Regulated	Not Regulated
<b>14.3. Transport hazard class(es)</b>	<b>DOT Hazard Class:</b> Not Applicable	<b>IMDG:</b> Not Applicable <b>Sub Class:</b> Not Applicable	<b>Air Class:</b> Not Applicable
<b>14.4. Packing group</b>	Not Applicable	Not Applicable	Not Applicable
<b>14.5. Environmental hazards</b>			
<b>IMDG</b>	Marine Pollutant: No		
<b>14.6. Special precautions for user</b>	No further information		

## 15. Regulatory information

<b>Regulatory Overview</b>	The regulatory data in Section 15 is not intended to be all-inclusive, only selected regulations are represented.
<b>Toxic Substance Control Act (TSCA)</b>	All components of this material are either listed or exempt from listing on the TSCA Inventory.
<b>WHMIS Classification</b>	D2B
<b>US EPA Tier II Hazards</b>	<b>Fire:</b> No <b>Sudden Release of Pressure:</b> No

# Safety Data Sheet

## SoloSite™ Conformable Wound Dressing

Reactive: No

Immediate (Acute): Yes

Delayed (Chronic): No

### **EPCRA 311/312 Chemicals and RQs:**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

### **EPCRA 302 Extremely Hazardous:**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

### **EPCRA 313 Toxic Chemicals:**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

### **Proposition 65 - Carcinogens (>0.0%):**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

### **Proposition 65 - Developmental Toxins (>0.0%):**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

### **Proposition 65 - Female Repro Toxins (>0.0%):**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

### **Proposition 65 - Male Repro Toxins (>0.0%):**

To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.

### **New Jersey RTK Substances (>1%):**

Glycerol

### **Pennsylvania RTK Substances (>1%):**

Glycerol

## 16. Other information

The information and recommendations contained herein are based upon data believed to be correct. However, no guarantee or warranty of any kind, expressed or implied, is made with respect to the information contained herein. We accept no responsibility and disclaim all liability for any harmful effects which may be caused by exposure to our products. Customers/users of this product must comply with all applicable health and safety laws, regulations, and orders.

The full text of the phrases appearing in section 3 is:

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

**This is the first version in the GHS SDS format. Listings of changes from previous versions in other formats are not applicable.**

End of Document