

+ PICO 7 sNPWT Sizing Chart

Smith+Nephew

PICO[◇]
Single Use Negative Pressure
Wound Therapy System



Code
66802006

Dressing size
15cm x 20cm

Incision length
15cm

Open wound
10cm x 15cm

Code
66802007

Dressing size
15cm x 30cm

Incision length
25cm

Open wound
10cm x 25cm

Code
66802002

Dressing size
10cm x 20cm

Incision length
15cm

Open wound
5cm x 15cm

Code
66802003

Dressing size
10cm x 30cm

Incision length
25cm

Open wound
5cm x 25cm


Code
66802004

Dressing size
10cm x 40cm

Incision length
35cm


Open wound
5cm x 30cm

+ PICO 7 sNPWT Sizing Chart



Smith+Nephew

PICO[®]
Single Use Negative Pressure
Wound Therapy System



Smith & Nephew Pty Ltd (Australia) smith-nephew.com/en-au. Smith & Nephew Ltd (New Zealand) smith-nephew.com/en-nz. [®]Trademark of Smith+Nephew. All trademarks acknowledged. © February 2025 Smith+Nephew. 46937-anz V1 02/25. This material is intended for healthcare professionals. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.