



Versatile dressing for minor acute wounds, designed to protect the wound against contamination

Low allergy adhesive dressing to prevent skin reactions in patients.<sup>1</sup>

 **smith&nephew**  
**CUTIPLAST<sup>◇</sup>**  
**STERIL**

Fabric Island Dressing

Supporting healthcare professionals

# Versatile dressing for minor acute wounds, designed to protect the wound against contamination

CUTIPLAST<sup>◊</sup> Steril and CUTIPLAST Plus Steril Dressings are absorbent wound dressings<sup>2</sup> designed to have a low-adherent pad which is intended to minimise patient and wound disturbance on removal.

## Features and benefits



### Protects the wound

CUTIPLAST Steril and CUTIPLAST Plus Steril Dressings are composed of an adhesive layer, which is designed to act as a barrier to help protect the wound against contamination and an absorbent pad which cushions the wound whilst managing exudate.<sup>2-3</sup>



### Ease of use

CUTIPLAST Steril and CUTIPLAST Plus Steril Dressings are both easy to apply to the patients wound.<sup>4</sup>



### Conformable<sup>2</sup>

CUTIPLAST Steril and CUTIPLAST Plus Steril Dressings are extensible, designed to adapt to patient movement and conform to highly contoured parts of the body.<sup>2,4</sup>



### Permeable to moisture vapour

CUTIPLAST Steril and CUTIPLAST Plus Steril Dressings are permeable to moisture vapour to prevent moisture build up underneath the dressing.<sup>2</sup>



### Low allergy adhesive

CUTIPLAST Steril and CUTIPLAST Plus Steril Dressings consist of a low allergy adhesive to help minimise skin reactions in patients.<sup>1</sup>

## Indications

To dress post-operative wounds, cuts, lacerations, grazes, or sutured wounds.

## Precautions

As with all adhesive products, apply and remove carefully.

## Warnings

Frequent dressing changes, particularly on patients with fragile skin or dermatological conditions, can result in skin reddening or damage.

| S&N Code                         | Size        | Carton |
|----------------------------------|-------------|--------|
| <b>CUTIPLAST Steril Dressing</b> |             |        |
| 66001473                         | 10cm x 8cm  | 50     |
| 66001474                         | 15cm x 8cm  | 50     |
| 66001475                         | 20cm x 10cm | 50     |
| 66001476                         | 25cm x 10cm | 50     |
| 66001477                         | 30cm x 10cm | 50     |
| 66001478                         | 7.2cm x 5cm | 100    |

| S&N Code                              | Size        | Carton |
|---------------------------------------|-------------|--------|
| <b>CUTIPLAST Plus Steril Dressing</b> |             |        |
| 66076825                              | 7.2cm x 5cm | 5      |
| 66176826                              | 10cm x 8cm  | 5      |
| 66176827                              | 15cm x 8cm  | 5      |
| 66076832                              | 20cm x 10cm | 5      |
| 66076833                              | 25cm x 10cm | 3      |
| 66176825                              | 7.2cm x 5cm | 5      |
| 66076826                              | 10cm x 8cm  | 5      |

Smith & Nephew  
Croxley Park  
Building 5, Lakeside,  
Hatters Lane, Watford,  
Hertfordshire WD18 8YE

[www.smith-nephew.com](http://www.smith-nephew.com)

T +44 (0) 1923 477100  
F +44 (0) 1923 477101

◊Trademark of Smith & Nephew  
All Trademarks acknowledged  
©March 2019 Smith & Nephew  
16062

Supporting healthcare professionals for over 150 years