Go beyond the repair

Optimized sharpness^{1, 2}

allo

CAPERY

CAP-FIX

An essential part of Hip Repair Solutions Learn more at **smith-nephew.com**

Smith
Nephew

CAP-FIX^O Capsular Management Family

CAP-FIX^o Suture Passer

Smith+Nephew now offers a family of products for all your Capsular Management needs. From open to close, start your successful capsulotomy with our range of CAP-FIX Blades before completing the repair using our range of CAP-FIX Suture Passers.

Sharpen your technique and go beyond the repair with the latest addition to Smith+Nephew's evolving hip repair solutions.

Optimized...

...sharpness, strength, suture throw and trajectory for capsule closure $^{\scriptscriptstyle 1,\,2}$

Compatible...

...with #2 Vicryl^{™2}, #2 ULTRABRAID°, MINITAPE° and MAGNUMWIRE° Sutures. Vicryl suture has been shown to be less damaged by the CAP-FIX Suture Passer than the Stryker SlingShot^{™3}



Streamlined...

...approach reduces the number of steps required for capsule closure compared to suture shuttling¹





#2 Vicryl Suture

Ability...

...to reach clinically relevant zones for closure of interportal and T-Capsulotomy¹

Available in 45° and 70° angles

CAP-FIX^{\$} Blade

A single piece design*...

- Available in curved and straight blades
- Textured handle to aid in gripping
- Working length able to fit down all standard Smith+Nephew hip length cannulas**

*applies only to pencil grip blade configurations **see IFU for compatible cannula 4.5 cannula compatibility Pointed tip Opinted tip Straight blade

Available with handle for increased leverage.

No breakage...

Optimized...

...sharpness to cut through hip capsule tissue to perform Interportal and T-Capsulotomy $^{\!\!\!\!\!^4}$

Ordering information

CAP-FIX° Blades	
Reference	Description
72205316	CAP-FIX Blade, Straight
72205317	CAP-FIX Blade, Curved
72205318	CAP-FIX Blade, Straight, with Hip Pac
72205319	CAP-FIX Blade, Curved, with Hip Pac
72205314	CAP-FIX Blade, Straight, with Handle
72205315	CAP-FIX Blade, Curved, with Handle
CAP-FIX Suture Passer	
Reference	Description
72205683	CAP-FIX Suture Passer, 45°
72205684	CAP-FIX Suture Passer, 70°

Learn more at smith-nephew.com

Sports Medicine Andover, MA 01810

www.smith-nephew.com Smith & Nephew, Inc. T +978 749 1000 150 Minuteman Road US Customer Service: +1 800 343 5717

°Trademark of Smith+Nephew. All trademarks acknowledged. ©2021 Smith+Nephew. All rights reserved. Printed in USA. 31451 V1 09/21

References

1. Smith+Nephew 2021. CAP-FIX Suture Passer. Internal Report 15011184 A 2. Smith+Nephew 2021. CAP-FIX Suture Passer. Internal Report 15011185 A 3. Smith+Nephew 2021. CAP-FIX Suture Passer. Internal Report 15011186 A 4. Smith+Nephew 2020. CAP-FIX Blade Claims Testing. Internal Report. 15009423 B. 5. Smith+Nephew 2020. CAP-FIX Blade Claims Testing, Extreme Bend Test of Blade. Internal Report. 15010366 A.